

Sanderson, D.C.W., Carmichael, L.A., Fisk, S., Key, P., Scott, E.M., Thompson, M., Alberti, A., Baden, M., Barth, S., Baumann, N., Bergman, M., Brockmann, R., Brown, P., Brunner, J., Buranasilp, S., Butz, B., Chan Sheot Harn, J., Das, R., Dittmar, F., Gemperle, C., Jung Heon, J., Iglovári, M., Kay, C., Kruspe, W., Lardner, C., MacFarlane, S., McKay, K., Muller, B., Niemi, E., Nuchter-Frangos, M., Onori, S., O'Sullivan, M., Pfordt, J., Poulima, I., Pugh, M., Sardone, N., Schleich, C., Song, I-S., Straub, I., Todoriki, S., Ward, A., Wiezorek, C. and Woo Suk, K. (2007) *Development of Proficiency Testing for Detection of Irradiated Food: Project E01068. Results of Second Round PSL and TL Trials, September 2006*. Project Report. Scottish Universities Environmental Research Centre, East Kilbride, UK.

<http://eprints.gla.ac.uk/44637>

Deposited on: 15 December 2010

**DEVELOPMENT OF PROFICIENCY TESTING
FOR DETECTION OF IRRADIATED FOOD PROJECT E01068**

FINAL REPORT NOVEMBER 2007

**VOLUME II: RESULTS OF SECOND ROUND PSL AND TL TRIALS
SEPTEMBER 2006**

**D.C.W. Sanderson, L.A. Carmichael, S. Fisk, P. Key, E.M. Scott and
M. Thompson**

DEVELOPMENT OF PROFICIENCY TESTING FOR DETECTION OF IRRADIATED FOOD

PROJECT E01068

RESULTS OF SECOND ROUND PSL AND TL TRIALS SEPTEMBER 2006

D.C.W. Sanderson¹, L.A. Carmichael¹, S. Fisk¹, P. Key², E.M. Scott³ and M. Thompson⁴
Angelo Alberti⁵, Michael Baden⁶, Stephan Barth⁷, Neil Baumann⁸, Meike Bergman⁹, Rainer Brockmann¹⁰, Peter Brown¹¹, Jürgen Brunner¹², Songphop Buranasilp¹³, Brigitte Bütz¹⁴, Joanne Chan Sheot Harn¹⁵, Rupa Das¹⁶, Frank Dittmar¹⁷, Claudius Gemperle¹⁸, Jo Jung Heon¹⁹, Maria Iglovári²⁰, Christina Kay²¹, Wolfgang Kruspe²², Caroline Lardner²³, Stuart MacFarlane²⁴, Keith McKay²⁵, Beate Müller²⁶, Esko Niemi²⁷, Magda Nüchter-Frangos⁹, Sandro Onori²⁸, Margaret O'Sullivan²⁹, Jürgen Pfordt³⁰, Irene Poulima³¹, Mike Pugh³², Nicola Sardone³³, Christine Schleich³⁴, In-Sang Song¹⁵, Irene Straub³⁵, Setsuko Todoriki³⁶, Andy Ward³⁷, Claus Wiezorek³⁸, Kang Woo Suk³⁹.

¹ Scottish Universities Environmental Research Centre (SUERC), UK

² Institute for Food Research (IFR), UK

³ University of Glasgow, UK

⁴ Birkbeck College, UK

⁵ Consiglio Nazionale delle Ricerche Istituto per la Sintesi Organica e la Fotoreattività (ISOF-CNR), Italy

⁶ Eurofins Scientific Group, Germany

⁷ Federal Research Centre of Nutrition & Food (BfEL), Germany

⁸ McCormicks (UK) Ltd, UK

⁹ Berliner Betrieb für Zentrale Gesundheitliche Aufgaben, Institute für Lebensmittel, Arzneimittel und Tierseuchen, Germany

¹⁰ Chemisches und Veterinaruntersuchungsamt Ostwestfalen-Lippe CVUA OWL, Germany

¹¹ Lincolne Sutton & Wood Ltd, UK

¹² LUA Sachsen, Germany

¹³ Thai President Foods Public Company Limited, Thailand

¹⁴ Bayerisches Landesamt für Gesundheit und Lebensmittelsicherheit, Germany

¹⁵ Centre for Analytical Science, Health Sciences Authority, Singapore

¹⁶ BI Nutraceuticals, USA

¹⁷ Landesbetrieb Hessisches Landeslabor, Germany

¹⁸ Kantonales Labor, Aarau, Switzerland

¹⁹ Nong Shim Co.Ltd, Korea

²⁰ Laboratory of Foodstuffs and Chemicals, Hungary

²¹ Manhattan Drug Company Inc, USA

²² Thüringer Landesamt für Lebensmittelsicherheit und Verbraucherschutz (TLLV), Germany

²³ Public Analyst Laboratory (Galway), Eire

²⁴ East Anglian Food Ingredients Ltd, UK

²⁵ Glasgow Scientific Services, UK

²⁶ Landeslabor Brandenburg, Standort Frankfurt (Oder), Germany

²⁷ Finnish Customs Laboratory, Finland

²⁸ Istituto Superiore di Sanità, Italy

²⁹ Public Analyst Laboratory (Cork), Eire

³⁰ LAVES - Lebensmittelinstitut Oldenburg

³¹ Food Division Laboratory, General Chemical State Laboratory of Greece, Greece

³² British Pepper and Spice Co, Ltd, UK

³³ Indena SpA, Italy

³⁴ Landesuntersuchungsamt, Institut für Lebensmittelchemie und Arzneimittelprüfung Mainz, Germany

³⁵ Tae Kyung Nong San Co Ltd, Korea

³⁶ Chemisches- und Veterinaruntersuchungsamt, Germany

³⁷ National Food Research Institute, Japan

³⁸ Unilever UK Foods, UK

³⁹ Chem. Landes Staatl. Vet. U. Amt., Germany

CONTENTS

| | |
|---|-------------|
| List of tables | iii |
| List of figures | vii |
| Summary | viii |
| Acknowledgments | ix |
| | |
| 1. Introduction | 1 |
| 2. Implementation of second round (PSL & TL) | 3 |
| 2.1 Outline of round 2 | 3 |
| 2.2 Participants | 3 |
| 2.3 Test material preparation and handling | 5 |
| 2.3.1 PSL Materials | 5 |
| 2.3.2 TL Materials | 6 |
| 2.3.3 Blending | 7 |
| 2.3.3.1 Mixing process | 7 |
| 3. PSL screening | 9 |
| 3.1 PSL Homogeneity testing (blends) | 9 |
| 3.2 PSL Homogeneity results (blends and paprika) | 9 |
| 3.3 Participants' results | 11 |
| 3.3.1 The paprika standard | 11 |
| 3.3.2 Screening results | 14 |
| 3.3.3 Summary statistics | 18 |
| 3.3.4 Qualitative results | 37 |
| 4. Calibrated PSL | 43 |
| 4.1 Calibrated PSL Homogeneity Testing | 43 |
| 4.2 Homogeneity Testing Results | 43 |
| 4.3 Participants' Results | 44 |
| 5. TL analysis | 50 |
| 5.1 TL Homogeneity Testing | 50 |
| 5.2 TL Homogeneity Testing Results | 50 |
| 5.3 Participants' Results | 59 |
| 5.3.1 Quantitative results | 59 |
| 5.3.2 Qualitative results | 75 |
| 6. Discussion and Conclusions | 77 |

APPENDICES

| | |
|---|------------|
| Appendix A: Study protocol | A.1 |
| Appendix B: Homogeneity testing raw data | B.1 |
| Appendix C: Participants' raw data | C.1 |
| Appendix D: Sample allocation | D.1 |
| Appendix E: Mixing the blends | E.1 |

LIST OF TABLES

| | | |
|-------------|--|-----|
| Table 2.1. | List of Round 2 Participants and their institutions | 4 |
| Table 2.2. | List of Round 2 Test Materials | 5 |
| Table 2.3. | PSL screening of 10 aliquots of each product prior to blending | 7 |
| Table 2.4. | Blend concentrations and predicted mean PSL estimate | 7 |
| Table 2.5. | Sample number allocations for blends | 8 |
| Table 2.6. | PSL screening results during paprika homogenisation | 8 |
| Table 3.1. | Reference values for the blended materials (screening) | 9 |
| Table 3.2. | Results of 10 replicate measurements from the paprika reference material | 11 |
| Table 3.3. | Participants' paprika standards – mean, standard deviation and standard error from 10 observations | 11 |
| Table 3.4 | Relativity instrumental sensitivities derived from measurements of the paprika standard (10 each for rounds 1 and 2) | 12 |
| Table 3.5 | Unirradiated Round 2 Test Materials: z-Scores for Participants 1-13 | 19 |
| Table 3.6 | Unirradiated Round 2 Test Materials: z-Scores for Participants 14-25 | 20 |
| Table 3.7 | Unirradiated Round 2 Test Materials: z-Scores for Participants 26-35 | 21 |
| Table 3.8 | Irradiated Round 2 Test Materials: z-Scores for Participants 1-13 | 22 |
| Table 3.9 | Irradiated Round 2 Test Materials: z-Scores for Participants 14-24 | 23 |
| Table 3.10 | Irradiated Round 2 Test Materials: z-Scores for Participants 25-35 | 24 |
| Table 3.11: | Blended Round 2 Test Materials: z-Scores for Participants 1-13 | 25 |
| Table 3.12: | Blended Round 2 Test Materials: z-Scores for Participants 14-24 | 25 |
| Table 3.13: | Blended Round 2 Test Materials: z-Scores for Participants 25-35 | 25 |
| Table 3.14 | Participants' qualitative results by laboratory | 39 |
| Table 3.15 | Participants' qualitative results by product | 40 |
| Table 3.16 | Homogeneity testing for blends (screening) with coefficient of variance | 41 |
| Table 3.17 | Comparison of participants' qualitative percentage with the reference set | 41 |
| Table 3.18 | PSL screening results from the blended samples by concentration and brightness | 42 |
| Table 4.1 | Type of source used for calibration by participants and whether they utilised the second data set to evaluate | 44 |
| Table 5.1 | LiF peak V temperatures and MDL for all laboratories | 59 |
| Table 5.2 | z-score for TL analyses – unirradiated samples | 64 |
| Table 5.3 | z-score for TL analyses – irradiated samples | 64 |
| Table 5.4 | z-score for TL analyses – blended samples | 65 |
| Table 5.5 | Qualitative results from TL participants | 75 |
| Table B.1 | Reference values based on log statistics for irradiated and unirradiated test materials – calibrated PSL | B.2 |
| Table B.2 | Reference data for blended samples | B.3 |
| Table B.3 | Raw data for unirradiated test materials – calibrated PSL | B.4 |
| Table B.4 | Raw data for irradiated test materials – calibrated PSL | B.5 |
| Table B.5 | Raw data for homogeneity testing of blended samples – screening | B.6 |
| Table B.6 | Raw data for homogeneity testing of blended samples – calibrated | B.6 |
| Table B.7 | Reference values for TL Homogeneity Testing based on Glow Ratio | B.6 |
| Table C.1 | Raw Data Returned By Laboratory 1 – screening of unirradiated samples | C.1 |
| Table C.2 | Raw Data Returned By Laboratory 1 – screening of irradiated samples and blends | C.2 |
| Table C.3 | Raw Data Returned By Laboratory 1 – calibrated PSL for unirradiated samples | C.3 |
| Table C.4 | Raw Data Returned By Laboratory 1 – calibrated PSL for irradiated samples and blends | C.4 |
| Table C.5 | Raw Data Returned By Laboratory 2 – screening of unirradiated samples | C.5 |

| | | |
|------------|---|------|
| Table C.6 | Raw Data Returned By Laboratory 2 – screening of irradiated samples and blends | C.6 |
| Table C.7 | Raw Data Returned By Laboratory 3 – screening of unirradiated samples | C.7 |
| Table C.8 | Raw Data Returned By Laboratory 3 – screening of irradiated samples and blends | C.8 |
| Table C.9 | Raw Data Returned By Laboratory 5 – screening of unirradiated samples | C.9 |
| Table C.10 | Raw Data Returned By Laboratory 5 – screening of irradiated samples and blends | C.10 |
| Table C.11 | Raw Data Returned By Laboratory 6 – screening of unirradiated samples | C.11 |
| Table C.12 | Raw Data Returned By Laboratory 6 – screening of irradiated samples and blends | C.12 |
| Table C.13 | Raw Data Returned By Laboratory 8 – screening of unirradiated samples | C.13 |
| Table C.14 | Raw Data Returned By Laboratory 8 – screening of irradiated samples and blends | C.14 |
| Table C.15 | Raw Data Returned By Laboratory 9 – screening of unirradiated samples | C.15 |
| Table C.16 | Raw Data Returned By Laboratory 9 – screening of irradiated samples and blends | C.16 |
| Table C.17 | Raw Data Returned By Laboratory 11 – screening of unirradiated samples | C.17 |
| Table C.18 | Raw Data Returned By Laboratory 11 – screening of irradiated samples and blends | C.18 |
| Table C.19 | Raw Data Returned By Laboratory 12 – screening of unirradiated samples | C.19 |
| Table C.20 | Raw Data Returned By Laboratory 12 – screening of irradiated samples and blends | C.20 |
| Table C.21 | Raw Data Returned By Laboratory 14 – screening of irradiated samples and blends | C.21 |
| Table C.22 | Raw Data Returned By Laboratory 14 – screening of irradiated samples and blends | C.22 |
| Table C.23 | Raw Data Returned By Laboratory 14 – calibrated PSL for unirradiated samples | C.23 |
| Table C.24 | Raw Data Returned By Laboratory 14 – calibrated PSL for irradiated samples and blends | C.24 |
| Table C.25 | Raw Data Returned By Laboratory 16 – screening of unirradiated samples | C.25 |
| Table C.26 | Raw Data Returned By Laboratory 16 – screening of irradiated samples and blends | C.26 |
| Table C.27 | Raw Data Returned By Laboratory 17 – screening of unirradiated samples | C.27 |
| Table C.28 | Raw Data Returned By Laboratory 17 – screening of irradiated samples and blends | C.28 |
| Table C.29 | Raw Data Returned By Laboratory 17 – calibrated PSL for unirradiated samples | C.29 |
| Table C.30 | Raw Data Returned By Laboratory 17 – calibrated PSL for irradiated samples and blends | C.30 |
| Table C.31 | Raw Data Returned By Laboratory 18 – screening of unirradiated samples | C.31 |
| Table C.32 | Raw Data Returned By Laboratory 18 – screening of irradiated samples and blends | C.32 |
| Table C.33 | Raw Data Returned By Laboratory 18 – calibrated PSL for unirradiated Samples | C.33 |
| Table C.34 | Raw Data Returned By Laboratory 18 – calibrated PSL for irradiated samples and blends | C.34 |
| Table C.35 | Raw Data Returned By Laboratory 19 – screening of unirradiated samples | C.35 |
| Table C.36 | Raw Data Returned By Laboratory 19 – screening of irradiated samples and blends | C.36 |
| Table C.37 | Raw Data Returned By Laboratory 19 – calibrated PSL for unirradiated samples | C.37 |

| | | |
|------------|---|------|
| Table C.38 | Raw Data Returned By Laboratory 19 – calibrated PSL for irradiated samples and blends | C.38 |
| Table C.39 | Raw Data Returned By Laboratory 20 – screening of unirradiated samples | C.39 |
| Table C.40 | Raw Data Returned By Laboratory 20 – screening of irradiated samples and blends | C.40 |
| Table C.41 | Raw Data Returned By Laboratory 21 – screening of unirradiated samples | C.41 |
| Table C.42 | Raw Data Returned By Laboratory 21 – screening of irradiated samples and blends | C.42 |
| Table C.43 | Raw Data Returned By Laboratory 22 – screening of unirradiated samples | C.43 |
| Table C.44 | Raw Data Returned By Laboratory 22 – screening of for irradiated samples and blends | C.44 |
| Table C.45 | Raw Data Returned By Laboratory 23 – screening of unirradiated samples | C.45 |
| Table C.46 | Raw Data Returned By Laboratory 23 – screening of irradiated samples and blends | C.46 |
| Table C.47 | Raw Data Returned By Laboratory 23 – calibrated PSL for unirradiated Samples | C.47 |
| Table C.48 | Raw Data Returned By Laboratory 23 – calibrated PSL for irradiated Samples and blends | C.48 |
| Table C.49 | Raw Data Returned By Laboratory 24 – screening of unirradiated samples | C.49 |
| Table C.50 | Raw Data Returned By Laboratory 24 – screening of irradiated samples and blends | C.50 |
| Table C.51 | Raw Data Returned By Laboratory 24 – calibrated PSL for unirradiated Samples | C.51 |
| Table C.52 | Raw Data Returned By Laboratory 24 – calibrated PSL for irradiated Samples and blends | C.52 |
| Table C.53 | Raw Data Returned By Laboratory 25 – screening of unirradiated samples | C.53 |
| Table C.54 | Raw Data Returned By Laboratory 25 – screening of irradiated samples and blends | C.54 |
| Table C.55 | Raw Data Returned By Laboratory 25 – calibrated PSL for unirradiated Samples | C.55 |
| Table C.56 | Raw Data Returned By Laboratory 25 – calibrated PSL for irradiated samples and blends | C.56 |
| Table C.57 | Raw Data Returned By Laboratory 26 – screening of unirradiated samples | C.57 |
| Table C.58 | Raw Data Returned By Laboratory 26 – screening of irradiated samples and blends | C.58 |
| Table C.59 | Raw Data Returned By Laboratory 26 – calibrated PSL for unirradiated samples | C.59 |
| Table C.60 | Raw Data Returned By Laboratory 26 – calibrated PSL for irradiated samples and blends | C.60 |
| Table C.61 | Raw Data Returned By Laboratory 27 – screening of unirradiated samples | C.61 |
| Table C.62 | Raw Data Returned By Laboratory 27 – screening of irradiated samples and blends | C.62 |
| Table C.63 | Raw Data Returned By Laboratory 27 – calibrated PSL for unirradiated Samples | C.63 |
| Table C.64 | Raw Data Returned By Laboratory 27 – calibrated PSL for irradiated samples and blends | C.64 |
| Table C.65 | Raw Data Returned By Laboratory 28 – screening of unirradiated samples | C.65 |
| Table C.66 | Raw Data Returned By Laboratory 28 – screening of irradiated samples and blends | C.66 |
| Table C.67 | Raw Data Returned By Laboratory 28 – calibrated PSL for unirradiated samples | C.67 |
| Table C.68 | Raw Data Returned By Laboratory 28 – calibrated PSL for irradiated Samples and blends | C.68 |

| | | |
|------------|---|------|
| Table C.69 | Raw Data Returned By Laboratory 29 – screening of unirradiated samples | C.69 |
| Table C.70 | Raw Data Returned By Laboratory 29 – screening of irradiated samples and blends | C.70 |
| Table C.71 | Raw Data Returned By Laboratory 30 – screening of unirradiated samples | C.71 |
| Table C.72 | Raw Data Returned By Laboratory 30 – screening of irradiated samples and blends | C.72 |
| Table C.73 | Raw Data Returned By Laboratory 31 – screening of unirradiated samples | C.73 |
| Table C.74 | Raw Data Returned By Laboratory 31 – screening of irradiated samples and blends | C.74 |
| Table C.75 | Raw Data Returned By Laboratory 33 – screening of unirradiated samples | C.75 |
| Table C.76 | Raw Data Returned By Laboratory 33 – screening of irradiated samples and blends | C.76 |
| Table C.77 | Raw Data Returned By Laboratory 34 – screening of unirradiated samples | C.77 |
| Table C.78 | Raw Data Returned By Laboratory 34 – screening of irradiated samples and blends | C.78 |
| Table C.79 | Raw Data Returned By Laboratory 35 – screening of unirradiated samples | C.79 |
| Table C.80 | Raw Data Returned By Laboratory 35 – screening of irradiated samples and blends | C.80 |
| Table C.81 | Raw Data Returned By Laboratory 35 – calibrated PSL for unirradiated samples | C.81 |
| Table C.82 | Raw Data Returned By Laboratory 35 – calibrated PSL for irradiated samples and blends | C.82 |
| Table C.83 | TL data for TL laboratory 1 | C.83 |
| Table C.84 | TL data for TL laboratory 2 | C.84 |
| Table C.85 | TL data for TL laboratory 3 | C.85 |
| Table C.86 | TL data for TL laboratory 4 | C.86 |
| Table C.87 | TL data for TL laboratory 5 | C.87 |
| Table C.88 | TL data for TL laboratory 6 | C.88 |
| Table C.89 | TL data for TL laboratory 7 | C.89 |
| Table C.90 | TL data for TL laboratory 8 | C.90 |
| Table C.91 | TL data for TL laboratory 10 | C.91 |
| Table C.92 | TL data for TL laboratory 11 | C.92 |
| Table C.93 | TL data for TL laboratory 13 | C.93 |
| Table C.94 | TL data for TL laboratory 14 | C.94 |
| Table C.95 | TL data for TL laboratory 15 | C.95 |
| Table C.96 | TL data for TL laboratory 18 | C.96 |
| Table C.97 | TL data for TL laboratory 24 | C.97 |
| Table D1. | Allocation of pots to laboratories for round 2 | D.1 |
| Table D2. | Correspondence between SP number and sample number as presented to Participants | D.2 |
| Table E.1 | PSL measurements during mixing for SP8578 Green Tea | E.1 |
| Table E.2 | PSL measurements during mixing for SP8535 Chilli Powder | E.1 |
| Table E.3 | PSL measurements during mixing for SP8531 Ginger 1% | E.2 |
| Table E.4 | PSL measurements during mixing for SP8531 Ginger 10% | E.3 |
| Table E.5 | PSL measurements during mixing for SP8577 Ginseng | E.3 |
| Table E.6 | PSL measurements during mixing for SP8533 Rosemary | E.4 |
| Table E.7 | PSL measurements during mixing for SP8539 Tarragon | E.4 |

LIST OF FIGURES

| | |
|---|----|
| Figure 3.1 Histogram of reference data to show position of blends | 10 |
| Figure 3.2 Mean paprika signal in photon counts obtained by laboratories in rounds 1 and 2 | 13 |
| Figure 3.3 Relative instrument sensitivities from measurements of paprika standards in rounds 1 and 2 standardised to SUERC results | 13 |
| Figure 3.4 Participants' screening data for each product | 14 |
| Figure 3.5 Participants' data for unirradiated samples by laboratory | 15 |
| Figure 3.6 Participants' data for irradiated samples by laboratory | 16 |
| Figure 3.7 Participants' data for blended samples by laboratory and concentration | 17 |
| Figure 3.8 Histogram of participants' screening results for all products | 18 |
| Figure 3.9 Participants' z-scores for unirradiated samples by sample | 27 |
| Figure 3.10 Participants' z-scores for unirradiated samples by laboratory | 28 |
| Figure 3.11 Participants' z-scores for irradiated samples by sample | 29 |
| Figure 3.12 Participants' z-scores for irradiated samples by laboratory | 30 |
| Figure 3.13 Time elapsed since irradiation vs z-score for all irradiated samples – participants' data | 31 |
| Figure 3.14 Participants' z-scores for blended samples by sample | 32 |
| Figure 3.15 Participants' z-scores for blended samples by laboratory | 33 |
| Figure 3.16 Contour plot for participants' data for unirradiated samples | 34 |
| Figure 3.17 Contour plot for participants' data for irradiated samples | 35 |
| Figure 3.18 Contour plot for participants' data for blended samples | 36 |
| Figure 4.1 Initial vs calibrated PSL for reference materials | 43 |
| Figure 4.2 Initial vs calibrated PSL for laboratories 1, 14 and 17 | 46 |
| Figure 4.3 Initial vs calibrated PSL for laboratories 18, 19 and 23 | 47 |
| Figure 4.4 Initial vs calibrated PSL for laboratories 24, 25 and 26 | 48 |
| Figure 4.5 Initial vs calibrated PSL for laboratories 27, 28 and 35 | 49 |
| Figure 5.1 Distribution of glow 1 peak shapes for 20 aliquots each of 3 blends | 52 |
| Figure 5.2 TL homogeneity testing glow 1 vs glow 2 scatter plot for SP8531 and SP8533 | 53 |
| Figure 5.3 TL homogeneity testing glow 1 vs glow 2 scatter plot for SP8535 and SP8539 | 54 |
| Figure 5.4 TL homogeneity testing glow 1 vs glow 2 scatter plot for SP8577 and SP8578 | 55 |
| Figure 5.5 TL homogeneity testing glow ratio plot for SP8531 and SP8533 | 56 |
| Figure 5.6 TL homogeneity testing glow ratio plot for SP8535 and SP8539 | 57 |
| Figure 5.7 TL homogeneity testing glow ratio plot for SP8577 and SP8578 | 58 |
| Figure 5.8 Participants' glow ratios for SP8531 and SP8533 | 61 |
| Figure 5.9 Participants' glow ratios for SP8535 and SP8539 | 62 |
| Figure 5.10 Participants' glow ratios for SP8577 and SP8578 | 63 |
| Figure 5.11 Participants' TL z-scores for unirradiated samples by laboratory | 66 |
| Figure 5.12 Participants' TL z-scores for irradiated samples by laboratory | 67 |
| Figure 5.13 Participants' TL z-scores for blended samples by laboratory | 68 |
| Figure 5.14 Participants' TL z-scores for unirradiated samples by sample | 69 |
| Figure 5.15 Participants' TL z-scores for irradiated samples by sample | 70 |
| Figure 5.16 Participants' TL z-scores for blended samples by sample | 71 |
| Figure 5.17 Contour plot for participants' TL data for unirradiated samples | 72 |
| Figure 5.18 Contour plot for participants' TL data for irradiated samples | 73 |
| Figure 5.19 Contour plot for participants' TL data for blended samples | 74 |

SUMMARY

This report summarises the results from Round 2 in which participants used the BS EN 13751:2002 photostimulated luminescence (PSL) method and the BS EN1788:2001 thermoluminescence (TL) method. Recognising the fundamentally qualitative nature of irradiation testing methods and the need to examine the relationship between underlying quantitative data and their associated qualitative outcomes led to a trial which attempts to utilise large numbers of test materials and large numbers of participants.

In this round, 33 test materials of herbs, spices, seasonings and dietary supplement ingredients were prepared in both irradiated and unirradiated forms and in addition, 6 products were selected for blending. The total number of samples distributed to each participant was 72. The selection aimed to include a range of sensitivities (from low to high), which together with 3 different blend concentrations would test proficiency at the limits of the method.

30 participating laboratories conducted PSL determinations and measurements from standard materials and returned results to the SUERC; 12 also carried out calibrated PSL. A further 15 laboratories participated in TL separations and measurements from a subset of 18 samples comprising 6 products presented in unirradiated, irradiated and blended forms.

Reference data for the PSL screening part of the study comprised values determined during round 1, supplemented with an additional 60 PSL determinations from the 6 blended products. Additional PSL screening from each product was also conducted during the blending process. For calibrated PSL analyses all 98 of the original samples used in round 1, plus the 6 blends were re-irradiated, with 10 fold replication, and then measured to determine the calibrated PSL response. Thus an additional 1040 PSL reference analyses were conducted in support of round 2. For the 18 TL samples 20 separations per sample were performed, followed by TL analysis, comprising 360 single aliquot TL determinations, equivalent to 180 EN1788 duplicated analyses. The reference analyses are documented in this report and used, where appropriate, to define assigned values for z score determination.

Participants returned results for PSL screening, calibrated PSL and TL analyses in a timely manner. The PSL screening results once again produced high quality results. Overall 99.9% of irradiated samples were successfully screened into intermediate and positive bands. Only 2.3% of the unirradiated samples screened into positive bands, representing a slight improvement in comparison with round 1. The blended mixtures produced intermediate or positive results in 95.8% of the cases containing 10% irradiated material, in 80.8% of samples with 1% irradiated material and in 37.5% of cases with 0.1% concentrations. Quantitative analyses based on z scores again were able to reveal both sample and laboratory specific differences in behaviour. In comparison with round 1 there has been marked improvement in sample handling with reduced levels of evidence for cross contamination of irradiated and unirradiated materials in the testing laboratories. There are still a few laboratories whose performance could be further improved. The identification and remediation of such problems is a significant achievement for PT work of this sort, and underlines its importance.

Calibrated PSL data were returned by 12 laboratories, each of whom succeeded in producing data sets that conformed to EN13751, despite having diverse access to irradiation facilities. There was some evidence, based on these data that x-ray irradiation and commercial gamma irradiation produced a greater spread in participants' data than the use of a LINAC for administering the calibration doses. Further work would be needed to establish whether such observations can be generalised, and to assess the sample-handling components of data scatter. For these samples

qualitative outcomes were not critically dependent on the calibrated response, however these data could now be used to design a more sensitive test of the calibrated system.

For the TL study 15 participants returned data in 241 out of 270 possible analyses, the remainder of determinations failing to produce sufficient minerals to satisfy EN1788. Of these some 92.6% of the unirradiated samples, 98.9% of the irradiated samples and 62% of the reported results from blends were correctly classified. A small proportion of unirradiated samples were identified as irradiated or containing irradiated material, from four laboratories. One irradiated sample was not successfully identified in one laboratory. For the blends the highest concentrations (10%) were all correctly identified in participants' results; while at lower concentrations laboratories detection rates for the irradiated component were again reduced. Quantitative analyses were also conducted using z scores on the recorded glow ratios. This appears to be informative and to identify laboratory to laboratory differences, although further work is needed to examine the distribution of the glow ratio variable and to explore other means of distinguishing between the detection conditions of each laboratory, for example by utilising absolute signal levels in comparison with detection limits.

ACKNOWLEDGMENTS

Support for this work by the Food Standard Agency, through project E01068, is gratefully acknowledged. We are also grateful to the following participants for the work and time they have contributed to the second round of this study:

Angelo Alberti, Michael Baden, Stephan Barth, Neil Baumann, Meike Bergman, Rainer Brockmann, Peter Brown, Jorgen Brunner, Songphop Buranasilp, Brigitte Butz, Joanne Chan Sheot Harn, Rupa Das, Frank Dittmar, Claudius Gemperle, Jo Jung Heon, Maria Iglovári, Christina Kay, Wolfgang Kruspe, Caroline Lardner, Stuart MacFarlane, Keith McKay, Beate Muller, Esko Niemi, Sandro Onori, Margaret O'Sullivan, Juergen Pfordt, Irene Poulima, Mike Pugh, Nicola Sardone, Christine Scleich, In-Sang Song, Irene Straub, Setsuko Todoriki, Andy Ward, Claus Wiezorek, Kang Woo Suk.

1. INTRODUCTION

Irradiation is used in many countries for the purposes of shelf life extension, reduction of spoilage and pathogen content, and retardation of ripening and sprouting processes in many different foods. UK^{1,2,3,4} and European Regulations^{5,6} require both licensing of plant and process and explicit product labelling at all stages of market presentation. However, in the absence of widespread consumer acceptance, there is little evidence of properly labelled products in the UK or in Europe. Several analytical methods for detection of food irradiation have been developed, of which the CEN international standards^{7,8,9,10,11,12,13,14,15,16} (i.e. within the UK the BS EN series) based on luminescence are in quite widespread use. Both photostimulated luminescence (PSL) screening and thermoluminescence (TL) analysis were used successfully in the 1996 MAFF survey of undeclared foods¹⁷ and also in the 2001 survey conducted by the Food Standards Agency¹⁸. These surveys have been successful in identifying undeclared irradiated spices and shellfish and particularly in drawing attention to the significant problems associated with dietary supplements¹⁹, which are the subject of current enforcement actions in the UK, and elsewhere in Europe. With this in mind the Food Standards Agency commissioned this project to assess the feasibility of developing a proficiency scheme appropriate for the detection of irradiated foods.

International Harmonised Protocols for conducting proficiency testing of analytical methods and laboratories are available^{20,21}, and form the basis for schemes such as those operated by

¹ Food (Control of Irradiation) Regulations, 1990, SI 2490

² The Food Irradiation Provisions (England) Regulations 2000, 2000, SI 2254

³ The Food Irradiation Provisions (Wales) Regulations 2001, 2001, WSI 1232 (W.66)

⁴ The Food Irradiation Provisions (Scotland) Regulations 2000, 2000, SSI 309

⁵ European Directive 1999/2/EC, On approximation of the laws of the Member States concerning foods and food ingredients treated with Ionising Radiation, OJEC, February 1999

⁶ European Directive 1999/3/EC, On the establishment of a community list of foods and food ingredients treated with ionising radiation, OJEC, February 1999

⁷ BS EN 1784:1996 Foodstuffs - Detection of Irradiated Food Containing fat – Gas Chromatographic Analysis of Hydrocarbons

⁸ BS EN 1785:1996 Foodstuffs - Detection of Irradiated Food Containing fat – Gas Chromatographic/Mass Spectrometric Analysis of Alkylcyclobutanones

⁹ BS EN 1786:1996 Foodstuffs - Detection of Irradiated Food Containing Bone – Method by ESR Spectroscopy

¹⁰ BS EN 1787:2000 Foodstuffs - Detection of Irradiated Food Containing Cellulose by ESR Spectroscopy

¹¹ BS EN 13708:2001 Foodstuffs - Detection of Irradiated Food Containing Crystalline Sugar by ESR Spectroscopy

¹² BS EN 13783:2001 Foodstuffs - Detection of Irradiated Using Direct Epifluorescent Filter technique/Aerobic Plate count (DEFT/APC) – Screening Method

¹³ BS EN 13784:2001 Foodstuffs - DNA Comet Assay for the Detection of Irradiated Foodstuffs – Screening method

¹⁴ BS EN 1788:1997 Foodstuffs - Detection of Irradiated Food From Which Silicate Minerals can be Isolated: Method by Thermoluminescence

¹⁵ BS EN 13751:2000 Foodstuffs - Detection of Irradiated Food Using Photostimulated Luminescence

¹⁶ BS EN 14569:2004 Foodstuffs – Microbiological Screening for Irradiated Food Using LAL/GNB Procedures.

¹⁷ MAFF, 1997, Undeclared Irradiation of Foodstuffs Surveillance Exercise, Food Surveillance Information Sheet, 102

¹⁸ Food Standards Agency 2002, Survey for Irradiated Foods – Herbs and Spices, Dietary Supplements and Prawns and Shrimps, Food Survey Information Sheet 25/02

¹⁹ European Commission, 2002, Report from the Commission on Food Irradiation for the period September 2000 to December 2001, OJEC, 255,2-12

²⁰ Thompson M., and Wood R, 1993, The International Harmonised Protocol for the Proficiency Testing of (Chemical) Analytical Laboratories, IUPAC/ISO/AOAC Protocol for Proficiency Testing, IUPAC, Geneva

FAPAS® and others for quantitative analysis of foods. However as noted in the earlier report on the first PT trial round ²² there are methodological questions concerning the adaptation of quantitative PT approaches to the qualitative criteria of irradiation tests.

In the first round²², 32 PSL screening laboratories examined 49 samples presented as irradiated and unirradiated pairs. Reference analyses were conducted at SUERC and used to define standard values for evaluation of z scores from participants' data. These were compared with qualitative classifications and it was shown that PSL screening outcomes were closely related to z score performance. Some performance differences between laboratories were noted, and suggestions made for enhancing sample handling integrity in the second round. In discussion with participants it was also agreed to use replicate samples in future rounds to bring the PT protocol into line with routine practice in the laboratories.

In this report the second PT trial round is discussed and results are presented. In addition to continuing to study PSL screening, the second round introduces calibrated PSL measurements and TL analysis. PSL screening was conducted by 29 laboratories, measuring 72 samples comprising 33 irradiated and unirradiated pairs, plus 6 blended mixtures of irradiated and unirradiated products. Calibrated PSL results were returned from the same sample sets by 11 participants. The TL analysis tasks involved measurement of 18 samples comprising irradiated, unirradiated and blended mixtures prepared from 6 products. TL analysis was undertaken by 16 laboratories.

This report outlines the design and preparation of round 2 together with details of the new reference analyses undertaken, mainly for the calibrated PSL and TL stages, and an analysis of participants' results.

²¹Thompson,M., et al, 2006, The International Harmonized Protocol For The Proficiency Testing Of Analytical Chemistry Laboratories, IUPAC Technical Report, in press

²² D.C.W. Sanderson, L.A. Carmichael, S. Fisk, P. Key, E.M. Scott And M. Thompson, 2005, Development of Proficiency Testing For Detection Of Irradiated Food. Project E01068: Results Of First Round PSL Trials, SUERC, East Kilbride

2. IMPLEMENTATION OF SECOND ROUND (PSL & TL)

In this section the work undertaken in preparation for round 2 is described. The design and protocols are presented, as are details of the sample preparation and presentation to participants. The samples were selected from materials used in round 1, thus permitting the retention of screening reference data from the earlier round. Blended mixtures were prepared and characterised. Additional reference analyses were undertaken and will be documented in the individual sections dealing with PSL screening, Calibrated PSL and TL.

2.1 Outline of round 2

Round 2 aimed to build on round 1 by repeating the PSL screening measurements on a reduced suite of materials now specifically including blends, with the introduction of duplicate aliquots, and by extending the study to calibrated PSL for those laboratories wishing to undertake such measurements. TL analysis was also introduced with a reduced set of samples, and with the participation of several new laboratories in addition to some round 1 participants who were able to conduct both methods.

2.2 Participants

All PSL round 1 participants indicated that they were willing to proceed to round 2 with a sub-set agreeing to conduct calibrated PSL. Subsequently one laboratory experienced technical problems and had to withdraw. Round 2 also included TL analysis; approximately 20 laboratories were approached and 16 agreed to participate of which 7 were already PSL participants. Two new PSL laboratories and one PSL / TL participant also joined for round 2. The full list of participants is given in Table 2.1 below.

Table 2.1. List of Round 2 Participants and their institutions

| Name | Address |
|------------------------|---|
| Dr Angelo Alberti | Consiglio Nazionale delle Ricerche Istituto per la Sintesi Organica e la Fotoreattività (ISOFCNR), Bologna, Italy |
| Dr Stephan Barth | Federal Research Centre of Nutrition & Food (BfEL), Karlsruhe, Germany |
| Neil Baumann | McCormicks (UK) Ltd, UK |
| Dr Rainer Brockmann | Chemisches und Veterinaruntersuchungsamt Ostwestfalen-Lippe CVUA OWL, Bielefeld, Germany |
| Peter Brown | Lincolne Sutton & Wood Ltd, UK |
| Songphop Buranasilp | Thai President Foods Public Company Limited, Thailand |
| Joanne Chan Sheot Harn | Centre for Analytical Science, Health Sciences Authority, Singapore |
| Rupa Das | BI Nutraceuticals, USA |
| Frank Dittmar | Landesbetrieb Hessisches Landeslabor, Kassel, Germany |
| Dr Claudius Gemperle | Kantonales Labor, Aarau, Switzerland |
| Jo Jung Heon | Nong Shim Co.Ltd, Korea |
| Christina Kay | Manhattan Drug Company Inc, USA |
| Dr Caroline Lardner | Public Analyst Laboratory (Galway), Eire |
| Stuart MacFarlane | East Anglian Food Ingredients Ltd, UK |
| Keith McKay | Glasgow Scientific Services, UK |
| Esko Niemi | Finnish Customs Laboratory, Finland |
| Dr Sandro Onori | Istituto Superiore di Sanita, Rome, Italy |
| Dr Juergen Pfordt | LAVES - Lebensmittelinstitut Oldenburg |
| Mike Pugh | British Pepper and Spice Co, Ltd, UK |
| Dr Nicola Sardone | Indena SpA, Settala, Italy |
| In-Sang Song | TaeKyung NongSan Co Ltd, Korea |
| Irene Straub | Chemisches- und Veterinaruntersuchungsamt, Karlsruhe, Germany |
| Setsuko Todoriki | National Food Research Institute, Japan |
| Dr Claus Wiezorek | Chem. Landes Staatl. Vet. U. Amt., Münster, Germany |
| Kang Woo Suk | Nong Shim Company, Korea |
| Meike Bergmann | Berliner Betrieb für Zentrale Gesundheitliche Aufgaben, Institute für Lebensmittel, Arzneimittel und Tierseuchen, Berlin, Germany |
| Andy Ward | Unilever UK Foods, Purfleet, Essex |
| Margaret O'Sullivan | Public Analyst Laboratory (Cork), Eire |
| Irene Poulima | Food Division Laboratory, General Chemical State Laboratory of Greece |
| Michael Baden | Eurofins Scientific Group, Hamburg, Germany |
| Wolfgang Kruspe | Thüringer Landesamt für Lebensmittelsicherheit und Verbraucherschutz (TLLV) |
| Christine Schleich | Landesuntersuchungsamt, Institut für Lebensmittelchemie und Arzneimittelprüfung Mainz |
| Jürgen Brunner | LUA Sachsen |
| Beate Müller | Landeslabor Brandenburg, Standort Frankfurt (Oder) Fachgebiet L4, Labor Bestrahlte Lebensmittel |
| Brigitte Butz | Bayerisches Landesamt für Gesundheit und Lebensmittelsicherheit |
| Maria Iglovári | Laboratory of Foodstuffs and Chemicals, Budapest, Hungary |

2.3 Test Material Preparation and Handling

2.3.1 PSL Materials

Having consulted participants from round 1 it was agreed that duplicated screening measurements from a reduced sample set would be utilised for the screening tasks in round 2 and that blended samples would be introduced. Therefore a subset of 33 of the original 49 samples was selected avoiding some of the replicated products and some of the products which had overflowed the instrument counter in irradiated form. 50 pots of each of these samples were sent to Isotron plc and irradiated to a mean dose of 4.6 (± 0.1) kGy on 27 March 2006.

The sample allocation comprised all 33 products in irradiated and unirradiated forms, plus 6 blended products (which were also used in the TL study under different labels), making a total of 72 samples. Since the products had already been measured in the presentational order of round 1, it was decided to present them in a random order. After considering alternatives it was decided to use a single random sequence for all laboratories (Table D.2 in Appendix D). Participants were informed that the samples were in random sequence, but were not informed which samples from round 1 had been selected or their order of presentation.

The 6 products selected for blending aimed to include a range of sensitivities (from low to high), products (2 herbs, 2 spices and 2 dietary supplement matrices were used), and concentrations levels (10%, 1% and 0.1% irradiated material). The blending process is described below.

Table 2.2 lists the products used in round 2. Blends were produced from SP8531 (Ginger 10% irradiated), SP8533 (Rosemary 1% irradiated), SP8535 (Chilli powder 0.1% irradiated), SP9539 (Tarragon 10% irradiated), SP8577 (Siberian Ginseng 1% irradiated) and SP8578 (Green Tea 0.1%).

Table 2.2. Round 2 Test Materials List

| Test Material No. | Product Description | Abbreviation | SUERC Ref No | Suppliers Code |
|-------------------|-------------------------|--------------|--------------|----------------|
| 1 | Basil | Basil | SP8512 | Not supplied |
| 2 | Curry Powder No 3 | CP3 | SP8513 | Not supplied |
| 3 | Barbeque Seasoning | BBQ | SP8514 | Not supplied |
| 4 | Ground Cinnamon | Cinn | SP8515 | Not supplied |
| 5 | Paprika | Pap | SP8516 | Not supplied |
| 6 | Ground Cumin | Cum | SP8517 | Not supplied |
| 7 | Oregano | Oreg | SP8518 | Not supplied |
| 8 | Ground Black Pepper | GBP | SP8519 | Not supplied |
| 9 | Chives | Chives | SP8520 | Not supplied |
| 10 | Mint | Mint | SP8521 | Not supplied |
| 11 | Sage | Sage | SP8523 | Not supplied |
| 12 | Ground Mixed Spice | GMS | SP8525 | Not supplied |
| 13 | Ground Coriander | Cori | SP8526 | Not supplied |
| 14 | Medium Curry Powder | MCP | SP8528 | Not supplied |
| 15 | Ground Nutmeg | Nut | SP8529 | Not supplied |
| 16 | Whole Black Pepper | WBP | SP8530 | Not supplied |
| 17 | Ground Ginger | Ginger | SP8531 | Not supplied |
| 18 | Thyme | Thyme | SP8532 | Not supplied |
| 19 | Rosemary | Rose | SP8533 | Not supplied |
| 20 | Turmeric | Turm | SP8534 | Not supplied |
| 21 | Chilli Powder | Chilli | SP8535 | Not supplied |
| 22 | Garlic Powder (Chinese) | Gar | SP8536 | Not supplied |
| 23 | Parsley | Pars | SP8538 | Not supplied |

| | | | | |
|----|--------------------------------|-------|--------|---------------------|
| 24 | Tarragon | Tar | SP8539 | Not supplied |
| 25 | Milk Thistle Seed Powder | MTSP | SP8571 | P13225 Batch 041104 |
| 26 | Alfalfa Herb Fine Powder | AHFP | SP8572 | P01102 Batch 4978 |
| 27 | Saw Palmetto Berry Fine Powder | SPBFP | SP8573 | P1906 Batch 031103 |
| 28 | Dong Quai Root Powder | DQRP | SP8574 | P04083 Batch 031209 |
| 29 | Guarana Seed Powder | GSP | SP8575 | P07385 Batch 7781 |
| 30 | Powdered Gingko biloba Leaves | GBL | SP8576 | P0711 Batch 1518 |
| 31 | Siberian Ginseng Root Powder | SGRP | SP8577 | P0714 Batch 040206 |
| 32 | Green Tea Powder | GTP | SP8578 | P07358 Batch 040419 |
| 33 | Echinacea Powder | EP | SP8579 | P05005 Batch 101673 |

To double-randomise the material for distribution, the PSL pots were arranged in numerical order according to the random number allocated in Round 1 (e.g. 208) and then paired with a new random number from a freshly generated list, placed in this new order and allocated a random laboratory destination on the basis of this number. The allocation is tabulated in Appendix D below.

2.3.2 TL Materials

TL laboratories were consulted at the outset to determine an approximate number of samples which could be considered practicable for the trial, bearing in mind that the sample preparation and handling for TL analysis is substantially more laborious than PSL analysis. On this basis it was decided to present 18 samples to the TL laboratories, comprising the 6 products described above, prepared in unirradiated, irradiated and blended forms. Once again a fixed allocation sequence was used (table 2.5 below) which was not communicated to participants in advance.

The TL samples were all re-potted so that the blends had the same appearance to participants as the end member samples for each product. For those laboratories participating both in PSL and TL stages the blends appeared twice, in both sample sets. Two PSL participants joined the TL trial after the products had been dispatched. Since very limited quantities of the TL blends remained at that stage it was decided to re-use remaining portions of the PSL samples as TL samples for these two participants. The correspondence between the PSL and TL sample codes was therefore revealed to these two laboratories after they had completed their PSL analysis.

TL laboratories also received two irradiated LiF chips, to provide temperature calibration checks on the TL readers used in the study. 50 LiF chips were subjected to the following pre-characterisation steps. Firstly, all the dosimeters were irradiated to a 100mGy dose, pre-heated for 5 minutes at 90°C and then read out in a Riso DA15 TL reader from 0-400°C at 5°Cs⁻¹ to precondition them and check for batch uniformity. The chips were then annealed at 400°C for 60 minutes, irradiated again to 100 mGy, annealed again at 90°C for 30 minutes and packed into laboratory sets for distribution. The dose was selected to remain within the linear range of output for LiF. Thermal processing was conducted in an external muffle furnace with temperature logging during annealing.

2.3.3 Blending

Initially, 10 aliquots of the recently irradiated portion of each of the 6 chosen products were dispensed into petri dishes and 60s PSL screening measurements conducted. These are

tabulated below, together with results for the same products from the Round 1 reference measurements, enabling comparison between the 2 irradiations:

| Product | Ginger SP8531 | Rosemary SP8533 | Chilli SP8535 | Tarragon SP8539 | Ginseng SP8577 | Green Tea SP8578 |
|---------|---------------|-----------------|---------------|-----------------|----------------|------------------|
| Aliquot | | | | | | |
| 1 | 297286 | 78706 | 133300 | 48723 | 1291448 | 117128 |
| 2 | 243906 | 71226 | 140675 | 36385 | 1126778 | 110876 |
| 3 | 407599 | 80194 | 134466 | 30577 | 1289684 | 125699 |
| 4 | 648877 | 78897 | 162891 | 26768 | 989735 | 112917 |
| 5 | 308605 | 80507 | 178017 | 73402 | 1262244 | 116855 |
| 6 | 385120 | 54319 | 145837 | 55013 | 1099444 | 118582 |
| 7 | 285810 | 71209 | 115018 | 94350 | 1135077 | 107604 |
| 8 | 225466 | 76635 | 124421 | 25738 | 1030373 | 94960 |
| 9 | 459122 | 85925 | 164007 | 22984 | 1159521 | 114507 |
| 10 | 384922 | 86673 | 156297 | 55672 | 958185 | 141958 |
| Mean | 364671 | 76429 | 145493 | 46961 | 1134249 | 116109 |
| SD | 124825 | 9312 | 19682 | 23338 | 120108 | 12145 |
| SE | 39473 | 2945 | 6224 | 7380 | 37982 | 3841 |
| CV | 34 | 12 | 14 | 50 | 11 | 10 |
| R1 mean | 319085 | 68466 | 132939 | 45325 | 917869 | 103468 |
| R1SD | 102844 | 4944 | 27094 | 68201 | 88866 | 28037 |
| R1 CV | 32 | 7 | 20 | 150 | 10 | 27 |

Table 2.3. PSL screening of 10 aliquots of each product prior to blending

On the basis of these data, the concentration of blend chosen for each product was as follows:

| Product | Concentration | Expected mean PSL |
|----------------------|---------------|-------------------|
| Ginger SP8531 | 10% | 36467 |
| Rosemary SP8533 | 1% | 764 |
| Chilli powder SP8535 | 0.1% | 145 |
| Tarragon SP8539 | 10% | 4696 |
| Ginseng SP8577 | 1% | 11342 |
| Green Tea SP8578 | 0.1% | 116 |

Table 2.4. Blend concentrations and predicted mean PSL estimate

Also tabulated are the expected mean PSL terminal counts, based on a simple multiplication of the mean of the 10 test screenings of the spike material and the concentration used.

2.3.3.1 Mixing process

Taking each product in turn, starting with the lowest (0.1%) spike concentration, and working with one product at a time, the mixing was performed as follows. To provide the unirradiated portion of the blend, pots 181-200 (as defined prior to round 1) of each product were mixed together and weighed. If insufficient material was available, further pots in the range 95-100, left over from Round 1, were included until the desired weight was reached (499.5g for 0.1%, 495g for 1% and 450g for 10%) to achieve a total weight of blended product of 500g. Full details of the mixing process, which involved sub-sampling the mixes in 10 fold replication 6

times during mixing and recording PSL data to judge the homogeneity and extent to which the expected values had been achieved, are included in Appendix E.

Following the blending, samples for TL analysis were allocated numbers cyclically as shown below, and pots were labelled and filled, with 7g per pot. It was noted that this was quite a small amount of material to work with, but it was not possible to provide more.

| | U | I | B |
|-----------|----|----|----|
| Ginger | 1 | 2 | 3 |
| Rosemary | 6 | 4 | 5 |
| Chilli | 8 | 9 | 7 |
| Tarragon | 10 | 11 | 12 |
| Ginseng | 15 | 13 | 14 |
| Green tea | 17 | 18 | 16 |

Table 2.5. Sample number allocations for blends

2.3.3.2 Homogenising the paprika standard

In Round 1 the spread of results from the participants' screening of the paprika standard was considerable and an attempt was made to reduce this for round 2. 600g of irradiated paprika were mixed together, then blended for 5 minutes with the same equipment as the blind blends, turned with a spoon then blended for a further 5 minutes before 10 aliquots were read out. This cycle was then repeated; the results are as follows:

| Cycle | 1 | 2 |
|---------|-------|-------|
| Aliquot | | |
| 1 | 64361 | 32076 |
| 2 | 57313 | 34681 |
| 3 | 43924 | 55685 |
| 4 | 30249 | 32568 |
| 5 | 42244 | 79514 |
| 6 | 37751 | 71847 |
| 7 | 38817 | 28091 |
| 8 | 29651 | 76939 |
| 9 | 35410 | 34671 |
| 10 | 77161 | 29606 |
| Mean | 45688 | 47568 |
| SD | 15637 | 21182 |
| CV(%) | 34 | 45 |

Table 2.6. PSL screening results during paprika homogenisation

It was concluded that mixing made little difference; 35 pots of standard were then dispensed. The last 5 pots also contained paprika which had not been part of the mixing experiment; PSL analysis confirmed that this had no discernible impact.

Sample packs were distributed between 8 and 10 May 2006. The protocol which participants were asked to follow is presented in Appendix A.

3. PSL SCREENING

3.1 PSL Homogeneity testing (blends)

Data from round 1 which had been used to calculate the reference values and z-scores were carried forward into round 2 to perform the same function with the new screening data from participants. It was however necessary to generate new data for the blends, and for the calibrated PSL section of the study described below.

On 20 July 2006, the new blends were screened using the procedure followed for the homogeneity testing performed at SUERC in July and August 2005, but without the intervening measurements of the paprika standard.

Ten aliquots of the irradiated paprika reference standard were also measured to repeat instrument comparisons undertaken in Round 1. Prior to measurement the material had been homogenised, so the measurements were also used to assess the efficiency of the mixing process.

3.2 PSL Homogeneity results (blends and paprika)

The full results from the blends are tabulated in Appendix B. Table 3.1 presents the reference values calculated from these data.

| | Linear | | | | Log | | | |
|--------|----------|----------|----------|---------|----------|----------|----------|--------|
| | Mean | SD | SE | CV | Mean | SD | SE | CV |
| SP8351 | 21015.4 | 15521.63 | 4908.371 | 73.86% | 4.241665 | 0.259966 | 0.082208 | 6.13% |
| SP8533 | 691.4 | 500.9658 | 158.4193 | 72.46% | 2.768233 | 0.246304 | 0.077888 | 8.90% |
| SP8535 | 931.6 | 1059.578 | 335.0679 | 113.74% | 2.843079 | 0.292019 | 0.092345 | 10.27% |
| SP8539 | 5148.143 | 2754.43 | 871.0273 | 53.50% | 3.453039 | 0.27816 | 0.087962 | 8.06% |
| SP8577 | 10073.7 | 3606.483 | 1140.47 | 35.80% | 3.979467 | 0.150172 | 0.047488 | 3.77% |
| SP8578 | 484.8 | 108.7114 | 34.37758 | 22.42% | 2.675532 | 0.099058 | 0.031325 | 3.70% |

Table 3.1. Reference values for the blended materials (screening)

Reference screening data Round 2




Figure 3.1 Histogram of reference data to show position of blends

Results for the round 2 measurements of the paprika standard are tabulated below.

| Paprika Test Material | Terminal Counts |
|------------------------------|------------------------|
| 1 | 64361 |
| 2 | 57313 |
| 3 | 43924 |
| 4 | 30249 |
| 5 | 42244 |
| 6 | 37751 |
| 7 | 38817 |
| 8 | 29651 |
| 9 | 35410 |
| 10 | 77161 |
| Mean | 45688 |
| SD | 15637 |
| CV(%) | 34 |

Table 3.2. Results of 10 replicate measurements from the paprika reference material, measured by SUERC.

3.3 Participants' Results

30 participants submitted PSL screening results. Participants were asked to return their PSL terminal counts and classifications for each of the pair of replicates used, plus an overall evaluation, and to return copies of summary and raw data files created by the instrument for possible audit purposes. In all cases the completed spreadsheets with terminal counts and individual classifications were returned. Not all participants supplied primary data however, and few utilised the duplicated readings to form a final evaluation.

3.3.1 The paprika standard

Participants also returned data from the irradiated paprika reference standard. The means, standard deviations and standard errors are tabulated in Table 3.3 and provide another assessment of the relative sensitivity of each instrument to this material and the variation from round 1.

| Lab | Mean | SD | SE |
|-----|-------|-------|-------|
| 1 | 42168 | 14416 | 4559 |
| 2 | 28778 | 21696 | 6861 |
| 3 | 62103 | 18360 | 5806 |
| 5 | 90382 | 26901 | 8507 |
| 6 | 61277 | 15860 | 5015 |
| 8 | 56797 | 11612 | 3672 |
| 9 | 62871 | 21934 | 6936 |
| 11 | 45860 | 23096 | 7304 |
| 13 | 57281 | 21704 | 6863 |
| 14 | 45295 | 20294 | 6418 |
| 16 | 50322 | 13450 | 4253 |
| 17 | 46629 | 16878 | 5337 |
| 18 | 53061 | 16733 | 5291 |
| 19 | 57171 | 14653 | 4634 |
| 20 | 35810 | 4820 | 1524 |
| 21 | 55819 | 17220 | 5446 |
| 22 | 40082 | 14358 | 4541 |
| 23 | 55323 | 13973 | 4419 |
| 24 | 65138 | 27344 | 8647 |
| 25 | 43101 | 11990 | 3792 |
| 26 | 37161 | 13042 | 4124 |
| 27 | 24407 | 7104 | 2247 |
| 28 | 70393 | 62679 | 19821 |
| 29 | 55652 | 21351 | 6752 |
| 30 | 34882 | 9459 | 2991 |
| 31 | 29275 | 14383 | 4548 |

Table 3.3. Participants' paprika standards – mean, standard deviation and standard error from 10 observations

The relative sensitivities of the instruments have been re-calculated using the round 2 data both for the participants and for the SUERC reference; the participants' instrumental sensitivities are shown in Table 3.4, together with the corresponding values from the first round. In round 2 the overall results appear to be quite similar to those obtained before. The data from round 2 range from 0.53 to 1.98; compared with 0.53 to 1.45, with similar

uncertainties (typically 10-20%) in both cases. The mean sensitivity between those laboratories participating in both rounds is very similar (0.94 ± 0.24 compared with 1.10 ± 0.31). The general similarity seems to imply that, as observed in homogeneity testing, the additional attempts to homogenise the paprika standard have made little difference to outcome, and also it seems that the spread in instrumental performance is comparable over the gap of approximately one year between the two rounds.

| | Round 1 | | Round 2 | |
|------|-----------|-------|-----------|-------|
| Lab | Rel sens. | Error | Rel sens. | Error |
| 1 | 0.79 | 0.15 | 0.92 | 0.14 |
| 2 | 0.7 | 0.08 | 0.63 | 0.16 |
| 3 | 0.95 | 0.12 | 1.36 | 0.19 |
| 5 | 1.22 | 0.15 | 1.98 | 0.28 |
| 6 | 1.18 | 0.21 | 1.34 | 0.18 |
| 8 | 1.28 | 0.22 | 1.24 | 0.16 |
| 9 | 0.7 | 0.07 | 1.38 | 0.21 |
| 11 | 1.17 | 0.16 | 1.00 | 0.19 |
| 13 | 1.02 | 0.12 | 1.25 | 0.20 |
| 14 | 1.06 | 0.16 | 0.99 | 0.18 |
| 16 | 1.01 | 0.12 | 1.10 | 0.15 |
| 17 | 1.12 | 0.14 | 1.02 | 0.16 |
| 18 | 1.08 | 0.16 | 1.16 | 0.17 |
| 19 | 0.97 | 0.11 | 1.25 | 0.17 |
| 20 | 0.7 | 0.07 | 0.78 | 0.09 |
| 21 | 0.65 | 0.08 | 1.22 | 0.18 |
| 22 | 0.59 | 0.08 | 0.88 | 0.14 |
| 23 | 1.45 | 0.2 | 1.21 | 0.16 |
| 24 | 0.81 | 0.09 | 1.43 | 0.24 |
| 25 | 0.8 | 0.08 | 0.94 | 0.13 |
| 26 | 1.03 | 0.13 | 0.81 | 0.13 |
| 27 | 0.53 | 0.06 | 0.53 | 0.08 |
| 28 | 1.25 | 0.27 | 1.54 | 0.47 |
| 29 | 0.83 | 0.09 | 1.22 | 0.20 |
| 30 | 0.75 | 0.09 | 0.76 | 0.11 |
| 31 | 0.72 | 0.16 | 0.64 | 0.12 |
| Mean | 0.94 | 0.13 | 1.10 | 0.18 |
| SD | 0.24 | 0.05 | 0.31 | 0.07 |

Table 3.4 Relativity instrumental sensitivities derived from measurements of the paprika standard (10 each for rounds 1 and 2)

Figures 3.1 and 3.2 show the relationship between the paprika standard results from the two rounds, firstly as raw data and then as relative sensitivities. In both cases there are positive correlations, implying that the laboratories have largely retained the same relative sensitivities between both rounds. Both figures also show the error bars corresponding to the standard errors on the mean counts, which appear to account for much of the dispersion about the trend lines in the data. Regression lines are marked on the figures. The correlation coefficients are equivalent between the two plots ($r=0.545$ cf $r=0.546$), confirming that a significant proportion of the variance in sensitivity is reproduced between rounds. It appears also from examination of the error bars that much of the deviation from either the trend lines or lines of

equality through the data is accounted for by the uncertainties of the data, which are driven by the heterogeneity of the paprika. At this stage it is not clear whether the relative sensitivities to all samples are well represented by the paprika data, although it was noted in round 1 that there was a marked correlation between paprika relative sensitivity and the pooled Z scores for irradiated samples.




Figure 3.2 Mean paprika signal in photon counts obtained by laboratories in rounds 1 and 2




Figure 3.3 Relative instrument sensitivities from measurements of paprika standards in rounds 1 and 2 standardised to SUERC results

3.3.2 Screening results

Raw data, comprising the PSL intensities for each test material and participant, are tabulated in Appendix C. These have been represented graphically in Figure 3.4, which shows all the participants' data by product. Unirradiated samples seem be mostly represented between negative and intermediate bands; irradiated samples again as expected straddle intermediate and positive bands and exhibit marked sample to sample sensitivity variations. The blended samples are generally located between unirradiated and irradiated results, again as expected.




Figure 3.4 Participants' screening data for each product.

Participants' data are also presented as scatter diagrams for each of the categories of sample (irradiated, unirradiated and blended) in Figures 3.4-3.6, arranged by laboratory.

From Figure 3.5 it can be seen that, while there are still some observations from unirradiated materials which are above the upper threshold, this is far less prevalent than in round one. This suggests that sample handling has greatly improved overall. Separation of irradiated and unirradiated samples has also improved as a result. This is reflected in the z scores and qualitative scores discussed below. Examination of the order of presentation of samples, which clearly contributed in round 1 to cross-contamination problems, explains two of the outliers; these two samples were run consecutively and the preceding sample was a high sensitivity irradiated product. It appear that cross-contamination has persisted to a limited extent, but judging by the proportion of unirradiated samples giving positive band results a only minority of laboratories now appear to be experiencing such problems.




Figure 3.5 Participants' data for unirradiated samples by laboratory




Figure 3.6 Participants' data for irradiated samples by laboratory

When arranged by laboratory the raw data from irradiated samples cover 4-5 orders of magnitude, representing the range of sample sensitivity variation. Laboratory to laboratory sensitivity variations (eg as seen in the paprika) clearly are very small in comparison with the range of variations from sample to sample.




Figure 3.7 Participants' data for blended samples by laboratory and concentration

From Figure 3.6 it can be seen that all but 2 of the results from irradiated samples exceed the lower threshold and that the great majority also exceed the upper threshold. This is consistent with the reference data produced at SUERC.

The results from the blended samples, shown in Figure 3.7, indicate significant dependence between concentration and screening outcome. Interestingly a significant proportion of results at 10% and 1% concentrations fall into the positive screening band, above the upper threshold, with the majority of the remainder in the intermediate band. By contrast none of the 0.1% blends produced terminal counts above the upper threshold, although a significant proportion would apparently be detected above the lower threshold. This confirms the importance of classifying intermediate results as requiring further attention, as recommended by EN13751, particularly for samples with low concentrations of irradiated materials.

Perhaps unsurprisingly the data also show that blends would be missed during routine analysis where concentrations of the irradiated component are low. The upper limit of the range for blends (apart from one outlier all results are below 10^5) is below that for pure irradiated samples ($10^7 - 10^8$).

Figure 3.8 presents all the participants' data as a histogram. It appears that the separation of irradiated and unirradiated samples has improved since round 1, with very few unirradiated examples falling into the positive band. General sample handling, randomised presentation rather than alternating irradiated and unirradiated materials, and corroboration from duplicate aliquots may all have contributed to the observed improvement. The blends, again as observed above, show reduced intensities in comparison with irradiated samples, with some examples reaching the negative bands.

Participants' screening data Round 2




Figure 3.8 Histogram of participants' screening results for all products

3.3.3 Summary Statistics

As described in the Round 1 report (section 2.3), z-scores were calculated for the screening measurements. The same calculations were applied in Round 2, with additional assigned values being derived for the blends.

Z-scores are tabulated below to allow each participant to assess their outcomes on a test material by test material basis. Tables 3.5, 3.6 and 3.7 show Z scores for each unirradiated sample for laboratories 1 to 13, 14 to 25 and 26 to 35 respectively. Similarly tables 3.8, 3.9 and 3.10, present results for irradiated samples and tables 3.11, 3.12, and 3.13 for the blends. Graphical representation, which reveals patterns in the data, follows.

Table 3.5 Unirradiated Round 2 Test Materials: z-Scores for Participants 1-13

| Test Material | Lab 1 | Lab 2 | Lab 3 | Lab 5 | Lab 6 | Lab 8 | Lab 9 | Lab 11 | Lab 12 | Lab 13 |
|---------------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|
| 1 | 0.84 | -0.17 | 0.53 | 2.21 | 1.57 | -0.19 | 0.66 | 1.37 | -1.60 | 0.56 |
| | 1.39 | 0.40 | -1.94 | 0.70 | -0.03 | 2.60 | 0.97 | 1.50 | -0.59 | 1.24 |
| 2 | 2.52 | -0.10 | 1.33 | 4.46 | 11.22 | 4.23 | -2.03 | 3.87 | -0.10 | 2.35 |
| | 2.65 | 2.53 | -1.55 | 4.21 | 10.10 | 2.40 | 4.45 | 3.38 | 3.00 | 3.03 |
| 3 | -0.75 | -1.39 | -0.49 | 0.56 | 1.37 | 0.15 | 1.62 | -0.23 | 0.50 | -0.20 |
| | -1.26 | -0.79 | 0.15 | 1.84 | 1.49 | 0.33 | 0.68 | 0.04 | 1.85 | 1.29 |
| 4 | 1.50 | 0.03 | 0.25 | 2.58 | 9.46 | -0.56 | -0.33 | 2.06 | 0.74 | -0.23 |
| | -1.54 | 0.71 | -2.44 | 1.83 | 8.60 | 1.65 | -0.10 | 2.25 | 0.55 | 1.13 |
| 5 | 0.48 | 1.08 | 2.21 | 4.14 | 5.43 | 3.82 | 2.17 | 0.94 | 2.48 | 5.30 |
| | -0.28 | 2.07 | -0.18 | 0.16 | 2.62 | 0.87 | 2.96 | 4.08 | 2.67 | 3.40 |
| 6 | 0.36 | 1.38 | 1.50 | -2.10 | 4.96 | 0.35 | | 0.72 | 0.96 | 0.92 |
| | 1.37 | 2.36 | 1.60 | 1.65 | 4.02 | 2.02 | -6.56 | 1.63 | -0.09 | 0.82 |
| 7 | -1.73 | 0.99 | -2.06 | 1.52 | 6.84 | 2.24 | 0.06 | 2.38 | 1.76 | 0.45 |
| | -0.46 | 1.26 | 0.06 | 1.88 | 5.70 | -0.01 | 1.46 | 1.45 | -2.28 | 2.91 |
| 8 | -3.07 | 2.34 | -2.66 | 1.37 | 12.59 | 4.01 | -2.46 | 2.73 | -1.00 | 0.05 |
| | 1.28 | 0.18 | 3.62 | -0.01 | 11.08 | 2.10 | 1.22 | 2.34 | -9.64 | 2.47 |
| 9 | -0.43 | 0.88 | 6.87 | -4.45 | 6.36 | 4.46 | 1.37 | 0.58 | 1.40 | 2.57 |
| | 4.69 | -0.57 | -5.51 | 1.03 | 5.10 | 1.37 | 2.54 | 5.80 | 1.59 | 3.98 |
| 10 | -1.93 | 0.57 | -0.36 | 2.09 | 5.62 | -0.42 | -2.60 | 1.96 | 1.22 | 2.57 |
| | 0.08 | -0.24 | 1.39 | 0.38 | 0.78 | 2.85 | 0.46 | 2.27 | -0.27 | 0.38 |
| 11 | 1.53 | 1.42 | 1.07 | 0.71 | 2.06 | 1.55 | 0.44 | 1.57 | 1.04 | 0.28 |
| | 0.22 | 1.25 | 1.42 | -1.41 | 2.95 | -0.24 | -0.81 | 1.57 | -0.81 | 1.22 |
| 12 | -0.50 | 0.32 | -1.72 | -0.34 | -0.74 | 0.44 | -0.62 | -0.71 | 1.40 | -0.50 |
| | -0.49 | -0.60 | -0.94 | 0.83 | 1.63 | -0.36 | -0.95 | -0.25 | -0.29 | -0.38 |
| 13 | 1.60 | 1.63 | 0.92 | -1.70 | 6.50 | -0.01 | 2.16 | 1.93 | 1.36 | 1.88 |
| | 0.25 | 2.80 | 1.12 | -0.47 | 4.99 | 0.22 | -0.52 | 1.06 | 4.10 | 0.44 |
| 14 | -0.71 | -0.32 | 0.55 | -1.38 | 1.10 | 0.81 | -0.68 | 0.44 | -0.20 | -0.10 |
| | -0.60 | 1.40 | 3.06 | -0.19 | 0.86 | 0.64 | 0.73 | 0.24 | -0.44 | 0.95 |
| 15 | 0.59 | 0.26 | 2.59 | 0.05 | 1.39 | 0.40 | -0.19 | 1.27 | -0.60 | 0.26 |
| | -0.71 | -0.01 | 0.60 | 1.40 | 1.12 | 0.72 | -0.75 | 0.31 | 1.36 | 0.64 |
| 16 | 0.69 | 1.35 | -4.63 | 1.38 | 3.43 | 1.46 | -1.29 | 1.70 | 1.40 | 1.41 |
| | -0.18 | 1.62 | 1.23 | 2.57 | 2.49 | 0.78 | 0.98 | 1.13 | 0.37 | 2.02 |
| 17 | -0.30 | 0.90 | -0.34 | 0.16 | 2.19 | 0.64 | -0.53 | 0.22 | -1.81 | 0.26 |
| | 0.31 | 0.60 | -0.48 | -0.62 | 3.21 | 0.49 | -0.93 | 0.27 | -0.94 | 1.00 |
| 18 | 0.31 | 2.48 | -4.71 | 2.07 | 12.45 | 3.58 | 5.22 | 2.48 | 0.50 | 2.78 |
| | 0.29 | 1.36 | 1.79 | 1.89 | 11.53 | 2.08 | 5.22 | 2.86 | 0.90 | 5.75 |
| 19 | 0.58 | 1.61 | 0.39 | 0.41 | 2.24 | -0.10 | 0.39 | 0.67 | 3.45 | 0.67 |
| | 0.73 | 1.18 | 1.45 | 5.53 | 2.63 | 2.22 | 1.77 | 1.86 | -1.44 | -1.30 |
| 20 | 0.87 | 0.08 | 1.70 | 0.52 | 2.29 | 1.04 | 1.79 | 0.52 | 0.45 | -1.33 |
| | -1.00 | 1.44 | 0.78 | 1.29 | 0.97 | 1.39 | 1.67 | 1.37 | -1.98 | -0.21 |
| 21 | -0.85 | -1.14 | 0.36 | 1.59 | 2.86 | 1.40 | -3.19 | 0.62 | 1.65 | 2.20 |
| | -0.21 | -0.62 | 1.92 | 3.15 | 2.31 | 1.73 | 0.21 | 1.25 | 3.09 | 1.83 |
| 22 | -0.89 | 1.01 | -1.53 | -0.37 | -0.42 | 0.11 | 0.76 | -0.29 | 1.64 | 0.16 |
| | -0.74 | -0.43 | 0.17 | 2.92 | -0.51 | 1.47 | -0.98 | 3.52 | 0.66 | -0.53 |
| 23 | -5.12 | 6.36 | -0.61 | 0.56 | 1.25 | -0.86 | -0.82 | -0.33 | -2.82 | 1.34 |
| | 1.20 | 7.05 | -1.25 | 1.59 | 3.52 | 0.80 | -1.23 | 0.86 | 0.94 | -0.31 |
| 24 | -2.09 | -0.64 | -0.72 | -1.04 | 2.08 | 0.21 | -4.84 | 0.87 | -0.31 | 1.20 |
| | 0.36 | -0.60 | 2.05 | 1.27 | 3.18 | 1.04 | 1.32 | -0.02 | 1.33 | -0.31 |
| 25 | 1.44 | -0.50 | -0.32 | 1.88 | 2.44 | -0.39 | 0.37 | 0.86 | 0.17 | -1.26 |
| | -0.55 | -0.30 | -0.12 | 1.01 | 1.04 | 1.37 | -2.15 | 0.51 | 0.33 | 0.25 |
| 26 | -0.44 | -0.32 | 2.33 | -0.26 | 3.45 | 3.18 | 0.46 | -0.44 | -2.32 | 1.44 |
| | -0.17 | 0.46 | -0.71 | 1.68 | 4.44 | 2.94 | -1.37 | 1.47 | -6.03 | 2.27 |
| 27 | -0.84 | 9.77 | 0.36 | 1.45 | 4.33 | 0.90 | 0.75 | 0.81 | -0.11 | 0.07 |
| | 0.49 | 10.61 | 2.03 | 0.28 | 2.46 | -1.23 | 0.67 | 2.07 | -0.57 | -1.63 |
| 28 | -0.15 | 0.73 | 0.21 | 2.07 | 1.47 | 2.60 | -1.21 | 1.28 | 0.56 | 2.71 |
| | -0.64 | -1.75 | 1.95 | 1.38 | 1.22 | 1.63 | -0.71 | 1.44 | 3.84 | 2.74 |
| 29 | 0.22 | -0.40 | 1.03 | 2.39 | 3.41 | -0.44 | 0.65 | 2.23 | 2.06 | 1.70 |
| | 1.89 | 0.52 | 0.38 | 2.01 | 2.33 | 0.54 | 1.61 | 1.86 | 1.58 | 1.07 |
| 30 | -2.56 | -0.94 | -0.09 | 0.02 | 0.18 | 1.63 | -0.89 | -0.17 | 0.85 | 0.06 |
| | -0.84 | 0.26 | 1.66 | 0.02 | 2.38 | 0.01 | -0.41 | 0.35 | 1.94 | -0.38 |
| 31 | -1.41 | -1.46 | 0.16 | 3.23 | 2.22 | 1.62 | -0.34 | 0.22 | 1.97 | 0.93 |
| | -0.92 | -0.76 | 0.53 | 1.50 | 0.27 | 3.22 | -0.90 | 0.39 | 0.61 | 0.25 |
| 32 | -3.21 | 2.25 | 1.84 | 0.10 | 16.65 | 2.68 | 1.13 | 4.03 | 2.43 | 0.63 |
| | -1.07 | 10.01 | 5.52 | 2.03 | 17.66 | 0.37 | -3.79 | 2.29 | 2.75 | 1.48 |
| 33 | -0.58 | 0.74 | -0.46 | 0.01 | 6.16 | 1.31 | -3.26 | 0.19 | 0.87 | 1.07 |
| | -0.33 | -0.20 | 0.74 | 2.45 | 6.14 | 2.22 | -0.07 | 0.67 | 2.22 | 3.35 |

Table 3.6 Unirradiated Round 2 Test Materials: z-Scores for Participants 14-25

| Test Material | Lah 14 | Lah 16 | Lah 17 | Lah 18 | Lah 19 | Lah 20 | Lah 21 | Lah 22 | Lah 23 | Lah 24 | Lah 25 |
|---------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 1 | -2.27 | 0.04 | 2.41 | 4.00 | 0.21 | 1.69 | -0.98 | 4.61 | 3.13 | 6.15 | -2.21 |
| | 1.16 | -0.70 | 1.97 | 4.38 | -1.17 | 0.81 | 1.59 | 21.20 | 1.37 | -0.03 | -0.59 |
| 2 | 3.18 | 2.92 | 3.01 | 13.93 | 5.19 | 3.58 | 0.63 | 1.83 | 3.71 | 1.83 | 2.64 |
| | 1.93 | 1.62 | 5.80 | 16.55 | 2.87 | 2.77 | -1.90 | 6.47 | 5.95 | 5.42 | 4.08 |
| 3 | -0.85 | 0.32 | 1.19 | 0.76 | 0.21 | 1.47 | -0.62 | -1.97 | 0.33 | 1.17 | 0.57 |
| | 0.89 | 0.67 | -0.19 | 1.64 | 0.05 | 0.38 | -0.50 | -0.86 | 1.64 | 2.61 | 3.28 |
| 4 | 0.11 | 0.63 | 1.42 | 6.81 | 2.54 | 1.22 | -0.17 | 2.54 | 2.41 | 1.58 | 1.99 |
| | 0.60 | 2.91 | -0.35 | 6.22 | 0.65 | 1.20 | -0.39 | 7.29 | 2.50 | 0.11 | 0.66 |
| 5 | 0.48 | 3.91 | 1.32 | 3.59 | 0.00 | -0.03 | 2.45 | -1.20 | 0.42 | 2.43 | 0.72 |
| | -0.80 | 5.11 | 1.23 | 0.75 | 1.80 | 0.31 | 1.29 | 1.24 | 0.37 | 0.57 | 4.24 |
| 6 | 0.93 | 2.17 | 2.06 | 3.73 | 2.08 | 1.66 | 1.09 | 0.65 | 3.06 | 2.15 | 1.83 |
| | -3.81 | 2.12 | -0.23 | 1.65 | 0.54 | 1.92 | 4.92 | 1.75 | 1.59 | 2.08 | 1.50 |
| 7 | 0.59 | 0.04 | 1.27 | 5.45 | 2.02 | 0.06 | 0.98 | 0.47 | 4.67 | 1.81 | 0.29 |
| | -0.15 | -1.56 | 0.64 | 3.75 | 1.77 | 1.29 | 0.26 | 1.16 | 2.11 | 1.37 | 2.08 |
| 8 | 2.79 | 0.98 | 1.64 | 9.50 | 2.47 | 0.42 | -0.06 | 1.03 | 7.47 | 1.24 | 2.08 |
| | -0.78 | 0.82 | 3.02 | 7.69 | 2.32 | 0.74 | -2.05 | 4.74 | 2.95 | 1.26 | 1.05 |
| 9 | 5.74 | 2.77 | 14.36 | 8.47 | 4.18 | 0.64 | 0.12 | 15.22 | 5.43 | 3.09 | 2.47 |
| | -0.50 | -0.40 | 15.07 | 9.74 | 2.94 | 1.92 | 6.61 | 21.64 | 5.16 | 2.84 | -0.80 |
| 10 | 3.91 | -2.68 | 2.64 | 1.66 | 2.36 | -1.79 | -0.24 | -1.58 | 1.89 | -0.85 | 0.57 |
| | 1.39 | -0.42 | 1.17 | 1.98 | 1.22 | -0.18 | 2.23 | -3.50 | 0.19 | -0.45 | -1.41 |
| 11 | 0.37 | 0.95 | 2.34 | 2.05 | 1.84 | 0.71 | 0.17 | 0.12 | 1.67 | 1.24 | 0.93 |
| | -1.02 | 1.44 | 1.44 | 2.06 | -0.03 | 0.64 | 2.00 | 0.60 | 1.30 | 1.11 | 1.38 |
| 12 | -0.40 | -0.57 | -0.61 | -0.60 | 0.49 | -0.97 | -1.74 | -0.29 | -0.58 | -1.34 | -1.29 |
| | -0.67 | -0.54 | -1.27 | -0.79 | 0.38 | -0.82 | -0.47 | 1.06 | -1.04 | -1.09 | -1.19 |
| 13 | 2.75 | 0.16 | 3.08 | 1.03 | 2.12 | 2.15 | -0.18 | 7.97 | 0.97 | 1.83 | 2.48 |
| | 0.58 | 1.32 | 1.77 | 4.24 | 1.11 | 0.84 | -1.70 | 1.21 | 2.03 | 2.56 | 0.14 |
| 14 | 2.40 | -0.09 | 1.27 | 8.54 | -0.66 | 0.35 | -0.32 | 0.51 | -0.15 | 0.47 | 0.05 |
| | -0.59 | -1.85 | -0.88 | 0.69 | 0.04 | -0.39 | -1.06 | -0.06 | -0.32 | -0.01 | -1.07 |
| 15 | 3.34 | 0.48 | -0.33 | 0.74 | 3.63 | 0.41 | 0.11 | 0.17 | 0.91 | 1.12 | -0.31 |
| | 0.13 | 0.99 | 1.32 | 0.54 | 1.38 | 0.46 | 1.36 | 0.27 | 0.82 | 0.40 | 0.81 |
| 16 | 2.09 | 0.33 | 2.16 | 9.15 | 1.50 | 1.53 | -1.01 | 1.33 | 1.48 | 1.02 | 1.73 |
| | 1.07 | -0.25 | 2.31 | 4.61 | 1.79 | 0.00 | 3.50 | 0.47 | 1.96 | -0.77 | 2.30 |
| 17 | -0.04 | -0.45 | 0.15 | 3.97 | 3.13 | -0.65 | 1.61 | 2.67 | 0.45 | -0.74 | -0.58 |
| | -0.45 | -0.47 | 0.33 | 3.30 | 0.54 | -1.13 | -0.47 | 3.32 | 0.85 | 0.06 | -0.23 |
| 18 | 2.35 | -2.69 | -0.56 | 22.04 | 0.03 | 1.69 | 0.46 | 21.83 | 2.38 | 2.52 | 1.95 |
| | 0.43 | 0.39 | 5.64 | 18.86 | 2.16 | 2.91 | 1.50 | 22.20 | 2.11 | 8.79 | 2.48 |
| 19 | -1.14 | -0.92 | 1.09 | 4.32 | 0.01 | 1.72 | 1.56 | 0.33 | -0.97 | 1.45 | 1.41 |
| | 1.15 | -0.59 | 1.09 | 1.86 | -0.03 | 2.70 | -3.89 | 0.48 | 1.16 | -0.08 | 0.37 |
| 20 | 1.26 | 0.10 | 1.16 | 2.08 | 1.67 | -0.06 | -0.52 | -1.10 | 1.27 | 1.01 | 0.73 |
| | 1.61 | 1.98 | 1.34 | 2.84 | 1.52 | -0.19 | -0.17 | 0.78 | 2.51 | -0.17 | 0.92 |
| 21 | -0.31 | 0.96 | -0.02 | 2.02 | 0.57 | -1.03 | -0.96 | 0.15 | 1.91 | 2.33 | 0.09 |
| | 0.15 | 2.10 | 0.19 | 2.24 | 0.96 | -1.66 | 0.25 | 0.79 | 5.62 | 1.38 | 0.43 |
| 22 | -0.93 | 1.23 | -0.56 | 1.86 | 2.95 | 1.95 | 0.08 | -0.93 | 0.13 | 1.15 | 1.29 |
| | -0.38 | 1.31 | -0.35 | 1.59 | 0.96 | 1.88 | -0.22 | 0.98 | 0.63 | 2.34 | 1.24 |
| 23 | 1.34 | -3.27 | 2.41 | 1.84 | 1.28 | -0.69 | -0.11 | -1.60 | 0.73 | 3.04 | 0.11 |
| | 0.17 | -0.53 | -0.04 | 5.29 | 1.43 | -1.19 | -3.18 | 2.94 | 1.55 | 1.38 | -0.53 |
| 24 | -2.77 | -2.14 | -0.14 | 2.41 | 0.31 | -0.92 | 0.22 | -1.37 | 0.99 | -1.68 | -2.39 |
| | 1.36 | -3.04 | 0.84 | 3.37 | 1.92 | -0.04 | -2.36 | -3.18 | -1.04 | -0.82 | 1.47 |
| 25 | 6.03 | -0.94 | 1.01 | 0.55 | 0.55 | -0.93 | -0.37 | -0.04 | 0.90 | 22.10 | 0.38 |
| | 1.06 | -0.84 | 0.27 | 0.45 | 0.87 | -0.94 | 0.82 | -1.53 | 0.72 | 22.62 | 0.98 |
| 26 | -2.45 | -2.75 | 0.27 | 3.64 | 3.09 | -0.92 | -1.72 | -2.58 | 0.71 | 0.57 | -0.28 |
| | 2.49 | -0.34 | -0.44 | 3.60 | 2.61 | 0.02 | -2.67 | -0.73 | 2.73 | 1.19 | 0.96 |
| 27 | 3.10 | -0.92 | 0.43 | 2.11 | -0.33 | 0.20 | 1.20 | 1.22 | 1.51 | 1.15 | 1.11 |
| | -0.08 | 1.09 | 2.09 | 0.74 | 1.01 | 1.07 | -1.75 | -2.45 | 1.74 | 1.02 | 0.86 |
| 28 | -0.48 | 1.08 | 1.11 | 1.34 | 1.53 | -0.13 | 0.71 | 4.23 | 1.58 | 2.25 | 1.94 |
| | 0.32 | 1.74 | -0.63 | 1.86 | 1.47 | 0.53 | 0.29 | 1.00 | 4.33 | 4.27 | 1.20 |
| 29 | -0.42 | 1.39 | 1.21 | 8.13 | 1.11 | 1.14 | 1.64 | 2.30 | 20.44 | -0.09 | 1.92 |
| | 1.11 | 1.30 | 1.59 | 2.51 | 0.89 | 0.42 | -0.80 | 1.21 | 19.31 | 1.86 | 2.29 |
| 30 | -0.50 | -0.91 | 3.96 | 1.18 | 1.63 | -0.54 | 0.56 | -1.39 | 0.03 | 0.35 | -1.05 |
| | -1.47 | 0.14 | -0.69 | 1.55 | 2.76 | -0.31 | 1.77 | -1.62 | 0.11 | 0.98 | -0.57 |
| 31 | 0.20 | 1.38 | -0.23 | 2.55 | 1.62 | -0.15 | -0.45 | -0.61 | 3.06 | 3.67 | 0.94 |
| | 0.64 | 1.72 | -0.46 | 3.02 | 1.85 | -0.56 | 1.60 | 3.86 | 1.19 | 1.83 | -0.20 |
| 32 | 2.27 | -2.63 | 2.93 | 24.81 | 2.68 | -3.00 | -2.81 | 12.39 | 3.28 | 1.61 | 2.07 |
| | -3.96 | -1.25 | 10.05 | 24.30 | 4.56 | -2.78 | 2.40 | 19.47 | 2.38 | -0.06 | -0.17 |
| 33 | 0.72 | 1.54 | 1.84 | 4.93 | 1.74 | -1.10 | 1.20 | 6.21 | 1.42 | 0.56 | 0.10 |
| | 0.19 | 1.77 | 0.27 | 5.61 | 2.06 | 0.13 | -0.31 | -1.32 | 0.11 | 1.51 | 0.77 |

Table 3.7 Unirradiated Round 2 Test Materials: z-Scores for Participants 26-35

| Test Material | Lab 26 | Lab 27 | Lab 28 | Lab 29 | Lab 30 | Lab 31 | Lab 33 | Lab 34 | Lab 35 |
|---------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 1 | 0.04 | 0.59 | 4.79 | 22.59 | -0.96 | -0.89 | 2.27 | 4.31 | 1.30 |
| | 0.65 | -0.59 | 3.94 | 18.72 | 0.59 | -0.05 | 0.53 | 3.06 | 2.22 |
| 2 | 2.64 | 0.98 | 0.81 | 2.81 | 1.27 | 8.86 | 5.60 | 4.17 | 4.00 |
| | 1.13 | 4.55 | 2.61 | 5.15 | 2.72 | 8.61 | 6.57 | 3.18 | 3.75 |
| 3 | 0.44 | -1.36 | -0.87 | 2.10 | -1.33 | -0.94 | 0.05 | 2.97 | 0.11 |
| | -0.90 | -1.57 | 1.65 | -1.49 | -1.60 | -1.12 | -0.63 | -1.28 | -0.18 |
| 4 | 1.59 | -0.33 | 1.86 | 5.21 | -0.84 | 2.06 | 1.07 | -0.19 | 4.76 |
| | 1.47 | 0.65 | 0.25 | 5.11 | 2.15 | 0.85 | 0.89 | -0.31 | 0.74 |
| 5 | 3.68 | -1.36 | 0.68 | 1.62 | 2.08 | 5.42 | -0.31 | -0.25 | 1.95 |
| | 2.07 | -1.55 | 1.85 | 1.61 | -0.83 | 5.75 | 0.35 | 0.16 | 6.45 |
| 6 | 2.17 | -0.17 | 1.68 | -5.16 | 0.73 | 1.99 | -3.78 | -0.22 | 2.64 |
| | 0.32 | -0.51 | 0.15 | 0.34 | 0.71 | -11.22 | 41.31 | 1.15 | 1.19 |
| 7 | 3.40 | -1.05 | 5.20 | 0.21 | 1.15 | 0.06 | 0.95 | 0.92 | 1.26 |
| | 1.91 | -0.62 | -0.13 | -0.16 | 1.03 | -6.69 | -2.17 | 1.42 | 1.33 |
| 8 | 0.15 | 1.81 | 1.35 | 7.16 | -0.49 | -0.88 | 3.65 | 0.91 | 3.98 |
| | 0.67 | -1.19 | 0.15 | 9.19 | 0.05 | -0.37 | 3.79 | 1.62 | 1.24 |
| 9 | -1.48 | -1.15 | 7.71 | 36.81 | -0.11 | 1.06 | 2.44 | 0.49 | 4.71 |
| | 3.46 | 0.12 | 9.29 | 33.84 | 1.97 | -1.55 | 0.91 | 3.83 | 3.58 |
| 10 | -0.30 | 1.34 | 1.70 | 0.44 | -0.27 | 3.38 | 0.65 | 1.12 | 4.00 |
| | -1.61 | 0.81 | 4.19 | 3.35 | -3.00 | -3.72 | 1.27 | 0.41 | 1.80 |
| 11 | 3.35 | 0.60 | 2.27 | 1.05 | -0.34 | 1.35 | -6.50 | -1.51 | 1.96 |
| | 2.37 | 0.46 | 0.15 | 1.73 | 2.45 | -0.90 | 2.31 | 0.32 | 1.66 |
| 12 | -1.11 | -1.90 | 0.93 | -0.83 | -1.33 | -0.94 | -1.40 | -1.49 | -0.51 |
| | -1.65 | -1.43 | -1.65 | -1.53 | -1.72 | -1.30 | -0.86 | -1.18 | 0.82 |
| 13 | 2.28 | 0.12 | 2.63 | 1.66 | 0.93 | 1.95 | 41.18 | 0.46 | 0.38 |
| | 0.69 | -1.40 | -1.30 | 0.52 | 2.22 | 1.41 | 1.32 | 0.45 | 2.40 |
| 14 | -0.05 | -1.21 | 0.15 | 0.13 | -0.66 | -1.17 | 0.83 | 0.93 | -0.64 |
| | -0.25 | 0.02 | 0.41 | -1.48 | -0.13 | -1.51 | -1.45 | 0.15 | -0.26 |
| 15 | 0.27 | 0.09 | 1.19 | 0.33 | 0.56 | -0.20 | 0.58 | 1.55 | 1.06 |
| | 0.48 | 0.33 | 1.09 | -0.17 | 0.23 | 0.24 | 0.93 | 0.78 | 1.10 |
| 16 | 1.20 | 2.23 | 3.10 | 9.38 | 0.91 | 1.45 | 2.57 | 1.00 | 0.72 |
| | 1.78 | 1.02 | 1.59 | 5.26 | 1.79 | 1.33 | 2.50 | 0.76 | 2.08 |
| 17 | -0.19 | -0.69 | -0.28 | 7.34 | -0.80 | 2.16 | 1.63 | -0.95 | 1.58 |
| | -0.70 | -1.35 | 0.63 | 6.51 | -0.60 | 0.66 | 0.66 | -0.88 | 1.14 |
| 18 | 11.93 | 9.33 | 9.49 | 26.36 | 0.69 | 4.49 | 15.69 | 2.37 | 0.63 |
| | 7.37 | 8.93 | 7.92 | 31.47 | 2.03 | 1.33 | 10.06 | 0.41 | 1.83 |
| 19 | 3.75 | -0.79 | 1.66 | 1.54 | 2.48 | 2.74 | -0.46 | 2.50 | 1.98 |
| | -1.11 | 1.35 | 2.03 | 1.65 | -0.36 | -0.84 | 3.18 | 0.83 | 1.75 |
| 20 | -0.11 | -2.56 | 2.42 | -0.34 | 0.39 | 0.14 | 6.37 | 0.18 | 2.07 |
| | 1.64 | 0.33 | 0.67 | -0.63 | 0.81 | 0.87 | 1.34 | -0.43 | 1.48 |
| 21 | -0.09 | -0.04 | -1.01 | -0.63 | -0.17 | -0.11 | -0.22 | 0.44 | 0.78 |
| | 0.22 | -1.55 | -1.60 | -1.08 | -1.11 | 0.26 | -0.85 | 0.78 | 3.43 |
| 22 | 0.38 | -0.77 | 0.10 | -0.60 | -0.85 | -1.86 | -1.20 | -3.63 | -0.20 |
| | -0.94 | -0.57 | -1.84 | 0.88 | -1.29 | -1.71 | 0.05 | -3.11 | 0.74 |
| 23 | -1.07 | -0.25 | 0.22 | 1.09 | -1.10 | -0.27 | -0.99 | 0.02 | -0.21 |
| | 1.04 | 0.02 | -5.53 | -0.77 | -0.53 | -0.09 | -2.04 | -0.09 | 3.20 |
| 24 | -0.31 | -1.02 | -0.47 | 3.49 | 0.34 | 0.79 | -0.76 | 1.20 | 3.46 |
| | -0.57 | 0.29 | 1.86 | 2.04 | -0.36 | 1.50 | -3.40 | -0.92 | 2.59 |
| 25 | -0.22 | -0.28 | 1.48 | -2.29 | 1.13 | -0.22 | 0.67 | 0.70 | 5.68 |
| | 0.35 | -0.65 | 1.87 | 0.10 | 0.18 | -2.84 | 0.59 | -0.10 | 2.17 |
| 26 | -1.03 | -1.18 | -0.96 | 1.46 | 0.59 | 0.94 | 1.49 | -2.29 | 2.46 |
| | 2.13 | 0.83 | 3.56 | 0.72 | 0.46 | 0.17 | -3.45 | 0.67 | 0.23 |
| 27 | 4.07 | -1.04 | 1.35 | -1.19 | -0.44 | 1.14 | 2.09 | 0.04 | 1.25 |
| | 2.91 | -1.49 | -1.25 | 0.04 | -0.13 | 0.20 | -1.49 | -0.09 | 0.64 |
| 28 | 0.70 | -1.43 | 2.38 | 6.91 | -1.07 | -3.07 | -1.75 | 3.67 | 2.69 |
| | 5.37 | -2.75 | 1.38 | 5.34 | -1.12 | -2.30 | -1.30 | -0.01 | 0.67 |
| 29 | 0.05 | 2.01 | 3.74 | 9.89 | 0.68 | 1.44 | 1.65 | 2.47 | 0.07 |
| | 1.73 | 0.82 | 6.68 | 8.44 | 0.19 | -10.03 | 1.32 | 1.88 | 0.84 |
| 30 | -1.95 | -2.24 | -2.19 | 1.18 | -2.84 | -1.36 | -3.28 | -0.77 | 10.02 |
| | -1.18 | -1.28 | -0.28 | 1.91 | 0.17 | -1.15 | -2.01 | 0.42 | 1.78 |
| 31 | -0.71 | -1.27 | -0.53 | 1.45 | -1.41 | -1.66 | -2.23 | 1.21 | -0.16 |
| | -0.21 | -2.29 | 1.08 | 0.61 | -2.56 | -2.59 | -2.85 | -1.63 | 0.39 |
| 32 | 0.91 | 0.10 | 3.60 | 44.68 | -2.06 | 20.06 | 9.80 | -0.97 | 2.49 |
| | 2.89 | -0.41 | 1.42 | 52.80 | -7.34 | 20.88 | -0.41 | 3.23 | 2.20 |
| 33 | 0.91 | -1.43 | -0.48 | 1.26 | 0.01 | 5.88 | 1.45 | 1.54 | 2.19 |
| | -0.29 | -0.99 | -0.85 | -0.70 | 1.34 | 3.84 | -0.97 | 0.98 | 2.92 |

Table 3.8 Irradiated Round 2 Test Materials: z-Scores for Participants 1-13

| Test Material | Lab 1 | Lab 2 | Lab 3 | Lab 5 | Lab 6 | Lab 8 | Lab 9 | Lab 11 | Lab 12 | Lab 13 |
|---------------|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 1 | -1.64 | -0.97 | -1.25 | 0.33 | 0.71 | 0.60 | -1.24 | 5.40 | 2.44 | 1.42 |
| | -1.45 | -2.50 | -1.52 | 0.35 | 1.24 | 3.62 | -0.95 | 0.03 | 2.05 | 0.13 |
| 2 | -4.83 | -5.27 | -7.40 | -17.43 | -6.69 | -12.71 | -2.58 | -45.32 | -4.03 | -4.05 |
| | -5.33 | -23.89 | -19.10 | -12.48 | -12.57 | -11.52 | -21.07 | -23.65 | -6.11 | -3.78 |
| 3 | -2.48 | -4.95 | -8.82 | -4.03 | -3.49 | -4.21 | -7.73 | -2.66 | -6.37 | |
| | -2.14 | -2.70 | -8.16 | -12.62 | -9.61 | -3.42 | -5.59 | -4.83 | -2.39 | |
| 4 | -2.64 | -2.16 | 1.44 | 1.71 | -0.90 | 2.66 | -1.25 | -0.59 | 0.37 | 0.79 |
| | -1.48 | -1.01 | 2.88 | -0.38 | 0.17 | 4.21 | -0.69 | -1.44 | 1.20 | 1.69 |
| 5 | -2.44 | -1.21 | 13.03 | 4.37 | 10.03 | 4.10 | 7.45 | 6.82 | 9.62 | 1.13 |
| | -2.19 | 2.56 | 13.69 | 4.86 | 9.60 | 4.84 | 11.00 | 0.39 | 7.45 | 0.73 |
| 6 | -1.14 | -2.53 | 4.39 | 2.37 | 0.85 | 1.94 | -1.62 | 1.36 | 1.47 | 0.21 |
| | -2.25 | -1.37 | 4.34 | 2.94 | 2.43 | 1.88 | -1.20 | 2.03 | 2.77 | 0.42 |
| 7 | -2.13 | -4.11 | 0.11 | -0.73 | 3.79 | 1.82 | -4.82 | 0.54 | -0.49 | 0.69 |
| | -1.39 | -2.31 | -2.44 | 0.66 | 3.65 | 0.04 | -3.22 | 1.86 | 3.89 | 0.86 |
| 8 | -0.31 | -2.21 | -0.32 | 0.19 | -0.78 | 1.44 | -2.27 | 0.81 | 1.71 | 0.00 |
| | -1.31 | -1.19 | -2.21 | 1.77 | 0.54 | 1.48 | -2.72 | -0.27 | 2.07 | -0.03 |
| 9 | -0.47 | -2.19 | 0.11 | 3.32 | 0.15 | 1.20 | 1.34 | -0.13 | 1.25 | 0.50 |
| | -1.03 | -0.10 | -0.38 | 1.33 | 1.69 | 0.73 | 0.87 | 2.58 | 2.49 | 0.56 |
| 10 | -3.02 | 4.20 | 0.30 | 2.24 | 2.13 | 1.46 | -0.87 | 0.71 | 0.37 | 4.19 |
| | -0.99 | -2.65 | -0.49 | 3.24 | 3.86 | 0.50 | -1.06 | 0.92 | 2.47 | 2.06 |
| 11 | -1.76 | -4.29 | -1.21 | 1.12 | 0.28 | -1.12 | -0.69 | -0.44 | 2.17 | 1.61 |
| | -0.21 | -2.87 | -1.34 | 0.04 | 2.85 | 2.55 | -2.80 | 1.37 | 0.40 | 2.63 |
| 12 | -1.59 | -0.14 | -1.66 | -0.98 | -1.27 | -1.03 | -1.20 | -0.37 | -0.68 | -1.09 |
| | -0.88 | -1.46 | -0.96 | -0.50 | -0.80 | -0.10 | -1.79 | 0.13 | -0.19 | -0.57 |
| 13 | -2.48 | -4.72 | 2.73 | 0.45 | 0.36 | 3.69 | 0.67 | 0.85 | 0.77 | -0.82 |
| | -3.57 | -0.89 | -1.03 | 3.19 | -0.40 | 1.18 | -1.32 | 1.16 | 3.08 | -0.37 |
| 14 | -6.43 | -24.34 | -1.65 | -16.19 | -10.57 | -26.12 | 2.43 | -6.36 | -23.97 | -2.69 |
| | -5.28 | 2.48 | -4.65 | -17.62 | -16.13 | -22.13 | -23.68 | -22.64 | -4.91 | -2.81 |
| 15 | -2.62 | -3.09 | -0.65 | 1.30 | 0.00 | -0.53 | -0.94 | -0.28 | 0.56 | -0.53 |
| | 0.66 | -2.35 | -1.62 | 0.14 | -2.11 | -0.55 | -1.45 | -0.89 | -0.18 | -1.95 |
| 16 | -0.40 | -1.86 | -0.36 | 0.53 | -0.37 | 0.54 | -0.18 | -0.13 | 2.05 | -0.17 |
| | -2.03 | -0.97 | -0.14 | 0.85 | -2.41 | 0.41 | -0.33 | -1.48 | -0.24 | -0.61 |
| 17 | -0.23 | -2.52 | -1.27 | 0.01 | 0.37 | 2.29 | -0.81 | -0.37 | 0.47 | 1.22 |
| | 0.86 | -1.89 | -0.39 | 1.82 | 1.30 | 0.98 | -1.23 | 0.75 | 1.04 | -0.67 |
| 18 | -3.38 | -3.49 | -1.45 | 2.05 | 0.77 | 0.30 | -2.66 | -0.62 | 7.59 | 2.02 |
| | 1.70 | -2.11 | -2.47 | 3.17 | 0.33 | 1.40 | -2.51 | 0.53 | 1.87 | 2.20 |
| 19 | -4.94 | -5.98 | -0.75 | 3.02 | 2.34 | 0.86 | -3.66 | 5.66 | 2.26 | 2.53 |
| | -4.89 | -5.20 | 1.24 | 1.66 | 3.84 | 2.32 | -4.03 | 1.81 | 3.21 | 1.88 |
| 20 | -3.15 | 1.32 | -2.00 | -2.11 | 0.70 | -0.05 | -1.12 | -0.74 | -0.05 | 1.70 |
| | -1.94 | 0.41 | -2.22 | -1.16 | -0.13 | 1.16 | -1.44 | 1.24 | 4.16 | 0.42 |
| 21 | -0.63 | -2.80 | -0.51 | 2.02 | 1.28 | 1.68 | -0.85 | -0.11 | 1.18 | 1.83 |
| | -2.27 | -1.53 | 2.56 | 2.51 | 2.39 | 1.07 | -0.33 | 0.98 | 1.16 | 1.16 |
| 22 | 0.89 | -3.57 | 2.02 | 4.73 | 3.08 | 2.89 | -3.03 | 1.84 | 1.62 | 0.28 |
| | 0.20 | -3.31 | -0.15 | 1.32 | 1.60 | 4.89 | -4.64 | 3.54 | 1.07 | -0.43 |
| 23 | -1.77 | -2.16 | 3.74 | -0.09 | 4.17 | 0.67 | -0.76 | 3.41 | 1.14 | 5.27 |
| | 0.23 | -0.64 | 1.79 | 0.95 | 1.64 | 1.95 | -0.89 | 4.99 | -0.29 | 4.70 |
| 24 | -0.93 | 0.40 | 1.14 | -0.74 | 0.54 | 0.85 | 0.48 | 0.37 | 0.63 | 0.31 |
| | -0.33 | -0.69 | 0.57 | -0.07 | 0.71 | 0.28 | -0.20 | 0.01 | 1.04 | 0.04 |
| 25 | -0.23 | -1.78 | -0.46 | 2.32 | 1.65 | 0.73 | -0.17 | -0.71 | 1.08 | 0.68 |
| | -1.53 | -1.11 | 0.51 | 0.35 | 0.47 | 3.56 | -0.02 | -0.05 | 2.14 | 0.79 |
| 26 | -0.88 | -1.17 | 0.57 | 1.11 | 2.74 | 2.77 | -2.10 | 0.21 | 1.09 | 2.28 |
| | 0.56 | -2.51 | 1.00 | -0.27 | 3.74 | 2.96 | -0.72 | 1.14 | 0.27 | 1.24 |
| 27 | -0.64 | -0.99 | 0.39 | 0.28 | 0.54 | 0.24 | 2.78 | 0.35 | -0.41 | 0.81 |
| | | -0.85 | -0.62 | 2.21 | -0.11 | -0.19 | 1.01 | -0.64 | -0.02 | -0.41 |
| 28 | -3.20 | -4.73 | -0.91 | 1.04 | 1.25 | 2.08 | -2.98 | -0.14 | 2.52 | 1.87 |
| | -3.92 | -3.19 | 0.19 | 0.72 | 3.20 | 2.91 | -2.78 | 1.71 | 1.68 | 1.09 |
| 29 | -1.70 | -0.83 | -0.87 | 1.60 | 3.65 | -0.37 | 0.63 | 0.64 | 0.55 | 1.06 |
| | -1.48 | 1.30 | -0.72 | 0.03 | 0.47 | 1.92 | -1.04 | 0.98 | 0.09 | 1.34 |
| 30 | -3.63 | -3.92 | 1.45 | 5.06 | 1.57 | 4.98 | -3.86 | -0.08 | 4.57 | 3.53 |
| | -2.97 | -5.20 | 2.12 | 2.21 | 1.38 | 2.70 | -2.41 | 0.97 | 5.32 | 1.58 |
| 31 | -4.44 | -4.11 | -5.16 | 3.97 | 0.48 | 3.67 | -4.35 | 3.39 | 2.15 | 1.30 |
| | -3.27 | -2.83 | -5.18 | 4.25 | 2.58 | 5.17 | -1.44 | 3.00 | 1.50 | 0.58 |
| 32 | 1.03 | -0.87 | 0.78 | 0.93 | 0.60 | 0.74 | -1.43 | 0.95 | 0.38 | -0.57 |
| | -1.79 | -1.43 | 0.55 | 1.52 | 1.17 | 2.11 | -0.96 | 1.55 | 0.18 | -0.04 |
| 33 | -2.27 | -2.62 | 1.42 | 0.24 | 2.01 | 8.06 | -0.77 | 5.28 | 1.56 | -0.77 |
| | -2.10 | -1.45 | 1.02 | 0.50 | 1.88 | 4.78 | -1.54 | 2.34 | 5.81 | -0.70 |

Table 3.9 Irradiated Round 2 Test Materials: z-Scores for Participants14-24

| Test Material | Lah 14 | Lah 16 | Lah 17 | Lah 18 | Lah 19 | Lah 20 | Lah 21 | Lah 22 | Lah 23 | Lah 24 |
|---------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 1 | -0.24 | 0.54 | 1.26 | 0.06 | 2.34 | -0.13 | -0.64 | -0.05 | 1.52 | 3.75 |
| | -0.49 | 0.91 | 0.63 | 1.67 | 0.53 | 0.19 | -1.03 | -2.91 | 0.43 | 2.79 |
| 2 | -4.48 | -5.25 | -3.57 | -18.36 | -2.63 | -3.58 | -12.88 | -40.85 | -3.02 | -4.58 |
| | -3.94 | -5.00 | -8.48 | -6.05 | -2.34 | -3.30 | -5.23 | -15.89 | -6.73 | -4.06 |
| 3 | -10.31 | -6.30 | -3.47 | -4.01 | -1.92 | -3.37 | -4.55 | -3.57 | -6.19 | -5.39 |
| | -5.23 | -5.05 | -3.19 | -9.27 | -1.61 | -3.83 | -6.67 | -4.34 | -6.76 | -6.40 |
| 4 | -1.58 | 3.54 | 0.01 | -0.27 | 1.67 | -2.09 | -0.60 | -2.11 | 1.11 | -0.40 |
| | -1.66 | 2.98 | 0.36 | -0.19 | 2.50 | -2.37 | -1.87 | 1.21 | -0.07 | -1.37 |
| 5 | 4.34 | -0.61 | 3.53 | 6.58 | 0.82 | -0.55 | 7.16 | -3.01 | 0.07 | 0.34 |
| | 4.49 | 0.16 | 0.95 | 3.98 | 0.82 | -1.09 | 6.92 | 0.76 | -1.06 | -0.06 |
| 6 | -0.67 | 1.56 | 1.76 | 1.08 | 2.70 | -0.30 | -1.81 | -1.68 | 1.28 | -0.56 |
| | -0.75 | 2.22 | 0.79 | 1.50 | 1.28 | 0.21 | -0.29 | -2.35 | 2.32 | 1.29 |
| 7 | -2.92 | 2.05 | -1.14 | 0.29 | 0.53 | -2.00 | 0.27 | -3.02 | 2.07 | 0.16 |
| | -1.42 | -0.02 | -0.34 | -0.69 | 1.46 | -2.46 | 0.66 | -2.87 | 1.84 | 1.42 |
| 8 | -0.41 | -1.04 | -1.49 | -1.24 | 0.37 | -2.34 | -1.92 | -2.71 | -0.93 | -1.40 |
| | -1.43 | -0.65 | -1.53 | 0.74 | -0.42 | -2.12 | -0.60 | -2.72 | -0.28 | -0.38 |
| 9 | -0.49 | 0.82 | 0.61 | 0.40 | 1.33 | -0.93 | 0.75 | -0.43 | 1.70 | 0.76 |
| | -0.53 | 1.46 | -0.12 | 0.46 | 0.73 | -0.10 | 0.27 | 0.74 | 5.49 | 0.55 |
| 10 | 0.76 | 0.64 | 0.20 | 0.16 | 0.62 | -1.79 | -1.15 | -2.64 | 3.10 | -0.77 |
| | -0.90 | 0.23 | -0.10 | -1.13 | 0.46 | -3.40 | -0.46 | -0.38 | -1.37 | -0.66 |
| 11 | 0.16 | -0.17 | 0.15 | -0.33 | 1.36 | -1.59 | -0.44 | -1.64 | 1.55 | 0.86 |
| | 1.93 | 1.28 | -1.80 | 0.31 | 0.87 | -2.15 | -0.95 | -2.61 | 0.60 | 2.13 |
| 12 | -1.59 | -0.20 | -1.76 | -0.96 | -0.32 | -0.91 | -1.12 | -1.25 | -0.67 | -1.11 |
| | -0.98 | -0.12 | 0.21 | -1.14 | -1.43 | -1.30 | -1.15 | -1.95 | 1.54 | -0.11 |
| 13 | 1.41 | 0.14 | -0.97 | -2.03 | 0.51 | 1.21 | -2.90 | 2.30 | 1.97 | 0.81 |
| | -0.41 | -0.60 | -1.37 | 6.48 | -0.46 | 1.57 | -5.12 | -0.39 | 0.60 | 1.74 |
| 14 | -18.82 | -3.13 | -19.09 | 1.83 | -3.43 | -2.72 | 0.11 | -28.60 | -1.97 | -21.28 |
| | -12.71 | -3.47 | -26.23 | -19.43 | -0.44 | -2.41 | -2.24 | -16.70 | -2.56 | -15.00 |
| 15 | -2.12 | -0.58 | 0.79 | 0.88 | 0.63 | -2.73 | -1.82 | -1.98 | -0.63 | 0.62 |
| | -1.20 | -0.69 | 0.60 | -1.21 | -2.49 | -2.28 | -1.09 | -1.37 | 0.73 | 2.96 |
| 16 | -0.78 | -1.16 | 0.30 | -1.19 | 1.44 | -2.04 | -0.56 | 0.12 | 0.09 | 0.82 |
| | -0.65 | -0.22 | -0.18 | -0.11 | -0.31 | -2.21 | -2.43 | -1.25 | -0.63 | -0.52 |
| 17 | -0.75 | 0.15 | 0.09 | -2.44 | -0.40 | -0.99 | -0.86 | -1.27 | -0.25 | 1.42 |
| | -0.93 | 1.57 | -0.74 | -2.75 | -0.25 | -1.12 | 0.39 | -0.83 | 1.90 | 2.98 |
| 18 | 0.33 | 1.21 | -0.27 | -0.50 | 0.43 | -2.67 | -0.83 | -1.56 | 0.48 | 2.21 |
| | -0.53 | 1.28 | -0.01 | 0.12 | 1.53 | -2.16 | 0.03 | -2.61 | 0.89 | 3.04 |
| 19 | -0.39 | -0.37 | 0.08 | -2.43 | 2.38 | -2.67 | -3.37 | -2.55 | 1.88 | 2.86 |
| | 0.37 | 3.25 | 0.54 | 0.15 | 0.20 | -2.20 | 2.84 | -1.89 | 1.04 | 2.04 |
| 20 | -2.04 | 0.96 | 0.22 | -0.21 | 0.14 | -1.31 | -0.70 | 2.86 | 5.84 | 1.05 |
| | -1.68 | 0.73 | -1.07 | 1.23 | -0.76 | -1.99 | -1.59 | 1.03 | -0.37 | 1.26 |
| 21 | 0.99 | 1.08 | -0.34 | -0.31 | 2.78 | -0.34 | -1.12 | -0.68 | 0.87 | 2.33 |
| | -0.90 | 1.01 | -1.72 | -0.09 | 1.11 | -0.35 | -0.41 | 0.15 | 0.89 | 2.22 |
| 22 | -1.21 | 2.06 | -0.18 | 0.19 | 1.32 | -0.03 | -5.34 | -2.31 | 2.66 | 4.36 |
| | -1.14 | 3.31 | -0.10 | 0.91 | -0.03 | 0.72 | -1.85 | -4.23 | 0.76 | 3.85 |
| 23 | 2.26 | 1.30 | 0.08 | 0.31 | -0.05 | -0.32 | -0.77 | 0.67 | 0.33 | 3.12 |
| | -0.86 | 2.35 | 2.87 | -0.06 | 2.25 | 0.16 | -0.12 | 0.43 | 1.99 | 4.57 |
| 24 | 0.74 | 0.08 | -0.21 | -0.31 | 0.69 | -0.59 | 0.20 | -0.40 | -0.11 | 0.70 |
| | 0.09 | 0.20 | -0.09 | 0.52 | -0.13 | -0.69 | -0.10 | -0.30 | -0.07 | 0.98 |
| 25 | 0.25 | -0.88 | 0.43 | 1.37 | 1.18 | 0.18 | -0.88 | -0.09 | 1.48 | 0.78 |
| | 0.15 | 0.34 | 0.47 | 2.22 | 1.20 | -0.32 | -0.20 | -1.57 | 1.56 | 1.95 |
| 26 | -3.97 | 0.81 | -1.75 | 0.87 | 1.59 | -0.09 | -1.46 | -1.73 | 1.39 | 2.86 |
| | -0.44 | 1.40 | -1.99 | -0.40 | 0.21 | -0.44 | -1.09 | -0.94 | 0.21 | 3.62 |
| 27 | -0.49 | -0.32 | -0.89 | -0.84 | 0.15 | -0.31 | -0.24 | -1.12 | -0.34 | 0.98 |
| | 0.91 | -0.53 | -0.83 | 0.30 | -0.77 | -0.25 | -0.47 | 0.63 | 0.11 | 0.92 |
| 28 | -1.53 | 1.86 | -1.67 | -9.43 | 0.61 | -1.63 | -2.33 | -4.27 | 1.78 | 2.97 |
| | -2.52 | 2.64 | -2.16 | -1.95 | -1.05 | -1.10 | -2.89 | -5.48 | 0.87 | 3.88 |
| 29 | -0.45 | -0.51 | -0.70 | 0.52 | 0.56 | -0.84 | 0.15 | 1.36 | 0.22 | 1.47 |
| | -0.09 | -0.82 | 0.04 | 0.59 | 0.09 | -0.85 | -0.13 | -0.21 | 0.24 | 0.62 |
| 30 | 1.37 | 3.17 | 2.35 | 1.42 | 2.70 | -2.30 | 0.11 | -1.97 | 4.38 | 5.79 |
| | 0.47 | 2.51 | 1.41 | -0.10 | 3.60 | -1.32 | -4.22 | -1.62 | 4.23 | 6.94 |
| 31 | -0.57 | -4.71 | -1.69 | -1.40 | 2.22 | -2.89 | -3.11 | -4.09 | 0.25 | -0.83 |
| | -3.42 | -3.37 | -0.91 | -2.28 | 2.32 | -2.16 | -2.51 | -4.21 | 1.51 | -0.09 |
| 32 | -1.11 | 0.68 | -0.95 | -0.17 | 2.56 | -1.23 | -1.13 | -1.94 | 2.35 | 0.13 |
| | 1.90 | -0.06 | -0.77 | -1.15 | -1.17 | -1.68 | -0.62 | -0.68 | 0.21 | 0.50 |
| 33 | 1.15 | 4.80 | 0.25 | 6.55 | 4.79 | 0.02 | 1.05 | -2.75 | 3.50 | 4.14 |
| | 0.10 | 4.27 | -1.50 | 9.21 | 0.37 | -0.78 | 4.02 | -1.76 | 2.16 | 2.49 |

Table 3.10 Irradiated Round 2 Test Materials: z-Scores for Participants 25-35

| Test Material | Lab 25 | Lab 26 | Lab 27 | Lab 28 | Lab 29 | Lab 30 | Lab 31 | Lab 33 | Lab 34 | Lab 35 |
|---------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 1 | -0.87 | -1.63 | -3.59 | 0.04 | -1.06 | -2.95 | -4.17 | -2.88 | 0.07 | 2.65 |
| | 0.41 | -2.08 | -2.58 | -0.38 | -1.13 | -3.56 | -5.19 | -2.76 | 0.07 | 0.63 |
| 2 | -3.70 | -2.67 | -7.09 | -11.12 | -15.65 | -3.35 | -4.95 | -2.87 | -8.99 | -10.29 |
| | -2.91 | -31.12 | -5.07 | -18.28 | -16.07 | -19.85 | -7.57 | -4.66 | -11.10 | -4.38 |
| 3 | -8.76 | -4.87 | -2.88 | -3.20 | -2.45 | -3.54 | -4.60 | -2.95 | -2.52 | -6.03 |
| | -3.91 | -2.90 | -3.26 | -3.51 | -3.90 | -10.66 | -4.14 | -7.63 | -2.68 | -6.13 |
| 4 | -1.00 | 3.75 | -4.09 | 1.12 | 0.34 | -2.04 | -4.91 | -1.06 | -3.26 | 2.68 |
| | -0.95 | -1.93 | -3.21 | -0.22 | 2.95 | -3.30 | -5.38 | -0.08 | -3.25 | -2.02 |
| 5 | -0.29 | 7.24 | -3.40 | 7.78 | 8.37 | 1.02 | 4.20 | 8.30 | 2.19 | 2.59 |
| | -0.59 | 9.18 | -2.61 | 8.98 | 7.68 | -0.01 | 1.05 | 6.26 | -1.83 | 1.89 |
| 6 | 1.09 | -1.83 | -2.26 | 0.49 | -0.09 | -3.17 | -3.48 | -2.45 | -6.21 | -0.01 |
| | -0.99 | -1.71 | -2.97 | 1.47 | -0.92 | -2.42 | -5.37 | -2.25 | -2.59 | -1.81 |
| 7 | 0.66 | -1.71 | -4.02 | -1.78 | 1.44 | -2.67 | -4.04 | -2.67 | -17.11 | -0.16 |
| | -3.23 | -2.26 | -4.07 | -0.91 | 1.58 | -2.24 | -5.20 | -2.87 | -16.45 | 0.60 |
| 8 | -2.11 | -2.33 | -4.38 | -1.28 | -0.76 | -2.89 | -2.74 | -3.00 | -3.19 | 2.76 |
| | -2.31 | -2.32 | -4.31 | -1.79 | -1.63 | -2.14 | -3.23 | -0.86 | -2.66 | -1.86 |
| 9 | 1.17 | 0.97 | 0.53 | -0.41 | 0.01 | -0.34 | -0.70 | -0.25 | -0.23 | 0.20 |
| | -0.78 | -0.84 | -1.00 | -0.63 | 0.95 | -0.16 | -0.44 | -1.62 | -0.43 | 0.28 |
| 10 | -1.94 | -1.54 | -4.76 | -1.42 | 0.87 | -1.51 | -4.14 | -1.59 | -4.69 | 0.52 |
| | -1.74 | -1.50 | -4.80 | -1.57 | 0.61 | -2.39 | -5.07 | 1.99 | -4.29 | -0.48 |
| 11 | -0.75 | -0.65 | -4.37 | -1.33 | 2.08 | -2.75 | -3.61 | -2.22 | -11.90 | 0.86 |
| | -0.06 | -3.55 | -2.55 | -0.94 | 1.35 | 0.51 | -4.68 | -1.76 | -1.68 | -1.08 |
| 12 | -1.22 | -1.28 | 0.50 | -0.17 | -2.14 | -1.90 | -2.58 | -1.69 | -1.33 | -0.98 |
| | -1.45 | -1.28 | -0.95 | -1.05 | -0.93 | 0.27 | -2.92 | 0.15 | -1.60 | -0.47 |
| 13 | 2.75 | -3.91 | -4.80 | 0.47 | -0.68 | -1.78 | -5.14 | 2.43 | -0.13 | -0.93 |
| | -3.25 | -4.14 | -2.81 | -0.62 | -0.14 | -2.90 | -6.57 | 0.56 | -0.37 | -2.97 |
| 14 | 0.77 | -11.14 | -9.86 | -10.58 | -4.32 | -12.28 | -6.65 | 0.81 | -9.00 | -24.76 |
| | 0.25 | -2.14 | -3.82 | -18.71 | -13.80 | -3.67 | -0.66 | -3.74 | 1.01 | -15.60 |
| 15 | -0.18 | -3.33 | -3.80 | -0.49 | -1.79 | -3.33 | -5.15 | -2.45 | -1.34 | 0.81 |
| | 2.88 | -2.85 | -5.46 | -1.53 | -3.33 | -0.42 | -5.94 | -3.07 | -1.34 | -1.35 |
| 16 | 0.33 | -1.28 | -0.71 | -1.57 | -0.40 | -1.15 | -0.66 | -1.34 | -3.87 | -0.77 |
| | 0.78 | -0.01 | -2.25 | 0.01 | -0.80 | -0.93 | -1.08 | -1.48 | -2.98 | -0.87 |
| 17 | 1.07 | -1.80 | -1.89 | 0.13 | -0.11 | -4.27 | -2.26 | -1.82 | -1.58 | 0.27 |
| | -0.54 | -1.21 | -2.30 | 0.39 | -1.28 | -1.98 | -2.79 | -1.64 | -5.86 | -0.41 |
| 18 | 0.44 | -3.18 | -4.91 | -0.99 | -0.86 | -3.19 | -6.32 | -3.29 | -2.35 | -0.17 |
| | 2.02 | 0.38 | -5.15 | -1.36 | -0.65 | -6.22 | -7.58 | -2.00 | -0.92 | -0.15 |
| 19 | -0.75 | -4.65 | -12.65 | -1.84 | 2.12 | -5.82 | -10.54 | -5.14 | -12.92 | 1.61 |
| | 0.03 | -5.70 | -10.71 | -3.76 | 2.13 | -10.36 | -14.45 | -7.82 | -22.88 | -0.52 |
| 20 | -2.50 | -1.44 | -4.99 | -0.38 | 2.03 | -0.43 | -2.13 | -1.77 | -6.05 | 0.88 |
| | -1.33 | -3.31 | -4.17 | -1.57 | 1.90 | -1.97 | -3.36 | -2.38 | -4.45 | -0.95 |
| 21 | -0.33 | -0.91 | -4.54 | 0.36 | -1.25 | -2.78 | -4.38 | -2.00 | -4.35 | 1.02 |
| | 0.74 | -2.60 | -4.10 | 0.19 | -1.56 | -1.70 | -5.45 | -1.90 | -4.07 | 1.27 |
| 22 | 0.34 | -0.97 | -7.09 | 0.60 | -1.45 | -2.99 | -3.65 | -3.79 | -3.27 | 0.09 |
| | 1.16 | -3.39 | -5.33 | -1.10 | -0.75 | -3.46 | -5.22 | -3.93 | -13.76 | 0.49 |
| 23 | -0.31 | 0.62 | -1.00 | 0.67 | 5.70 | 1.58 | -1.35 | -0.48 | -1.81 | -1.37 |
| | -0.67 | 0.90 | -1.87 | -0.30 | 0.47 | -1.53 | -1.46 | 0.76 | -2.72 | 0.08 |
| 24 | 0.54 | -0.20 | -0.11 | 0.23 | 0.20 | -0.54 | -1.11 | -0.06 | -1.40 | 0.08 |
| | -0.61 | -0.31 | -1.09 | 0.09 | -0.30 | -0.96 | -1.42 | -0.87 | -0.69 | 0.17 |
| 25 | 1.13 | 0.23 | -1.78 | -0.32 | -0.54 | -0.86 | -1.74 | 0.50 | -2.42 | 1.13 |
| | -0.02 | -0.37 | -0.39 | -0.40 | -1.00 | 0.57 | -2.73 | -0.88 | -1.94 | -0.29 |
| 26 | -0.87 | -1.21 | -4.80 | 2.03 | 0.33 | -0.44 | -3.52 | -3.17 | -6.49 | -1.11 |
| | -0.85 | -2.65 | -4.37 | -0.99 | -1.20 | -3.26 | -4.68 | -1.02 | -5.40 | 0.14 |
| 27 | -0.39 | -0.79 | -4.47 | -0.55 | -1.29 | -1.41 | -1.36 | -0.83 | -0.97 | -0.73 |
| | -0.46 | -0.82 | -0.94 | 0.82 | -0.77 | -1.12 | -1.57 | -0.47 | -0.48 | -0.15 |
| 28 | -0.10 | -5.38 | -7.82 | -1.48 | -0.64 | -1.56 | -7.37 | -2.63 | -11.17 | -1.21 |
| | -0.77 | -3.19 | -4.81 | -1.16 | 0.19 | -5.38 | -9.48 | -4.63 | -12.20 | -1.18 |
| 29 | -0.21 | -1.29 | -1.60 | 0.27 | 3.43 | -0.96 | -1.46 | -1.61 | -1.27 | 0.21 |
| | -0.71 | -1.10 | -0.49 | -0.36 | 3.93 | -1.01 | -1.72 | -2.85 | -0.55 | 0.62 |
| 30 | -0.04 | -2.64 | -8.19 | 0.58 | -0.91 | -2.49 | -8.48 | -4.11 | 2.41 | 1.20 |
| | -1.41 | -4.06 | -7.43 | 1.79 | -0.44 | -4.29 | -11.06 | -5.26 | 4.50 | -1.27 |
| 31 | -0.97 | -5.27 | -10.46 | -0.84 | -0.19 | -7.66 | -11.71 | -6.20 | -8.05 | -1.76 |
| | -1.86 | -3.01 | -9.57 | -1.32 | 1.06 | -7.56 | -14.23 | -6.66 | -9.70 | -1.87 |
| 32 | 2.41 | -2.45 | -3.44 | -0.74 | 1.11 | -2.09 | -3.08 | -2.04 | -2.44 | -0.42 |
| | 1.07 | -2.30 | -3.21 | -0.71 | -1.38 | -1.40 | -3.98 | -1.47 | -3.52 | -1.00 |
| 33 | -0.98 | -5.62 | -7.13 | 1.39 | 0.32 | -3.72 | -4.37 | 0.33 | -8.74 | 1.93 |
| | 5.83 | -3.42 | -3.76 | -1.60 | 12.76 | -3.95 | -6.36 | 1.63 | -5.13 | 4.08 |

Table 3.11: Blended Round 2 Test Materials: z-Scores for Participants1-13

| Test Material | Lab 1 | Lab 2 | Lab 3 | Lab 5 | Lab 6 | Lab 8 | Lab 9 | Lab 11 | Lab 12 | Lab 13 |
|---------------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|
| 17 | 1.54 | -0.88 | 0.35 | 0.77 | 4.54 | 0.43 | -0.52 | 1.30 | 0.21 | 0.87 |
| | 0.55 | -1.06 | 2.38 | 1.34 | 0.97 | 1.45 | 0.10 | 0.24 | 2.05 | 1.18 |
| 19 | -0.30 | -0.08 | -0.01 | 2.10 | 2.30 | 0.23 | 0.07 | 1.52 | 0.53 | 0.34 |
| | -1.58 | -0.87 | 3.03 | -0.80 | 0.67 | 0.13 | -2.19 | -0.67 | 0.47 | 0.05 |
| 21 | -1.35 | 2.88 | -0.95 | 1.44 | 1.56 | 2.07 | -0.30 | 1.74 | 1.86 | 2.34 |
| | -1.19 | -0.14 | -2.35 | 0.37 | 5.40 | 1.66 | 0.65 | 1.49 | -1.10 | -0.91 |
| 24 | -0.41 | -0.81 | -0.59 | 0.00 | 0.20 | 0.43 | -0.66 | 0.27 | 0.25 | 1.80 |
| | -0.84 | -0.70 | 1.36 | 0.38 | 1.90 | 1.56 | 0.50 | 0.54 | 0.06 | 1.48 |
| 31 | -0.29 | 0.02 | 0.07 | 0.16 | 1.82 | 1.09 | -0.57 | 0.60 | 0.29 | -0.49 |
| | -1.30 | -1.36 | -0.44 | -0.52 | 2.74 | -1.00 | -1.94 | -0.77 | -0.22 | 0.02 |
| 32 | -0.62 | 0.00 | 0.47 | 1.11 | 1.11 | 3.59 | 0.77 | 3.06 | 0.92 | 0.80 |
| | -0.99 | 1.10 | 0.13 | 1.28 | 4.79 | 1.71 | 0.74 | 9.03 | 0.12 | -0.82 |

Table 3.12: Blended Round 2 Test Materials: z-Scores for Participants 14-24

| Test Material | Lab 14 | Lab 15 | Lab 16 | Lab 17 | Lab 18 | Lab 19 | Lab 20 | Lab 22 | Lab 23 | Lab 24 |
|---------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 17 | 1.81 | 0.81 | 0.71 | 0.09 | 1.81 | -0.78 | 0.34 | 1.33 | -5.79 | 2.01 |
| | 1.04 | -0.23 | 2.26 | 0.58 | 0.66 | -0.11 | 0.55 | -0.27 | 1.58 | 2.42 |
| 19 | -0.25 | 0.23 | -0.34 | 2.26 | -0.53 | -0.03 | -0.13 | 1.04 | 1.33 | 0.82 |
| | 0.72 | -0.84 | -0.84 | 1.39 | -0.13 | -1.79 | -1.01 | -0.64 | -0.60 | 0.23 |
| 21 | -0.56 | 5.43 | 1.79 | -0.20 | 1.46 | -0.35 | -0.14 | -1.98 | 0.72 | 2.17 |
| | -0.71 | 0.42 | 0.74 | 3.81 | 1.50 | -0.03 | -0.52 | -0.02 | 0.19 | 2.29 |
| 24 | 1.15 | 0.30 | 1.96 | 0.05 | 1.27 | -0.71 | 0.55 | -0.19 | -6.38 | 2.33 |
| | 1.99 | 0.29 | 0.98 | 0.68 | 1.71 | 0.21 | 0.25 | 0.47 | 2.47 | 1.00 |
| 31 | -0.55 | 0.02 | -0.17 | 2.46 | 0.64 | -0.29 | -0.07 | 2.64 | 0.53 | 0.29 |
| | 1.13 | -0.57 | 0.16 | -0.09 | -0.69 | -1.73 | -0.81 | -0.94 | -0.15 | -0.09 |
| 32 | 3.02 | 5.41 | 0.43 | -0.65 | 1.48 | -1.26 | 1.41 | -0.40 | 0.50 | 3.69 |
| | 0.25 | -0.10 | 0.42 | 4.30 | 0.52 | 0.17 | -0.68 | 0.50 | 0.48 | 2.19 |

Table 3.13: Blended Round 2 Test Materials: z-Scores for Participants 25-35

| Test Material | Lab 25 | Lab 26 | Lab 27 | Lab 28 | Lab 29 | Lab 30 | Lab 31 | Lab 33 | Lab 34 | Lab 35 |
|---------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 17 | 1.66 | 0.07 | -0.80 | 0.90 | 0.39 | 0.61 | -1.08 | -1.09 | -0.59 | 2.29 |
| | 0.68 | 2.12 | -0.05 | 1.78 | 0.40 | 0.20 | -0.23 | 0.46 | -0.49 | 0.91 |
| 19 | 1.02 | -0.28 | 0.07 | 0.09 | 0.28 | 0.85 | 1.48 | 1.30 | -0.74 | 0.62 |
| | -0.58 | -1.41 | -1.03 | -0.85 | -0.49 | -0.96 | -1.86 | -1.23 | -1.74 | 1.49 |
| 21 | -0.02 | -0.46 | -2.18 | -0.11 | 0.45 | -0.65 | -2.81 | -0.29 | 0.13 | 0.64 |
| | -1.22 | 0.94 | -0.29 | 1.55 | 4.96 | -1.23 | -0.11 | -1.05 | 0.61 | 0.47 |
| 24 | 1.38 | -0.30 | -0.85 | 2.72 | 0.36 | -0.80 | -1.42 | 0.19 | -1.40 | 0.64 |
| | 0.33 | 1.15 | -1.16 | 0.35 | 0.40 | -0.06 | -0.15 | 0.27 | -0.26 | 0.70 |
| 31 | 0.44 | -0.42 | 0.06 | 0.46 | -0.03 | -0.69 | 0.75 | -0.06 | -0.15 | 0.58 |
| | -0.38 | -1.36 | -1.27 | -0.13 | -0.66 | -1.67 | -1.86 | -1.30 | -1.93 | 0.39 |
| 32 | 0.16 | 0.94 | -2.06 | -0.20 | -0.62 | -1.28 | -3.37 | -0.31 | -0.08 | -1.09 |
| | 0.39 | -1.56 | -0.88 | 1.55 | 4.29 | -0.61 | -0.05 | -0.55 | -1.14 | 0.26 |

Figures 3.9-3.12 show the z-scores for irradiated, unirradiated, first by sample and then by laboratory.

In figure 3.9 it appears that for the unirradiated samples there are still a number of products with a wider range of z-scores than the majority. This suggests that these may be more problematic than the other materials. Samples 6 (Ground Cumin SP8517), 9 (Chives SP8520), 13 (Ground Coriander SP8526), 18 (Thyme SP8532), 25 (Milk Thistle SP8571), 29 (Guarana SP8575) and 32 (Green Tea SP8578) have produced outlying results. For these products with dispersed z-scores it is clear that several laboratories have produced the results appearing as outliers; with the use of duplicate aliquots a pair of outlying observations may represent a single laboratory's results from one product. It is not at this stage clear why these materials should have presented difficulties to some laboratories; they are not, for instance, multi-ingredient or salt-containing, nor is the order of presentation an explanatory factor for all of them.

When analysed by laboratory, as in figure 3.10 it is again apparent that there are more positive than negative outlying z scores, suggesting that contamination of unirradiated samples may be occurring occasionally, but that it affects only a minority of laboratories in these data. In comparison with round 1 this represents a considerable improvement in results, although it appears that a small number of laboratories could still find room for improvement. Other than for these laboratory-sample combinations the majority of unirradiated test materials show closely zero centred z scores.

For the irradiated samples, shown in figures 3.11 and 3.12, there is a greater spread of z-scores. When analysed by sample there are three samples with a larger number of outliers. Samples 2 (Curry powder no 3 SP8513), 14 (Medium Curry Powder SP8528) and 19 (Rosemary SP8533) show systematically lower z scores than the other samples, implying that participants were unable to retrieve the same signal levels as the reference laboratory despite the similarities of response to the paprika standard. The first two of these are salt-containing products and it is possible that participants have still not corrected for counter overflows. Without the full set of PSL data files it is not possible to confirm this. For the other sample another explanation is needed. When arranged by laboratory, it appears that all participants have at least one or two samples showing low-z score values, implying some sort of signal loss or under-estimation relative to the reference data. Some laboratories show more extreme values or a greater proportion of such values, for reasons which are not at this stage clear.

To assess the possibility that some fading has occurred between the irradiation data and the time at which participants had measured their samples the results were examined as a function of post-irradiation delay. In the case of the reference data the values had been determined during round 1, using material that had been irradiated in 2005 and measured during the study with a post irradiation delay of approximately 40 days. Participants have reported their measurement dates, from which it was possible to prepare figure 3.15. There is no obvious relationship between Z scores from irradiated samples and post irradiation delays ranging from 45 to 112 days. Apparently fading is unlikely to explain the negative z-scores. A plot of only those samples with negative outliers also supports this. Considering other potential explanations it is noted that laboratory bleaching might be a factor. Also it is noted that the irradiation dose received in round 1 was close to the requested value of 10 ± 1 kGy (8.5-9.9 kGy), whereas in round 2 a lower dose was delivered. However from table 2.3 it appears that the signal levels from the Rosemary sample using in blending were very close to the original reference values despite the dose difference, making it unlikely that dose response effects are responsible.

Round 2 Z scores unirradiated samples




Figure 3.9 Participants' z-scores for unirradiated samples by sample

Round 2
Zscores for Participants data from unirradiated samples




Figure 3.10 Participants' z-scores for unirradiated samples by laboratory




Figure 3.11 Participants' z-scores for irradiated samples by sample

Round 2
Zscores for Participants data for irradiated samples
corrected




Figure 3.12 Participants' z-scores for irradiated samples by laboratory




Figure 3.13 Time elapsed since irradiation vs z-score for all irradiated samples – participants' data

For the blends, displayed in figures 3.14 and 3.15, the data sets are quite well zero centred with only a very small number of outlying observations, which seem to occur randomly as far as dependence on sample or laboratory is concerned. There is no significant difference between samples, nor any marked dependence on the concentrations of the blends. From this it appears that the performance with the 6 blended samples was well matched between reference and participants' data sets.




Figure 3.14 Participants' z-scores for blended samples by sample

Round 2
Zscores for Participants data for blended samples




Figure 3.15 Participants' z-scores for blended samples by laboratory

Figures 3.16, 3.17 and 3.18 present participants z-scores in the form of 2-d colour contour plots. This forms of presentation permits visual identification of the interaction between test materials and participants with outlying results. Figure 3.16 confirms the observations concerning test material handling of unirradiated materials in the participants identified above and but does not indicate that any particular sample uniformly presented difficulties. In Figure 3.17 the irradiated test materials from three samples tended to be underestimated and two to

be overestimated. Figure 3.18, which shows the blended samples, indicates a more even distribution of results.




Figure 3.16 Contour plot for participants' data for unirradiated samples

Z-scores for irradiated samples




Figure 3.17 Contour plot for participants' data for irradiated samples

Z-scores for blended samples




Figure 3.18 Contour plot for participants' data for blended samples

3.3.4 Qualitative Results

As with round 1, participants' results were classified into the 3 conventional screening bands (negative, intermediate and positive) using the two arbitrary threshold values (set at 700 and 5000 counts in keeping with BS EN 13751:2002) established during development of the method and in the validation studies. These threshold values do not take account of individual test material sensitivities, and therefore may be of limited applicability to matrices for which research or validation data are missing. Neither do they take account of individual instrument sensitivities and one potential opportunity furnished by this study would be to assess the benefits of standardising data sets prior to classification to results from the paprika reference standard or test material mean values prior to classification. At this stage however the qualitative outcomes are presented in descriptive form based on the data as received.

Table 3.14 shows these data tabulated by participant. Table 3.15 shows them by test material. Table 3.16 presents the qualitative outcomes of the reference data for the blends. Table 3.17 shows the qualitative outcomes from round 2 in comparison with those obtained from round 1 and reference data. Table 3.18 examines the results from blended samples as functions of concentration and the brightness of the material.

Overall the performance is again quite encouraging, with remarkably consistent performance from laboratories between the two rounds. As can be seen from table 3.17 in the present round 99.9% of irradiated samples fell into the intermediate or positive bands, compared with 99.7% in round 1. The proportions of unirradiated samples giving positive results has fallen very slightly to 2.3% from 2.8%. Both of these results represent slight improvements in performance between the two rounds, and again imply that the overall qualitative outcomes are relatively insensitive to the quantitative differences observed in the previous sections. For additional detail it is necessary to examine individual laboratory and sample performance.

Table 3.14 indicates that 2 laboratories (18 and 29) have still not overcome the cross-contamination problems they experienced in round 1 despite the randomised order of presentation, with a large number of positive results from the unirradiated samples. Laboratory 8 has greatly improved performance.

From Table 3.15 it is clear that there are some products (SP8514, SP8516, SP8536, SP8574 and SP8577) where the majority of observations for the unirradiated material fall into the intermediate band. This was also the case for these products in round 1, suggesting that the result is product-associated. There are also samples (SP8525, SP8531) where the negative/intermediate split is approximately even; again this repeats the round 1 observations. Two other products (SP8576, SP8578) have a significant number of intermediate observations, and this was also the case in round 1. Overall, however, 75% of observations from unirradiated samples are negative and 98% are either negative or intermediate. For the irradiated materials, four products (SP8530, SP8538, SP8573 and SP8575) are wholly or predominantly intermediate rather than positive. This was the same in round 1 for these products, reflecting their relatively low sensitivity to radiation. Overall, 99.9% of irradiated samples were identified as positive or intermediate.

For the blends it is not possible to compare with round 1. One of the 10% blends (SP8531) produced 95% positive outcomes, but the remainder were negative not intermediate. The other 10% blend (SP8539) produced 88% intermediate results. Neither of these products appears to be problematic in pure irradiated form, implying that the results for the blends are associated with the presence of a mixture of irradiated and unirradiated material. For the 1% blends, SP8533 produced 62% intermediate and 38% negative results and SP8577, which had elevated screening results in its unirradiated form, was 90% positive. With the concentration

of irradiated material at 0.1%, SP8535 (which again had almost as many intermediate results as negative for the unirradiated samples) has more intermediate than negative outcomes, whereas SP8578 is 87% negative. Table 3.18 examines the extent to which the outcomes can be related to the concentrations and brightnesses of the irradiated component in these samples. Once the expected signal levels, taken from table 2.4 are considered, the differences between the two examples in the 10% and 1% concentration blends become easier to understand. Thus the majority of results from the high sensitivity ginger sample at 10% concentration fall into the positive band, whereas the same concentration in tarragon leaves many results in the intermediate category. At 1% concentrations the high sensitivity ginseng produces many positive determinations, unlike the less sensitive rosemary. When all pairs are considered together however there is a noticeable trend in detection rate with increasing concentration. At 0.1% concentration 37.5% of blends would be selected for further investigation, on the basis of intermediate or positive screening, rising to 80.8% at 1% concentration and 95.8% at 10% concentration. On this basis it seems that both sensitivity and concentration are relevant in determining screening outcomes from blends. These blends mixed irradiated and unirradiated portions of the same product. In commercial applications where the brightness of minor irradiated components may differ from that of the matrix the situation may be more complex. Such scenarios were however studied under an earlier FSA/MAFF investigation of blending at SURRC.

Table 3.16 presents the homogeneity testing data for the same blends (10 aliquots of each). These show a similar pattern to the participants' results, with some evidence of a correlation between "unexpected" outcome and the coefficient of variance of the 10 results. This supports the conclusion that sample heterogeneity is a major factor.

Table 3.17 compares the participants' results with the round 1 data (from a larger number of samples measured as single aliquots) and with round 1 reference data. Again as in round 1, the correspondence between participants' and reference performance suggests that the slight quantitative underestimation implied by z-scores makes little difference in practice to qualitative outcomes at least for these test materials.

The unirradiated test materials still show a small proportion (c. 2%) of false positive band results in participant data sets that might be attributable to test material handling problems in certain laboratories (indeed the same laboratories in both rounds) and test materials. The main detrimental effect of elevated readings in unirradiated samples would be to select a greater proportion for other investigations than necessary which clearly represents a cost but does not necessarily lead to errors.

It should be noted that the qualitative outcomes tabulated below are based comparison between the reported terminal counts and the thresholds, as in round 1. A major difference in round 2 was the use of duplicate aliquots, with potential for conflict between classifications. Participants were asked to evaluate the outcomes of their screenings, but this was not done consistently. Those participants who did evaluate used a variety of approaches both to the evaluation and to dealing with different outcomes between aliquots; some returned a positive, intermediate or negative classification, some returned a decision as to whether the material was irradiated, some regarded samples with conflicting results for the two aliquots as inconclusive, some always chose the higher category, others chose a mixture of higher and lower, particularly for the subset of laboratories who used calibrated PSL to inform their decisions. The tables below treat the two aliquots from each sample as separate outcomes because of this variation in participant response.

| Lab | Unirradiated | | | Irradiated | | | Blends | | |
|------------|--------------|--------|-------|------------|-------|--------|--------|--------|--------|
| | N | I | P | N | I | P | N | I | P |
| 1 | 53 | 13 | 0 | 0 | 7 | 59 | 5 | 3 | 4 |
| 2 | 47 | 19 | 0 | 0 | 7 | 59 | 5 | 3 | 4 |
| 3 | 51 | 15 | 0 | 0 | 5 | 61 | 3 | 4 | 5 |
| 5 | 51 | 13 | 2 | 0 | 4 | 62 | 2 | 6 | 4 |
| 6 | 47 | 19 | 0 | 0 | 3 | 63 | 6 | 3 | 3 |
| 8 | 49 | 16 | 1 | 0 | 5 | 61 | 2 | 6 | 4 |
| 9 | 52 | 13 | 1 | 0 | 5 | 61 | 5 | 3 | 4 |
| 11 | 50 | 15 | 1 | 0 | 3 | 63 | 2 | 6 | 4 |
| 12 | 50 | 14 | 2 | 0 | 5 | 61 | 3 | 5 | 4 |
| 13 | 49 | 17 | 0 | 0 | 3 | 63 | 3 | 5 | 4 |
| 14 | 50 | 16 | 0 | 0 | 6 | 60 | 4 | 3 | 5 |
| 16 | 52 | 14 | 0 | 0 | 6 | 60 | 4 | 4 | 4 |
| 17 | 48 | 18 | 0 | 0 | 6 | 60 | 4 | 4 | 4 |
| 18 | 30 | 27 | 9 | 0 | 5 | 61 | 0 | 7 | 5 |
| 19 | 47 | 17 | 2 | 0 | 5 | 61 | 3 | 5 | 4 |
| 20 | 56 | 8 | 2 | 0 | 8 | 58 | 6 | 2 | 4 |
| 21 | 52 | 14 | 0 | 0 | 7 | 59 | 5 | 3 | 4 |
| 22 | 41 | 21 | 4 | 0 | 6 | 60 | 3 | 6 | 3 |
| 23 | 45 | 18 | 3 | 0 | 5 | 61 | 4 | 6 | 2 |
| 24 | 49 | 12 | 5 | 0 | 2 | 64 | 0 | 8 | 4 |
| 25 | 56 | 10 | 0 | 0 | 8 | 58 | 2 | 6 | 4 |
| 26 | 54 | 12 | 0 | 0 | 8 | 58 | 4 | 4 | 4 |
| 27 | 56 | 10 | 0 | 0 | 10 | 56 | 4 | 6 | 2 |
| 28 | 49 | 17 | 0 | 0 | 6 | 60 | 2 | 6 | 4 |
| 29 | 38 | 17 | 11 | 0 | 5 | 61 | 1 | 7 | 4 |
| 30 | 56 | 10 | 0 | 0 | 7 | 59 | 5 | 3 | 4 |
| 31 | 50 | 16 | 0 | 0 | 10 | 56 | 4 | 6 | 2 |
| 33 | 47 | 17 | 2 | 0 | 8 | 58 | 4 | 4 | 4 |
| 34 | 56 | 10 | 0 | 2 | 8 | 56 | 6 | 2 | 4 |
| 35 | 47 | 19 | 0 | 0 | 6 | 60 | 2 | 5 | 5 |
| Total | 1478 | 457 | 45 | 2 | 179 | 1799 | 103 | 141 | 116 |
| Percentage | 74.65% | 23.08% | 2.27% | 0.10% | 9.40% | 89.24% | 28.61% | 39.17% | 32.22% |

Table 3.14 Participants' qualitative results by laboratory

| | Unirradiated | | | Irradiated | | | Blends | | | |
|-----------|--------------|-------|------|------------|------|-------|-------------|-------|-------|-------|
| SP Number | N | I | P | N | I | P | SP Number | N | I | P |
| 8512 | 56 | 1 | 3 | 0 | 0 | 60 | 8531(10%) | 3 | 0 | 57 |
| 8513 | 52 | 7 | 1 | 0 | 0 | 60 | 8533 (1%) | 22 | 38 | 0 |
| 8514 | 1 | 59 | 0 | 0 | 0 | 60 | 8535 (0.1%) | 23 | 37 | 0 |
| 8515 | 54 | 6 | 0 | 0 | 2 | 58 | 8539 (10%) | 2 | 53 | 5 |
| 8516 | 4 | 56 | 0 | 0 | 0 | 60 | 8577 (1%) | 1 | 5 | 54 |
| 8517 | 54 | 4 | 2 | 0 | 0 | 60 | 8578 (0.1%) | 52 | 8 | 0 |
| 8518 | 56 | 4 | 0 | 0 | 0 | 60 | Total | 103 | 141 | 116 |
| 8519 | 54 | 6 | 0 | 0 | 0 | 60 | Percent | 28.61 | 39.17 | 32.22 |
| 8520 | 51 | 7 | 2 | 0 | 0 | 60 | | | | |
| 8521 | 60 | 0 | 0 | 0 | 0 | 60 | | | | |
| 8523 | 60 | 0 | 0 | 0 | 0 | 60 | | | | |
| 8525 | 29 | 31 | 0 | 0 | 0 | 60 | | | | |
| 8526 | 55 | 4 | 1 | 0 | 0 | 60 | | | | |
| 8528 | 57 | 3 | 0 | 0 | 0 | 60 | | | | |
| 8529 | 57 | 3 | 0 | 0 | 0 | 60 | | | | |
| 8530 | 58 | 2 | 0 | 2 | 58 | 0 | | | | |
| 8531 | 27 | 30 | 3 | 0 | 0 | 60 | | | | |
| 8532 | 40 | 14 | 6 | 0 | 0 | 60 | | | | |
| 8533 | 60 | 0 | 0 | 0 | 0 | 60 | | | | |
| 8534 | 59 | 1 | 0 | 0 | 4 | 60 | | | | |
| 8535 | 33 | 27 | 0 | 0 | 0 | 60 | | | | |
| 8536 | 3 | 48 | 9 | 0 | 0 | 60 | | | | |
| 8538 | 57 | 3 | 0 | 0 | 40 | 20 | | | | |
| 8539 | 60 | 0 | 0 | 0 | 0 | 60 | | | | |
| 8571 | 56 | 2 | 2 | 0 | 0 | 60 | | | | |
| 8572 | 59 | 1 | 0 | 0 | 0 | 60 | | | | |
| 8573 | 58 | 2 | 0 | 0 | 44 | 16 | | | | |
| 8574 | 2 | 57 | 1 | 0 | 0 | 60 | | | | |
| 8575 | 54 | 4 | 2 | 0 | 31 | 29 | | | | |
| 8576 | 49 | 11 | 0 | 0 | 0 | 60 | | | | |
| 8577 | 1 | 49 | 10 | 0 | 0 | 60 | | | | |
| 8578 | 47 | 10 | 3 | 0 | 0 | 60 | | | | |
| 8579 | 55 | 5 | 0 | 0 | 0 | 60 | | | | |
| Total | 1478 | 457 | 45 | 2 | 179 | 1799 | | | | |
| Percent | 74.65 | 23.08 | 2.27 | 0.10 | 9.04 | 90.86 | | | | |

Table 3.15 Participants' qualitative results by product

| SP Number | Concentration | N | I | P | CV |
|-----------|---------------|-------|-------|-------|-----|
| 8531 | 10% | 0 | 0 | 10 | 6% |
| 8533 | 1% | 7 | 3 | 0 | 9% |
| 8535 | 0.1% | 6 | 4 | 0 | 10% |
| 8539 | 10% | 0 | 8 | 2 | 8% |
| 8577 | 1% | 0 | 0 | 10 | 4% |
| 8578 | 0.1% | 10 | 0 | 0 | 4% |
| Total | | 23 | 15 | 22 | |
| Percent | | 38.33 | 25.00 | 36.67 | |

Table 3.16 Homogeneity testing for blends (screening) with coefficient of variance

| | Irradiated Test Materials | | Unirradiated Test Materials | | Blended Test Materials | |
|---------------|---------------------------|-----------|-----------------------------|-----|------------------------|-----------|
| | Neg | Pos / Int | Neg / Int | Pos | Neg | Pos / Int |
| Participants | | | | | | |
| Round 2 | | | | | | |
| Total | 2 | 1978 | 1935 | 45 | 103 | 257 |
| % | 0.1 | 99.9 | 97.7 | 2.3 | 28.6 | 71.4 |
| | | | | | | |
| Round 1 | 5 | 1514 | 1477 | 42 | | |
| Total | | | | | | |
| % | 0.3 | 99.7 | 97.2 | 2.8 | | |
| | | | | | | |
| Reference Set | | | | | | |
| % | 0 | 490 | 489 | 1 | 23 | 37 |
| Percent | 0 | 100 | 99.8 | 0.2 | 38.3 | 61.7 |

Table 3.17 Comparison of participants' qualitative percentage with the reference set

| Concentration | Sample | Product | Expected PSL ¹ | Negative | Intermediate | Positive | Selected for investigation ² |
|---------------|--------|-----------|---------------------------|----------|--------------|----------|---|
| 10% | SP8531 | Ginger | 36467 | 3 | 0 | 57 | 57 |
| | SP8539 | Tarragon | 4696 | 2 | 53 | 5 | 58 |
| | | | Total | 5 | 53 | 62 | 115 |
| | | | % | 4.2 | 44.2 | 51.7 | 95.8 |
| 1% | SP8533 | Rosemary | 764 | 22 | 38 | 0 | 38 |
| | SP8577 | Ginseng | 11342 | 1 | 5 | 54 | 59 |
| | | | Total | 23 | 43 | 54 | 97 |
| | | | % | 19.2 | 35.8 | 45.0 | 80.8 |
| 0.10% | SP8535 | Chilli | 145 | 23 | 37 | 0 | 37 |
| | SP8578 | Green Tea | 116 | 52 | 8 | 0 | 8 |
| | | | | 75 | 45 | 0 | 45 |
| | | | | 62.5 | 37.5 | 0 | 37.5 |
| All blends | | | Total | 103 | 141 | 116 | 257 |
| | | | % | 28.6 | 39.2 | 32.2 | 71.4 |

Table 3.18 PSL screening results from the blended samples by concentration and brightness. 1. Expected PSL signals in photon counts from table 2.4. 2. The numbers which would be selected for further investigation under EN13751, ie the sum of intermediate and positive outcomes.

4. CALIBRATED PSL

4.1 Calibrated PSL Homogeneity Testing

Round 2 introduced calibrated PSL measurements for those laboratories able to irradiate their material. 12 of the screening laboratories participated in this part of the trial. To provide reference data for comparison with these extra results, calibrated PSL homogeneity testing was carried out at SUERC in July and August 2006.

All the petri dishes from the 2005 homogeneity testing, including those products not sent to participants for round 2, plus those used for the blends, were given a dose of 1kGy at Isotron and re-read. This provided reference data comparable with participants' calibrated measurements and with the homogeneity testing for the screening.

4.2 Homogeneity Testing Results

The PSL terminal counts for the calibrated measurements are fully tabulated in Appendix B. Figure 4. 1 presents the data as a scatter plot with initial PSL against calibrated PSL; the data show good separation of irradiated and unirradiated products, with the blends overlapping both categories as expected.




Figure 4.1 Initial vs calibrated PSL for reference materials

4.3 Participants' Results

Participants were asked to return their classification data, PSL terminal counts, summary and raw data files, as for the screening. Some participants indicated that they had re-dispensed material to be irradiated, rather than sealing the petri-dishes used for screening and re-measuring the same material after irradiation. For some laboratories it is not clear which procedure was followed (the calibrated homogeneity testing was performed on the screened material) and neither is it clear whether the thick or thin methods of preparation were used. For each of the 12 laboratories which returned calibrated data, a plot equivalent to Figure 4.1 was produced. These are displayed below as Figures 4.2 – 4.5. Full data are in Appendix C.

Table 4.1 shows the type of source used by participants, where this was stated in their returns. It also shows whether they used the information obtained from the second measurement to inform their judgment about the samples.

| Laboratory | Type of source | Whether used calibrated results to evaluate | |
|------------|----------------------------|---|---|
| | | Y | N |
| 1 | ^{60}Co | Y | |
| 14 | X-ray 60kV 0.2mm Al filter | Y | |
| 17 | 10 MeV LINAC | Y | |
| 18 | ^{60}Co | | N |
| 19 | ^{60}Co | Y | |
| 23 | ^{60}Co | Y | |
| 24 | ^{60}Co | | N |
| 25 | ^{60}Co | | N |
| 26 | ^{60}Co | | N |
| 27 | Not stated | | N |
| 28 | Not stated | Y | |
| 35 | ^{60}Co | Y | |

Table 4.1 Type of source used for calibration by participants and whether they utilised the second data set to evaluate

From figures 4.2 to 4.5 it can be seen that the overall scatter distribution for the participants is very similar to that for the reference data. The small cluster of results with high initial PSL (salt-containing materials) has been distinguished by almost all the participants, and most have the irradiated samples clustering along a gradient. Laboratory 18 is an exception to this, with a much more diffuse pattern; this laboratory has already been seen to have a sample-handling issue. Laboratory 14 also shows a slightly more random distribution for the irradiated samples, and laboratory 17 appears to have the least dispersed data set from irradiated samples of all laboratories, including the reference laboratory. When the source types used are considered in combination with these data it appears that the quality of the calibrated relationship may depend on the irradiation facility used. For laboratories utilising ^{60}Co sources the data appear to be very similar to the reference set, which was re-irradiated remotely using the Isotron commercial ^{60}Co facility. Laboratory 17 used an in-house LINAC, which apparently results in a very well controlled dose distribution from sample to sample, reflecting in high quality calibrated PSL data. The x-ray system used by laboratory 14 appears to have generated additional scatter. It would be useful to assess the extent to which these findings are generic to their respective technologies in future work. It is also noted that the evaluation methods differed between participants, and that for these samples there were rather few cases where screening results from low-sensitivity irradiated samples might have been re-interpreted following the calibrated response.

The use of z-scores was not considered appropriate here. Each participant was responsible for the necessary irradiation; dose and the time elapsed between irradiation and read-out have not been assessed. Figure 3.9 above implies that time elapsed will not be important within the time-scale of the study. Since it is unclear to what extent the availability of calibrated results affected participants' evaluation of outcome (some indicated that they had taken it into account; in some other cases this is implicit), qualitative tables have not been prepared at this stage. Nevertheless, it is clear that overall performance was good, irradiations having been performed and repeat measurements conducted successfully. Further analyses of these data may prove informative in the future.






Figure 4.2 Initial vs calibrated PSL for laboratories 1, 14 and 17






Figure 4.3 Initial vs calibrated PSL for laboratories 18, 19 and 23






Figure 4.4 Initial vs calibrated PSL for laboratories 24, 25 and 26






Figure 4.5 Initial vs calibrated PSL for laboratories 27, 28 and 35

5. THERMOLUMINESCENCE ANALYSIS

5.1 TL Homogeneity Testing

Homogeneity testing was conducted at SUERC on randomly selected excess material from that prepared for distribution to participants, to characterise the TL response and assess the intrinsic variability of all samples within one laboratory. For irradiated and unirradiated products sent to participants, 2 discs were prepared from each of 10 pots; 20 discs were prepared from each of the blended materials. (Previously, 2 discs of each of the unirradiated products as originally purchased had been analysed to determine whether there was any evidence that they contained irradiated components.) To ensure that there was no cross-contamination between irradiated, unirradiated and blended samples, the unirradiated products were separated and read out first, then the irradiated and finally the blended products; the procedure was carried out with the usual quality assurance steps using process and glassware blanks.

5.2 TL Homogeneity Testing Results

EN1788 calls for TL quantification of glow 1 signals, glow 2 signals, the glow ratio and identification of whether or not a peak can be observed in the 150-250°C region. Where glow 2 intensities are less than 10 times the laboratory minimum detectable level, defined by the mean and 3 standard deviations of full-process blanks the mineral yield is considered unsatisfactory. An additional criterion for use in UK surveys, which has yet to be added to EN1788, was introduced in 2001 which requires any sample where the presence of irradiated material has been identified on the basis of a low temperature peak in glow 1 to have intensity which exceeds the minimum detectable level.

The reference laboratory routinely checks temperature calibration using encapsulated natural CaF_2 grains and monitoring the positions of 3 well characterised TL peaks. This material presents the grains in a manner which more closely resembles the thermal contact of the minerals extracted from food samples than the 0.8 mm LiF chips suggested in the EN1788 annexe, and does not require the complex annealing cycle of LiF. In this study however LiF chips were prepared and distributed, since participants would be more familiar with this approach. The reference laboratory obtained peak V temperature estimates of 249, 250, 250, and 251 °C from 4 retained LiF chips from the same batch as those distributed in pairs to participants. An additional 6 chips drawn from a new batch of LiF obtained from Thermo-electron for possible use in later trial rounds was also measured at SUERC, giving values of 250, 251, 249, 248, 251 and 249 °C, in July 2006. It was concluded that the temperature control of the SURRC reader used gave a mean value of $250 \pm 1^\circ\text{C}$ for peak V during this phase of the study. In earlier interlaboratory trials, between 1991 and 1995, the SURRC reader had given results between 245°C and 258°C with paired reproducibility of $<5^\circ\text{C}$, which compared with other laboratories peak V temperature range estimates from 190-328°C, with paired reproducibility of up to 30°C ²³. Given the importance of identifying the presence of peaks in the 150-250°C region, particularly for blended mixtures, temperature control accuracy and precision remains of concern to TL analysis.

23. Sanderson D.C.W, Carmichael L.A., Naylor J.D., 1996, Establishing Luminescence Detection Methods for Irradiated Fruits, Vegetables and Shellfish, MAFF Project 1B073, Third Report, SURRC, East Kilbride.p8-12



A further distinction which has been used in the reference laboratory to help to describe the nature of low temperature peak evidence is to note the general shape of the glow 1 peak, with respect to signals in the 150-250°C region. In these terms a “dominant” peak is the main signal in glow 1, a “major” peak describes the largest signal, although other similar signals may also occur at higher temperatures due to geological residuals, a “minor” peak means that there are larger higher temperature signals in the glow curve than the low temperature peak, and a “shoulder” is a small inflected signal component. If no signal is present the term “absent” is used. These descriptions have been used in UK surveys, although they are not in general international use. The peak shape descriptors of the reference analysis are discussed below for information.

All the irradiated materials showed “dominant” G1 peaks in the 150-250°C range. For the unirradiated materials with the exception of SP8539 Tarragon, peaks in this range were “absent”. In the case of SP8539 a low temperature signal in glow 1 was noted in 11 out of the 20 aliquots, which implies, according to EN1788, the presence of irradiated material before treatment. For the blends, results were more mixed, as might be expected. The 10% blends (SP8531 Ginger and SP8539) and one of the 1% blends (SP8533 Rosemary) showed “dominant” G1 peaks in all aliquots. The other blends (SP8577 Ginseng at 1% and 0.1% blends SP8535 Chilli powder and SP8578 Green tea) produced a range of G1 glow shapes. The distribution of these shapes for 20 aliquots is shown in Figure 5.1. As can be seen the 1% Ginseng sample shows both “major” and “minor” peaks in equal proportions. At 0.1% the Green Tea samples shows a significant proportion of “absent” outcomes; more than 50% of analyses have produced “shoulders” with a number of aliquots giving “minor” signals. The other 0.1% concentration blend produces a mixture of “shoulder”, “minor”, “major” and “dominant” in most cases.



Figures 5.2-5.4 are scatter plots of glow 1 against glow 2 log intensities (integrated for 220°-240°C) for each product, with the three irradiation categories indicated by different symbols. These plots show complete separation of irradiated and unirradiated materials. For the blends, the degree of separation varies. The 10% and 1% blends the blended results occupy a discrete area, whereas for the 0.1% the blended and unirradiated samples overlap. This is consistent with the glow shapes described above, but underlines the importance of utilising glow shape information at low concentrations. Figures 5.5-5.7 plot glow ratio for the same temperature interval, again showing very good discrimination between irradiated and unirradiated samples, with concentration dependence in the glow ratios from blends. Table B.7 in Appendix B summarises the reference values of glow ratio.

Figure 5.1 Distribution of glow 1 peak shapes for 20 aliquots each of 3 blends



SP8535 Chilli Powder 0.1% Blend



SP8578 Green Tea 0.1% Blend



SP8577 Ginseng 1% Blend







Figure 5.2 TL homogeneity testing glow 1 vs glow 2 scatter plot for SP8531 and SP8533





Figure 5.3 TL homogeneity testing glow 1 vs glow 2 scatter plot for SP8535 and SP8539




Figure 5.4 TL homogeneity testing glow 1 vs glow 2 scatter plot for SP8577 and SP8578





Figure 5.5 TL homogeneity testing glow ratio plot for SP8531 and SP8533




Figure 5.6 TL homogeneity testing glow ratio plot for SP8535 and SP8539




Figure 5.7 TL homogeneity testing glow ratio plot for SP8577 and SP8578

5.3 Participants' Results

5.3.1 Quantitative results

TL results were returned by 16 of the 18 laboratories which received samples. TL participants were asked to return LiF temperature control data, first and second glow data, glow ratio, comments on glow curve shape and classification.

Table 5.1 presents participants' LiF data as a mean for peak V for each laboratory, and the minimum detectable level (MDL). All raw data appear in Appendix C. The units for MDL depend upon the equipment, which is also tabulated. The value for lab 5 LiF appears to be an error; it is suspected that the wrong peak was identified since this magnitude of temperature discrepancy could not have allowed them to perform as well as they have. The range of peak V estimates from the other laboratories varies from 220-293 °C, which although a slight reduction from that encountered in some of the early trials in the 1990's still represents a very significant variation from laboratory to laboratory. This is of potential concern when dealing with blends classified by the presence or absence of peaks in the EN1788 temperature band of 150-250°C.

| Labs | TL Reader | Mean Peak V | MDL |
|------|--|-------------|----------|
| 1 | TLD 3500 | 234 | 0,37 nC |
| 2 | Harshaw TLD 4000 | 275 | 7.5nC |
| 3 | TLD 3500 | 234 | 3,12nC |
| 4 | Harshaw TLD-Reader 3500 | 245 | 20 nC |
| 5 | Harshaw 3500 | 90.5 | 0.047nC |
| 6 | HARSHAW TLD 3500 | 235 | 0,098 nC |
| 7 | Harshaw M3500 | 255 | 0.13nC |
| 8 | Risö TL/OSL-DA15 | 254 | 502 |
| 10 | Harshaw 3500 | 246 | 0.54 |
| 11 | TL-DA-15; RISO | 257 | 358.49 |
| 13 | Risoe TL-DA-15 Series No 194-05/2005-b | 263 | 2492 |
| 14 | Harshaw 3500 | 220 | 0.354nC |
| 15 | TL-DA-10 | 293 | 636 |
| 18 | Harshow TL 3500 | 232 | 2.56nC |
| 24 | Harshaw TLD 3500 | 228 | 1.9nC |

Table 5.1 LiF peak V temperatures and MDL for all laboratories

Figures 5.8-5.10 present glow ratios for each laboratory and each product, with the 3 irradiation categories given different symbols. Intensities could not be plotted analogously to Figures 5.2-4 because participants used a variety of units for their measurements.

Figures 5.11-5.13 display *z* scores by laboratory; figures 5.14-16 display them by sample. The data are tabulated in Table 5.2. Contour plots have also been produced (Figures 5.17-19).

From the figures it can be seen that participants' data distinguish the irradiated samples but are markedly less well separated by glow ratio than the homogeneity testing data for the 10% and 1% blends. For the 0.1% blends where the unirradiated and blended results overlap in the homogeneity testing, the same is true for the participants but with more scatter. Some of the scatter may be explained by the use of different instruments with several different types of equipment with different irradiation sources, but there is also likely to be a component

associated with efficiency of mineral separation. In some instances the samples could have been mislabelled leading to a lower glow ratio for the irradiated material than for either the unirradiated or the blend. This is also apparent in the z -score plots which have corresponding outliers. Even without these outliers, the z -scores show a wider distribution than the equivalent for PSL. This is not unexpected since the distribution of a ratio is likely to follow a Cauchy distribution rather than a Gaussian, and thus to have broader tails.


Figures 5.11-13 show clearly that some laboratories have more closely spaced z -scores than others; for unirradiated samples labs 1, 13 and 18 had the smallest spread. These 3 laboratories also performed well with irradiated samples, if an outlier for lab 1 is discounted, but labs 6, 8, 11, and 14 also did well. With the blends, only lab 18 has a tight dispersion of z -score.

Plots of z -score by sample show that samples 6 and 9 (unirradiated and irradiated Rosemary), 10, 11 and 12 (all 3 Tarragon samples) and 7 (0.1% Chilli blend) produced tighter groupings. This suggests that higher mineral yields from some products may affect results. There is no overall correlation between blend concentration and spread of z -score.

The contour plots reveal that SP8533 Rosemary consistently produces z -scores which differ from zero, despite the tighter grouping. These plots also show that labs 2, 3 and 4 have higher scores for unirradiated materials; higher scores for irradiated materials are seen for labs 3, 5, and 7 with slightly elevated scores for lab 11 but lower scores for lab 14, and labs 2 and 4 also have higher scores for the blends, which are otherwise more varied.

Given the limit amount of sample available (some participants were unable to obtain sufficient minerals from some sample to be able to return a result) the overall performance is quite good, although there are some issues still to be addressed.

SP8531 Ginger



SP8533 Rosemary





Figure 5.8 Participants' glow ratios for SP8531 and SP8533

SP8535 Chilli



SP8539 Tarragon





Figure 5.9 Participants' glow ratios for SP8535 and SP8539

SP8577 Ginseng



SP8578 Green Tea




Figure 5.10 Participants' glow ratios for SP8577 and SP8578

Unirradiated

| | Sample 1 Ginger | | Sample 6 Rosemary | | Sample 8 Chilli | | Sample10 Tarragon | | Sample15 Ginseng | | Sample17 Green Tea | |
|--------|--------------------|-------|----------------------|------|--------------------|-------|----------------------|-------|---------------------|-------|-----------------------|-------|
| Lab 1 | -2.61 | -1.19 | 1.03 | 1.92 | -3.49 | -1.85 | 0.09 | -0.12 | 0.12 | -0.76 | -4.95 | -2.41 |
| Lab 2 | -0.36 | -2.61 | 5.51 | 5.51 | 11.81 | 11.81 | 4.19 | 4.19 | 21.03 | 21.03 | 20.34 | 20.34 |
| Lab 3 | 0.69 | 7.37 | 3.96 | 5.51 | 4.96 | 5.93 | 3.32 | 1.73 | 5.43 | 7.15 | 2.54 | 0.37 |
| Lab 4 | 2.09 | 11.56 | 4.95 | 5.51 | 11.81 | 11.81 | 4.45 | 4.19 | 9.61 | 21.03 | 8.48 | 20.34 |
| Lab 5 | -0.36 | -0.36 | 4.07 | 5.51 | -4.74 | 11.81 | 2.46 | 4.19 | 2.44 | 5.53 | 1.61 | 20.34 |
| Lab 6 | -1.19 | -1.19 | 5.51 | 5.51 | -1.35 | -1.42 | -0.05 | 4.19 | -0.71 | 2.23 | -4.45 | -1.83 |
| Lab 7 | 5.94 | 5.35 | 2.53 | 3.68 | 5.75 | 5.13 | 1.68 | 0.82 | 0.21 | -0.06 | 0.86 | 0.66 |
| Lab 8 | 0.23 | 2.11 | 3.08 | 5.51 | -1.52 | -1.85 | -2.65 | -1.45 | 2.44 | 1.45 | 0.12 | -0.93 |
| Lab 10 | -0.90 | -1.65 | 1.05 | 5.51 | -0.43 | -1.31 | 1.20 | -0.62 | -1.14 | -1.70 | -8.49 | -4.17 |
| Lab 11 | -0.36 | 0.69 | 0.56 | 5.51 | 1.09 | -0.28 | -1.85 | -0.54 | 2.63 | 4.35 | 0.38 | 0.74 |
| Lab 13 | -0.36 | 0.69 | 2.58 | 1.03 | 0.64 | 1.37 | -0.42 | 0.09 | 1.97 | 1.45 | 0.12 | 0.12 |
| Lab 14 | -2.05 | 0.91 | 2.22 | 5.51 | -2.11 | -3.41 | 0.87 | 0.30 | -6.77 | -8.94 | -5.60 | -5.58 |
| Lab 15 | -5.01 | 11.56 | 1.78 | 5.51 | -4.70 | -3.66 | 1.13 | 4.19 | -5.41 | -4.79 | -6.57 | -9.21 |
| Lab 18 | 4.50 | 5.01 | 1.59 | 1.59 | 2.62 | 2.29 | 0.97 | 0.60 | 6.14 | 5.78 | 5.20 | 4.96 |
| Lab 24 | 1.38 | -0.36 | 0.90 | 5.51 | 2.16 | 1.34 | 1.38 | 0.60 | 2.81 | -0.08 | 1.32 | 4.24 |

Table 5.2 z-score for TL analyses – unirradiated samples. Laboratory numbers are those allocated for TL apart from 18 and 24 which are PSL identifiers (these labs joined later)

Irradiated

| | Sample2 Ginger | | Sample 4 Rosemary | | Sample 9 Chilli | | Sample11 Tarragon | | Sample13 Ginseng | | Sample18 Green Tea | |
|--------|-------------------|-------|----------------------|--------|--------------------|--------|----------------------|--------|---------------------|-------|-----------------------|-------|
| Lab 1 | -4.29 | -3.44 | -5.41 | -22.02 | -2.11 | -4.16 | -0.64 | -4.67 | -1.62 | -7.50 | -2.58 | -5.71 |
| Lab 2 | -1.59 | -1.78 | -5.66 | -5.66 | -3.39 | -3.39 | -3.67 | -3.67 | -5.42 | -5.42 | -6.60 | -6.60 |
| Lab 3 | 6.27 | 2.94 | -0.20 | -5.66 | 1.34 | 4.15 | -1.48 | -0.23 | 0.95 | 2.36 | 7.47 | 8.24 |
| Lab 4 | -6.70 | -5.24 | -4.47 | -5.66 | -0.13 | -3.99 | -0.58 | -3.67 | -0.97 | -5.42 | -54.43 | -6.60 |
| Lab 5 | -13.15 | -5.24 | 2.17 | -5.66 | 1.38 | 2.59 | 0.88 | -3.67 | 4.56 | 7.50 | 4.33 | -6.60 |
| Lab 6 | -2.38 | -1.14 | -5.66 | -5.66 | -1.40 | -3.39 | -3.61 | -3.67 | -3.46 | -3.10 | -2.94 | -1.16 |
| Lab 7 | 8.95 | 5.12 | -2.73 | -4.48 | 5.14 | 2.33 | -1.86 | -2.15 | 7.23 | 2.85 | 13.53 | 9.22 |
| Lab 8 | -3.47 | -2.96 | -4.94 | -5.66 | -2.75 | -2.34 | -2.20 | -2.71 | -4.56 | -5.27 | -3.53 | -3.34 |
| Lab 10 | 4.19 | -0.71 | -19.29 | -5.66 | -3.38 | -5.39 | -3.94 | -15.51 | -2.39 | -2.60 | -6.93 | -6.64 |
| Lab 11 | -2.42 | -1.73 | -0.75 | -5.66 | 1.98 | 0.90 | 0.59 | -0.62 | -1.90 | -2.91 | -0.96 | -2.18 |
| Lab 13 | -0.98 | -1.20 | 0.77 | -0.01 | -0.83 | -0.70 | -3.00 | -0.54 | -2.09 | -2.31 | -1.86 | -2.25 |
| Lab 14 | -3.74 | -4.63 | -6.55 | -5.66 | -3.75 | -3.86 | -2.31 | -6.63 | -2.78 | -3.71 | -5.03 | -5.17 |
| Lab 15 | 8.08 | -2.54 | -1.25 | -5.66 | -6.92 | -6.19 | -3.68 | -3.67 | -6.70 | -8.25 | -7.13 | -9.67 |
| Lab 18 | -1.98 | -2.29 | 0.01 | 0.42 | -1.31 | -1.25 | 0.43 | 0.21 | 2.48 | 1.22 | -1.46 | -1.57 |
| Lab 24 | -3.28 | -4.34 | -45.52 | -5.66 | -0.08 | -12.29 | -2.59 | -2.81 | 15.73 | -2.58 | -4.93 | -4.58 |

Table 5.3 z-score for TL analyses – irradiated samples. Laboratory numbers are those allocated for TL apart from 18 and 24 which are PSL identifiers (these labs joined later)

Blends

| | Sample 3 Ginger | | Sample 5 Rosemary | | Sample 7 Chilli | | Sample12 Tarragon | | Sample14 Ginseng | | Sample16 Green Tea | |
|--------|--------------------|-------|----------------------|-------|--------------------|-------|----------------------|-------|---------------------|--------|-----------------------|-------|
| Lab 1 | -1.48 | -5.32 | -0.56 | 0.19 | 0.61 | -1.04 | -1.92 | -2.24 | -2.11 | -6.03 | -1.99 | -0.54 |
| Lab 2 | -0.57 | -1.95 | 5.01 | 5.01 | 3.30 | 3.30 | 1.44 | 1.44 | 13.94 | 13.94 | 7.13 | 7.13 |
| Lab 3 | 4.38 | 0.90 | 1.50 | 5.01 | 1.66 | 0.61 | 1.43 | 0.29 | -1.42 | -4.03 | 3.00 | 1.18 |
| Lab 4 | -0.70 | 1.85 | 5.01 | 5.01 | 0.98 | 3.30 | 2.72 | 1.44 | -0.56 | 13.94 | 4.58 | 7.13 |
| Lab 5 | -0.49 | 0.40 | 2.77 | 5.01 | -0.03 | 1.25 | 1.93 | 1.44 | 1.74 | 2.04 | -0.16 | 7.13 |
| Lab 6 | 0.64 | -1.26 | 5.01 | 5.01 | -1.00 | 3.30 | -1.23 | 1.44 | -6.03 | -11.06 | -0.54 | -1.26 |
| Lab 7 | 4.16 | 5.09 | 1.11 | 1.56 | 2.11 | 1.20 | 0.92 | -0.39 | -5.39 | -1.99 | 2.02 | -0.08 |
| Lab 8 | -0.17 | 1.21 | -4.91 | -5.58 | -1.14 | -1.44 | -0.70 | -0.60 | -0.53 | -12.72 | -1.07 | -1.07 |
| Lab 10 | -0.57 | -0.63 | 0.37 | 5.01 | -0.33 | -0.31 | -1.23 | -0.98 | -7.53 | -6.56 | -1.16 | -1.05 |
| Lab 11 | -0.73 | 0.07 | 5.01 | 5.01 | -1.11 | -0.54 | 0.49 | -2.48 | -3.34 | -5.01 | 3.93 | -0.40 |
| Lab 13 | -0.30 | 0.14 | 0.63 | -0.47 | -0.10 | -0.10 | 0.66 | -0.59 | -4.49 | -5.20 | 0.38 | 1.29 |
| Lab 14 | -1.17 | -0.84 | 0.06 | 5.01 | -0.80 | -0.92 | 0.05 | 0.74 | -7.36 | -5.44 | -1.95 | -0.29 |
| Lab 15 | 0.81 | -4.11 | -2.65 | 5.01 | -1.83 | -2.13 | -0.70 | 1.44 | -2.02 | -14.01 | -1.78 | -3.88 |
| Lab 18 | -0.12 | -0.44 | -0.66 | -0.31 | 0.48 | 0.40 | 0.29 | -0.52 | -1.37 | -1.04 | 2.09 | 1.83 |
| Lab 24 | -0.77 | 1.55 | 0.65 | 5.01 | 0.42 | 0.13 | -1.11 | -1.21 | -2.06 | -1.11 | 0.85 | 0.92 |

Table 5.4 z-score for TL analyses – blended samples. Laboratory numbers are those allocated for TL apart from 18 and 24 which are PSL identifiers (these labs joined later)

Round 2
Zscores for Participants data from unirradiated samples




Figure 5.11 Participants' TL z-scores for unirradiated samples by laboratory

Round 2
Zscores for Participants data from irradiated samples




Figure 5.12 Participants' TL z-scores for irradiated samples by laboratory

Round 2
Zscores for Participants data from blended samples




Figure 5.13 Participants' TL z -scores for blended samples by laboratory

Round 2
Zscores for Participants data from unirradiated samples




Figure 5.14 Participants' TL z-scores for unirradiated samples by sample
Sample 1 = Ginger, sample 6 = Rosemary, sample 8 = Chilli, sample 10= Tarragon,
sample 15 = Ginseng, sample 17 = Green tea

Round 2
Zscores for Participants data from irradiated samples




Figure 5.15 Participants' TL z-scores for irradiated samples by sample
Sample 2 = Ginger, sample 4 = Rosemary, sample 9 = Chilli, sample 11= Tarragon,
sample 13 = Ginseng, sample 18 = Green tea

Round 2
Zscores for Participants data from blended samples




Figure 5.16 Participants' TL z -scores for blended samples by sample
Sample 3 = Ginger, sample 5 = Rosemary, sample 7 = Chilli, sample 12= Tarragon,
sample 14 = Ginseng, sample 16 = Green tea

Z-scores for unirradiated TL Samples




Figure 5.17 Contour plot for participants' TL data for unirradiated samples

Z-scores for irradiated TL Samples




Figure 5.18 Contour plot for participants' TL data for irradiated samples

Z-scores for blended TL Samples




Figure 5.19 Contour plot for participants' TL data for blended samples

5.3.2 Qualitative results

Participants were asked to return classifications of the samples as well as descriptions of the glow curves. These can both be used to assess whether they correctly identified the irradiation status of the materials. Table 5.5 details their classifications compared with the known status of the material, both by laboratory and by product.

| Lab | Unirradiated | | | Irradiated | | | Blends | | |
|--------------------------------|--------------|-----------|-------|------------|-----------|-------|---------|-----------|--------|
| | Correct | Incorrect | Other | Correct | Incorrect | Other | Correct | Incorrect | Other |
| 1 | 6 | 0 | 0 | 6 | 0 | 0 | 3 | 3 | 0 |
| 2 | 1 | 0 | 5 | 1 | 0 | 5 | 1 | 0 | 5 |
| 3 | 6 | 0 | 0 | 6 | 0 | 0 | 3 | 3 | 0 |
| 4 | 3 | 3 | 0 | 6 | 0 | 0 | 2 | 3 | 1 |
| 5 | 4 | 1 | 1 | 6 | 0 | 0 | 4 | 2 | 0 |
| 6 | 5 | 0 | 1 | 6 | 0 | 0 | 2 | 2 | 2 |
| 7 | 6 | 0 | 0 | 6 | 0 | 0 | 2 | 2 | 2 |
| 8 | 6 | 0 | 0 | 6 | 0 | 0 | 3 | 3 | 0 |
| 10 | 6 | 0 | 0 | 6 | 0 | 0 | 5 | 0 | 1 |
| 11 | 6 | 0 | 0 | 5 | 1 | 0 | 5 | 1 | 0 |
| 13 | 6 | 0 | 0 | 6 | 0 | 0 | 3 | 3 | 0 |
| 14 | 4 | 2 | 0 | 6 | 0 | 0 | 5 | 1 | 0 |
| 15 | 4 | 0 | 2 | 4 | 0 | 2 | 2 | 2 | 2 |
| 18 | 6 | 0 | 0 | 6 | 0 | 0 | 6 | 0 | 0 |
| 24 | 6 | 0 | 0 | 6 | 0 | 0 | 2 | 4 | 0 |
| Total | 75 | 6 | 9 | 82 | 1 | 7 | 48 | 29 | 13 |
| Percentage | 83.33% | 6.67% | 10% | 91.11% | 1.11% | 7.78% | 53.34% | 32.22% | 14.44% |
| Percentage of results returned | 92.59% | 7.40% | | 98.80% | 1.20% | | 62.34% | 37.66% | |

| Product | Unirradiated | | | Irradiated | | | Blends | | |
|--------------------------------|--------------|-----------|-------|------------|-----------|-------|---------|-----------|--------|
| | Correct | Incorrect | Other | Correct | Incorrect | Other | Correct | Incorrect | Other |
| SP8531 | 14 | 1 | 0 | 14 | 1 | 0 | 15 | 0 | 0 |
| SP8533 | 10 | 1 | 4 | 13 | 0 | 2 | 3 | 8 | 4 |
| SP8535 | 14 | 0 | 1 | 14 | 0 | 1 | 5 | 7 | 3 |
| SP8539 | 11 | 2 | 2 | 13 | 0 | 2 | 12 | 0 | 3 |
| SP8577 | 13 | 1 | 1 | 14 | 0 | 1 | 9 | 4 | 2 |
| SP8578 | 13 | 1 | 1 | 14 | 0 | 1 | 4 | 10 | 1 |
| Total | 75 | 6 | 9 | 82 | 1 | 7 | 48 | 29 | 13 |
| Percentage | 83.33% | 6.67% | 10% | 91.11% | 1.11% | 7.78% | 53.34% | 32.22% | 14.44% |
| Percentage of results returned | 92.59% | 7.40% | | 98.80% | 1.20% | | 62.34% | 37.66% | |

Other : Participants have either stated that they were unable to obtain sufficient minerals, presented G2 results less than 10 x MDL or supplied no data for that sample.

Table 5.5 Qualitative results from TL participants

Overall performance for unirradiated and irradiated materials is good, with correct identification in 83% and 91% of cases respectively. If those cases where participants were unable to recover sufficient minerals are excluded the correct identification rate rises to 93% and 99% respectively. The incorrect results from unirradiated samples were obtained in laboratories 4 (3 cases), 14 (2 cases) and 5 (1 case), and for the irradiated samples laboratory 11 (1 case). It was noted from reference analyses that small low temperature peaks were observed in the tarragon sample in 11 out of 20 aliquots prepared. This may account for two of the “incorrect” classifications from unirradiated samples in participants’ data, however the other 4 are not represented in reference data, and are thus unexplained.

The success rate for blends is lower, as to be expected both from the nature of the materials and from the quantitative data described above. From table 5.5 it is apparent that the 10% blends are clearly easier to identify, with all reported analyses giving correct outcomes, and 3 cases where insufficient mineral yields were encountered. At the lower concentrations detection rates are poorer, and the concentration dependence is less apparent.

The TL results from this round of the proficiency testing study have confirmed that participants are capable of exceeding the IUPAC criteria for collaborative trials (of 80% success rates from a quorum of at least 5-10 laboratories studying a minimum of 5 blind samples). The performance overall therefore can be considered to affirm the validity of the EN1788 standard, with a new generation of laboratory practitioners in comparison with the original blind studies of the 1990’s. That itself is a useful result at this stage, albeit incidental to the aims of the PT development. However the observations of 7-15% failures to recover sufficient minerals from the samples, which had all produced adequate recoveries in reference analyses, and a small proportion (7%) of incorrect classifications from unirradiated materials does suggest varying levels of proficiency in different laboratories. For blended mixtures the results at 10% concentration were very satisfactory. However at both lower concentration levels detection rates were reduced. Reference analyses suggested that the performance at 1% concentration might have been better, although at 0.1% concentration a spread of outcomes was expected. At this stage it is unclear whether individual laboratory performance differences drive the outcomes for dilute blends, or whether the results are mainly reflecting the statistical fluctuations in mineral yields and brightnesses, although these features may be identifiable in future work of this sort.

6. DISCUSSION AND CONCLUSIONS

As shown in the previous sections the second round of the trial was successfully implemented, with test material distribution, additional homogeneity testing and the second measurement round taking place within a few months.

The calibrated PSL reference data sets obtained in the organising laboratory resemble earlier published data from the original PSL research and inter-laboratory trial data sets, as does TL homogeneity testing.

PSL screening results from 30 participating laboratories were returned in a timely manner and generally show equivalence to the SUERC reference data. For some participants the round 1 difficulties with cross contamination have been overcome or reduced; for others there are still problems. The participants that still experienced these difficulties tended to have greater proportions of unirradiated test materials in higher qualitative classification categories than the reference data set or data from the other participants. The 3 new laboratories who had not taken part in round 1 all performed well.

The data sets have been explored using descriptive statistics and in tabular and graphical forms. There is still evidence of additional inter-laboratory variation in comparison with the reference set even when outliers have been taken into account. There remains scope for further examination of this, in particular to examine the effect of standardisation of instrument response to the paprika standard, although it is noted that such procedures are not currently within the scope of the EN13751 standard. As shown in Sections 3.3.4 99.9% of irradiated samples fell into the intermediate or positive bands, compared with 99.7% in round 1. The proportions of unirradiated samples giving positive results has fallen very slightly to 2.3% from 2.8%. Both of these results represent slight improvements in qualitative outcome between the two rounds and suggest that a stable level of performance has been reached. As shown in the previous report there is a correlation between the qualitative performance level and pooled z scores implying that for PSL screening at least it may be possible to combine quantitative indicators with predictive qualitative performance. The data from round 2 are very encouraging – since there has been improvement, readily seen in z score plots, in sample handling in many laboratories, and it has been shown that qualitative outcomes are well reproduced using these, relatively large sets of samples. It may therefore be possible, utilising these results to develop a design for a “sparse” study utilising smaller numbers of samples, and thus capable of being sustained beyond the timescale of this developmental project.

Calibrated PSL measurements were undertaken successfully by 12 laboratories using a diverse set of irradiation facilities. The results can be shown to be broadly comparable to the reference analyses, with good discrimination between irradiated and unirradiated samples in most cases, and blends being distributed in the calibrated PSL plots in the loci between unirradiated and irradiated samples. In the present study the proportion of samples with very low sensitivity, where calibrated PSL measurements can very usefully call into question the status of negative screening results, was extremely low. Therefore it has not been possible to demonstrate a large impact on qualitative outcomes from screening. Nonetheless interlaboratory differences, potentially attributed to the use of different sources, and perhaps also reflecting variations in sample handling within the guidelines of EN13751, were observed. It would be of interest to explore these differences further in future calibrated PSL rounds, with a sample set that is more critically selected to reveal the impact of such measurements on difficult cases. Against these objectives the results from this round are an encouraging starting point.

For the TL study this round included 18 samples from 6 products, comprising unirradiated, irradiated and blended conditions. The reference analyses showed that it was possible to recover minerals from all products, and that the data sets recovered showed a high degree of internal consistency in both separating irradiated and unirradiated examples and in placing the blends in a concentration ordered sequence. For one sample the reference data identified the presence of small low temperature peaks in the untreated material. Reference analyses comprised 20 preparations per product producing 360 individual TL determinations, equivalent to 180 conventional duplicated TL analyses. All samples in the reference set satisfied the EN1788 requirement for glow 2 sensitivity. Fifteen participants returned data in 241 out of 270 possible TL analyses, the remainder of determinations failing to produce sufficient minerals to satisfy EN1788. Of these some 92.6% of the unirradiated samples 98.9% of the irradiated samples and 62% of the reported results from blends were correctly classified. A small proportion of unirradiated samples were identified as irradiated or containing irradiated material, from four laboratories. One irradiated sample was not successfully identified in one laboratory. For the blends the highest concentrations (10%) were all correctly identified in participants' results; while at lower concentrations laboratories detection rates for the irradiated component were reduced. The qualitative performance is encouraging overall, bearing in mind the sample sizes distributed and the differing equipment used by participants. Temperature control, however, could be improved. Quantitative analyses were also conducted using z scores on the recorded glow ratios. This appears to be informative and to identify some laboratory to laboratory differences, although further work is needed to examine the distribution of the glow ratio variable and to explore other means of distinguishing between the detection conditions of each laboratory, for example by utilising absolute signal levels in comparison with detection limits.

The second round of proficiency testing analysis under this project has again resulted in a high return rate of extremely good and useful data from laboratories engaged in PSL and TL analysis of irradiated foods. Participation in such work represents an important commitment by participants and by the study organisers to the ongoing task of ensuring high analytical quality in routine determinations. It is to be hoped that this work will lead to sustained activities in support of food irradiation analyses.

APPENDIX A: TEST SAMPLE DISTRIBUTION AND STUDY PROTOCOL

**PLEASE FAX YOUR REPLY TO LORNA CARMICHAEL / SAFFRON FISK AT
44 1355229898**

LABORATORY :

We have received the Proficiency test materials on _____ (date)

Were the test materials in good condition on arrival ? _____

Any other comments :

Food Standards Agency Sponsored Development of Proficiency Testing for Detection of Irradiated Foods. Round 2 (PSL and TL)

1. ROUND 2

The second round of the Food Standards Agency project to develop PT schemes for detection of irradiated foods is now underway. Many thanks for your continued participation in this project. Round 2 involves both PSL measurements, following the very interesting first round trial, and also for some laboratories introduces TL measurements. The data returns should be made electronically using the EXCEL spreadsheets issued. Please enter your **lab number as marked on the samples** (it may be different for the PSL and TL parts) in the data sheets.

The protocols for PSL and TL analysis are summarised here.

2. PSL ANALYSIS

All PSL participants are asked to return screening results from round 2, as outlined below in section 2.1. Some participants have also indicated a desire to continue to perform calibrated PSL measurements as well; this is outlined in section 2.2.

This time you will receive a new sample of irradiated paprika standard material, plus a set of 72 numbered samples for PSL analysis, and some petri dishes.

2.1 PSL Screening

(i) Please set your system up in your normal manner, ready for use in accordance with EN13751.

(ii) Before working on the test samples please dispense and measure **10 different portions of the new irradiated paprika standard**. This is to provide us with a measurement of the sensitivity of your instrument to this new material, and does not form part of the blind trial. For round 2, rather than dispensing these standard samples by volume, please weigh **2 g.** of the standard into the petri dish, shake or spread it to cover the base, and record the counting data. The material may be discarded after measurement. If you wish to re-use the sample petri dish for all 10 paprika measurements you may. But please be careful not to spread irradiated paprika to the outsides of the dishes or to the instrument chamber.

(iii) Then measure the samples using a new pair of petri dishes for each of the numbered test samples. You should dispense and measure these in the numbered sequence using your normal procedure.

Measure each sample in duplicate, and record the PSL screening results in the EXCEL spreadsheet provided by disk and email.

Please be careful when dispensing, handling and measuring the samples to avoid cross-contamination, or contamination of the sample chamber. Some of the samples are irradiated, and may have high sensitivities.

PSL screening results should be returned to the organisers by email (to D.Sanderson@suerc.gla.ac.uk &/or L.Carmichael@suerc.gla.ac.uk) both using the EXCEL spreadsheet and also with a copy of the PSL summary files for the runs used.

2.2 Calibrated PSL

For those laboratories continuing with calibrated PSL measurements, seal your petri dishes and bag them after screening measurements. They should then be irradiated using a gamma or electron beam facility to a dose of 1 kGy. Then after a suitable post-irradiation delay (at least 24 hours) they should be re-measured and evaluated.

Return these results also using the combination of EXCEL spreadsheet and back-up PSL summary file.

3. TL laboratories

Participating TL laboratories are supplied with two LiF TLD:100 chips for temperature scale calibration, and a set of 18 numbered test samples.

The two LiF chips should be measured for TL using the same system as used to read the trial samples, and the full glow curve returned to the organisers (in its normal format) plus an evaluation of the temperature of the main (LiF peak V) band.

The test samples should be prepared in duplicate following EN1788 compatible procedures, including appropriate quality control measures (particularly blanks). Prepared samples should then be read out by TL (typically from 0-400°C at 5°C s⁻¹) to produce **first glow** TL data. The samples should then be irradiated using a ⁶⁰Co, ⁹⁰Sr or equivalent radiation source to a dose of 1 kGy, stored or preheated (eg at 50°C for one hour) and then remeasured to obtain **second glow** TL data. Return data for each sample comprising the first glow and second glow intensities, the TL glow ratio, observations on glow1 peak shapes and your evaluation. The samples may be irradiated, untreated, or mixtures containing irradiated materials. Please evaluate your data in these terms.

The TL results should be returned both in the supplied EXCEL spreadsheet (in the TL results page) and with additional supporting data files for the LiF glow curves plus any additional observations that you wish to have considered for evaluation purposes.

4. Return dates

Please return **PSL screening results no later than 2nd June**; for **calibrated PSL by 16th June**, and for **TL results no later than 30th June 2006**.

APPENDIX B – RAW DATA

Table B.1 Reference values based on log statistics for irradiated and unirradiated test materials – calibrated PSL

| SP Number | Description | Irradiated cal | | | Unirradiated cal | | |
|-----------|----------------------|------------------------|------|-------|------------------------|------|-------|
| | | Log ₁₀ mean | SD | CV | Log ₁₀ mean | SD | CV |
| 8511 | Italian seasoning | 5.29 | 0.07 | 1.39% | 5.19 | 0.06 | 1.20% |
| 8512 | Basil | 4.61 | 0.21 | 4.53% | 4.30 | 0.06 | 1.39% |
| 8513 | Curry no 3 | 7.23 | 0.03 | 0.42% | 6.94 | 0.04 | 0.61% |
| 8514 | BBQ seasoning | 7.54 | 0.10 | 1.27% | 7.28 | 0.07 | 0.94% |
| 8515 | Cinnamon | 4.25 | 0.13 | 2.97% | 3.88 | 0.12 | 3.03% |
| 8516 | Paprika | 5.58 | 0.06 | 1.07% | 5.11 | 0.05 | 1.06% |
| 8517 | Cumin | 5.63 | 0.07 | 1.25% | 5.48 | 0.07 | 1.21% |
| 8518 | Oregano | 5.46 | 0.10 | 1.77% | 5.33 | 0.04 | 0.75% |
| 8519 | Ground black pepper | 4.75 | 0.21 | 4.46% | 4.47 | 0.21 | 4.60% |
| 8520 | Chives | 4.22 | 0.18 | 4.22% | 4.22 | 0.33 | 7.80% |
| 8521 | Mint | 4.48 | 0.07 | 1.54% | 4.21 | 0.04 | 0.96% |
| 8522 | Ground white pepper | 3.17 | 0.28 | 8.98% | 2.96 | 0.19 | 6.35% |
| 8523 | Sage | 5.43 | 0.13 | 2.37% | 5.25 | 0.08 | 1.49% |
| 8524 | Thai seasoning | 7.71 | 0.16 | 2.07% | 7.64 | 0.06 | 0.72% |
| 8525 | Ground mixed spice | 4.53 | 0.17 | 3.86% | 4.24 | 0.19 | 4.43% |
| 8526 | Ground coriander | 4.38 | 0.16 | 3.75% | 4.18 | 0.09 | 2.27% |
| 8527 | Steak seasoning | 7.87 | 0.06 | 0.72% | 7.41 | 0.14 | 1.94% |
| 8528 | Medium curry | 7.15 | 0.02 | 0.27% | 6.98 | 0.05 | 0.69% |
| 8529 | Nutmeg | 4.93 | 0.08 | 1.69% | 4.71 | 0.17 | 3.67% |
| 8530 | Whole black pepper | 3.29 | 0.16 | 4.79% | 3.16 | 0.16 | 4.93% |
| 8531 | Ginger | 5.38 | 0.13 | 2.39% | 5.11 | 0.08 | 1.59% |
| 8532 | Thyme | 4.90 | 0.07 | 1.49% | 4.71 | 0.03 | 0.67% |
| 8533 | Rosemary | 4.77 | 0.06 | 1.17% | 4.50 | 0.03 | 0.65% |
| 8534 | Turmeric | 4.15 | 0.02 | 0.59% | 3.93 | 0.09 | 2.29% |
| 8535 | Chilli powder | 5.18 | 0.05 | 1.03% | 4.96 | 0.05 | 1.09% |
| 8536 | Garlic powder | 5.90 | 0.08 | 1.33% | 5.82 | 0.07 | 1.15% |
| 8537 | Mixed herbs | 5.22 | 0.07 | 1.42% | 5.04 | 0.08 | 1.60% |
| 8538 | Parsley | 3.76 | 0.23 | 6.12% | 3.52 | 0.18 | 5.13% |
| 8539 | Tarragon | 4.49 | 0.29 | 6.45% | 4.19 | 0.13 | 3.04% |
| 8540 | Basil | 4.42 | 0.05 | 1.04% | 4.20 | 0.08 | 1.88% |
| 8541 | Mint | 4.67 | 0.05 | 1.17% | 4.41 | 0.12 | 2.61% |
| 8542 | Parsley | 4.80 | 0.05 | 1.14% | 4.62 | 0.07 | 1.51% |
| 8543 | Tarragon | 3.83 | 0.11 | 2.89% | 3.75 | 0.19 | 4.99% |
| 8544 | Thyme | 6.22 | 0.02 | 0.36% | 6.01 | 0.05 | 0.85% |
| 8545 | Cinnamon | 4.06 | 0.13 | 3.23% | 3.79 | 0.08 | 2.11% |
| 8546 | Ground coriander | 3.72 | 0.13 | 3.49% | 3.42 | 0.09 | 2.67% |
| 8547 | Cumin | 4.26 | 0.06 | 1.40% | 4.04 | 0.06 | 1.40% |
| 8548 | Paprika | 5.11 | 0.06 | 1.27% | 4.79 | 0.12 | 2.48% |
| 8549 | Turmeric | 4.96 | 0.09 | 1.82% | 4.65 | 0.12 | 2.64% |
| 8571 | Milk thistle seed | 5.14 | 0.06 | 1.23% | 4.96 | 0.08 | 1.64% |
| 8572 | Alfalfa herb | 4.54 | 0.09 | 2.03% | 4.30 | 0.07 | 1.52% |
| 8573 | Saw palmetto berry | 3.55 | 0.18 | 4.99% | 3.31 | 0.22 | 6.60% |
| 8574 | Dong quai root | 6.03 | 0.04 | 0.70% | 5.92 | 0.04 | 0.68% |
| 8575 | Guarana seed | 3.67 | 0.20 | 5.33% | 3.32 | 0.12 | 3.57% |
| 8576 | Ginkgo biloba leaves | 5.33 | 0.04 | 0.72% | 5.19 | 0.07 | 1.25% |
| 8577 | Siberian ginseng | 5.87 | 0.05 | 0.90% | 5.68 | 0.05 | 0.86% |
| 8578 | Green tea | 4.98 | 0.09 | 1.82% | 4.75 | 0.13 | 2.68% |
| 8579 | Echinacea | 4.99 | 0.07 | 1.43% | 4.86 | 0.13 | 2.70% |
| 8580 | Garlic powder | 5.90 | 0.05 | 0.87% | 5.70 | 0.05 | 0.82% |

Table B.2 Reference data for blended samples

| | Linear | | | | Log | | | |
|--------|----------|----------|---------|---------|------|------|------|--------|
| | Mean | SD | SE | CV | Mean | SD | SE | CV |
| SP8351 | 21015.40 | 15521.63 | 4908.37 | 73.86% | 4.24 | 0.26 | 0.08 | 6.13% |
| SP8533 | 691.40 | 500.97 | 158.42 | 72.46% | 2.77 | 0.25 | 0.08 | 8.90% |
| SP8535 | 931.60 | 1059.58 | 335.07 | 113.74% | 2.84 | 0.29 | 0.09 | 10.27% |
| SP8539 | 5148.14 | 2754.43 | 871.03 | 53.50% | 3.45 | 0.28 | 0.09 | 8.06% |
| SP8577 | 10073.70 | 3606.48 | 1140.47 | 35.80% | 3.98 | 0.15 | 0.05 | 3.77% |
| SP8578 | 484.80 | 108.71 | 34.38 | 22.42% | 2.68 | 0.10 | 0.03 | 3.70% |

a) screening

| | Linear | | | | Log | | | |
|--------|-----------|----------|----------|--------|------|------|------|-------|
| | Mean | SD | SE | CV | Mean | SD | SE | CV |
| SP8351 | 165260.20 | 31390.18 | 9926.45 | 18.99% | 5.21 | 0.08 | 0.03 | 1.62% |
| SP8533 | 40215.30 | 3157.07 | 998.35 | 7.85% | 4.60 | 0.03 | 0.01 | 0.73% |
| SP8535 | 106267.40 | 37812.79 | 11957.45 | 35.58% | 5.01 | 0.13 | 0.04 | 2.54% |
| SP8539 | 23122.70 | 10679.65 | 3377.20 | 46.19% | 4.33 | 0.18 | 0.06 | 4.19% |
| SP8577 | 559935.40 | 52882.23 | 16722.83 | 9.44% | 5.75 | 0.04 | 0.01 | 0.72% |
| SP8578 | 61180.70 | 13827.13 | 4372.52 | 22.60% | 4.78 | 0.10 | 0.03 | 2.11% |

b) calibrated

Table B.3 Raw data for unirradiated test materials – calibrated PSL

| SP | Pot 1 | Pot 2 | Pot 3 | Pot 4 | Pot 5 | Pot 6 | Pot 7 | Pot 8 | Pot 9 | Pot 10 | Mean | SD | CV |
|------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|---------|---------|
| 8511 | 124311 | 174636 | 172345 | 138274 | 195529 | 136257 | 158792 | 148507 | 172136 | 138622 | 155941 | 22389 | 14.36% |
| 8512 | 16132 | 24527 | 20025 | 17319 | 23798 | 16826 | 20080 | 20625 | 19930 | 20623 | 19989 | 2747 | 13.74% |
| 8513 | 8027697 | 9531145 | 7595899 | 8792558 | 10295791 | 9618417 | 9439472 | 8437471 | 7949878 | 8547162 | 8823549 | 869839 | 9.86% |
| 8514 | 15029686 | 17761159 | 16596225 | 18435764 | 21324919 | 19868612 | 26472466 | 20821701 | 18883158 | 16972465 | 19216616 | 3203575 | 16.67% |
| 8515 | 7660 | 7001 | 5922 | 6766 | 8724 | 5644 | 5603 | 8183 | 13797 | 8563 | 7786 | 2410 | 30.95% |
| 8516 | 164891 | 128701 | 126253 | 149987 | 120886 | 127646 | 121181 | 126010 | 105271 | 139896 | 131072 | 16709 | 12.75% |
| 8517 | 317713 | 289945 | 315485 | 426044 | 337408 | 254494 | 313305 | 291928 | 247751 | 289703 | 308378 | 49876 | 16.17% |
| 8518 | 223361 | 180899 | 211856 | 250231 | 233291 | 200282 | 213769 | 215075 | 208589 | 194836 | 213219 | 19586 | 9.19% |
| 8519 | 19778 | 74041 | 36757 | 16855 | 21550 | 27256 | 30711 | 17580 | 32240 | 49489 | 32626 | 17677 | 54.18% |
| 8520 | 15384 | 19390 | 7809 | 21671 | 12093 | 23581 | 12291 | 7511 | 11333 | 105097 | 23616 | 29154 | 123.45% |
| 8521 | 18087 | 18090 | 15275 | 16624 | 14426 | 16499 | 14670 | 15071 | 15511 | 18550 | 16280 | 1526 | 9.37% |
| 8522 | 628 | 633 | 885 | 811 | 617 | 839 | 1950 | 853 | 1989 | 741 | 995 | 523 | 52.60% |
| 8523 | 145429 | 253022 | 138421 | 192922 | 159090 | 171456 | 152341 | 188885 | 195457 | 193705 | 179073 | 33658 | 18.80% |
| 8524 | 32793765 | 40019983 | 43555176 | 45856292 | 44042568 | 43382528 | 48548840 | 43388938 | 50013171 | 51167164 | 44276843 | 5304845 | 11.98% |
| 8525 | 12446 | 13570 | 29077 | 12098 | 11673 | 16579 | 27123 | 35912 | 19605 | 10963 | 18905 | 8807 | 46.59% |
| 8526 | 16278 | 12418 | 14607 | 18475 | 10888 | 21062 | 14528 | 19365 | 15198 | 11734 | 15455 | 3368 | 21.79% |
| 8527 | 36635438 | 18376016 | 30018868 | 32679996 | 41455927 | 13457906 | 25453478 | 26997806 | 23208559 | 23734635 | 27201863 | 8353141 | 30.71% |
| 8528 | 7457131 | 9729985 | 9470072 | 9012609 | 10150807 | 10695855 | 9678367 | 9728232 | 9858988 | 11303502 | 9708555 | 1019189 | 10.50% |
| 8529 | 67145 | 47572 | 41695 | 102914 | 37116 | 72495 | 43703 | 69252 | 27079 | 39122 | 54809 | 22763 | 41.53% |
| 8530 | 896 | 1049 | 1254 | 1937 | 1726 | 1057 | 1838 | 2017 | 1016 | 2494 | 1528 | 544 | 35.59% |
| 8531 | 97108 | 119508 | 145583 | 163155 | 134736 | 150802 | 142284 | 92927 | 146342 | 124628 | 131707 | 23049 | 17.50% |
| 8532 | 48463 | 52082 | 53367 | 57135 | 44638 | 52486 | 54640 | 53221 | 47641 | 51496 | 51517 | 3660 | 7.10% |
| 8533 | 32396 | 28314 | 30600 | 30632 | 34551 | 32080 | 35019 | 31314 | 29208 | 30439 | 31455 | 2133 | 6.78% |
| 8534 | 12440 | 10591 | 8385 | 7255 | 10228 | 8521 | 9336 | 6818 | 7080 | 6874 | 8753 | 1885 | 21.54% |
| 8535 | 97975 | 87569 | 96507 | 82548 | 105819 | 104967 | 98528 | 71314 | 81938 | 93997 | 92116 | 11054 | 12.00% |
| 8536 | 711968 | 672927 | 737937 | 656307 | 708094 | 635172 | 462190 | 561803 | 775055 | 730130 | 665158 | 93420 | 14.04% |
| 8537 | 104280 | 166959 | 126955 | 103555 | 108178 | 102796 | 109003 | 102544 | 80538 | 111785 | 111659 | 22513 | 20.16% |
| 8538 | 3067 | 4354 | 5131 | 5000 | 1698 | 3948 | 2724 | 2177 | 2038 | 5093 | 3523 | 1346 | 38.20% |
| 8539 | 24728 | 15381 | 12146 | 14151 | 14704 | 10902 | 14550 | 17689 | 11794 | 26024 | 16207 | 5220 | 32.21% |
| 8540 | 16867 | 13927 | 15001 | 16051 | 11121 | 15800 | 15765 | 23051 | 17283 | 15398 | 16026 | 3012 | 18.80% |
| 8541 | 33755 | 21131 | 24790 | 32832 | 22170 | 26179 | 20982 | 44762 | 19953 | 21906 | 26846 | 7959 | 29.65% |
| 8542 | 36526 | 47135 | 31425 | 55561 | 38599 | 47783 | 40647 | 43809 | 41754 | 37895 | 42113 | 6823 | 16.20% |
| 8543 | 5064 | 4323 | 7390 | 11625 | 9852 | 4884 | 6218 | 5515 | 2660 | 4181 | 6171 | 2745 | 44.48% |
| 8544 | 1316975 | 1139076 | 990324 | 1044169 | 1082905 | 1024919 | 917660 | 975545 | 904996 | 899760 | 1029633 | 127813 | 12.41% |
| 8545 | 5503 | 5531 | 4465 | 7140 | 8200 | 6361 | 6864 | 5101 | 6941 | 6880 | 6299 | 1126 | 17.88% |
| 8546 | 2543 | 2485 | 2257 | 3327 | 3399 | 2773 | 3240 | 2727 | 1716 | 2294 | 2676 | 535 | 19.98% |
| 8547 | 12175 | 10784 | 9745 | 11472 | 10429 | 10630 | 14715 | 10373 | 11214 | 9179 | 11072 | 1535 | 13.86% |
| 8548 | 63863 | 61574 | 62022 | 42968 | 64923 | 102073 | 59892 | 52478 | 42657 | 87431 | 63988 | 18429 | 28.80% |
| 8549 | 62861 | 44309 | 48757 | 33702 | 34810 | 78848 | 36108 | 43696 | 33598 | 46479 | 46317 | 14554 | 31.42% |
| 8571 | 116238 | 110391 | 78923 | 74572 | 76921 | 89288 | 112818 | 73695 | 91127 | 111566 | 93554 | 17522 | 18.73% |
| 8572 | 19754 | 21092 | 15996 | 24569 | 23732 | 19682 | 20217 | 16126 | 22163 | 17258 | 20059 | 2970 | 14.81% |
| 8573 | 2799 | 1486 | 1378 | 4241 | 3920 | 1735 | 3356 | 1113 | 1575 | 1239 | 2284 | 1185 | 51.89% |
| 8574 | 820184 | 738431 | 867873 | 864771 | 1000367 | 883548 | 729644 | 836450 | 772735 | 820145 | 833415 | 79155 | 9.50% |
| 8575 | 1777 | 1782 | 1906 | 1875 | 1950 | 2584 | 2716 | 3735 | 1443 | 1934 | 2170 | 668 | 30.77% |
| 8576 | 165002 | 150555 | 129431 | 201677 | 141517 | 149142 | 174504 | 181116 | 151989 | 124632 | 156957 | 23860 | 15.20% |
| 8577 | 435986 | 431089 | 412854 | 473832 | 502088 | 493176 | 539856 | 582503 | 422430 | 479106 | 477292 | 54639 | 11.45% |
| 8578 | 41693 | 45999 | 61686 | 80830 | 60196 | 70724 | 94483 | 48400 | 45544 | 40315 | 58987 | 18235 | 30.91% |
| 8579 | 104369 | 108365 | 58761 | 42956 | 67785 | 54471 | 77439 | 78608 | 95777 | 60437 | 74897 | 22117 | 29.53% |
| 8580 | 570045 | 478819 | 571735 | 509948 | 543500 | 509107 | 542497 | 500422 | 404174 | 456626 | 508687 | 52377 | 10.30% |

Table B.4 Raw data for irradiated test materials – calibrated PSL

| SP | Pot 1 | Pot 2 | Pot 3 | Pot 4 | Pot 5 | Pot 6 | Pot 7 | Pot 8 | Pot 9 | Pot 10 | Mean | SD | CV |
|------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|---------|
| 8511 | 205012 | 133563 | 219983 | 222758 | 187541 | 198042 | 185207 | 210055 | 241356 | 166984 | 197050 | 30820 | 15.64% |
| 8512 | 36950 | 35885 | 156351 | 33391 | 36952 | 37758 | 31493 | 37326 | 31321 | 32618 | 47005 | 38500 | 81.91% |
| 8513 | 15940295 | 17454502 | 15954032 | 16321188 | 16432017 | 18495317 | 17214216 | 16004401 | 18980855 | 18833104 | 17162993 | 1224675 | 7.14% |
| 8514 | 30511729 | 22979030 | 25790752 | 33689888 | 35409912 | 40965138 | 46446285 | 39403028 | 35437912 | 41488329 | 35212200 | 7296649 | 20.72% |
| 8515 | 15945 | 15663 | 16898 | 17102 | 14765 | 18146 | 21633 | 14783 | 38359 | 14839 | 18813 | 7177 | 38.15% |
| 8516 | 423074 | 296422 | 394453 | 357865 | 386735 | 390009 | 439086 | 334595 | 343968 | 466795 | 383300 | 51597 | 13.46% |
| 8517 | 459632 | 493724 | 561377 | 443819 | 397695 | 333192 | 375511 | 447807 | 348913 | 461538 | 432321 | 69608 | 16.10% |
| 8518 | 294686 | 235946 | 501377 | 306226 | 331260 | 254512 | 281540 | 247562 | 289073 | 237693 | 297988 | 78046 | 26.19% |
| 8519 | 99216 | 54381 | 38136 | 65622 | 164537 | 42712 | 37531 | 51174 | 37685 | 42341 | 63334 | 40259 | 63.57% |
| 8520 | 13805 | 25207 | 20134 | 12246 | 10910 | 27739 | 10089 | 14175 | 13658 | 31898 | 17986 | 7757 | 43.13% |
| 8521 | 43866 | 30922 | 32631 | 24904 | 25888 | 27673 | 31083 | 27239 | 29989 | 30981 | 30518 | 5327 | 17.46% |
| 8522 | 2105 | 1166 | 1304 | 1110 | 1442 | 8552 | 1154 | 925 | 1130 | 995 | 1988 | 2330 | 117.19% |
| 8523 | 170225 | 267709 | 339267 | 286942 | 172395 | 265976 | 358465 | 213080 | 328792 | 396191 | 279904 | 77498 | 27.69% |
| 8524 | 34697103 | 28929879 | 33729531 | 38753155 | 62642226 | 67833960 | 76377140 | 63789246 | 66209421 | 68778642 | 54174030 | 17870766 | 32.99% |
| 8525 | 30345 | 42252 | 61127 | 46435 | 29394 | 22544 | 27471 | 58205 | 24415 | 19269 | 36146 | 14966 | 41.40% |
| 8526 | 20737 | 61332 | 34861 | 20442 | 24639 | 20850 | 17567 | 20073 | 18547 | 21958 | 26101 | 13299 | 50.95% |
| 8527 | 62859702 | 79048552 | 65244265 | 65935412 | 78088765 | 75960545 | 75871798 | 76847128 | 97474796 | 66238607 | 74356957 | 10184395 | 13.70% |
| 8528 | 13714686 | 14710803 | 13500559 | 13803206 | 13946229 | 13475184 | 15299674 | 14978058 | 14132903 | 14236428 | 14179773 | 627983 | 4.43% |
| 8529 | 92647 | 90515 | 69013 | 104811 | 115615 | 99213 | 61485 | 80791 | 77382 | 83563 | 87504 | 16480 | 18.83% |
| 8530 | 1321 | 2188 | 2164 | 1497 | 1812 | 1460 | 4397 | 1595 | 1693 | 2795 | 2092 | 922 | 44.08% |
| 8531 | 309141 | 264624 | 237577 | 241497 | 196963 | 448462 | 250712 | 164018 | 181971 | 186088 | 248105 | 83139 | 33.51% |
| 8532 | 76978 | 78579 | 76000 | 64025 | 78543 | 80745 | 117998 | 76509 | 67305 | 91295 | 80798 | 14994 | 18.56% |
| 8533 | 75863 | 51555 | 56422 | 60876 | 54742 | 53208 | 50745 | 67276 | 62215 | 54833 | 58774 | 7926 | 13.49% |
| 8534 | 14282 | 13654 | 15181 | 14648 | 14823 | 12920 | 13747 | 15646 | 14407 | 13770 | 14308 | 810 | 5.66% |
| 8535 | 165649 | 158439 | 144594 | 182215 | 134819 | 163733 | 153205 | 134490 | 119652 | 153930 | 151073 | 18190 | 12.04% |
| 8536 | 663803 | 954660 | 632301 | 1119455 | 802072 | 834005 | 630645 | 806464 | 752535 | 817118 | 801306 | 150944 | 18.84% |
| 8537 | 131628 | 163098 | 213688 | 136943 | 186619 | 193846 | 138559 | 167776 | 153987 | 197092 | 168324 | 28523 | 16.95% |
| 8538 | 6486 | 4198 | 2725 | 4958 | 3994 | 20527 | 6871 | 6530 | 6126 | 4881 | 6730 | 5026 | 74.68% |
| 8539 | 21296 | 21930 | 24484 | 22308 | 24476 | 38455 | 32744 | 192770 | 23319 | 24257 | 42604 | 53044 | 124.51% |
| 8540 | 23759 | 25601 | 24694 | 24960 | 30842 | 24896 | 28904 | 28961 | 27221 | 21756 | 26159 | 2771 | 10.59% |
| 8541 | 41298 | 48859 | 61036 | 44537 | 48744 | 45028 | 44358 | 51046 | 41847 | 39869 | 46662 | 6197 | 13.28% |
| 8542 | 60206 | 51148 | 73894 | 72176 | 64167 | 62968 | 68380 | 52286 | 69279 | 60426 | 63493 | 7745 | 12.20% |
| 8543 | 7378 | 5889 | 6836 | 7694 | 5572 | 9357 | 5419 | 4262 | 9056 | 8525 | 6999 | 1701 | 24.30% |
| 8544 | 1681549 | 1714240 | 1604792 | 1554809 | 1801888 | 1701555 | 1704925 | 1803930 | 1571401 | 1642304 | 1678139 | 86272 | 5.14% |
| 8545 | 10845 | 13067 | 10181 | 8298 | 10586 | 17496 | 20534 | 7860 | 10695 | | 12174 | 4229 | 34.74% |
| 8546 | 10303 | 8203 | 6366 | 3294 | 5060 | 4989 | 5037 | 7231 | 5107 | 6052 | 3286 | 1992 | 60.61% |
| 8547 | 19889 | 18599 | 18668 | 17661 | 22505 | 15548 | 18699 | 15718 | 14473 | 20556 | 18232 | 2470 | 13.55% |
| 8548 | 112485 | 128866 | 111231 | 135256 | 128072 | 168499 | 145477 | 142765 | 101446 | 122285 | 129638 | 19579 | 15.10% |
| 8549 | 74809 | 86536 | 82440 | 87361 | 99000 | 82916 | 131758 | 92327 | 125180 | 68104 | 93043 | 20580 | 22.12% |
| 8571 | 150936 | 156871 | 152879 | 178261 | 112120 | 132837 | 131899 | 117584 | 122694 | 146859 | 140294 | 20484 | 14.60% |
| 8572 | 40157 | 39564 | 30496 | 52932 | 26932 | 37761 | 29039 | 29226 | 29321 | 38212 | 35364 | 7973 | 22.55% |
| 8573 | 2753 | 3757 | 5573 | 3585 | 2725 | 2884 | 2101 | 8330 | 2620 | 4021 | 3835 | 1859 | 48.49% |
| 8574 | 1171395 | 1164631 | 1057568 | 999988 | 1063220 | 964566 | 1085723 | 949237 | 1274135 | 974717 | 1070518 | 105952 | 9.90% |
| 8575 | 4055 | 4391 | 4724 | 4022 | 14378 | 5508 | 4833 | 3156 | 4600 | 2631 | 5230 | 3320 | 63.47% |
| 8576 | 189752 | 242044 | 211731 | 208093 | 220578 | 237932 | 207004 | 241089 | 214728 | 190608 | 216356 | 19164 | 8.86% |
| 8577 | 903049 | 819138 | 746424 | 772496 | 755818 | 711804 | 724061 | 751250 | 794219 | 561518 | 753978 | 87059 | 11.55% |
| 8578 | 134763 | 129319 | 97620 | 86958 | 77496 | 76716 | 81962 | 101404 | 76716 | 102775 | 96573 | 21219 | 21.97% |
| 8579 | 134763 | 129319 | 97620 | 86958 | 77496 | 76716 | 81962 | 101404 | 76716 | 102775 | 96573 | 21219 | 21.97% |
| 8580 | 739471 | 781331 | 1041743 | 801144 | 827719 | 726787 | 834829 | 781923 | 686038 | 709231 | 793022 | 100477 | 12.67% |

Table B.5 Raw data for homogeneity testing of blended samples - screening

| SP | Pot 1 | Pot 2 | Pot 3 | Pot 4 | Pot 5 | Pot 6 | Pot 7 | Pot 8 | Pot 9 | Pot 10 | Mean | SD | CV |
|------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|-------|-------|---------|
| 8531 | 8989 | 45122 | 11516 | 12373 | 17718 | 12305 | 11552 | 53863 | 16332 | 20384 | 21015 | 15522 | 73.86% |
| 8533 | 2021 | 242 | 786 | 394 | 551 | 833 | 392 | 522 | 662 | 511 | 691 | 501 | 72.46% |
| 8535 | 645 | 381 | 593 | 489 | 3898 | 560 | 456 | 836 | 1006 | 452 | 932 | 1060 | 113.74% |
| 8539 | 3194 | 6254 | 10240 | 2068 | 1247 | 3283 | 2334 | 2430 | 1367 | 2526 | 3494 | 2754 | 78.83% |
| 8577 | 7080 | 17855 | 11459 | 12239 | 9616 | 7144 | 5819 | 10736 | 6876 | 11913 | 10074 | 3606 | 35.80% |
| 8578 | 339 | 604 | 491 | 670 | 491 | 352 | 564 | 379 | 505 | 453 | 485 | 109 | 22.42% |

Table B.6 Raw data for homogeneity testing of blended samples – calibrated

| SP | Pot 1 | Pot 2 | Pot 3 | Pot 4 | Pot 5 | Pot 6 | Pot 7 | Pot 8 | Pot 9 | Pot 10 | Mean | SD | CV |
|------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|----------|----------|--------|
| 8531 | 193375 | 151100 | 212606 | 149594 | 140099 | 190742 | 113002 | 159045 | 196617 | 146422 | 165260.2 | 31390.18 | 18.99% |
| 8533 | 46740 | 40186 | 40409 | 43853 | 35456 | 39001 | 40635 | 39033 | 39346 | 37494 | 40215.3 | 3157.065 | 7.85% |
| 8535 | 115051 | 74743 | 77814 | 115621 | 83172 | 86666 | 104694 | 205160 | 92172 | 107581 | 106267.4 | 37812.79 | 35.58% |
| 8539 | 26057 | 31705 | 47843 | 14624 | 13660 | 14195 | 25237 | 13989 | 21232 | 22685 | 23122.7 | 10679.65 | 46.19% |
| 8577 | 552352 | 615177 | 554081 | 492794 | 502941 | 602737 | 616370 | 626743 | 492731 | 543428 | 559935.4 | 52882.23 | 9.44% |
| 8578 | 39329 | 65200 | 69402 | 80236 | 48451 | 53277 | 51719 | 65173 | 56711 | 82309 | 61180.7 | 13827.13 | 22.60% |

Table B.7 Reference values for TL Homogeneity Testing based on Glow Ratio

| | Mean | SD |
|-------------|----------|----------|
| Ginger U | -2.44737 | 0.211783 |
| Rosemary U | -2.94587 | 0.534393 |
| Chilli U | -2.85467 | 0.241645 |
| Tarragon U | -2.45061 | 0.584921 |
| Ginseng U | -2.31409 | 0.110042 |
| Green Tea U | -2.4126 | 0.118619 |
| Ginger I | 0.49417 | 0.094234 |
| Rosemary I | 0.295098 | 0.052178 |
| Chilli I | 0.220145 | 0.065008 |
| Tarragon I | 0.357058 | 0.097278 |
| Ginseng I | 0.31036 | 0.05731 |
| Green Tea I | 0.308111 | 0.046717 |
| Ginger B | -0.5351 | 0.288699 |
| Rosemary B | -1.66684 | 0.332742 |
| Chilli B | -2.4516 | 0.742994 |
| Tarragon B | -0.73895 | 0.512392 |
| Ginseng B | -1.25363 | 0.08994 |
| Green Tea B | -2.34651 | 0.328888 |

APPENDIX C – PARTICIPANTS’ RAW DATA

Table C.1 Raw Data Returned By Laboratory 1 – screening of unirradiated samples

| SP number | Status | Code number | Lab 1 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 345 | 380 |
| 8513 | U | 31 | 367 | 377 |
| 8514 | U | 56 | 1074 | 908 |
| 8515 | U | 42 | 350 | 186 |
| 8516 | U | 25 | 977 | 813 |
| 8517 | U | 58 | 373 | 485 |
| 8518 | U | 54 | 238 | 295 |
| 8519 | U | 39 | 166 | 309 |
| 8520 | U | 68 | 327 | 524 |
| 8521 | U | 3 | 239 | 302 |
| 8523 | U | 12 | 397 | 292 |
| 8525 | U | 1 | 818 | 820 |
| 8526 | U | 59 | 385 | 267 |
| 8528 | U | 10 | 335 | 346 |
| 8529 | U | 11 | 416 | 254 |
| 8530 | U | 23 | 246 | 205 |
| 8531 | U | 29 | 622 | 876 |
| 8532 | U | 27 | 321 | 320 |
| 8533 | U | 14 | 309 | 316 |
| 8534 | U | 4 | 329 | 238 |
| 8535 | U | 53 | 482 | 580 |
| 8536 | U | 19 | 1374 | 1492 |
| 8538 | U | 61 | 167 | 428 |
| 8539 | U | 20 | 286 | 410 |
| 8571 | U | 9 | 589 | 371 |
| 8572 | U | 6 | 357 | 371 |
| 8573 | U | 60 | 235 | 313 |
| 8574 | U | 70 | 1285 | 1150 |
| 8575 | U | 71 | 257 | 381 |
| 8576 | U | 13 | 319 | 441 |
| 8577 | U | 55 | 2314 | 2600 |
| 8578 | U | 28 | 322 | 398 |
| 8579 | U | 57 | 308 | 321 |

Non-negative results shaded

Table C.2 Raw Data Returned By Laboratory 1 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 1 | |
|-----------|--------|-------------|----------|-------------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 25370 | 26239 |
| 8513 | I | 26 | 12812648 | 12073490 |
| 8514 | I | 64 | 16211998 | 18186510 |
| 8515 | I | 35 | 8267 | 11678 |
| 8516 | I | 18 | 184264 | 190536 |
| 8517 | I | 5 | 384688 | 324003 |
| 8518 | I | 33 | 234603 | 266577 |
| 8519 | I | 24 | 57810 | 43263 |
| 8520 | I | 69 | 13056 | 10325 |
| 8521 | I | 38 | 17858 | 25178 |
| 8523 | I | 32 | 246477 | 344781 |
| 8525 | I | 8 | 22605 | 36455 |
| 8526 | I | 50 | 15856 | 12942 |
| 8528 | I | 62 | 8938396 | 9686189 |
| 8529 | I | 51 | 66732 | 154530 |
| 8530 | I | 21 | 1742 | 947 |
| 8531 | I | 45 | 283706 | 399343 |
| 8532 | I | 49 | 62853 | 142591 |
| 8533 | I | 15 | 48223 | 48418 |
| 8534 | I | 65 | 5735 | 7590 |
| 8535 | I | 66 | 115857 | 84560 |
| 8536 | I | 44 | 1149989 | 1035595 |
| 8538 | I | 16 | 2075 | 3864 |
| 8539 | I | 48 | 14193 | 22554 |
| 8571 | I | 34 | 593782 | 665066 |
| 8572 | I | 41 | 32733 | 44585 |
| 8573 | I | 7 | 2534 | 6833 / 1684 |
| 8574 | I | 40 | 1023518 | 941434 |
| 8575 | I | 63 | 1877 | 2146 |
| 8576 | I | 52 | 191167 | 203506 |
| 8577 | I | 22 | 157651 | 114192 |
| 8578 | I | 37 | 130081 | 63914 |
| 8579 | I | 43 | 97608 | 99827 |
| 8531 | B | 72 | 43805 | 13634 |
| 8533 | B | 17 | 802 | 364 |
| 8535 | B | 30 | 571 | 572 |
| 8539 | B | 36 | 1033 | 1237 |
| 8577 | B | 2 | 5976 | 7695 |
| 8578 | B | 46 | 361 | 378 |

Table C.3 Raw Data Returned By Laboratory 1 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 1 | |
|-----------|--------|-------------|-----------|-----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 15023 | 16009 |
| 8513 | U | 31 | 9931242 | 11409380 |
| 8514 | U | 56 | 56912073* | 63534172* |
| 8515 | U | 42 | 7185 | 10605 |
| 8516 | U | 25 | 701000 | 358609 |
| 8517 | U | 58 | 266870 | 213119 |
| 8518 | U | 54 | 154833 | 149483 |
| 8519 | U | 39 | 16730 | 15109 |
| 8520 | U | 68 | 11567 | 20614 |
| 8521 | U | 3 | 12439 | 13845 |
| 8523 | U | 12 | 126552 | 109817 |
| 8525 | U | 1 | 13264 | 29187 |
| 8526 | U | 59 | 10136 | 11195 |
| 8528 | U | 10 | 37391569* | 35153118* |
| 8529 | U | 11 | 41986 | 29839 |
| 8530 | U | 23 | 1098 | 875 |
| 8531 | U | 29 | 119518 | 156964 |
| 8532 | U | 27 | 37817 | 28825 |
| 8533 | U | 14 | 28043 | 31168 |
| 8534 | U | 4 | 3454 | 4729 |
| 8535 | U | 53 | 46042 | 72711 |
| 8536 | U | 19 | 457399 | 464508 |
| 8538 | U | 61 | 2669 | 5397 |
| 8539 | U | 20 | 16538 | 10474 |
| 8571 | U | 9 | 62632 | 67571 |
| 8572 | U | 6 | 11870 | 15930 |
| 8573 | U | 60 | 1745 | 2281 |
| 8574 | U | 70 | 618897 | 696057 |
| 8575 | U | 71 | 2870 | 2656 |
| 8576 | U | 13 | 130508 | 125174 |
| 8577 | U | 55 | 342980 | 319177 |
| 8578 | U | 28 | 47965 | 41214 |
| 8579 | U | 57 | 47285 | 31911 |

* participant has identified but not corrected overflows

Table C.4 Raw Data Returned By Laboratory 1 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 1 | |
|-----------|--------|-------------|-----------|------------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 23954 | 23175 |
| 8513 | I | 26 | 40961101* | 35408813* |
| 8514 | I | 64 | 88375230* | 82523750** |
| 8515 | I | 35 | 9680 | 8226 |
| 8516 | I | 18 | 619880 | 765705 |
| 8517 | I | 5 | 327568 | 313269 |
| 8518 | I | 33 | 200509 | 207198 |
| 8519 | I | 24 | 28836 | 35852 |
| 8520 | I | 69 | 13881 | 23372 |
| 8521 | I | 38 | 19575 | 17509 |
| 8523 | I | 32 | 216782 | 164077 |
| 8525 | I | 8 | 23883 | 43121 |
| 8526 | I | 50 | 12961 | 15644 |
| 8528 | I | 62 | 36146369* | 39914149* |
| 8529 | I | 51 | 86807 | 200646 |
| 8530 | I | 21 | 1737 | 1112 |
| 8531 | I | 45 | 216683 | 313006 |
| 8532 | I | 49 | 79474 | 104691 |
| 8533 | I | 15 | 45818 | 40928 |
| 8534 | I | 65 | 6641 | 8054 |
| 8535 | I | 66 | 84442 | 89264 |
| 8536 | I | 44 | 611404 | 870954 |
| 8538 | I | 16 | 3332 | 3129 |
| 8539 | I | 48 | 23912 | 23328 |
| 8571 | I | 34 | 485860 | 517646 |
| 8572 | I | 41 | 23980 | 34659 |
| 8573 | I | 7 | 1183 | 1718 |
| 8574 | I | 40 | 1091575 | 1084291 |
| 8575 | I | 63 | 3050 | 5430 |
| 8576 | I | 52 | 164187 | 173440 |
| 8577 | I | 22 | 155204 | 118657 |
| 8578 | I | 37 | 114151 | 69101 |
| 8579 | I | 43 | 85701 | 111550 |
| 8531 | B | 72 | 149224 | 168478 |
| 8533 | B | 17 | 28150 | 27213 |
| 8535 | B | 30 | 53318 | 61713 |
| 8539 | B | 36 | 8569 | 10435 |
| 8577 | B | 2 | 368419 | 437021 |
| 8578 | B | 46 | 46881 | 103791 |

* participant has identified but not corrected overflows

Table C.5 Raw Data Returned By Laboratory 2 – screening of unirradiated samples

| SP number | Status | Code number | Lab 2 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 288 | 319 |
| 8513 | U | 31 | 220 | 368 |
| 8514 | U | 56 | 873 | 1060 |
| 8515 | U | 42 | 258 | 297 |
| 8516 | U | 25 | 1129 | 1434 |
| 8517 | U | 58 | 487 | 629 |
| 8518 | U | 54 | 377 | 394 |
| 8519 | U | 39 | 359 | 264 |
| 8520 | U | 68 | 369 | 323 |
| 8521 | U | 3 | 320 | 291 |
| 8523 | U | 12 | 387 | 372 |
| 8525 | U | 1 | 1397 | 763 |
| 8526 | U | 59 | 388 | 532 |
| 8528 | U | 10 | 374 | 608 |
| 8529 | U | 11 | 367 | 331 |
| 8530 | U | 23 | 283 | 300 |
| 8531 | U | 29 | 1216 | 1031 |
| 8532 | U | 27 | 455 | 380 |
| 8533 | U | 14 | 362 | 339 |
| 8534 | U | 4 | 287 | 363 |
| 8535 | U | 53 | 442 | 515 |
| 8536 | U | 19 | 3829 | 1759 |
| 8538 | U | 61 | 924 | 1024 |
| 8539 | U | 20 | 354 | 356 |
| 8571 | U | 9 | 376 | 394 |
| 8572 | U | 6 | 363 | 405 |
| 8573 | U | 60 | 2319 | 2777 |
| 8574 | U | 70 | 1568 | 894 |
| 8575 | U | 71 | 222 | 276 |
| 8576 | U | 13 | 433 | 542 |
| 8577 | U | 55 | 2285 | 2701 |
| 8578 | U | 28 | 553 | 1191 |
| 8579 | U | 57 | 382 | 328 |

Non-negative results shaded

Table C.6 Raw Data Returned By Laboratory 2 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 2 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 28655 | 21669 |
| 8513 | I | 26 | 12153397 | 1323394 |
| 8514 | I | 64 | 6992832 | 15041014 |
| 8515 | I | 35 | 9536 | 13448 |
| 8516 | I | 18 | 217089 | 359230 |
| 8517 | I | 5 | 310442 | 370979 |
| 8518 | I | 33 | 167288 | 227632 |
| 8519 | I | 24 | 33212 | 44696 |
| 8520 | I | 69 | 6338 | 15280 |
| 8521 | I | 38 | 60398 | 19003 |
| 8523 | I | 32 | 142492 | 194092 |
| 8525 | I | 8 | 60332 | 24770 |
| 8526 | I | 50 | 10415 | 21349 |
| 8528 | I | 62 | 2554108 | 16658787 |
| 8529 | I | 51 | 59218 | 71593 |
| 8530 | I | 21 | 1009 | 1409 |
| 8531 | I | 45 | 138189 | 168069 |
| 8532 | I | 49 | 61799 | 77150 |
| 8533 | I | 15 | 44840 | 47362 |
| 8534 | I | 65 | 16186 | 13087 |
| 8535 | I | 66 | 76481 | 97465 |
| 8536 | I | 44 | 583379 | 606412 |
| 8538 | I | 16 | 1839 | 2947 |
| 8539 | I | 48 | 40015 | 17109 |
| 8571 | I | 34 | 612883 | 694439 |
| 8572 | I | 41 | 30806 | 23109 |
| 8573 | I | 7 | 1926 | 2143 |
| 8574 | I | 40 | 856570 | 1024047 |
| 8575 | I | 63 | 3161 | 11340 |
| 8576 | I | 52 | 185977 | 164907 |
| 8577 | I | 22 | 107139 | 126581 |
| 8578 | I | 37 | 80555 | 69986 |
| 8579 | I | 43 | 93250 | 108730 |
| 8531 | B | 72 | 10308 | 10771 |
| 8533 | B | 17 | 321 | 395 |
| 8535 | B | 30 | 659 | 704 |
| 8539 | B | 36 | 1627 | 1189 |
| 8577 | B | 2 | 25846 | 9547 |
| 8578 | B | 46 | 459 | 609 |

Table C.7 Raw Data Returned By Laboratory 3 – screening of unirradiated samples

| SP number | Status | Code number | Lab 3 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 326 | 210 |
| 8513 | U | 31 | 291 | 166 |
| 8514 | U | 56 | 1169 | 1437 |
| 8515 | U | 42 | 270 | 154 |
| 8516 | U | 25 | 1484 | 833 |
| 8517 | U | 58 | 502 | 516 |
| 8518 | U | 54 | 225 | 322 |
| 8519 | U | 39 | 176 | 431 |
| 8520 | U | 68 | 640 | 205 |
| 8521 | U | 3 | 287 | 352 |
| 8523 | U | 12 | 356 | 387 |
| 8525 | U | 1 | 366 | 611 |
| 8526 | U | 59 | 320 | 338 |
| 8528 | U | 10 | 478 | 971 |
| 8529 | U | 11 | 885 | 417 |
| 8530 | U | 23 | 80 | 276 |
| 8531 | U | 29 | 610 | 562 |
| 8532 | U | 27 | 143 | 407 |
| 8533 | U | 14 | 300 | 353 |
| 8534 | U | 4 | 380 | 324 |
| 8535 | U | 53 | 685 | 1079 |
| 8536 | U | 19 | 974 | 2435 |
| 8538 | U | 61 | 327 | 297 |
| 8539 | U | 20 | 350 | 526 |
| 8571 | U | 9 | 392 | 410 |
| 8572 | U | 6 | 525 | 344 |
| 8573 | U | 60 | 304 | 436 |
| 8574 | U | 70 | 1395 | 2066 |
| 8575 | U | 71 | 311 | 267 |
| 8576 | U | 13 | 508 | 705 |
| 8577 | U | 55 | 3361 | 3671 |
| 8578 | U | 28 | 531 | 764 |
| 8579 | U | 57 | 314 | 382 |

Non-negative results shaded

Table C.8 Raw Data Returned By Laboratory 3 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 3 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 27255 | 25912 |
| 8513 | I | 26 | 9425926 | 2342188 |
| 8514 | I | 64 | 1873048 | 2351920 |
| 8515 | I | 35 | 27921 | 42858 |
| 8516 | I | 18 | 1454923 | 1589263 |
| 8517 | I | 5 | 901051 | 894510 |
| 8518 | I | 33 | 343884 | 222553 |
| 8519 | I | 24 | 57661 | 33253 |
| 8520 | I | 69 | 16675 | 13576 |
| 8521 | I | 38 | 31290 | 27368 |
| 8523 | I | 32 | 277585 | 270206 |
| 8525 | I | 8 | 21630 | 34545 |
| 8526 | I | 50 | 42076 | 20805 |
| 8528 | I | 62 | 12481875 | 10124928 |
| 8529 | I | 51 | 110584 | 86368 |
| 8530 | I | 21 | 1768 | 1922 |
| 8531 | I | 45 | 204514 | 270074 |
| 8532 | I | 49 | 85808 | 72769 |
| 8533 | I | 15 | 64793 | 74555 |
| 8534 | I | 65 | 7487 | 7110 |
| 8535 | I | 66 | 118457 | 213128 |
| 8536 | I | 44 | 1365643 | 981635 |
| 8538 | I | 16 | 11466 | 6276 |
| 8539 | I | 48 | 70740 | 45639 |
| 8571 | I | 34 | 553473 | 552686 |
| 8572 | I | 41 | 44636 | 49004 |
| 8573 | I | 7 | 5693 | 2569 |
| 8574 | I | 40 | 1334573 | 1516933 |
| 8575 | I | 63 | 3081 | 3370 |
| 8576 | I | 52 | 308715 | 328936 |
| 8577 | I | 22 | 148935 | 189186 |
| 8578 | I | 37 | 122117 | 115421 |
| 8579 | I | 43 | 158702 | 150525 |
| 8531 | B | 72 | 21479 | 12221 |
| 8533 | B | 17 | 2257 | 1271 |
| 8535 | B | 30 | 693 | 731 |
| 8539 | B | 36 | 19730 | 2141 |
| 8577 | B | 2 | 6866 | 11234 |
| 8578 | B | 46 | 277 | 488 |

Table C.9 Raw Data Returned By Laboratory 5 – screening of unirradiated samples

| SP number | Status | Code number | Lab 5 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 440 | 336 |
| 8513 | U | 31 | 536 | 511 |
| 8514 | U | 56 | 1641 | 2490 |
| 8515 | U | 42 | 438 | 375 |
| 8516 | U | 25 | 2364 | 903 |
| 8517 | U | 58 | 196 | 522 |
| 8518 | U | 54 | 412 | 438 |
| 8519 | U | 39 | 313 | 257 |
| 8520 | U | 68 | 226 | 374 |
| 8521 | U | 3 | 382 | 313 |
| 8523 | U | 12 | 327 | 199 |
| 8525 | U | 1 | 903 | 1944 |
| 8526 | U | 59 | 158 | 220 |
| 8528 | U | 10 | 278 | 388 |
| 8529 | U | 11 | 338 | 565 |
| 8530 | U | 23 | 285 | 366 |
| 8531 | U | 29 | 804 | 522 |
| 8532 | U | 27 | 426 | 414 |
| 8533 | U | 14 | 301 | 663 |
| 8534 | U | 4 | 310 | 354 |
| 8535 | U | 53 | 979 | 1543 |
| 8536 | U | 19 | 1816 | 10777 |
| 8538 | U | 61 | 389 | 454 |
| 8539 | U | 20 | 334 | 469 |
| 8571 | U | 9 | 653 | 534 |
| 8572 | U | 6 | 366 | 480 |
| 8573 | U | 60 | 385 | 299 |
| 8574 | U | 70 | 2123 | 1815 |
| 8575 | U | 71 | 428 | 392 |
| 8576 | U | 13 | 518 | 518 |
| 8577 | U | 55 | 7001 | 4630 |
| 8578 | U | 28 | 447 | 541 |
| 8579 | U | 57 | 339 | 505 |

Table C.10 Raw Data Returned By Laboratory 5 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 5 | |
|-----------|--------|-------------|---------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 36351 | 36529 |
| 8513 | I | 26 | 2857873 | 5152301 |
| 8514 | I | 64 | 9583298 | 514910 |
| 8515 | I | 35 | 30233 | 16228 |
| 8516 | I | 18 | 457332 | 488801 |
| 8517 | I | 5 | 660384 | 720439 |
| 8518 | I | 33 | 298382 | 377938 |
| 8519 | I | 24 | 67050 | 106411 |
| 8520 | I | 69 | 64881 | 27921 |
| 8521 | I | 38 | 43376 | 51329 |
| 8523 | I | 32 | 459919 | 364415 |
| 8525 | I | 8 | 34091 | 47218 |
| 8526 | I | 50 | 27492 | 45917 |
| 8528 | I | 62 | 4517074 | 4086563 |
| 8529 | I | 51 | 181758 | 135146 |
| 8530 | I | 21 | 2470 | 2790 |
| 8531 | I | 45 | 306196 | 541091 |
| 8532 | I | 49 | 150867 | 180809 |
| 8533 | I | 15 | 84503 | 76808 |
| 8534 | I | 65 | 7290 | 9101 |
| 8535 | I | 66 | 192335 | 211131 |
| 8536 | I | 44 | 2063093 | 1227166 |
| 8538 | I | 16 | 3502 | 4822 |
| 8539 | I | 48 | 16466 | 27692 |
| 8571 | I | 34 | 1343951 | 1381383 |
| 8572 | I | 41 | 50109 | 37345 |
| 8573 | I | 7 | 5194 | 23739 |
| 8574 | I | 40 | 1673935 | 1612547 |
| 8575 | I | 63 | 13577 | 5299 |
| 8576 | I | 52 | 434031 | 331636 |
| 8577 | I | 22 | 296764 | 181895 |
| 8578 | I | 37 | 127073 | 147275 |
| 8579 | I | 43 | 135809 | 140487 |
| 8531 | B | 72 | 27648 | 17454 |
| 8533 | B | 17 | 1255 | 727 |
| 8535 | B | 30 | 2858 | 774 |
| 8539 | B | 36 | 1705 | 2038 |
| 8577 | B | 2 | 15699 | 14003 |
| 8578 | B | 46 | 515 | 634 |

Table C.11 Raw Data Returned By Laboratory 6 – screening of unirradiated samples

| SP number | Status | Code number | Lab 6 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 393 | 295 |
| 8513 | U | 31 | 2007 | 1612 |
| 8514 | U | 56 | 2137 | 2219 |
| 8515 | U | 42 | 1835 | 1535 |
| 8516 | U | 25 | 3231 | 1640 |
| 8517 | U | 58 | 1242 | 970 |
| 8518 | U | 54 | 1012 | 835 |
| 8519 | U | 39 | 1551 | 1249 |
| 8520 | U | 68 | 611 | 544 |
| 8521 | U | 3 | 576 | 328 |
| 8523 | U | 12 | 450 | 555 |
| 8525 | U | 1 | 697 | 3297 |
| 8526 | U | 59 | 1442 | 960 |
| 8528 | U | 10 | 558 | 522 |
| 8529 | U | 11 | 563 | 508 |
| 8530 | U | 23 | 439 | 360 |
| 8531 | U | 29 | 2512 | 4446 |
| 8532 | U | 27 | 2265 | 1954 |
| 8533 | U | 14 | 399 | 424 |
| 8534 | U | 4 | 421 | 335 |
| 8535 | U | 53 | 1418 | 1206 |
| 8536 | U | 19 | 1773 | 1687 |
| 8538 | U | 61 | 431 | 605 |
| 8539 | U | 20 | 528 | 621 |
| 8571 | U | 9 | 744 | 537 |
| 8572 | U | 6 | 614 | 704 |
| 8573 | U | 60 | 716 | 479 |
| 8574 | U | 70 | 1853 | 1750 |
| 8575 | U | 71 | 544 | 422 |
| 8576 | U | 13 | 534 | 807 |
| 8577 | U | 55 | 5499 | 3451 |
| 8578 | U | 28 | 2297 | 2539 |
| 8579 | U | 57 | 925 | 921 |

Non-negative results shaded

Table C.12 Raw Data Returned By Laboratory 6 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 6 | |
|-----------|--------|-------------|----------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 39024 | 42940 |
| 8513 | I | 26 | 10258283 | 5097008 |
| 8514 | I | 64 | 11523933 | 1435797 |
| 8515 | I | 35 | 13905 | 19107 |
| 8516 | I | 18 | 974182 | 920214 |
| 8517 | I | 5 | 522103 | 665723 |
| 8518 | I | 33 | 645083 | 630362 |
| 8519 | I | 24 | 50417 | 74159 |
| 8520 | I | 69 | 16983 | 32513 |
| 8521 | I | 38 | 42587 | 57013 |
| 8523 | I | 32 | 383866 | 669992 |
| 8525 | I | 8 | 28146 | 38664 |
| 8526 | I | 50 | 27026 | 23443 |
| 8528 | I | 62 | 6690026 | 4537385 |
| 8529 | I | 51 | 130386 | 76019 |
| 8530 | I | 21 | 1766 | 821 |
| 8531 | I | 45 | 342272 | 459322 |
| 8532 | I | 49 | 122730 | 114389 |
| 8533 | I | 15 | 80574 | 89546 |
| 8534 | I | 65 | 14004 | 11546 |
| 8535 | I | 66 | 166956 | 206351 |
| 8536 | I | 44 | 1606337 | 1282496 |
| 8538 | I | 16 | 13138 | 5985 |
| 8539 | I | 48 | 44503 | 50762 |
| 8571 | I | 34 | 957256 | 1174846 |
| 8572 | I | 41 | 71131 | 88111 |
| 8573 | I | 7 | 6400 | 3842 |
| 8574 | I | 40 | 1715013 | 2153031 |
| 8575 | I | 63 | 46430 | 6879 |
| 8576 | I | 52 | 312142 | 306764 |
| 8577 | I | 22 | 251074 | 187240 |
| 8578 | I | 37 | 116852 | 134786 |
| 8579 | I | 43 | 171389 | 168463 |
| 8531 | B | 72 | 263840 | 19647 |
| 8533 | B | 17 | 1018 | 1721 |
| 8535 | B | 30 | 3271 | 2365 |
| 8539 | B | 36 | 4349 | 16388 |
| 8577 | B | 2 | 16377 | 13989 |
| 8578 | B | 46 | 1622 | 1412 |

Table C.13 Raw Data Returned By Laboratory 8 – screening of unirradiated samples

| SP number | Status | Code number | Lab 8 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 287 | 472 |
| 8513 | U | 31 | 513 | 359 |
| 8514 | U | 56 | 1439 | 1526 |
| 8515 | U | 42 | 228 | 361 |
| 8516 | U | 25 | 2188 | 1074 |
| 8517 | U | 58 | 372 | 575 |
| 8518 | U | 54 | 465 | 318 |
| 8519 | U | 39 | 456 | 347 |
| 8520 | U | 68 | 513 | 386 |
| 8521 | U | 3 | 285 | 417 |
| 8523 | U | 12 | 399 | 262 |
| 8525 | U | 1 | 1505 | 895 |
| 8526 | U | 59 | 249 | 265 |
| 8528 | U | 10 | 515 | 491 |
| 8529 | U | 11 | 387 | 437 |
| 8530 | U | 23 | 290 | 251 |
| 8531 | U | 29 | 1055 | 971 |
| 8532 | U | 27 | 543 | 427 |
| 8533 | U | 14 | 278 | 398 |
| 8534 | U | 4 | 339 | 360 |
| 8535 | U | 53 | 926 | 1021 |
| 8536 | U | 19 | 2356 | 4926 |
| 8538 | U | 61 | 315 | 403 |
| 8539 | U | 20 | 401 | 453 |
| 8571 | U | 9 | 385 | 580 |
| 8572 | U | 6 | 591 | 572 |
| 8573 | U | 60 | 342 | 216 |
| 8574 | U | 70 | 2394 | 1922 |
| 8575 | U | 71 | 220 | 277 |
| 8576 | U | 13 | 701 | 517 |
| 8577 | U | 55 | 4762 | 6975 |
| 8578 | U | 28 | 577 | 459 |
| 8579 | U | 57 | 419 | 486 |

Table C.14 Raw Data Returned By Laboratory 8 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 8 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 38255 | 66437 |
| 8513 | I | 26 | 5013003 | 5777319 |
| 8514 | I | 64 | 9005901 | 11765444 |
| 8515 | I | 35 | 40121 | 63660 |
| 8516 | I | 18 | 441489 | 487474 |
| 8517 | I | 5 | 617732 | 612078 |
| 8518 | I | 33 | 460973 | 340141 |
| 8519 | I | 24 | 96630 | 97581 |
| 8520 | I | 69 | 26495 | 21684 |
| 8521 | I | 38 | 38005 | 32330 |
| 8523 | I | 32 | 283332 | 627800 |
| 8525 | I | 8 | 33015 | 61945 |
| 8526 | I | 50 | 50372 | 31481 |
| 8528 | I | 62 | 2256229 | 2982021 |
| 8529 | I | 51 | 113874 | 113438 |
| 8530 | I | 21 | 2481 | 2360 |
| 8531 | I | 45 | 627485 | 415256 |
| 8532 | I | 49 | 113881 | 135969 |
| 8533 | I | 15 | 72593 | 80428 |
| 8534 | I | 65 | 11781 | 15579 |
| 8535 | I | 66 | 180232 | 160485 |
| 8536 | I | 44 | 1560909 | 2115277 |
| 8538 | I | 16 | 4432 | 6578 |
| 8539 | I | 48 | 56541 | 36352 |
| 8571 | I | 34 | 1305640 | 1510217 |
| 8572 | I | 41 | 71578 | 74609 |
| 8573 | I | 7 | 5032 | 3614 |
| 8574 | I | 40 | 1888955 | 2080817 |
| 8575 | I | 63 | 4168 | 16497 |
| 8576 | I | 52 | 430785 | 347337 |
| 8577 | I | 22 | 199962 | 404080 |
| 8578 | I | 37 | 121127 | 171051 |
| 8579 | I | 43 | 379977 | 246661 |
| 8531 | B | 72 | 22529 | 22544 |
| 8533 | B | 17 | 1333 | 1417 |
| 8535 | B | 30 | 811 | 1453 |
| 8539 | B | 36 | 3088 | 1500 |
| 8577 | B | 2 | 19498 | 33047 |
| 8578 | B | 46 | 691 | 699 |

Table C.15 Raw Data Returned By Laboratory 9 – screening of unirradiated samples

| SP number | Status | Code number | Lab 9 | |
|-----------|--------|-------------|-------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 334 | 353 |
| 8513 | U | 31 | 151 | 535 |
| 8514 | U | 56 | 2320 | 1706 |
| 8515 | U | 42 | 239 | 251 |
| 8516 | U | 25 | 1468 | 1780 |
| 8517 | U | 58 | <0 | 61 |
| 8518 | U | 54 | 322 | 408 |
| 8519 | U | 39 | 181 | 306 |
| 8520 | U | 68 | 386 | 430 |
| 8521 | U | 3 | 221 | 316 |
| 8523 | U | 12 | 307 | 229 |
| 8525 | U | 1 | 752 | 607 |
| 8526 | U | 59 | 447 | 217 |
| 8528 | U | 10 | 338 | 503 |
| 8529 | U | 11 | 309 | 250 |
| 8530 | U | 23 | 162 | 262 |
| 8531 | U | 29 | 547 | 438 |
| 8532 | U | 27 | 708 | 708 |
| 8533 | U | 14 | 300 | 371 |
| 8534 | U | 4 | 386 | 378 |
| 8535 | U | 53 | 244 | 656 |
| 8536 | U | 19 | 3357 | 1307 |
| 8538 | U | 61 | 317 | 298 |
| 8539 | U | 20 | 191 | 472 |
| 8571 | U | 9 | 460 | 256 |
| 8572 | U | 6 | 405 | 314 |
| 8573 | U | 60 | 331 | 325 |
| 8574 | U | 70 | 1010 | 1132 |
| 8575 | U | 71 | 284 | 356 |
| 8576 | U | 13 | 437 | 478 |
| 8577 | U | 55 | 2983 | 2609 |
| 8578 | U | 28 | 495 | 304 |
| 8579 | U | 57 | 199 | 335 |

Non-negative results shaded

Table C.16 Raw Data Returned By Laboratory 9 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 9 | |
|-----------|--------|-------------|----------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 27287 | 28794 |
| 8513 | I | 26 | 16746256 | 1851754 |
| 8514 | I | 64 | 2715956 | 5639558 |
| 8515 | I | 35 | 12510 | 14795 |
| 8516 | I | 18 | 690769 | 1109477 |
| 8517 | I | 5 | 357292 | 380636 |
| 8518 | I | 33 | 148278 | 194996 |
| 8519 | I | 24 | 32684 | 28622 |
| 8520 | I | 69 | 28114 | 23031 |
| 8521 | I | 38 | 25680 | 24876 |
| 8523 | I | 32 | 310872 | 196710 |
| 8525 | I | 8 | 29494 | 19700 |
| 8526 | I | 50 | 28601 | 19714 |
| 8528 | I | 62 | 16602262 | 2675379 |
| 8529 | I | 51 | 102620 | 90064 |
| 8530 | I | 21 | 1892 | 1787 |
| 8531 | I | 45 | 236624 | 207396 |
| 8532 | I | 49 | 70549 | 72271 |
| 8533 | I | 15 | 52791 | 51415 |
| 8534 | I | 65 | 9187 | 8531 |
| 8535 | I | 66 | 111036 | 122731 |
| 8536 | I | 44 | 633165 | 495479 |
| 8538 | I | 16 | 2840 | 2725 |
| 8539 | I | 48 | 42484 | 24962 |
| 8571 | I | 34 | 599227 | 794649 |
| 8572 | I | 41 | 25237 | 33911 |
| 8573 | I | 7 | 36889 | 9241 |
| 8574 | I | 40 | 1050054 | 1074282 |
| 8575 | I | 63 | 7602 | 2783 |
| 8576 | I | 52 | 187041 | 214453 |
| 8577 | I | 22 | 159955 | 165979 |
| 8578 | I | 37 | 69955 | 78913 |
| 8579 | I | 43 | 118905 | 107451 |
| 8531 | B | 72 | 12773 | 11780 |
| 8533 | B | 17 | 621 | 780 |
| 8535 | B | 30 | 731 | 475 |
| 8539 | B | 36 | 699 | 818 |
| 8577 | B | 2 | 8604 | 12463 |
| 8578 | B | 46 | 549 | 561 |

Table C.17 Raw Data Returned By Laboratory 11 – screening of un irradiated samples

| SP number | Status | Code number | Lab 11 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 379 | 388 |
| 8513 | U | 31 | 478 | 434 |
| 8514 | U | 56 | 1271 | 1387 |
| 8515 | U | 42 | 393 | 409 |
| 8516 | U | 25 | 1092 | 2330 |
| 8517 | U | 58 | 409 | 520 |
| 8518 | U | 54 | 476 | 407 |
| 8519 | U | 39 | 380 | 359 |
| 8520 | U | 68 | 359 | 580 |
| 8521 | U | 3 | 376 | 390 |
| 8523 | U | 12 | 401 | 401 |
| 8525 | U | 1 | 713 | 961 |
| 8526 | U | 59 | 420 | 332 |
| 8528 | U | 10 | 464 | 438 |
| 8529 | U | 11 | 538 | 374 |
| 8530 | U | 23 | 305 | 270 |
| 8531 | U | 29 | 833 | 858 |
| 8532 | U | 27 | 455 | 484 |
| 8533 | U | 14 | 313 | 376 |
| 8534 | U | 4 | 310 | 359 |
| 8535 | U | 53 | 738 | 887 |
| 8536 | U | 19 | 1898 | 14830 |
| 8538 | U | 61 | 341 | 407 |
| 8539 | U | 20 | 442 | 388 |
| 8571 | U | 9 | 515 | 475 |
| 8572 | U | 6 | 357 | 466 |
| 8573 | U | 60 | 335 | 440 |
| 8574 | U | 70 | 1777 | 1843 |
| 8575 | U | 71 | 412 | 378 |
| 8576 | U | 13 | 500 | 551 |
| 8577 | U | 55 | 3409 | 3555 |
| 8578 | U | 28 | 659 | 555 |
| 8579 | U | 57 | 349 | 378 |

Non-negative results shaded

Table C.18 Raw Data Returned By Laboratory 11 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 11 | |
|-----------|--------|-------------|----------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 92107 | 34431 |
| 8513 | I | 26 | 103113 | 1362074 |
| 8514 | I | 64 | 15247364 | 7288572 |
| 8515 | I | 35 | 15213 | 11837 |
| 8516 | I | 18 | 634327 | 268979 |
| 8517 | I | 5 | 564823 | 626304 |
| 8518 | I | 33 | 370651 | 464098 |
| 8519 | I | 24 | 80293 | 58550 |
| 8520 | I | 69 | 15089 | 47461 |
| 8521 | I | 38 | 33531 | 34743 |
| 8523 | I | 32 | 328068 | 486060 |
| 8525 | I | 8 | 51606 | 72291 |
| 8526 | I | 50 | 29625 | 31377 |
| 8528 | I | 62 | 8978821 | 2877814 |
| 8529 | I | 51 | 121618 | 103858 |
| 8530 | I | 21 | 1928 | 1163 |
| 8531 | I | 45 | 271178 | 386885 |
| 8532 | I | 49 | 98112 | 118126 |
| 8533 | I | 15 | 101792 | 77585 |
| 8534 | I | 65 | 10036 | 15885 |
| 8535 | I | 66 | 127918 | 157657 |
| 8536 | I | 44 | 1330046 | 1721133 |
| 8538 | I | 16 | 10367 | 16936 |
| 8539 | I | 48 | 38849 | 29391 |
| 8571 | I | 34 | 1270579 | 1222728 |
| 8572 | I | 41 | 41337 | 50489 |
| 8573 | I | 7 | 5489 | 2528 |
| 8574 | I | 40 | 1459218 | 1810176 |
| 8575 | I | 63 | 7655 | 9380 |
| 8576 | I | 52 | 267254 | 295124 |
| 8577 | I | 22 | 139744 | 164785 |
| 8578 | I | 37 | 127468 | 148546 |
| 8579 | I | 43 | 263425 | 178994 |
| 8531 | B | 72 | 37949 | 20481 |
| 8533 | B | 17 | 672 | 797 |
| 8535 | B | 30 | 1933 | 1042 |
| 8539 | B | 36 | 1848 | 1733 |
| 8577 | B | 2 | 17399 | 27486 |
| 8578 | B | 46 | 666 | 3719 |

Table C.19 Raw Data Returned By Laboratory 12 – screening of unirradiated samples

| SP number | Status | Code number | Lab 12 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 223 | 267 |
| 8513 | U | 31 | 220 | 403 |
| 8514 | U | 56 | 1610 | 2494 |
| 8515 | U | 42 | 299 | 287 |
| 8516 | U | 25 | 1582 | 1658 |
| 8517 | U | 58 | 436 | 331 |
| 8518 | U | 54 | 429 | 217 |
| 8519 | U | 39 | 223 | 65 |
| 8520 | U | 68 | 387 | 394 |
| 8521 | U | 3 | 345 | 290 |
| 8523 | U | 12 | 354 | 229 |
| 8525 | U | 1 | 2837 | 933 |
| 8526 | U | 59 | 360 | 755 |
| 8528 | U | 10 | 387 | 362 |
| 8529 | U | 11 | 265 | 556 |
| 8530 | U | 23 | 286 | 230 |
| 8531 | U | 29 | 268 | 436 |
| 8532 | U | 27 | 331 | 353 |
| 8533 | U | 14 | 481 | 226 |
| 8534 | U | 4 | 306 | 201 |
| 8535 | U | 53 | 996 | 1513 |
| 8536 | U | 19 | 5394 | 3171 |
| 8538 | U | 61 | 235 | 412 |
| 8539 | U | 20 | 372 | 473 |
| 8571 | U | 9 | 439 | 456 |
| 8572 | U | 6 | 275 | 164 |
| 8573 | U | 60 | 275 | 249 |
| 8574 | U | 70 | 1509 | 3166 |
| 8575 | U | 71 | 396 | 354 |
| 8576 | U | 13 | 606 | 744 |
| 8577 | U | 55 | 5175 | 3746 |
| 8578 | U | 28 | 563 | 581 |
| 8579 | U | 57 | 390 | 486 |

Non-negative results shaded

Table C.20 Raw Data Returned By Laboratory 12 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 12 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 53557 | 49887 |
| 8513 | I | 26 | 14089009 | 11001777 |
| 8514 | I | 64 | 4319319 | 16715966 |
| 8515 | I | 35 | 20271 | 25944 |
| 8516 | I | 18 | 922293 | 690070 |
| 8517 | I | 5 | 574654 | 702196 |
| 8518 | I | 33 | 310569 | 656315 |
| 8519 | I | 24 | 104397 | 116189 |
| 8520 | I | 69 | 27043 | 45739 |
| 8521 | I | 38 | 31634 | 45112 |
| 8523 | I | 32 | 577206 | 393358 |
| 8525 | I | 8 | 41980 | 58165 |
| 8526 | I | 50 | 29174 | 44932 |
| 8528 | I | 62 | 2622708 | 9937559 |
| 8529 | I | 51 | 150619 | 124756 |
| 8530 | I | 21 | 4368 | 1852 |
| 8531 | I | 45 | 354004 | 422850 |
| 8532 | I | 49 | 368979 | 146646 |
| 8533 | I | 15 | 80117 | 85660 |
| 8534 | I | 65 | 11777 | 31280 |
| 8535 | I | 66 | 163716 | 163240 |
| 8536 | I | 44 | 1284868 | 1181836 |
| 8538 | I | 16 | 5127 | 3291 |
| 8539 | I | 48 | 47738 | 65751 |
| 8571 | I | 34 | 1125758 | 1057626 |
| 8572 | I | 41 | 49909 | 41860 |
| 8573 | I | 7 | 3035 | 4102 |
| 8574 | I | 40 | 1988114 | 1804061 |
| 8575 | I | 63 | 7248 | 5506 |
| 8576 | I | 52 | 414450 | 444712 |
| 8577 | I | 22 | 218021 | 283802 |
| 8578 | I | 37 | 110446 | 104964 |
| 8579 | I | 43 | 161514 | 282455 |
| 8531 | B | 72 | 19799 | 20209 |
| 8533 | B | 17 | 1871 | 606 |
| 8535 | B | 30 | 998 | 849 |
| 8539 | B | 36 | 3840 | 2467 |
| 8577 | B | 2 | 18118 | 13101 |
| 8578 | B | 46 | 369 | 487 |

Table C.21 Raw Data Returned By Laboratory 14 – screening of unirradiated samples

| SP number | Status | Code number | Lab 14 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 198 | 365 |
| 8513 | U | 31 | 418 | 327 |
| 8514 | U | 56 | 1040 | 1828 |
| 8515 | U | 42 | 262 | 290 |
| 8516 | U | 25 | 976 | 717 |
| 8517 | U | 58 | 433 | 125 |
| 8518 | U | 54 | 352 | 311 |
| 8519 | U | 39 | 383 | 230 |
| 8520 | U | 68 | 577 | 325 |
| 8521 | U | 3 | 472 | 352 |
| 8523 | U | 12 | 302 | 218 |
| 8525 | U | 1 | 868 | 729 |
| 8526 | U | 59 | 524 | 292 |
| 8528 | U | 10 | 806 | 347 |
| 8529 | U | 11 | 1174 | 349 |
| 8530 | U | 23 | 331 | 267 |
| 8531 | U | 29 | 719 | 572 |
| 8532 | U | 27 | 446 | 327 |
| 8533 | U | 14 | 237 | 337 |
| 8534 | U | 4 | 352 | 374 |
| 8535 | U | 53 | 564 | 645 |
| 8536 | U | 19 | 1342 | 1807 |
| 8538 | U | 61 | 437 | 367 |
| 8539 | U | 20 | 259 | 475 |
| 8571 | U | 9 | 1712 | 540 |
| 8572 | U | 6 | 270 | 537 |
| 8573 | U | 60 | 550 | 277 |
| 8574 | U | 70 | 1193 | 1430 |
| 8575 | U | 71 | 221 | 317 |
| 8576 | U | 13 | 470 | 392 |
| 8577 | U | 55 | 3395 | 3770 |
| 8578 | U | 28 | 554 | 299 |
| 8579 | U | 57 | 381 | 349 |

Non-negative results shaded

Table C.22 Raw Data Returned By Laboratory 14 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 14 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 32773 | 31281 |
| 8513 | I | 26 | 13358561 | 14234501 |
| 8514 | I | 64 | 1130402 | 6365294 |
| 8515 | I | 35 | 11347 | 11061 |
| 8516 | I | 18 | 455956 | 465202 |
| 8517 | I | 5 | 413321 | 408508 |
| 8518 | I | 33 | 205260 | 265050 |
| 8519 | I | 24 | 56235 | 41743 |
| 8520 | I | 69 | 12950 | 12741 |
| 8521 | I | 38 | 33784 | 25559 |
| 8523 | I | 32 | 373461 | 548651 |
| 8525 | I | 8 | 22546 | 34149 |
| 8526 | I | 50 | 32906 | 23379 |
| 8528 | I | 62 | 3757791 | 5760608 |
| 8529 | I | 51 | 75991 | 96123 |
| 8530 | I | 21 | 1513 | 1585 |
| 8531 | I | 45 | 240597 | 227576 |
| 8532 | I | 49 | 114305 | 99591 |
| 8533 | I | 15 | 66449 | 70117 |
| 8534 | I | 65 | 7413 | 8067 |
| 8535 | I | 66 | 157987 | 110003 |
| 8536 | I | 44 | 835039 | 844666 |
| 8538 | I | 16 | 7243 | 2753 |
| 8539 | I | 48 | 51948 | 31314 |
| 8571 | I | 34 | 864916 | 655695 |
| 8572 | I | 41 | 16911 | 35983 |
| 8573 | I | 7 | 2842 | 8538 |
| 8574 | I | 40 | 1242739 | 1106701 |
| 8575 | I | 63 | 3973 | 4928 |
| 8576 | I | 52 | 306522 | 281536 |
| 8577 | I | 22 | 177387 | 172988 |
| 8578 | I | 37 | 75930 | 162278 |
| 8579 | I | 43 | 153095 | 133370 |
| 8531 | B | 72 | 51402 | 34795 |
| 8533 | B | 17 | 1055 | 1816 |
| 8535 | B | 30 | 589 | 481 |
| 8539 | B | 36 | 4500 | 5860 |
| 8577 | B | 2 | 7851 | 27144 |
| 8578 | B | 46 | 403 | 502 |

Table C.23 Raw Data Returned By Laboratory 14 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 14 | |
|-----------|--------|-------------|----------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 3414 | 4554 |
| 8513 | U | 31 | 3751559 | 3600757 |
| 8514 | U | 56 | 12955295 | 6848282 |
| 8515 | U | 42 | 1321 | 823 |
| 8516 | U | 25 | 73654 | 69123 |
| 8517 | U | 58 | 44709 | 35881 |
| 8518 | U | 54 | 26934 | 22509 |
| 8519 | U | 39 | 4916 | 3830 |
| 8520 | U | 68 | 1015 | 3081 |
| 8521 | U | 3 | 2470 | 3288 |
| 8523 | U | 12 | 23280 | 14410 |
| 8525 | U | 1 | 2125 | 4111 |
| 8526 | U | 59 | 1544 | 3689 |
| 8528 | U | 10 | 5385095 | 3631067 |
| 8529 | U | 11 | 6899 | 6788 |
| 8530 | U | 23 | 454 | 441 |
| 8531 | U | 29 | 20114 | 33118 |
| 8532 | U | 27 | 7290 | 5013 |
| 8533 | U | 14 | 5539 | 3735 |
| 8534 | U | 4 | 1194 | 806 |
| 8535 | U | 53 | 10981 | 8047 |
| 8536 | U | 19 | 76491 | 75512 |
| 8538 | U | 61 | 657 | 643 |
| 8539 | U | 20 | 1439 | 2916 |
| 8571 | U | 9 | 14178 | 13237 |
| 8572 | U | 6 | 2352 | 2082 |
| 8573 | U | 60 | 810 | 323 |
| 8574 | U | 70 | 83300 | 81123 |
| 8575 | U | 71 | 560 | 584 |
| 8576 | U | 13 | 18566 | 23641 |
| 8577 | U | 55 | 52547 | 43508 |
| 8578 | U | 28 | 11761 | 6354 |
| 8579 | U | 57 | 7756 | 6503 |

Table C.24 Raw Data Returned By Laboratory 14 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 14 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 14865 | 15272 |
| 8513 | I | 26 | 8878492 | 7655782 |
| 8514 | I | 64 | 13467789 | 12423490 |
| 8515 | I | 35 | 3717 | 6571 |
| 8516 | I | 18 | 747328 | 266928 |
| 8517 | I | 5 | 59854 | 60715 |
| 8518 | I | 33 | 62871 | 59281 |
| 8519 | I | 24 | 27239 | 23461 |
| 8520 | I | 69 | 16711 | 10827 |
| 8521 | I | 38 | 11941 | 19098 |
| 8523 | I | 32 | 94996 | 147651 |
| 8525 | I | 8 | 3571 | 4723 |
| 8526 | I | 50 | 8348 | 7972 |
| 8528 | I | 62 | 12684636 | 6109259 |
| 8529 | I | 51 | 25273 | 35103 |
| 8530 | I | 21 | 712 | 991 |
| 8531 | I | 45 | 17267 | 24804 |
| 8532 | I | 49 | 44324 | 41338 |
| 8533 | I | 15 | 16842 | 19668 |
| 8534 | I | 65 | 1738 | 3063 |
| 8535 | I | 66 | 55848 | 58278 |
| 8536 | I | 44 | 121044 | 29730 |
| 8538 | I | 16 | 1682 | 1182 |
| 8539 | I | 48 | 14011 | 15447 |
| 8571 | I | 34 | 192636 | 145546 |
| 8572 | I | 41 | 3801 | 4575 |
| 8573 | I | 7 | 428 | 797 |
| 8574 | I | 40 | 275377 | 237089 |
| 8575 | I | 63 | 1470 | 1462 |
| 8576 | I | 52 | 116371 | 107371 |
| 8577 | I | 22 | 38033 | 32022 |
| 8578 | I | 37 | 26743 | 37787 |
| 8579 | I | 43 | 12617 | 15067 |
| 8531 | B | 72 | 30791 | 30446 |
| 8533 | B | 17 | 3581 | 4275 |
| 8535 | B | 30 | 11568 | 5904 |
| 8539 | B | 36 | 2741 | 2453 |
| 8577 | B | 2 | 97359 | 48740 |
| 8578 | B | 46 | 16699 | 7144 |

Table C.25 Raw Data Returned By Laboratory 16 –screening of unirradiated samples

| SP number | Status | Code number | Lab 16 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 299 | 262 |
| 8513 | U | 31 | 397 | 308 |
| 8514 | U | 56 | 1517 | 1700 |
| 8515 | U | 42 | 292 | 469 |
| 8516 | U | 25 | 2240 | 2995 |
| 8517 | U | 58 | 598 | 590 |
| 8518 | U | 54 | 321 | 245 |
| 8519 | U | 39 | 296 | 289 |
| 8520 | U | 68 | 439 | 328 |
| 8521 | U | 3 | 219 | 285 |
| 8523 | U | 12 | 346 | 389 |
| 8525 | U | 1 | 778 | 792 |
| 8526 | U | 59 | 261 | 357 |
| 8528 | U | 10 | 399 | 243 |
| 8529 | U | 11 | 398 | 484 |
| 8530 | U | 23 | 228 | 202 |
| 8531 | U | 29 | 573 | 566 |
| 8532 | U | 27 | 198 | 325 |
| 8533 | U | 14 | 245 | 258 |
| 8534 | U | 4 | 288 | 399 |
| 8535 | U | 53 | 815 | 1135 |
| 8536 | U | 19 | 4321 | 4514 |
| 8538 | U | 61 | 220 | 331 |
| 8539 | U | 20 | 284 | 249 |
| 8571 | U | 9 | 339 | 347 |
| 8572 | U | 6 | 259 | 362 |
| 8573 | U | 60 | 231 | 356 |
| 8574 | U | 70 | 1696 | 1969 |
| 8575 | U | 71 | 338 | 331 |
| 8576 | U | 13 | 435 | 530 |
| 8577 | U | 55 | 4502 | 4886 |
| 8578 | U | 28 | 341 | 391 |
| 8579 | U | 57 | 435 | 452 |

Table C.26 Raw Data Returned By Laboratory 16 –screening of irradiated samples and blends

| SP number | Status | Code number | Lab 16 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 37832 | 40453 |
| 8513 | I | 26 | 12190108 | 12558760 |
| 8514 | I | 64 | 4429161 | 6772705 |
| 8515 | I | 35 | 52134 | 44052 |
| 8516 | I | 18 | 235279 | 260766 |
| 8517 | I | 5 | 582754 | 645251 |
| 8518 | I | 33 | 479017 | 336745 |
| 8519 | I | 24 | 46825 | 52423 |
| 8520 | I | 69 | 22546 | 29499 |
| 8521 | I | 38 | 33141 | 30898 |
| 8523 | I | 32 | 348041 | 475984 |
| 8525 | I | 8 | 58070 | 61115 |
| 8526 | I | 50 | 25911 | 22566 |
| 8528 | I | 62 | 11256910 | 10990372 |
| 8529 | I | 51 | 112475 | 109397 |
| 8530 | I | 21 | 1312 | 1862 |
| 8531 | I | 45 | 319774 | 500269 |
| 8532 | I | 49 | 131856 | 133355 |
| 8533 | I | 15 | 66531 | 85899 |
| 8534 | I | 65 | 14888 | 14115 |
| 8535 | I | 66 | 160744 | 158637 |
| 8536 | I | 44 | 1375426 | 1663245 |
| 8538 | I | 16 | 5383 | 7453 |
| 8539 | I | 48 | 31024 | 34063 |
| 8571 | I | 34 | 578207 | 658764 |
| 8572 | I | 41 | 47035 | 53367 |
| 8573 | I | 7 | 3253 | 2768 |
| 8574 | I | 40 | 1842047 | 2015902 |
| 8575 | I | 63 | 3840 | 3172 |
| 8576 | I | 52 | 363010 | 341079 |
| 8577 | I | 22 | 133893 | 181461 |
| 8578 | I | 37 | 119224 | 98923 |
| 8579 | I | 43 | 247390 | 230868 |
| 8531 | B | 72 | 28279 | 20920 |
| 8533 | B | 17 | 516 | 691 |
| 8535 | B | 30 | 816 | 708 |
| 8539 | B | 36 | 1659 | 1969 |
| 8577 | B | 2 | 62417 | 62015 |
| 8578 | B | 46 | 521 | 463 |

Table C.27 Raw Data Returned By Laboratory 17 –screening of unirradiated samples

| SP number | Status | Code number | Lab 17 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 456 | 422 |
| 8513 | U | 31 | 404 | 697 |
| 8514 | U | 56 | 2014 | 1289 |
| 8515 | U | 42 | 344 | 238 |
| 8516 | U | 25 | 1197 | 1170 |
| 8517 | U | 58 | 582 | 319 |
| 8518 | U | 54 | 395 | 355 |
| 8519 | U | 39 | 325 | 396 |
| 8520 | U | 68 | 1275 | 1361 |
| 8521 | U | 3 | 407 | 343 |
| 8523 | U | 12 | 480 | 389 |
| 8525 | U | 1 | 759 | 492 |
| 8526 | U | 59 | 573 | 403 |
| 8528 | U | 10 | 585 | 320 |
| 8529 | U | 11 | 293 | 547 |
| 8530 | U | 23 | 336 | 347 |
| 8531 | U | 29 | 803 | 885 |
| 8532 | U | 27 | 279 | 757 |
| 8533 | U | 14 | 334 | 334 |
| 8534 | U | 4 | 346 | 357 |
| 8535 | U | 53 | 613 | 651 |
| 8536 | U | 19 | 1644 | 1837 |
| 8538 | U | 61 | 513 | 356 |
| 8539 | U | 20 | 381 | 440 |
| 8571 | U | 9 | 534 | 449 |
| 8572 | U | 6 | 394 | 357 |
| 8573 | U | 60 | 309 | 442 |
| 8574 | U | 70 | 1707 | 1152 |
| 8575 | U | 71 | 324 | 355 |
| 8576 | U | 13 | 1086 | 454 |
| 8577 | U | 55 | 3062 | 2901 |
| 8578 | U | 28 | 591 | 1196 |
| 8579 | U | 57 | 457 | 354 |

Non-negative results shaded

Table C.28 Raw Data Returned By Laboratory 17 –screening of irradiated samples and blends

| SP number | Status | Code number | Lab 17 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 43101 | 38425 |
| 8513 | I | 26 | 14885775 | 8297428 |
| 8514 | I | 64 | 11573165 | 12763174 |
| 8515 | I | 35 | 18182 | 20216 |
| 8516 | I | 18 | 409244 | 289787 |
| 8517 | I | 5 | 600805 | 517517 |
| 8518 | I | 33 | 278004 | 318696 |
| 8519 | I | 24 | 40968 | 40522 |
| 8520 | I | 69 | 20615 | 15170 |
| 8521 | I | 38 | 30767 | 29255 |
| 8523 | I | 32 | 373203 | 244611 |
| 8525 | I | 8 | 20178 | 76299 |
| 8526 | I | 50 | 21039 | 19537 |
| 8528 | I | 62 | 3687049 | 2238305 |
| 8529 | I | 51 | 159666 | 152249 |
| 8530 | I | 21 | 2264 | 1894 |
| 8531 | I | 45 | 314308 | 241702 |
| 8532 | I | 49 | 103749 | 108166 |
| 8533 | I | 15 | 68709 | 70954 |
| 8534 | I | 65 | 12539 | 9279 |
| 8535 | I | 66 | 122454 | 94076 |
| 8536 | I | 44 | 977140 | 989948 |
| 8538 | I | 16 | 3686 | 8761 |
| 8539 | I | 48 | 24737 | 27337 |
| 8571 | I | 34 | 775803 | 836543 |
| 8572 | I | 41 | 27164 | 25818 |
| 8573 | I | 7 | 2083 | 2171 |
| 8574 | I | 40 | 1222017 | 1153774 |
| 8575 | I | 63 | 3418 | 5340 |
| 8576 | I | 52 | 336217 | 307503 |
| 8577 | I | 22 | 185691 | 187410 |
| 8578 | I | 37 | 78928 | 82772 |
| 8579 | I | 43 | 136005 | 108082 |
| 8531 | B | 72 | 26734 | 56354 |
| 8533 | B | 17 | 2115 | 1024 |
| 8535 | B | 30 | 556 | 620 |
| 8539 | B | 36 | 1655 | 3153 |
| 8577 | B | 2 | 17701 | 11070 |
| 8578 | B | 46 | 561 | 521 |

Table C.29 Raw Data Returned By Laboratory 17 –calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 17 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 21922 | 21886 |
| 8513 | U | 31 | 355723 | 1037613 |
| 8514 | U | 56 | 14279475 | 7673281 |
| 8515 | U | 42 | 7996 | 9973 |
| 8516 | U | 25 | 250528 | 223736 |
| 8517 | U | 58 | 258126 | 293771 |
| 8518 | U | 54 | 234452 | 219798 |
| 8519 | U | 39 | 32830 | 21612 |
| 8520 | U | 68 | 10532 | 7423 |
| 8521 | U | 3 | 20087 | 16994 |
| 8523 | U | 12 | 230725 | 179810 |
| 8525 | U | 1 | 21543 | 22531 |
| 8526 | U | 59 | 11190 | 12728 |
| 8528 | U | 10 | 13741888 | 14121970 |
| 8529 | U | 11 | 51858 | 42678 |
| 8530 | U | 23 | 1224 | 1120 |
| 8531 | U | 29 | 139241 | 118772 |
| 8532 | U | 27 | 57114 | 50246 |
| 8533 | U | 14 | 42794 | 39440 |
| 8534 | U | 4 | 4722 | 4989 |
| 8535 | U | 53 | 70962 | 160183 |
| 8536 | U | 19 | 818606 | 843894 |
| 8538 | U | 61 | 2429 | 2692 |
| 8539 | U | 20 | 12634 | 16200 |
| 8571 | U | 9 | 96656 | 110085 |
| 8572 | U | 6 | 15373 | 15582 |
| 8573 | U | 60 | 2465 | 4077 |
| 8574 | U | 70 | 690094 | 710691 |
| 8575 | U | 71 | 2696 | 3072 |
| 8576 | U | 13 | 154812 | 196081 |
| 8577 | U | 55 | 443187 | 490339 |
| 8578 | U | 28 | 51312 | 66654 |
| 8579 | U | 57 | 90770 | 131240 |

Table C.30 Raw Data Returned By Laboratory 17 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 17 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 32591 | 35645 |
| 8513 | I | 26 | 14364557 | 11298692 |
| 8514 | I | 64 | 3660960 | 13620196 |
| 8515 | I | 35 | 18037 | 22421 |
| 8516 | I | 18 | 703010 | 368718 |
| 8517 | I | 5 | 508825 | 491621 |
| 8518 | I | 33 | 349158 | 314843 |
| 8519 | I | 24 | 43392 | 39630 |
| 8520 | I | 69 | 17130 | 24037 |
| 8521 | I | 38 | 30501 | 41027 |
| 8523 | I | 32 | 285549 | 224930 |
| 8525 | I | 8 | 23133 | 59476 |
| 8526 | I | 50 | 20609 | 21604 |
| 8528 | I | 62 | 10180002 | 6041691 |
| 8529 | I | 51 | 111331 | 157626 |
| 8530 | I | 21 | 3659 | 2068 |
| 8531 | I | 45 | 251614 | 201030 |
| 8532 | I | 49 | 103292 | 100880 |
| 8533 | I | 15 | 64306 | 59023 |
| 8534 | I | 65 | 9858 | 8609 |
| 8535 | I | 66 | 114560 | 88641 |
| 8536 | I | 44 | 842700 | 969748 |
| 8538 | I | 16 | 4785 | 4637 |
| 8539 | I | 48 | 21666 | 29064 |
| 8571 | I | 34 | 682208 | 815835 |
| 8572 | I | 41 | 23117 | 22941 |
| 8573 | I | 7 | 1426 | 2931 |
| 8574 | I | 40 | 1018235 | 1013676 |
| 8575 | I | 63 | 3025 | 4684 |
| 8576 | I | 52 | 324406 | 317028 |
| 8577 | I | 22 | 159233 | 180862 |
| 8578 | I | 37 | 107418 | 96582 |
| 8579 | I | 43 | 118545 | 100403 |
| 8531 | B | 72 | 146856 | 174606 |
| 8533 | B | 17 | 37008 | 41469 |
| 8535 | B | 30 | 65898 | 139647 |
| 8539 | B | 36 | 13346 | 12077 |
| 8577 | B | 2 | 449981 | 509025 |
| 8578 | B | 46 | 78125 | 43503 |

Table C.31 Raw Data Returned By Laboratory 18 – screening of unirradiated samples

| SP number | Status | Code number | Lab 18 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 606 | 648 |
| 8513 | U | 31 | 3406 | 5680 |
| 8514 | U | 56 | 1754 | 2334 |
| 8515 | U | 42 | 1057 | 934 |
| 8516 | U | 25 | 2073 | 1041 |
| 8517 | U | 58 | 900 | 522 |
| 8518 | U | 54 | 800 | 600 |
| 8519 | U | 39 | 998 | 770 |
| 8520 | U | 68 | 742 | 834 |
| 8521 | U | 3 | 363 | 377 |
| 8523 | U | 12 | 449 | 450 |
| 8525 | U | 1 | 763 | 676 |
| 8526 | U | 59 | 330 | 785 |
| 8528 | U | 10 | 4547 | 498 |
| 8529 | U | 11 | 439 | 408 |
| 8530 | U | 23 | 1470 | 564 |
| 8531 | U | 29 | 6798 | 4669 |
| 8532 | U | 27 | 10614 | 6362 |
| 8533 | U | 14 | 550 | 376 |
| 8534 | U | 4 | 406 | 463 |
| 8535 | U | 53 | 1111 | 1183 |
| 8536 | U | 19 | 6062 | 5235 |
| 8538 | U | 61 | 471 | 787 |
| 8539 | U | 20 | 554 | 638 |
| 8571 | U | 9 | 480 | 468 |
| 8572 | U | 6 | 630 | 627 |
| 8573 | U | 60 | 444 | 330 |
| 8574 | U | 70 | 1802 | 2026 |
| 8575 | U | 71 | 1656 | 441 |
| 8576 | U | 13 | 645 | 691 |
| 8577 | U | 55 | 5954 | 6659 |
| 8578 | U | 28 | 5151 | 4899 |
| 8579 | U | 57 | 757 | 845 |

Non-negative results shaded

Table C.32 Raw Data Returned By Laboratory 18 –screening of irradiated samples and blends

| SP number | Status | Code number | Lab 18 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 34651 | 46529 |
| 8513 | I | 26 | 2556449 | 11080339 |
| 8514 | I | 64 | 9644356 | 1612226 |
| 8515 | I | 35 | 16772 | 17139 |
| 8516 | I | 18 | 614815 | 434441 |
| 8517 | I | 5 | 541114 | 577184 |
| 8518 | I | 33 | 355128 | 300246 |
| 8519 | I | 24 | 44105 | 78681 |
| 8520 | I | 69 | 18883 | 19376 |
| 8521 | I | 38 | 30554 | 24573 |
| 8523 | I | 32 | 336180 | 386396 |
| 8525 | I | 8 | 34557 | 30699 |
| 8526 | I | 50 | 17269 | 84973 |
| 8528 | I | 62 | 15924188 | 3602589 |
| 8529 | I | 51 | 163417 | 95879 |
| 8530 | I | 21 | 1298 | 1947 |
| 8531 | I | 45 | 141332 | 128409 |
| 8532 | I | 49 | 99972 | 110511 |
| 8533 | I | 15 | 57574 | 69023 |
| 8534 | I | 65 | 11333 | 15834 |
| 8535 | I | 66 | 123102 | 128365 |
| 8536 | I | 44 | 1033955 | 1153979 |
| 8538 | I | 16 | 3965 | 3532 |
| 8539 | I | 48 | 22874 | 43711 |
| 8571 | I | 34 | 797903 | 732433 |
| 8572 | I | 41 | 47694 | 36289 |
| 8573 | I | 7 | 2169 | 5270 |
| 8574 | I | 40 | 496058 | 1183363 |
| 8575 | I | 63 | 7105 | 7396 |
| 8576 | I | 52 | 307923 | 266704 |
| 8577 | I | 22 | 234614 | 289431 |
| 8578 | I | 37 | 96189 | 75055 |
| 8579 | I | 43 | 311452 | 442033 |
| 8531 | B | 72 | 18400 | 17949 |
| 8533 | B | 17 | 814 | 860 |
| 8535 | B | 30 | 3176 | 3642 |
| 8539 | B | 36 | 6920 | 2672 |
| 8577 | B | 2 | 8891 | 7624 |
| 8578 | B | 46 | 1129 | 1263 |

Table C.33 Raw Data Returned By Laboratory 18 –calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 18 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 30949 | 34642 |
| 8513 | U | 31 | 7499685 | 1953902 |
| 8514 | U | 56 | 2766261 | 7730987 |
| 8515 | U | 42 | 16368 | 18455 |
| 8516 | U | 25 | 827848 | 766659 |
| 8517 | U | 58 | 309302 | 421435 |
| 8518 | U | 54 | 184209 | 257628 |
| 8519 | U | 39 | 41760 | 48480 |
| 8520 | U | 68 | 11814 | 18582 |
| 8521 | U | 3 | 20465 | 23531 |
| 8523 | U | 12 | 324630 | 193046 |
| 8525 | U | 1 | 30657 | 20642 |
| 8526 | U | 59 | 19061 | 19172 |
| 8528 | U | 10 | 5572252 | 15462266 |
| 8529 | U | 11 | 125778 | 153102 |
| 8530 | U | 23 | 3778 | 5143 |
| 8531 | U | 29 | 312012 | 139239 |
| 8532 | U | 27 | 152351 | 94043 |
| 8533 | U | 14 | 51571 | 48277 |
| 8534 | U | 4 | 8458 | 7466 |
| 8535 | U | 53 | 69545 | 83882 |
| 8536 | U | 19 | 764495 | 651559 |
| 8538 | U | 61 | 6162 | 3789 |
| 8539 | U | 20 | 28266 | 18319 |
| 8571 | U | 9 | 115950 | 119771 |
| 8572 | U | 6 | 26481 | 24957 |
| 8573 | U | 60 | 3603 | 5565 |
| 8574 | U | 70 | 809630 | 841360 |
| 8575 | U | 71 | 5304 | 6335 |
| 8576 | U | 13 | 198513 | 203921 |
| 8577 | U | 55 | 561080 | 573338 |
| 8578 | U | 28 | 92993 | 70998 |
| 8579 | U | 57 | 69862 | 68238 |

Table C.34 Raw Data Returned By Laboratory 18 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 18 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 53309 | 40590 |
| 8513 | I | 26 | 2520711 | 15000999 |
| 8514 | I | 64 | 8080384 | 9929371 |
| 8515 | I | 35 | 50446 | 40241 |
| 8516 | I | 18 | 11224436 | 1107243 |
| 8517 | I | 5 | 479267 | 530441 |
| 8518 | I | 33 | 402348 | 338781 |
| 8519 | I | 24 | 72426 | 53176 |
| 8520 | I | 69 | 63653 | 15846 |
| 8521 | I | 38 | 41827 | 48583 |
| 8523 | I | 32 | 499723 | 486317 |
| 8525 | I | 8 | 50406 | 82007 |
| 8526 | I | 50 | 22360 | 22223 |
| 8528 | I | 62 | 1231703 | 1412590 |
| 8529 | I | 51 | 88313 | 108259 |
| 8530 | I | 21 | 3763 | 35677 |
| 8531 | I | 45 | 236673 | 267181 |
| 8532 | I | 49 | 81887 | 95789 |
| 8533 | I | 15 | 65932 | 69839 |
| 8534 | I | 65 | 35145 | 33306 |
| 8535 | I | 66 | 216240 | 135457 |
| 8536 | I | 44 | 988980 | 1114127 |
| 8538 | I | 16 | 14540 | 18063 |
| 8539 | I | 48 | 34803 | 43978 |
| 8571 | I | 34 | 862857 | 1023492 |
| 8572 | I | 41 | 41386 | 55128 |
| 8573 | I | 7 | 53336 | 23223 |
| 8574 | I | 40 | 1385570 | 1466601 |
| 8575 | I | 63 | 54361 | 15287 |
| 8576 | I | 52 | 430079 | 288509 |
| 8577 | I | 22 | 201699 | 226910 |
| 8578 | I | 37 | 97082 | 108374 |
| 8579 | I | 43 | 135001 | 114013 |
| 8531 | B | 72 | 148885 | 116898 |
| 8533 | B | 17 | 44560 | 604150 |
| 8535 | B | 30 | 127813 | 96466 |
| 8539 | B | 36 | 48809 | 46300 |
| 8577 | B | 2 | 519170 | 513987 |
| 8578 | B | 46 | 66130 | 56621 |

Table C.35 Raw Data Returned By Laboratory 19 – screening of unirradiated samples

| SP number | Status | Code number | Lab 19 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 308 | 241 |
| 8513 | U | 31 | 619 | 393 |
| 8514 | U | 56 | 1465 | 1390 |
| 8515 | U | 42 | 435 | 293 |
| 8516 | U | 25 | 870 | 1342 |
| 8517 | U | 58 | 584 | 391 |
| 8518 | U | 54 | 448 | 430 |
| 8519 | U | 39 | 366 | 358 |
| 8520 | U | 68 | 500 | 446 |
| 8521 | U | 3 | 394 | 345 |
| 8523 | U | 12 | 427 | 275 |
| 8525 | U | 1 | 1560 | 1449 |
| 8526 | U | 59 | 442 | 337 |
| 8528 | U | 10 | 340 | 414 |
| 8529 | U | 11 | 1310 | 559 |
| 8530 | U | 23 | 292 | 311 |
| 8531 | U | 29 | 4247 | 997 |
| 8532 | U | 27 | 307 | 432 |
| 8533 | U | 14 | 283 | 281 |
| 8534 | U | 4 | 378 | 368 |
| 8535 | U | 53 | 727 | 816 |
| 8536 | U | 19 | 10909 | 3723 |
| 8538 | U | 61 | 433 | 443 |
| 8539 | U | 20 | 407 | 516 |
| 8571 | U | 9 | 479 | 516 |
| 8572 | U | 6 | 584 | 546 |
| 8573 | U | 60 | 262 | 350 |
| 8574 | U | 70 | 1880 | 1853 |
| 8575 | U | 71 | 317 | 301 |
| 8576 | U | 13 | 702 | 867 |
| 8577 | U | 55 | 4762 | 5029 |
| 8578 | U | 28 | 577 | 695 |
| 8579 | U | 57 | 450 | 474 |

Non-negative results shaded

Table C.36 Raw Data Returned By Laboratory 19 –screening of irradiated samples and blends

| SP number | Status | Code number | Lab 19 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 52539 | 37709 |
| 8513 | I | 26 | 16637143 | 17226362 |
| 8514 | I | 64 | 19647258 | 21815714 |
| 8515 | I | 35 | 29840 | 38275 |
| 8516 | I | 18 | 284797 | 284714 |
| 8517 | I | 5 | 694776 | 558473 |
| 8518 | I | 33 | 369839 | 433608 |
| 8519 | I | 24 | 70660 | 55980 |
| 8520 | I | 69 | 27975 | 21677 |
| 8521 | I | 38 | 32985 | 32117 |
| 8523 | I | 32 | 484311 | 435721 |
| 8525 | I | 8 | 53286 | 25129 |
| 8526 | I | 50 | 27766 | 23141 |
| 8528 | I | 62 | 11022426 | 13586093 |
| 8529 | I | 51 | 153515 | 69011 |
| 8530 | I | 21 | 3476 | 1804 |
| 8531 | I | 45 | 269072 | 281963 |
| 8532 | I | 49 | 116150 | 138858 |
| 8533 | I | 15 | 80793 | 69270 |
| 8534 | I | 65 | 12308 | 9988 |
| 8535 | I | 66 | 222613 | 161575 |
| 8536 | I | 44 | 1228023 | 1000500 |
| 8538 | I | 16 | 3544 | 7224 |
| 8539 | I | 48 | 49876 | 26336 |
| 8571 | I | 34 | 1134014 | 1145009 |
| 8572 | I | 41 | 55597 | 41331 |
| 8573 | I | 7 | 4692 | 2284 |
| 8574 | I | 40 | 1592395 | 1313447 |
| 8575 | I | 63 | 7287 | 5503 |
| 8576 | I | 52 | 347521 | 378101 |
| 8577 | I | 22 | 223357 | 224667 |
| 8578 | I | 37 | 191809 | 74759 |
| 8579 | I | 43 | 247137 | 138215 |
| 8531 | B | 72 | 51400 | 37287 |
| 8533 | B | 17 | 852 | 1544 |
| 8535 | B | 30 | 488 | 1069 |
| 8539 | B | 36 | 2604 | 1821 |
| 8577 | B | 2 | 15779 | 15935 |
| 8578 | B | 46 | 667 | 533 |

Table C.37 Raw Data Returned By Laboratory 19 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 19 | |
|-----------|--------|-------------|---------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 29910 | 20067 |
| 8513 | U | 31 | 6631594 | 5535342 |
| 8514 | U | 56 | 9992404 | 8408140 |
| 8515 | U | 42 | 16608 | 9045 |
| 8516 | U | 25 | 123928 | 131153 |
| 8517 | U | 58 | 256487 | 210764 |
| 8518 | U | 54 | 183336 | 241370 |
| 8519 | U | 39 | 19704 | 18714 |
| 8520 | U | 68 | 16210 | 6189 |
| 8521 | U | 3 | 16747 | 20614 |
| 8523 | U | 12 | 217303 | 195918 |
| 8525 | U | 1 | 51295 | 18440 |
| 8526 | U | 59 | 37325 | 8121 |
| 8528 | U | 10 | 7040103 | 6075185 |
| 8529 | U | 11 | 74783 | 81973 |
| 8530 | U | 23 | 922 | 957 |
| 8531 | U | 29 | 1097759 | 217329 |
| 8532 | U | 27 | 49780 | 45589 |
| 8533 | U | 14 | 41107 | 30378 |
| 8534 | U | 4 | 5476 | 5150 |
| 8535 | U | 53 | 71974 | 63919 |
| 8536 | U | 19 | 898811 | 599040 |
| 8538 | U | 61 | 3037 | 2542 |
| 8539 | U | 20 | 14254 | 12957 |
| 8571 | U | 9 | 91521 | 64153 |
| 8572 | U | 6 | 11944 | 18854 |
| 8573 | U | 60 | 1371 | 2071 |
| 8574 | U | 70 | 589556 | 719899 |
| 8575 | U | 71 | 3637 | 3555 |
| 8576 | U | 13 | 240896 | 174644 |
| 8577 | U | 55 | 555057 | 496069 |
| 8578 | U | 28 | 65733 | 68470 |
| 8579 | U | 57 | 53635 | 76942 |

Table C.38 Raw Data Returned By Laboratory 19 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 19 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 33331 | 35588 |
| 8513 | I | 26 | 13702222 | 10956695 |
| 8514 | I | 64 | 6012367 | 5257530 |
| 8515 | I | 35 | 15496 | 15039 |
| 8516 | I | 18 | 273258 | 254668 |
| 8517 | I | 5 | 367894 | 312092 |
| 8518 | I | 33 | 249887 | 186239 |
| 8519 | I | 24 | 46078 | 63232 |
| 8520 | I | 69 | 12498 | 24975 |
| 8521 | I | 38 | 22054 | 17798 |
| 8523 | I | 32 | 203170 | 170348 |
| 8525 | I | 8 | 18899 | 27940 |
| 8526 | I | 50 | 14846 | 21037 |
| 8528 | I | 62 | 10096556 | 10514290 |
| 8529 | I | 51 | 60112 | 104768 |
| 8530 | I | 21 | 2101 | 1428 |
| 8531 | I | 45 | 304807 | 225705 |
| 8532 | I | 49 | 90340 | 88443 |
| 8533 | I | 15 | 49569 | 47809 |
| 8534 | I | 65 | 9176 | 10762 |
| 8535 | I | 66 | 268069 | 97307 |
| 8536 | I | 44 | 832417 | 1500341 |
| 8538 | I | 16 | 4808 | 9056 |
| 8539 | I | 48 | 15216 | 17996 |
| 8571 | I | 34 | 946070 | 592174 |
| 8572 | I | 41 | 25732 | 37775 |
| 8573 | I | 7 | 3821 | 12065 |
| 8574 | I | 40 | 1308998 | 1107425 |
| 8575 | I | 63 | 10594 | 5306 |
| 8576 | I | 52 | 289368 | 266744 |
| 8577 | I | 22 | 176204 | 126784 |
| 8578 | I | 37 | 75248 | 94695 |
| 8579 | I | 43 | 110847 | 84905 |
| 8531 | B | 72 | 166594 | 109436 |
| 8533 | B | 17 | 34767 | 34968 |
| 8535 | B | 30 | 71564 | 74025 |
| 8539 | B | 36 | 12577 | 6352 |
| 8577 | B | 2 | 535661 | 514532 |
| 8578 | B | 46 | 92645 | 77801 |

Table C.39 Raw Data Returned By Laboratory 20 – screening of unirradiated samples

| SP number | Status | Code number | Lab 20 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 401 | 343 |
| 8513 | U | 31 | 452 | 386 |
| 8514 | U | 56 | 2210 | 1549 |
| 8515 | U | 42 | 330 | 329 |
| 8516 | U | 25 | 864 | 937 |
| 8517 | U | 58 | 523 | 560 |
| 8518 | U | 54 | 322 | 396 |
| 8519 | U | 39 | 273 | 286 |
| 8520 | U | 68 | 361 | 406 |
| 8521 | U | 3 | 243 | 293 |
| 8523 | U | 12 | 327 | 322 |
| 8525 | U | 1 | 599 | 663 |
| 8526 | U | 59 | 446 | 313 |
| 8528 | U | 10 | 452 | 367 |
| 8529 | U | 11 | 388 | 395 |
| 8530 | U | 23 | 294 | 213 |
| 8531 | U | 29 | 513 | 391 |
| 8532 | U | 27 | 401 | 488 |
| 8533 | U | 14 | 368 | 428 |
| 8534 | U | 4 | 280 | 274 |
| 8535 | U | 53 | 457 | 381 |
| 8536 | U | 19 | 6352 | 6125 |
| 8538 | U | 61 | 323 | 300 |
| 8539 | U | 20 | 340 | 387 |
| 8571 | U | 9 | 340 | 339 |
| 8572 | U | 6 | 334 | 381 |
| 8573 | U | 60 | 294 | 355 |
| 8574 | U | 70 | 1291 | 1500 |
| 8575 | U | 71 | 319 | 269 |
| 8576 | U | 13 | 467 | 487 |
| 8577 | U | 55 | 3126 | 2835 |
| 8578 | U | 28 | 329 | 336 |
| 8579 | U | 57 | 283 | 346 |

Non-negative results shaded

Table C.40 Raw Data Returned By Laboratory 20 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 20 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 33452 | 35478 |
| 8513 | I | 26 | 14858885 | 15377677 |
| 8514 | I | 64 | 11980728 | 10253402 |
| 8515 | I | 35 | 9739 | 8967 |
| 8516 | I | 18 | 237260 | 220833 |
| 8517 | I | 5 | 437423 | 473200 |
| 8518 | I | 33 | 239935 | 221992 |
| 8519 | I | 24 | 31968 | 34153 |
| 8520 | I | 69 | 10768 | 15284 |
| 8521 | I | 38 | 21979 | 16758 |
| 8523 | I | 32 | 255761 | 226620 |
| 8525 | I | 8 | 35856 | 27433 |
| 8526 | I | 50 | 31661 | 33866 |
| 8528 | I | 62 | 11584632 | 11834355 |
| 8529 | I | 51 | 65011 | 72857 |
| 8530 | I | 21 | 943 | 884 |
| 8531 | I | 45 | 223537 | 214373 |
| 8532 | I | 49 | 70452 | 76525 |
| 8533 | I | 15 | 56582 | 58526 |
| 8534 | I | 65 | 8788 | 7507 |
| 8535 | I | 66 | 122532 | 122175 |
| 8536 | I | 44 | 999716 | 1120749 |
| 8538 | I | 16 | 3255 | 3778 |
| 8539 | I | 48 | 18450 | 17116 |
| 8571 | I | 34 | 690187 | 740668 |
| 8572 | I | 41 | 38781 | 36009 |
| 8573 | I | 7 | 3269 | 3448 |
| 8574 | I | 40 | 1228297 | 1306115 |
| 8575 | I | 63 | 3146 | 3126 |
| 8576 | I | 52 | 216786 | 237672 |
| 8577 | I | 22 | 174435 | 154015 |
| 8578 | I | 37 | 73684 | 65787 |
| 8579 | I | 43 | 131867 | 118807 |
| 8531 | B | 72 | 10967 | 11406 |
| 8533 | B | 17 | 551 | 660 |
| 8535 | B | 30 | 681 | 573 |
| 8539 | B | 36 | 899 | 935 |
| 8577 | B | 2 | 8454 | 6161 |
| 8578 | B | 46 | 470 | 492 |

Table C.41 Raw Data Returned By Laboratory 21 – screening of unirradiated samples

| SP number | Status | Code number | Lab 21 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 249 | 394 |
| 8513 | U | 31 | 254 | 155 |
| 8514 | U | 56 | 1121 | 1164 |
| 8515 | U | 42 | 247 | 236 |
| 8516 | U | 25 | 1574 | 1188 |
| 8517 | U | 58 | 451 | 1230 |
| 8518 | U | 54 | 376 | 333 |
| 8519 | U | 39 | 255 | 192 |
| 8520 | U | 68 | 344 | 625 |
| 8521 | U | 3 | 291 | 388 |
| 8523 | U | 12 | 288 | 444 |
| 8525 | U | 1 | 362 | 831 |
| 8526 | U | 59 | 238 | 158 |
| 8528 | U | 10 | 374 | 304 |
| 8529 | U | 11 | 346 | 555 |
| 8530 | U | 23 | 172 | 446 |
| 8531 | U | 29 | 1816 | 566 |
| 8532 | U | 27 | 329 | 389 |
| 8533 | U | 14 | 359 | 155 |
| 8534 | U | 4 | 259 | 275 |
| 8535 | U | 53 | 467 | 663 |
| 8536 | U | 19 | 2316 | 1975 |
| 8538 | U | 61 | 352 | 223 |
| 8539 | U | 20 | 402 | 275 |
| 8571 | U | 9 | 387 | 510 |
| 8572 | U | 6 | 299 | 262 |
| 8573 | U | 60 | 365 | 193 |
| 8574 | U | 70 | 1561 | 1420 |
| 8575 | U | 71 | 359 | 202 |
| 8576 | U | 13 | 574 | 720 |
| 8577 | U | 55 | 2905 | 4744 |
| 8578 | U | 28 | 335 | 561 |
| 8579 | U | 57 | 412 | 322 |

Non-negative results shaded

Table C.42 Raw Data Returned By Laboratory 21 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 21 | |
|-----------|--------|-------------|----------|----------------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 30450 | 28338 |
| 8513 | I | 26 | 4908542 | 12208098 |
| 8514 | I | 64 | 8032934 | 3895753 |
| 8515 | I | 35 | 15186 | 10408 |
| 8516 | I | 18 | 663722 | 642956 |
| 8517 | I | 5 | 346549 | 437913 |
| 8518 | I | 33 | 353983 | 378227 |
| 8519 | I | 24 | 36107 | 53234 |
| 8520 | I | 69 | 21916 | 17856 |
| 8521 | I | 38 | 24478 | 27500 |
| 8523 | I | 32 | 328564 | 293756 |
| 8525 | I | 8 | 31118 | 30553 |
| 8526 | I | 50 | 14653 | 9666 |
| 8528 | I | 62 | 14120491 | 11980180 |
| 8529 | I | 51 | 81957 | 98700 |
| 8530 | I | 21 | 1640 | 815 |
| 8531 | I | 45 | 232648 | 345110 |
| 8532 | I | 49 | 94843 | 108897 |
| 8533 | I | 15 | 53871 | 83420 |
| 8534 | I | 65 | 10113 | 8233 |
| 8535 | I | 66 | 105484 | 120737 |
| 8536 | I | 44 | 445418 | 758151 |
| 8538 | I | 16 | 2828 | 3467 |
| 8539 | I | 48 | 34026 | 26960 |
| 8571 | I | 34 | 675978 | 715956 |
| 8572 | I | 41 | 28907 | 31315 |
| 8573 | I | 7 | 3459 | 2898 |
| 8574 | I | 40 | 1131817 | 1061134 |
| 8575 | I | 63 | 5688 | 4801 |
| 8576 | I | 52 | 272157 | 180950 |
| 8577 | I | 22 | 133891 | 158767 |
| 8578 | I | 37 | 75450 | 85802 |
| 8579 | I | 43 | 151155 | 223303 |
| 8531 | B | 72 | 21392 | 24240 |
| 8533 | B | 17 | 801 | 675 |
| 8535 | B | 30 | 639 | 664 |
| 8539 | B | 36 | 1486 | 1691 |
| 8577 | B | 2 | 9081 | 15508 |
| 8578 | B | 46 | 421 | 406 |

Table C.43 Raw Data Returned By Laboratory 22 – screening of unirradiated samples

| SP number | Status | Code number | Lab 22 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 675 | 13024 |
| 8513 | U | 31 | 321 | 794 |
| 8514 | U | 56 | 723 | 1035 |
| 8515 | U | 42 | 435 | 1169 |
| 8516 | U | 25 | 650 | 1174 |
| 8517 | U | 58 | 402 | 536 |
| 8518 | U | 54 | 345 | 388 |
| 8519 | U | 39 | 298 | 506 |
| 8520 | U | 68 | 1380 | 2493 |
| 8521 | U | 3 | 249 | 199 |
| 8523 | U | 12 | 285 | 319 |
| 8525 | U | 1 | 938 | 2261 |
| 8526 | U | 59 | 2144 | 346 |
| 8528 | U | 10 | 473 | 403 |
| 8529 | U | 11 | 354 | 368 |
| 8530 | U | 23 | 282 | 235 |
| 8531 | U | 29 | 3283 | 4704 |
| 8532 | U | 27 | 10258 | 10888 |
| 8533 | U | 14 | 297 | 304 |
| 8534 | U | 4 | 234 | 324 |
| 8535 | U | 53 | 645 | 777 |
| 8536 | U | 19 | 1345 | 3772 |
| 8538 | U | 61 | 282 | 555 |
| 8539 | U | 20 | 318 | 244 |
| 8571 | U | 9 | 418 | 296 |
| 8572 | U | 6 | 265 | 343 |
| 8573 | U | 60 | 366 | 166 |
| 8574 | U | 70 | 3458 | 1668 |
| 8575 | U | 71 | 419 | 324 |
| 8576 | U | 13 | 398 | 381 |
| 8577 | U | 55 | 2795 | 8135 |
| 8578 | U | 28 | 1508 | 3037 |
| 8579 | U | 57 | 932 | 273 |

Non-negative results shaded

Table C.44 Raw Data Returned By Laboratory 22 – screening of for irradiated samples and blends

| SP number | Status | Code number | Lab 22 | |
|-----------|--------|-------------|----------|---------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 33926 | 20086 |
| 8513 | I | 26 | 175500 | 3432284 |
| 8514 | I | 64 | 11215921 | 8620105 |
| 8515 | I | 35 | 9678 | 26023 |
| 8516 | I | 18 | 170810 | 282709 |
| 8517 | I | 5 | 353685 | 319248 |
| 8518 | I | 33 | 201582 | 206797 |
| 8519 | I | 24 | 28666 | 28581 |
| 8520 | I | 69 | 13319 | 21781 |
| 8521 | I | 38 | 19051 | 27895 |
| 8523 | I | 32 | 253076 | 205132 |
| 8525 | I | 8 | 28411 | 17715 |
| 8526 | I | 50 | 38824 | 23461 |
| 8528 | I | 62 | 1896464 | 4357782 |
| 8529 | I | 51 | 78621 | 91919 |
| 8530 | I | 21 | 2116 | 1269 |
| 8531 | I | 45 | 204294 | 234938 |
| 8532 | I | 49 | 84275 | 71122 |
| 8533 | I | 15 | 57089 | 59818 |
| 8534 | I | 65 | 23167 | 15118 |
| 8535 | I | 66 | 114687 | 134509 |
| 8536 | I | 44 | 707049 | 527493 |
| 8538 | I | 16 | 4432 | 4114 |
| 8539 | I | 48 | 21420 | 23216 |
| 8571 | I | 34 | 614400 | 607210 |
| 8572 | I | 41 | 27308 | 32317 |
| 8573 | I | 7 | 1740 | 6877 |
| 8574 | I | 40 | 903721 | 784854 |
| 8575 | I | 63 | 11770 | 4572 |
| 8576 | I | 52 | 223658 | 231192 |
| 8577 | I | 22 | 163068 | 112965 |
| 8578 | I | 37 | 61532 | 84639 |
| 8579 | I | 43 | 91713 | 104471 |
| 8531 | B | 72 | 38620 | 15526 |
| 8533 | B | 17 | 503 | 767 |
| 8535 | B | 30 | 1405 | 4110 |
| 8539 | B | 36 | 1887 | 1559 |
| 8577 | B | 2 | 4807 | 8295 |
| 8578 | B | 46 | 472 | 531 |

Table C.45 Raw Data Returned By Laboratory 23 – screening of unirradiated samples

| SP number | Status | Code number | Lab 23 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 519 | 379 |
| 8513 | U | 31 | 463 | 718 |
| 8514 | U | 56 | 1526 | 2330 |
| 8515 | U | 42 | 423 | 431 |
| 8516 | U | 25 | 962 | 950 |
| 8517 | U | 58 | 755 | 514 |
| 8518 | U | 54 | 701 | 455 |
| 8519 | U | 39 | 747 | 392 |
| 8520 | U | 68 | 561 | 547 |
| 8521 | U | 3 | 373 | 306 |
| 8523 | U | 12 | 410 | 376 |
| 8525 | U | 1 | 776 | 574 |
| 8526 | U | 59 | 324 | 432 |
| 8528 | U | 10 | 393 | 374 |
| 8529 | U | 11 | 469 | 453 |
| 8530 | U | 23 | 291 | 322 |
| 8531 | U | 29 | 949 | 1182 |
| 8532 | U | 27 | 448 | 429 |
| 8533 | U | 14 | 243 | 338 |
| 8534 | U | 4 | 353 | 437 |
| 8535 | U | 53 | 1076 | 3164 |
| 8536 | U | 19 | 2377 | 3129 |
| 8538 | U | 61 | 399 | 451 |
| 8539 | U | 20 | 450 | 334 |
| 8571 | U | 9 | 520 | 499 |
| 8572 | U | 6 | 419 | 555 |
| 8573 | U | 60 | 390 | 410 |
| 8574 | U | 70 | 1899 | 3537 |
| 8575 | U | 71 | 30006 | 23013 |
| 8576 | U | 13 | 519 | 527 |
| 8577 | U | 55 | 6726 | 4301 |
| 8578 | U | 28 | 612 | 560 |
| 8579 | U | 57 | 427 | 345 |

Non-negative results shaded

Table C.46 Raw Data Returned By Laboratory 23 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 23 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 45246 | 37063 |
| 8513 | I | 26 | 15896915 | 10211898 |
| 8514 | I | 64 | 4588391 | 3786666 |
| 8515 | I | 35 | 25275 | 17774 |
| 8516 | I | 18 | 257610 | 221475 |
| 8517 | I | 5 | 557752 | 655278 |
| 8518 | I | 33 | 481186 | 462532 |
| 8519 | I | 24 | 48345 | 58470 |
| 8520 | I | 69 | 32739 | 161977 |
| 8521 | I | 38 | 50165 | 23591 |
| 8523 | I | 32 | 505684 | 411124 |
| 8525 | I | 8 | 42012 | 188388 |
| 8526 | I | 50 | 36498 | 28250 |
| 8528 | I | 62 | 12205124 | 11717020 |
| 8529 | I | 51 | 111113 | 157203 |
| 8530 | I | 21 | 2099 | 1599 |
| 8531 | I | 45 | 281605 | 554185 |
| 8532 | I | 49 | 117186 | 125239 |
| 8533 | I | 15 | 77956 | 73484 |
| 8534 | I | 65 | 46254 | 10930 |
| 8535 | I | 66 | 154340 | 155002 |
| 8536 | I | 44 | 1506284 | 1128405 |
| 8538 | I | 16 | 3978 | 6675 |
| 8539 | I | 48 | 26863 | 27758 |
| 8571 | I | 34 | 936202 | 1058850 |
| 8572 | I | 41 | 53316 | 41328 |
| 8573 | I | 7 | 3212 | 4559 |
| 8574 | I | 40 | 1824842 | 1641073 |
| 8575 | I | 63 | 5945 | 6000 |
| 8576 | I | 52 | 406881 | 401368 |
| 8577 | I | 22 | 240863 | 245679 |
| 8578 | I | 37 | 181815 | 105907 |
| 8579 | I | 43 | 208448 | 174872 |
| 8531 | B | 72 | 545 | 382 |
| 8533 | B | 17 | 1438 | 2379 |
| 8535 | B | 30 | 1706 | 996 |
| 8539 | B | 36 | 1933 | 2578 |
| 8577 | B | 2 | 12254 | 11356 |
| 8578 | B | 46 | 495 | 528 |

Table C.47 Raw Data Returned By Laboratory 23 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 23 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 27126 | 30702 |
| 8513 | U | 31 | 12843016 | 11370888 |
| 8514 | U | 56 | 11310263 | 13722745 |
| 8515 | U | 42 | 27518 | 13978 |
| 8516 | U | 25 | 227545 | 182518 |
| 8517 | U | 58 | 384708 | 375937 |
| 8518 | U | 54 | 324582 | 375462 |
| 8519 | U | 39 | 35896 | 42675 |
| 8520 | U | 68 | 11790 | 18583 |
| 8521 | U | 3 | 32173 | 39145 |
| 8523 | U | 12 | 235602 | 298993 |
| 8525 | U | 1 | 19730 | 30394 |
| 8526 | U | 59 | 18221 | 19675 |
| 8528 | U | 10 | 8898895 | 8207120 |
| 8529 | U | 11 | 78299 | 88134 |
| 8530 | U | 23 | 2545 | 3770 |
| 8531 | U | 29 | 151393 | 194283 |
| 8532 | U | 27 | 97669 | 94046 |
| 8533 | U | 14 | 77856 | 79257 |
| 8534 | U | 4 | 13468 | 16453 |
| 8535 | U | 53 | 130416 | 177003 |
| 8536 | U | 19 | 787179 | 722801 |
| 8538 | U | 61 | 10955 | 4272 |
| 8539 | U | 20 | 21746 | 29762 |
| 8571 | U | 9 | 218646 | 121024 |
| 8572 | U | 6 | 38890 | 36596 |
| 8573 | U | 60 | 5220 | 3230 |
| 8574 | U | 70 | 1250682 | 1336745 |
| 8575 | U | 71 | 171240 | 164152 |
| 8576 | U | 13 | 232214 | 276429 |
| 8577 | U | 55 | 804372 | 699137 |
| 8578 | U | 28 | 72516 | 59597 |
| 8579 | U | 57 | 103985 | 122123 |

Table C.48 Raw Data Returned By Laboratory 23 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 23 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 504421 | 58521 |
| 8513 | I | 26 | 27017713 | 25182247 |
| 8514 | I | 64 | 30236335 | 25610617 |
| 8515 | I | 35 | 18761 | 31518 |
| 8516 | I | 18 | 424749 | 479657 |
| 8517 | I | 5 | 604036 | 580568 |
| 8518 | I | 33 | 492562 | 520034 |
| 8519 | I | 24 | 132360 | 45453 |
| 8520 | I | 69 | 31401 | 101567 |
| 8521 | I | 38 | 47971 | 52836 |
| 8523 | I | 32 | 649531 | 587723 |
| 8525 | I | 8 | 35339 | 50470 |
| 8526 | I | 50 | 23419 | 22656 |
| 8528 | I | 62 | 20302293 | 19018463 |
| 8529 | I | 51 | 108602 | 106855 |
| 8530 | I | 21 | 2674 | 3369 |
| 8531 | I | 45 | 324485 | 261227 |
| 8532 | I | 49 | 124995 | 130054 |
| 8533 | I | 15 | 72620 | 84671 |
| 8534 | I | 65 | 17547 | 16753 |
| 8535 | I | 66 | 251896 | 321500 |
| 8536 | I | 44 | 903064 | 1723331 |
| 8538 | I | 16 | 8657 | 6510 |
| 8539 | I | 48 | 33254 | 30441 |
| 8571 | I | 34 | 1079961 | 1102512 |
| 8572 | I | 41 | 58928 | 51631 |
| 8573 | I | 7 | 3755 | 4770 |
| 8574 | I | 40 | 1985661 | 2278841 |
| 8575 | I | 63 | 8751 | 10166 |
| 8576 | I | 52 | 306748 | 288243 |
| 8577 | I | 22 | 212946 | 211089 |
| 8578 | I | 37 | 95855 | 129065 |
| 8579 | I | 43 | 162973 | 218089 |
| 8531 | B | 72 | 5489 | 3325 |
| 8533 | B | 17 | 50186 | 61315 |
| 8535 | B | 30 | 131323 | 136991 |
| 8539 | B | 36 | 25605 | 28374 |
| 8577 | B | 2 | 804184 | 709120 |
| 8578 | B | 46 | 65388 | 71690 |

Table C.49 Raw Data Returned By Laboratory 24 – screening of unirradiated samples

| SP number | Status | Code number | Lab 24 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 889 | 295 |
| 8513 | U | 31 | 321 | 647 |
| 8514 | U | 56 | 2001 | 3197 |
| 8515 | U | 42 | 356 | 262 |
| 8516 | U | 25 | 1563 | 999 |
| 8517 | U | 58 | 595 | 584 |
| 8518 | U | 54 | 433 | 402 |
| 8519 | U | 39 | 307 | 308 |
| 8520 | U | 68 | 452 | 442 |
| 8521 | U | 3 | 271 | 284 |
| 8523 | U | 12 | 371 | 360 |
| 8525 | U | 1 | 472 | 555 |
| 8526 | U | 59 | 409 | 498 |
| 8528 | U | 10 | 467 | 408 |
| 8529 | U | 11 | 508 | 387 |
| 8530 | U | 23 | 264 | 181 |
| 8531 | U | 29 | 486 | 760 |
| 8532 | U | 27 | 458 | 1257 |
| 8533 | U | 14 | 353 | 279 |
| 8534 | U | 4 | 337 | 275 |
| 8535 | U | 53 | 1216 | 922 |
| 8536 | U | 19 | 4144 | 7877 |
| 8538 | U | 61 | 563 | 440 |
| 8539 | U | 20 | 304 | 345 |
| 8571 | U | 9 | 71742 | 80906 |
| 8572 | U | 6 | 411 | 448 |
| 8573 | U | 60 | 361 | 351 |
| 8574 | U | 70 | 2212 | 3491 |
| 8575 | U | 71 | 239 | 378 |
| 8576 | U | 13 | 552 | 621 |
| 8577 | U | 55 | 7781 | 5009 |
| 8578 | U | 28 | 519 | 440 |
| 8579 | U | 57 | 371 | 433 |

Table C.50 Raw Data Returned By Laboratory 24 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 24 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 68073 | 57082 |
| 8513 | I | 26 | 13188717 | 14036104 |
| 8514 | I | 64 | 6029965 | 4281496 |
| 8515 | I | 35 | 16111 | 12053 |
| 8516 | I | 18 | 267268 | 253218 |
| 8517 | I | 5 | 420056 | 559096 |
| 8518 | I | 33 | 347250 | 430311 |
| 8519 | I | 24 | 42122 | 56764 |
| 8520 | I | 69 | 22019 | 20094 |
| 8521 | I | 38 | 26129 | 26618 |
| 8523 | I | 32 | 434720 | 572184 |
| 8525 | I | 8 | 31213 | 61576 |
| 8526 | I | 50 | 29379 | 34955 |
| 8528 | I | 62 | 3163677 | 4909169 |
| 8529 | I | 51 | 152932 | 278455 |
| 8530 | I | 21 | 2751 | 1667 |
| 8531 | I | 45 | 477079 | 779454 |
| 8532 | I | 49 | 154759 | 177001 |
| 8533 | I | 15 | 83581 | 78880 |
| 8534 | I | 65 | 15188 | 15973 |
| 8535 | I | 66 | 204070 | 199679 |
| 8536 | I | 44 | 1951431 | 1805891 |
| 8538 | I | 16 | 9474 | 14847 |
| 8539 | I | 48 | 50254 | 62748 |
| 8571 | I | 34 | 843559 | 906310 |
| 8572 | I | 41 | 72996 | 85907 |
| 8573 | I | 7 | 8981 | 8610 |
| 8574 | I | 40 | 2094402 | 2328608 |
| 8575 | I | 63 | 12552 | 7546 |
| 8576 | I | 52 | 464995 | 518246 |
| 8577 | I | 22 | 202436 | 270595 |
| 8578 | I | 37 | 103837 | 113930 |
| 8579 | I | 43 | 226762 | 182473 |
| 8531 | B | 72 | 58172 | 70512 |
| 8533 | B | 17 | 2309 | 1034 |
| 8535 | B | 30 | 1210 | 849 |
| 8539 | B | 36 | 3285 | 2672 |
| 8577 | B | 2 | 20234 | 34121 |
| 8578 | B | 46 | 798 | 780 |

Table C.51 Raw Data Returned By Laboratory 24 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 24 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 28072 | 29786 |
| 8513 | U | 31 | 7215163 | 7322224 |
| 8514 | U | 56 | 3300452 | 15224325 |
| 8515 | U | 42 | 11294 | 12144 |
| 8516 | U | 25 | 419918 | 222887 |
| 8517 | U | 58 | 338561 | 331646 |
| 8518 | U | 54 | 180908 | 237446 |
| 8519 | U | 39 | 60163 | 37850 |
| 8520 | U | 68 | 8489 | 18507 |
| 8521 | U | 3 | 28153 | 26921 |
| 8523 | U | 12 | 168843 | 198062 |
| 8525 | U | 1 | 23264 | 19699 |
| 8526 | U | 59 | 14232 | 20890 |
| 8528 | U | 10 | 6701364 | 5934251 |
| 8529 | U | 11 | 58782 | 56797 |
| 8530 | U | 23 | 730 | 1881 |
| 8531 | U | 29 | 184594 | 172743 |
| 8532 | U | 27 | 56592 | 60078 |
| 8533 | U | 14 | 27131 | 34118 |
| 8534 | U | 4 | 6155 | 6689 |
| 8535 | U | 53 | 111238 | 78712 |
| 8536 | U | 19 | 988672 | 810026 |
| 8538 | U | 61 | 2447 | 2393 |
| 8539 | U | 20 | 17905 | 28610 |
| 8571 | U | 9 | 55868 | 54423 |
| 8572 | U | 6 | 23451 | 21197 |
| 8573 | U | 60 | 2788 | 2278 |
| 8574 | U | 70 | 992615 | 844343 |
| 8575 | U | 71 | 3130 | 2923 |
| 8576 | U | 13 | 292025 | 233056 |
| 8577 | U | 55 | 625338 | 701776 |
| 8578 | U | 28 | 59992 | 66165 |
| 8579 | U | 57 | 108282 | 189162 |

Table C.52 Raw Data Returned By Laboratory 24 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 24 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 53105 | 42231 |
| 8513 | I | 26 | 15070049 | 16120632 |
| 8514 | I | 64 | 2136081 | 16201835 |
| 8515 | I | 35 | 23149 | 17031 |
| 8516 | I | 18 | 503773 | 424515 |
| 8517 | I | 5 | 335117 | 376897 |
| 8518 | I | 33 | 420331 | 400253 |
| 8519 | I | 24 | 73899 | 43259 |
| 8520 | I | 69 | 22872 | 31494 |
| 8521 | I | 38 | 29386 | 30361 |
| 8523 | I | 32 | 387850 | 374088 |
| 8525 | I | 8 | 45811 | 36817 |
| 8526 | I | 50 | 23850 | 19137 |
| 8528 | I | 62 | 3679937 | 1488252 |
| 8529 | I | 51 | 103851 | 92981 |
| 8530 | I | 21 | 2867 | 1265 |
| 8531 | I | 45 | 256626 | 455275 |
| 8532 | I | 49 | 98842 | 99641 |
| 8533 | I | 15 | 48685 | 64334 |
| 8534 | I | 65 | 9554 | 10203 |
| 8535 | I | 66 | 173397 | 144224 |
| 8536 | I | 44 | 1159459 | 974808 |
| 8538 | I | 16 | 5349 | 2175 |
| 8539 | I | 48 | 36361 | 36438 |
| 8571 | I | 34 | 1084074 | 897816 |
| 8572 | I | 41 | 42478 | 50558 |
| 8573 | I | 7 | 3942 | 5174 |
| 8574 | I | 40 | 1631181 | 1404564 |
| 8575 | I | 63 | 8764 | 7794 |
| 8576 | I | 52 | 373326 | 406720 |
| 8577 | I | 22 | 152163 | 238487 |
| 8578 | I | 37 | 113047 | 133386 |
| 8579 | I | 43 | 94092 | 239583 |
| 8531 | B | 72 | 143690 | 303328 |
| 8533 | B | 17 | 35147 | 36846 |
| 8535 | B | 30 | 79547 | 65724 |
| 8539 | B | 36 | 15141 | 10599 |
| 8577 | B | 2 | 788876 | 725733 |
| 8578 | B | 46 | 92927 | 111189 |

Table C.53 Raw Data Returned By Laboratory 25 – screening of unirradiated samples

| SP number | Status | Code number | Lab 25 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 200 | 267 |
| 8513 | U | 31 | 376 | 498 |
| 8514 | U | 56 | 1650 | 3975 |
| 8515 | U | 42 | 388 | 294 |
| 8516 | U | 25 | 1034 | 2424 |
| 8517 | U | 58 | 547 | 502 |
| 8518 | U | 54 | 335 | 453 |
| 8519 | U | 39 | 346 | 299 |
| 8520 | U | 68 | 427 | 316 |
| 8521 | U | 3 | 320 | 254 |
| 8523 | U | 12 | 345 | 383 |
| 8525 | U | 1 | 486 | 519 |
| 8526 | U | 59 | 488 | 259 |
| 8528 | U | 10 | 416 | 303 |
| 8529 | U | 11 | 295 | 452 |
| 8530 | U | 23 | 307 | 346 |
| 8531 | U | 29 | 533 | 648 |
| 8532 | U | 27 | 418 | 455 |
| 8533 | U | 14 | 351 | 299 |
| 8534 | U | 4 | 321 | 332 |
| 8535 | U | 53 | 633 | 698 |
| 8536 | U | 19 | 4467 | 4338 |
| 8538 | U | 61 | 364 | 331 |
| 8539 | U | 20 | 274 | 483 |
| 8571 | U | 9 | 461 | 530 |
| 8572 | U | 6 | 365 | 434 |
| 8573 | U | 60 | 358 | 339 |
| 8574 | U | 70 | 2064 | 1743 |
| 8575 | U | 71 | 383 | 418 |
| 8576 | U | 13 | 424 | 464 |
| 8577 | U | 55 | 4056 | 3089 |
| 8578 | U | 28 | 543 | 435 |
| 8579 | U | 57 | 344 | 384 |

Non-negative results shaded

Table C.54 Raw Data Returned By Laboratory 25 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 25 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 29214 | 36917 |
| 8513 | I | 26 | 14653081 | 16094050 |
| 8514 | I | 64 | 1916374 | 9972190 |
| 8515 | I | 35 | 13495 | 13672 |
| 8516 | I | 18 | 245643 | 235948 |
| 8517 | I | 5 | 542001 | 393641 |
| 8518 | I | 33 | 378162 | 194684 |
| 8519 | I | 24 | 34252 | 32276 |
| 8520 | I | 69 | 26156 | 11492 |
| 8521 | I | 38 | 21430 | 22150 |
| 8523 | I | 32 | 306643 | 356629 |
| 8525 | I | 8 | 28998 | 24835 |
| 8526 | I | 50 | 42238 | 13722 |
| 8528 | I | 62 | 14782595 | 14257131 |
| 8529 | I | 51 | 124675 | 272642 |
| 8530 | I | 21 | 2292 | 2713 |
| 8531 | I | 45 | 426703 | 257506 |
| 8532 | I | 49 | 116332 | 150228 |
| 8533 | I | 15 | 64815 | 68459 |
| 8534 | I | 65 | 6667 | 8754 |
| 8535 | I | 66 | 122774 | 150600 |
| 8536 | I | 44 | 1058034 | 1199078 |
| 8538 | I | 16 | 3270 | 2918 |
| 8539 | I | 48 | 44393 | 18191 |
| 8571 | I | 34 | 831817 | 762982 |
| 8572 | I | 41 | 32813 | 32961 |
| 8573 | I | 7 | 3070 | 2902 |
| 8574 | I | 40 | 1466454 | 1357369 |
| 8575 | I | 63 | 4586 | 3394 |
| 8576 | I | 52 | 268313 | 235802 |
| 8577 | I | 22 | 221093 | 166053 |
| 8578 | I | 37 | 184267 | 131590 |
| 8579 | I | 43 | 115655 | 283417 |
| 8531 | B | 72 | 47173 | 39754 |
| 8533 | B | 17 | 861 | 706 |
| 8535 | B | 30 | 1382 | 937 |
| 8539 | B | 36 | 1961 | 2219 |
| 8577 | B | 2 | 9462 | 10068 |
| 8578 | B | 46 | 359 | 518 |

Table C.55 Raw Data Returned By Laboratory 25 –calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 25 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 20722 | 19467 |
| 8513 | U | 31 | 4949295 | 15483075 |
| 8514 | U | 56 | 1561909 | 6885856 |
| 8515 | U | 42 | 10262 | 46242 |
| 8516 | U | 25 | 216214 | 310170 |
| 8517 | U | 58 | 250737 | 299545 |
| 8518 | U | 54 | 219527 | 247300 |
| 8519 | U | 39 | 33393 | 22196 |
| 8520 | U | 68 | 9120 | 30445 |
| 8521 | U | 3 | 20363 | 24259 |
| 8523 | U | 12 | 280020 | 233123 |
| 8525 | U | 1 | 21319 | 15307 |
| 8526 | U | 59 | 11863 | 12818 |
| 8528 | U | 10 | 2765768 | 15303573 |
| 8529 | U | 11 | 166965 | 57115 |
| 8530 | U | 23 | 1925 | 1129 |
| 8531 | U | 29 | 132154 | 144437 |
| 8532 | U | 27 | 62912 | 52774 |
| 8533 | U | 14 | 39799 | 42736 |
| 8534 | U | 4 | 8506 | 7720 |
| 8535 | U | 53 | 81897 | 74912 |
| 8536 | U | 19 | 818134 | 897340 |
| 8538 | U | 61 | 31581 | 21682 |
| 8539 | U | 20 | 14893 | 23452 |
| 8571 | U | 9 | 114067 | 130362 |
| 8572 | U | 6 | 32710 | 27749 |
| 8573 | U | 60 | 2731 | 3133 |
| 8574 | U | 70 | 962765 | 943902 |
| 8575 | U | 71 | 2686 | 3995 |
| 8576 | U | 13 | 172218 | 153058 |
| 8577 | U | 55 | 488270 | 509935 |
| 8578 | U | 28 | 44827 | 52558 |
| 8579 | U | 57 | 78604 | 88978 |

Table C.56 Raw Data Returned By Laboratory 25 –calibrated PSL for irradiated samples

| SP number | Status | Code number | Lab 25 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 30193 | 30498 |
| 8513 | I | 26 | 1108494 | 13503502 |
| 8514 | I | 64 | 1838472 | 2257048 |
| 8515 | I | 35 | 15776 | 20003 |
| 8516 | I | 18 | 1010070 | 641343 |
| 8517 | I | 5 | 487732 | 513865 |
| 8518 | I | 33 | 322461 | 243261 |
| 8519 | I | 24 | 48468 | 86335 |
| 8520 | I | 69 | 24894 | 24949 |
| 8521 | I | 38 | 22181 | 29020 |
| 8523 | I | 32 | 254489 | 247665 |
| 8525 | I | 8 | 32982 | 30871 |
| 8526 | I | 50 | 30262 | 23250 |
| 8528 | I | 62 | 14607130 | 12687175 |
| 8529 | I | 51 | 132348 | 141299 |
| 8530 | I | 21 | 2000 | 2381 |
| 8531 | I | 45 | 310014 | 187733 |
| 8532 | I | 49 | 115090 | 96452 |
| 8533 | I | 15 | 60087 | 66517 |
| 8534 | I | 65 | 9170 | 11202 |
| 8535 | I | 66 | 142478 | 143540 |
| 8536 | I | 44 | 1282496 | 1478708 |
| 8538 | I | 16 | 9669 | 7520 |
| 8539 | I | 48 | 33377 | 26874 |
| 8571 | I | 34 | 752683 | 828586 |
| 8572 | I | 41 | 19319 | 32076 |
| 8573 | I | 7 | 4113 | 3127 |
| 8574 | I | 40 | 1409525 | 1365439 |
| 8575 | I | 63 | 18467 | 10513 |
| 8576 | I | 52 | 249669 | 324279 |
| 8577 | I | 22 | 176398 | 129438 |
| 8578 | I | 37 | 81046 | 119961 |
| 8579 | I | 43 | 146900 | 444267 |
| 8531 | B | 72 | 158914 | 203422 |
| 8533 | B | 17 | 43588 | 37107 |
| 8535 | B | 30 | 86856 | 68608 |
| 8539 | B | 36 | 15337 | 42143 |
| 8577 | B | 2 | 421424 | 458753 |
| 8578 | B | 46 | 74752 | 46971 |

Table C.57 Raw Data Returned By Laboratory 26 – screening of unirradiated samples

| SP number | Status | Code number | Lab 26 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 299 | 333 |
| 8513 | U | 31 | 376 | 280 |
| 8514 | U | 56 | 1580 | 1022 |
| 8515 | U | 42 | 357 | 348 |
| 8516 | U | 25 | 2116 | 1435 |
| 8517 | U | 58 | 599 | 369 |
| 8518 | U | 54 | 566 | 440 |
| 8519 | U | 39 | 263 | 283 |
| 8520 | U | 68 | 297 | 468 |
| 8521 | U | 3 | 289 | 248 |
| 8523 | U | 12 | 609 | 484 |
| 8525 | U | 1 | 548 | 385 |
| 8526 | U | 59 | 462 | 301 |
| 8528 | U | 10 | 404 | 382 |
| 8529 | U | 11 | 368 | 399 |
| 8530 | U | 23 | 274 | 310 |
| 8531 | U | 29 | 661 | 498 |
| 8532 | U | 27 | 2084 | 1000 |
| 8533 | U | 14 | 504 | 238 |
| 8534 | U | 4 | 278 | 376 |
| 8535 | U | 53 | 600 | 658 |
| 8536 | U | 19 | 2734 | 1338 |
| 8538 | U | 61 | 305 | 418 |
| 8539 | U | 20 | 372 | 358 |
| 8571 | U | 9 | 401 | 458 |
| 8572 | U | 6 | 329 | 511 |
| 8573 | U | 60 | 677 | 528 |
| 8574 | U | 70 | 1559 | 4477 |
| 8575 | U | 71 | 247 | 367 |
| 8576 | U | 13 | 358 | 414 |
| 8577 | U | 55 | 2734 | 3075 |
| 8578 | U | 28 | 484 | 589 |
| 8579 | U | 57 | 393 | 323 |

Non-negative results shaded

Table C.58 Raw Data Returned By Laboratory 26 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 26 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 25419 | 23414 |
| 8513 | I | 26 | 16559751 | 559458 |
| 8514 | I | 64 | 7192896 | 14047800 |
| 8515 | I | 35 | 55476 | 10228 |
| 8516 | I | 18 | 670984 | 869742 |
| 8517 | I | 5 | 345413 | 352157 |
| 8518 | I | 33 | 252406 | 229573 |
| 8519 | I | 24 | 32061 | 32190 |
| 8520 | I | 69 | 24041 | 11203 |
| 8521 | I | 38 | 22909 | 23079 |
| 8523 | I | 32 | 313808 | 167297 |
| 8525 | I | 8 | 27924 | 27976 |
| 8526 | I | 50 | 12125 | 11619 |
| 8528 | I | 62 | 6429549 | 12062249 |
| 8529 | I | 51 | 55725 | 62914 |
| 8530 | I | 21 | 1251 | 2017 |
| 8531 | I | 45 | 173274 | 208672 |
| 8532 | I | 49 | 64915 | 115261 |
| 8533 | I | 15 | 49224 | 45715 |
| 8534 | I | 65 | 8526 | 5515 |
| 8535 | I | 66 | 109885 | 79495 |
| 8536 | I | 44 | 866877 | 599006 |
| 8538 | I | 16 | 4360 | 4751 |
| 8539 | I | 48 | 25083 | 22986 |
| 8571 | I | 34 | 548036 | 682493 |
| 8572 | I | 41 | 30514 | 22423 |
| 8573 | I | 7 | 2246 | 2199 |
| 8574 | I | 40 | 793854 | 1024138 |
| 8575 | I | 63 | 2394 | 2682 |
| 8576 | I | 52 | 209878 | 183694 |
| 8577 | I | 22 | 176795 | 151985 |
| 8578 | I | 37 | 54087 | 56220 |
| 8579 | I | 43 | 62883 | 83922 |
| 8531 | B | 72 | 18219 | 14555 |
| 8533 | B | 17 | 1953 | 1127 |
| 8535 | B | 30 | 579 | 527 |
| 8539 | B | 36 | 1149 | 1188 |
| 8577 | B | 2 | 8149 | 13203 |
| 8578 | B | 46 | 587 | 332 |

Table C.59 Raw Data Returned By Laboratory 26 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 26 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 15310 | 21613 |
| 8513 | U | 31 | 2129713 | 16378840 |
| 8514 | U | 56 | 7768798 | 1600473 |
| 8515 | U | 42 | 6521 | 5133 |
| 8516 | U | 25 | 266702 | 314550 |
| 8517 | U | 58 | 234983 | 206775 |
| 8518 | U | 54 | 189129 | 54597 |
| 8519 | U | 39 | 24156 | 17205 |
| 8520 | U | 68 | 5869 | 7831 |
| 8521 | U | 3 | 31953 | 18983 |
| 8523 | U | 12 | 127250 | 165050 |
| 8525 | U | 1 | 14327 | 12180 |
| 8526 | U | 59 | 32538 | 33714 |
| 8528 | U | 10 | 7433412 | 13388576 |
| 8529 | U | 11 | 49926 | 59052 |
| 8530 | U | 23 | 1179 | 1386 |
| 8531 | U | 29 | 127828 | 137471 |
| 8532 | U | 27 | 64647 | 74173 |
| 8533 | U | 14 | 35994 | 33338 |
| 8534 | U | 4 | 8322 | 4930 |
| 8535 | U | 53 | 51645 | 138912 |
| 8536 | U | 19 | 807699 | 562688 |
| 8538 | U | 61 | 4117 | 3955 |
| 8539 | U | 20 | 10241 | 11359 |
| 8571 | U | 9 | 79796 | 96872 |
| 8572 | U | 6 | 18811 | 19654 |
| 8573 | U | 60 | 4973 | 1546 |
| 8574 | U | 70 | 714857 | 642996 |
| 8575 | U | 71 | 2858 | 2070 |
| 8576 | U | 13 | 142537 | 143351 |
| 8577 | U | 55 | 364815 | 321459 |
| 8578 | U | 28 | 45394 | 38528 |
| 8579 | U | 57 | 45989 | 198385 |

Table C.60 Raw Data Returned By Laboratory 26 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 26 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 31854 | 22213 |
| 8513 | I | 26 | 4390921 | 12317747 |
| 8514 | I | 64 | 9127918 | 16348653 |
| 8515 | I | 35 | 44587 | 11409 |
| 8516 | I | 18 | 543342 | 667190 |
| 8517 | I | 5 | 418990 | 314295 |
| 8518 | I | 33 | 206819 | 250137 |
| 8519 | I | 24 | 43126 | 30881 |
| 8520 | I | 69 | 20383 | 13495 |
| 8521 | I | 38 | 24793 | 40103 |
| 8523 | I | 32 | 278812 | 138647 |
| 8525 | I | 8 | 16262 | 39742 |
| 8526 | I | 50 | 24230 | 12410 |
| 8528 | I | 62 | 14703655 | 9232486 |
| 8529 | I | 51 | 77897 | 71235 |
| 8530 | I | 21 | 1623 | 2149 |
| 8531 | I | 45 | 155144 | 247579 |
| 8532 | I | 49 | 61320 | 65505 |
| 8533 | I | 15 | 55912 | 55872 |
| 8534 | I | 65 | 10639 | 9642 |
| 8535 | I | 66 | 174029 | 147468 |
| 8536 | I | 44 | 975020 | 666993 |
| 8538 | I | 16 | 4349 | 9652 |
| 8539 | I | 48 | 22737 | 24318 |
| 8571 | I | 34 | 573623 | 606758 |
| 8572 | I | 41 | 29138 | 25545 |
| 8573 | I | 7 | 7587 | 1989 |
| 8574 | I | 40 | 784677 | 1090956 |
| 8575 | I | 63 | 6277 | 3895 |
| 8576 | I | 52 | 190257 | 290478 |
| 8577 | I | 22 | 155336 | 129428 |
| 8578 | I | 37 | 73647 | 80672 |
| 8579 | I | 43 | 88303 | 51407 |
| 8531 | B | 72 | 1747980 | 86628 |
| 8533 | B | 17 | 30876 | 30276 |
| 8535 | B | 30 | 76750 | 64145 |
| 8539 | B | 36 | 9850 | 12807 |
| 8577 | B | 2 | 374151 | 361286 |
| 8578 | B | 46 | 53382 | 36909 |

Table C.61 Raw Data Returned By Laboratory 27 – screening of unirradiated samples

| SP number | Status | Code number | Lab 27 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 330 | 267 |
| 8513 | U | 31 | 272 | 546 |
| 8514 | U | 56 | 880 | 821 |
| 8515 | U | 42 | 239 | 293 |
| 8516 | U | 25 | 625 | 598 |
| 8517 | U | 58 | 324 | 297 |
| 8518 | U | 54 | 267 | 287 |
| 8519 | U | 39 | 333 | 217 |
| 8520 | U | 68 | 306 | 344 |
| 8521 | U | 3 | 350 | 329 |
| 8523 | U | 12 | 319 | 309 |
| 8525 | U | 1 | 327 | 443 |
| 8526 | U | 59 | 258 | 171 |
| 8528 | U | 10 | 291 | 412 |
| 8529 | U | 11 | 344 | 376 |
| 8530 | U | 23 | 341 | 264 |
| 8531 | U | 29 | 501 | 346 |
| 8532 | U | 27 | 1371 | 1285 |
| 8533 | U | 14 | 250 | 348 |
| 8534 | U | 4 | 182 | 300 |
| 8535 | U | 53 | 609 | 393 |
| 8536 | U | 19 | 1468 | 1629 |
| 8538 | U | 61 | 345 | 359 |
| 8539 | U | 20 | 335 | 406 |
| 8571 | U | 9 | 395 | 363 |
| 8572 | U | 6 | 322 | 426 |
| 8573 | U | 60 | 225 | 204 |
| 8574 | U | 70 | 963 | 713 |
| 8575 | U | 71 | 392 | 296 |
| 8576 | U | 13 | 339 | 406 |
| 8577 | U | 55 | 2393 | 1873 |
| 8578 | U | 28 | 447 | 425 |
| 8579 | U | 57 | 268 | 288 |

Non-negative results shaded

Table C.62 Raw Data Returned By Laboratory 27 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 27 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 17748 | 21353 |
| 8513 | I | 26 | 9781532 | 12446254 |
| 8514 | I | 64 | 14178201 | 12438359 |
| 8515 | I | 35 | 5376 | 6982 |
| 8516 | I | 18 | 162206 | 180251 |
| 8517 | I | 5 | 323287 | 290101 |
| 8518 | I | 33 | 169886 | 168575 |
| 8519 | I | 24 | 17595 | 17982 |
| 8520 | I | 69 | 19973 | 10475 |
| 8521 | I | 38 | 13313 | 13224 |
| 8523 | I | 32 | 139971 | 207579 |
| 8525 | I | 8 | 93256 | 34780 |
| 8526 | I | 50 | 10275 | 14904 |
| 8528 | I | 62 | 7033536 | 10726789 |
| 8529 | I | 51 | 49435 | 32341 |
| 8530 | I | 21 | 1549 | 871 |
| 8531 | I | 45 | 168518 | 147755 |
| 8532 | I | 49 | 49121 | 47231 |
| 8533 | I | 15 | 28033 | 32120 |
| 8534 | I | 65 | 3735 | 4515 |
| 8535 | I | 66 | 54807 | 59590 |
| 8536 | I | 44 | 341192 | 445812 |
| 8538 | I | 16 | 2635 | 2011 |
| 8539 | I | 48 | 26832 | 12517 |
| 8571 | I | 34 | 330842 | 360796 |
| 8572 | I | 41 | 14158 | 15521 |
| 8573 | I | 7 | 2900 | 2006 |
| 8574 | I | 40 | 598362 | 848532 |
| 8575 | I | 63 | 1992 | 3877 |
| 8576 | I | 52 | 124406 | 133610 |
| 8577 | I | 22 | 107184 | 151422 |
| 8578 | I | 37 | 42183 | 44700 |
| 8579 | I | 43 | 51499 | 80283 |
| 8531 | B | 72 | 10832 | 10473 |
| 8533 | B | 17 | 569 | 303 |
| 8535 | B | 30 | 729 | 726 |
| 8539 | B | 36 | 1467 | 1261 |
| 8577 | B | 2 | 4487 | 4680 |
| 8578 | B | 46 | 443 | 388 |

Table C.63 Raw Data Returned By Laboratory 27 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 27 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 11600 | 9067 |
| 8513 | U | 31 | 16000206 | 15730733 |
| 8514 | U | 56 | 5627783 | 9432607 |
| 8515 | U | 42 | 4232 | 4900 |
| 8516 | U | 25 | 376925 | 403688 |
| 8517 | U | 58 | 134257 | 119890 |
| 8518 | U | 54 | 99347 | 108499 |
| 8519 | U | 39 | 17116 | 16398 |
| 8520 | U | 68 | 4299 | 4764 |
| 8521 | U | 3 | 9814 | 12021 |
| 8523 | U | 12 | 73494 | 97374 |
| 8525 | U | 1 | 8216 | 13853 |
| 8526 | U | 59 | 6649 | 5674 |
| 8528 | U | 10 | 14646318 | 14737471 |
| 8529 | U | 11 | 44203 | 25408 |
| 8530 | U | 23 | 1232 | 671 |
| 8531 | U | 29 | 68540 | 72863 |
| 8532 | U | 27 | 30822 | 23222 |
| 8533 | U | 14 | 18461 | 17718 |
| 8534 | U | 4 | 3297 | 3085 |
| 8535 | U | 53 | 36321 | 32635 |
| 8536 | U | 19 | 386610 | 343622 |
| 8538 | U | 61 | 3574 | 7439 |
| 8539 | U | 20 | 25319 | 7594 |
| 8571 | U | 9 | 69007 | 40059 |
| 8572 | U | 6 | 11023 | 9545 |
| 8573 | U | 60 | 1563 | 1329 |
| 8574 | U | 70 | 385299 | 373003 |
| 8575 | U | 71 | 1332 | 1826 |
| 8576 | U | 13 | 77212 | 76844 |
| 8577 | U | 55 | 209376 | 227951 |
| 8578 | U | 28 | 30095 | 29954 |
| 8579 | U | 57 | 43222 | 37506 |

Table C.64 Raw Data Returned By Laboratory 27 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 27 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 18592 | 13374 |
| 8513 | I | 26 | 10267835 | 11961542 |
| 8514 | I | 64 | 4697571 | 9058260 |
| 8515 | I | 35 | 7892 | 9850 |
| 8516 | I | 18 | 911440 | 417701 |
| 8517 | I | 5 | 208407 | 283216 |
| 8518 | I | 33 | 162248 | 177075 |
| 8519 | I | 24 | 33444 | 24919 |
| 8520 | I | 69 | 14039 | 5926 |
| 8521 | I | 38 | 15334 | 12735 |
| 8523 | I | 32 | 134858 | 262085 |
| 8525 | I | 8 | 22992 | 8493 |
| 8526 | I | 50 | 10250 | 10817 |
| 8528 | I | 62 | 12147665 | 10603278 |
| 8529 | I | 51 | 62001 | 36012 |
| 8530 | I | 21 | 2452 | 5958 |
| 8531 | I | 45 | 127917 | 152187 |
| 8532 | I | 49 | 46667 | 36620 |
| 8533 | I | 15 | 25478 | 27767 |
| 8534 | I | 65 | 7435 | 9605 |
| 8535 | I | 66 | 52416 | 87467 |
| 8536 | I | 44 | 481115 | 336968 |
| 8538 | I | 16 | 4478 | 3327 |
| 8539 | I | 48 | 19280 | 17527 |
| 8571 | I | 34 | 419362 | 372177 |
| 8572 | I | 41 | 15399 | 19795 |
| 8573 | I | 7 | 2233 | 1432 |
| 8574 | I | 40 | 792196 | 629634 |
| 8575 | I | 63 | 3280 | 1856 |
| 8576 | I | 52 | 106570 | 106426 |
| 8577 | I | 22 | 81070 | 104830 |
| 8578 | I | 37 | 37597 | 41842 |
| 8579 | I | 43 | 47667 | 58239 |
| 8531 | B | 72 | 50833 | 78412 |
| 8533 | B | 17 | 21120 | 20396 |
| 8535 | B | 30 | 32068 | 36240 |
| 8539 | B | 36 | 5641 | 11812 |
| 8577 | B | 2 | 217115 | 204438 |
| 8578 | B | 46 | 18840 | 16442 |

Table C.65 Raw Data Returned By Laboratory 28 – screening of unirradiated samples

| SP number | Status | Code number | Lab 28 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 697 | 599 |
| 8513 | U | 31 | 263 | 374 |
| 8514 | U | 56 | 1031 | 2340 |
| 8515 | U | 42 | 377 | 270 |
| 8516 | U | 25 | 1026 | 1360 |
| 8517 | U | 58 | 527 | 353 |
| 8518 | U | 54 | 767 | 312 |
| 8519 | U | 39 | 312 | 263 |
| 8520 | U | 68 | 692 | 800 |
| 8521 | U | 3 | 365 | 488 |
| 8523 | U | 12 | 473 | 287 |
| 8525 | U | 1 | 2084 | 385 |
| 8526 | U | 59 | 508 | 176 |
| 8528 | U | 10 | 427 | 460 |
| 8529 | U | 11 | 522 | 501 |
| 8530 | U | 23 | 410 | 298 |
| 8531 | U | 29 | 631 | 1049 |
| 8532 | U | 27 | 1407 | 1093 |
| 8533 | U | 14 | 365 | 386 |
| 8534 | U | 4 | 430 | 318 |
| 8535 | U | 53 | 460 | 387 |
| 8536 | U | 19 | 2346 | 824 |
| 8538 | U | 61 | 370 | 157 |
| 8539 | U | 20 | 363 | 511 |
| 8571 | U | 9 | 595 | 651 |
| 8572 | U | 6 | 332 | 623 |
| 8573 | U | 60 | 377 | 215 |
| 8574 | U | 70 | 2278 | 1818 |
| 8575 | U | 71 | 589 | 1176 |
| 8576 | U | 13 | 342 | 490 |
| 8577 | U | 55 | 2851 | 4193 |
| 8578 | U | 28 | 632 | 509 |
| 8579 | U | 57 | 313 | 295 |

Non-negative results shaded

Table C.66 Raw Data Returned By Laboratory 28 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 28 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 34514 | 31954 |
| 8513 | I | 26 | 6058085 | 2581105 |
| 8514 | I | 64 | 12702482 | 11417309 |
| 8515 | I | 35 | 25336 | 17021 |
| 8516 | I | 18 | 721293 | 846620 |
| 8517 | I | 5 | 493770 | 574905 |
| 8518 | I | 33 | 249094 | 289214 |
| 8519 | I | 24 | 43563 | 37591 |
| 8520 | I | 69 | 13424 | 12220 |
| 8521 | I | 38 | 23381 | 22817 |
| 8523 | I | 32 | 270452 | 294572 |
| 8525 | I | 8 | 59002 | 32631 |
| 8526 | I | 50 | 27572 | 22464 |
| 8528 | I | 62 | 6684711 | 3786616 |
| 8529 | I | 51 | 115030 | 88303 |
| 8530 | I | 21 | 1124 | 2031 |
| 8531 | I | 45 | 317801 | 345380 |
| 8532 | I | 49 | 92440 | 87119 |
| 8533 | I | 15 | 59988 | 52430 |
| 8534 | I | 65 | 10907 | 8264 |
| 8535 | I | 66 | 139954 | 135493 |
| 8536 | I | 44 | 1100775 | 849671 |
| 8538 | I | 16 | 4432 | 3273 |
| 8539 | I | 48 | 34955 | 31447 |
| 8571 | I | 34 | 842672 | 804175 |
| 8572 | I | 41 | 61052 | 32013 |
| 8573 | I | 7 | 2724 | 7944 |
| 8574 | I | 40 | 1249000 | 1296348 |
| 8575 | I | 63 | 6105 | 4181 |
| 8576 | I | 52 | 284462 | 318662 |
| 8577 | I | 22 | 153910 | 151195 |
| 8578 | I | 37 | 83264 | 83890 |
| 8579 | I | 43 | 158020 | 106559 |
| 8531 | B | 72 | 29890 | 88722 |
| 8533 | B | 17 | 1613 | 715 |
| 8535 | B | 30 | 741 | 949 |
| 8539 | B | 36 | 1649 | 2605 |
| 8577 | B | 2 | 9181 | 8900 |
| 8578 | B | 46 | 675 | 675 |

Table C.67 Raw Data Returned By Laboratory 28 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 28 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 14618 | 13865 |
| 8513 | U | 31 | 9407648 | 3979723 |
| 8514 | U | 56 | 9203401 | 12948678 |
| 8515 | U | 42 | 8874 | 11065 |
| 8516 | U | 25 | 297535 | 326027 |
| 8517 | U | 58 | 217881 | 179946 |
| 8518 | U | 54 | 134110 | 142555 |
| 8519 | U | 39 | 29316 | 9800 |
| 8520 | U | 68 | 6834 | 5403 |
| 8521 | U | 3 | 11519 | 11993 |
| 8523 | U | 12 | 92443 | 109171 |
| 8525 | U | 1 | 21324 | 13289 |
| 8526 | U | 59 | 7496 | 16734 |
| 8528 | U | 10 | 7256934 | 7465576 |
| 8529 | U | 11 | 27216 | 29152 |
| 8530 | U | 23 | 2193 | 800 |
| 8531 | U | 29 | 64947 | 96916 |
| 8532 | U | 27 | 33080 | 36471 |
| 8533 | U | 14 | 21942 | 22668 |
| 8534 | U | 4 | 3525 | 3934 |
| 8535 | U | 53 | 61319 | 46895 |
| 8536 | U | 19 | 176580 | 144961 |
| 8538 | U | 61 | 1797 | 2941 |
| 8539 | U | 20 | 7852 | 7480 |
| 8571 | U | 9 | 56120 | 59988 |
| 8572 | U | 6 | 18550 | 15193 |
| 8573 | U | 60 | 1337 | 1000 |
| 8574 | U | 70 | 434201 | 422043 |
| 8575 | U | 71 | 2128 | 1368 |
| 8576 | U | 13 | 87906 | 84666 |
| 8577 | U | 55 | 270492 | 288981 |
| 8578 | U | 28 | 31229 | 23462 |
| 8579 | U | 57 | 84443 | 46584 |

Table C.68 Raw Data Returned By Laboratory 28 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 28 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 16537 | 17132 |
| 8513 | I | 26 | 16159920 | 16101535 |
| 8514 | I | 64 | 9977615 | 12462322 |
| 8515 | I | 35 | 16425 | 8843 |
| 8516 | I | 18 | 412536 | 426159 |
| 8517 | I | 5 | 191597 | 224690 |
| 8518 | I | 33 | 161349 | 165810 |
| 8519 | I | 24 | 29670 | 26908 |
| 8520 | I | 69 | 9380 | 12095 |
| 8521 | I | 38 | 13630 | 13468 |
| 8523 | I | 32 | 137885 | 183532 |
| 8525 | I | 8 | 23262 | 17001 |
| 8526 | I | 50 | 10657 | 11424 |
| 8528 | I | 62 | 15897800 | 15702167 |
| 8529 | I | 51 | 75594 | 60978 |
| 8530 | I | 21 | 1286 | 1855 |
| 8531 | I | 45 | 134319 | 148807 |
| 8532 | I | 49 | 37122 | 35979 |
| 8533 | I | 15 | 31235 | 27021 |
| 8534 | I | 65 | 4152 | 3768 |
| 8535 | I | 66 | 78065 | 74444 |
| 8536 | I | 44 | 429927 | 384415 |
| 8538 | I | 16 | 2882 | 2573 |
| 8539 | I | 48 | 12902 | 14123 |
| 8571 | I | 34 | 412730 | 352615 |
| 8572 | I | 41 | 19162 | 29087 |
| 8573 | I | 7 | 1963 | 3321 |
| 8574 | I | 40 | 347727 | 373901 |
| 8575 | I | 63 | 2597 | 1634 |
| 8576 | I | 52 | 172478 | 227121 |
| 8577 | I | 22 | 77808 | 82069 |
| 8578 | I | 37 | 43518 | 46422 |
| 8579 | I | 43 | 74403 | 49514 |
| 8531 | B | 72 | 91023 | 111084 |
| 8533 | B | 17 | 14509 | 16178 |
| 8535 | B | 30 | 40177 | 38616 |
| 8539 | B | 36 | 6190 | 7259 |
| 8577 | B | 2 | 293107 | 242503 |
| 8578 | B | 46 | 37265 | 36529 |

Table C.69 Raw Data Returned By Laboratory 29 – screening of unirradiated samples

| SP number | Status | Code number | Lab 29 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 16712 | 8375 |
| 8513 | U | 31 | 389 | 614 |
| 8514 | U | 56 | 2711 | 844 |
| 8515 | U | 42 | 758 | 742 |
| 8516 | U | 25 | 1286 | 1283 |
| 8517 | U | 58 | 88 | 371 |
| 8518 | U | 54 | 330 | 310 |
| 8519 | U | 39 | 714 | 955 |
| 8520 | U | 68 | 10066 | 7657 |
| 8521 | U | 3 | 315 | 442 |
| 8523 | U | 12 | 355 | 416 |
| 8525 | U | 1 | 658 | 416 |
| 8526 | U | 59 | 391 | 287 |
| 8528 | U | 10 | 425 | 270 |
| 8529 | U | 11 | 376 | 312 |
| 8530 | U | 23 | 1545 | 647 |
| 8531 | U | 29 | 44580 | 27987 |
| 8532 | U | 27 | 21277 | 48399 |
| 8533 | U | 14 | 358 | 364 |
| 8534 | U | 4 | 267 | 254 |
| 8535 | U | 53 | 514 | 451 |
| 8536 | U | 19 | 1607 | 3577 |
| 8538 | U | 61 | 421 | 319 |
| 8539 | U | 20 | 650 | 525 |
| 8571 | U | 9 | 248 | 432 |
| 8572 | U | 6 | 465 | 420 |
| 8573 | U | 60 | 218 | 284 |
| 8574 | U | 70 | 6348 | 4446 |
| 8575 | U | 71 | 2503 | 1780 |
| 8576 | U | 13 | 644 | 739 |
| 8577 | U | 55 | 4577 | 3747 |
| 8578 | U | 28 | 36794 | 82173 |
| 8579 | U | 57 | 416 | 302 |

Non-negative results shaded

Table C.70 Raw Data Returned By Laboratory 29 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 29 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 28205 | 27841 |
| 8513 | I | 26 | 3529971 | 3359649 |
| 8514 | I | 64 | 16387379 | 10000048 |
| 8515 | I | 35 | 20094 | 43690 |
| 8516 | I | 18 | 780543 | 712315 |
| 8517 | I | 5 | 452033 | 397588 |
| 8518 | I | 33 | 432165 | 442444 |
| 8519 | I | 24 | 50746 | 39327 |
| 8520 | I | 69 | 16046 | 23797 |
| 8521 | I | 38 | 34442 | 32939 |
| 8523 | I | 32 | 567007 | 484013 |
| 8525 | I | 8 | 15558 | 35266 |
| 8526 | I | 50 | 22205 | 24615 |
| 8528 | I | 62 | 10357120 | 5337156 |
| 8529 | I | 51 | 82544 | 55783 |
| 8530 | I | 21 | 1746 | 1502 |
| 8531 | I | 45 | 294775 | 203923 |
| 8532 | I | 49 | 94328 | 97675 |
| 8533 | I | 15 | 79323 | 79358 |
| 8534 | I | 65 | 19081 | 18543 |
| 8535 | I | 66 | 102872 | 97026 |
| 8536 | I | 44 | 805865 | 896548 |
| 8538 | I | 16 | 21111 | 4159 |
| 8539 | I | 48 | 34136 | 23117 |
| 8571 | I | 34 | 897698 | 1012810 |
| 8572 | I | 41 | 42489 | 30617 |
| 8573 | I | 7 | 1514 | 2280 |
| 8574 | I | 40 | 1377648 | 1516052 |
| 8575 | I | 63 | 40872 | 54925 |
| 8576 | I | 52 | 247135 | 258343 |
| 8577 | I | 22 | 145944 | 130159 |
| 8578 | I | 37 | 132950 | 70849 |
| 8579 | I | 43 | 137262 | 704778 |
| 8531 | B | 72 | 22056 | 21588 |
| 8533 | B | 17 | 736 | 735 |
| 8535 | B | 30 | 840 | 682 |
| 8539 | B | 36 | 2074 | 1865 |
| 8577 | B | 2 | 11159 | 7698 |
| 8578 | B | 46 | 1469 | 1260 |

Table C.71 Raw Data Returned By Laboratory 30 – screening of unirradiated samples

| SP number | Status | Code number | Lab 30 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 250 | 330 |
| 8513 | U | 31 | 288 | 382 |
| 8514 | U | 56 | 888 | 815 |
| 8515 | U | 42 | 215 | 401 |
| 8516 | U | 25 | 1436 | 711 |
| 8517 | U | 58 | 411 | 408 |
| 8518 | U | 54 | 387 | 379 |
| 8519 | U | 39 | 240 | 259 |
| 8520 | U | 68 | 337 | 408 |
| 8521 | U | 3 | 290 | 211 |
| 8523 | U | 12 | 256 | 493 |
| 8525 | U | 1 | 475 | 368 |
| 8526 | U | 59 | 321 | 455 |
| 8528 | U | 10 | 340 | 395 |
| 8529 | U | 11 | 411 | 363 |
| 8530 | U | 23 | 258 | 311 |
| 8531 | U | 29 | 472 | 527 |
| 8532 | U | 27 | 341 | 423 |
| 8533 | U | 14 | 414 | 267 |
| 8534 | U | 4 | 303 | 326 |
| 8535 | U | 53 | 587 | 447 |
| 8536 | U | 19 | 1400 | 1105 |
| 8538 | U | 61 | 304 | 331 |
| 8539 | U | 20 | 409 | 369 |
| 8571 | U | 9 | 549 | 440 |
| 8572 | U | 6 | 412 | 405 |
| 8573 | U | 60 | 256 | 274 |
| 8574 | U | 70 | 1044 | 1033 |
| 8575 | U | 71 | 286 | 255 |
| 8576 | U | 13 | 303 | 533 |
| 8577 | U | 55 | 2313 | 1757 |
| 8578 | U | 28 | 361 | 214 |
| 8579 | U | 57 | 339 | 421 |

Non-negative results shaded

Table C.72 Raw Data Returned By Laboratory 30 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 30 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 19956 | 17853 |
| 8513 | I | 26 | 15269333 | 2140690 |
| 8514 | I | 64 | 11317068 | 1001799 |
| 8515 | I | 35 | 9897 | 6792 |
| 8516 | I | 18 | 292597 | 254947 |
| 8517 | I | 5 | 281399 | 315492 |
| 8518 | I | 33 | 214204 | 230262 |
| 8519 | I | 24 | 27260 | 33955 |
| 8520 | I | 69 | 13830 | 14929 |
| 8521 | I | 38 | 23039 | 19880 |
| 8523 | I | 32 | 198913 | 403476 |
| 8525 | I | 8 | 18382 | 79661 |
| 8526 | I | 50 | 18085 | 14649 |
| 8528 | I | 62 | 5938571 | 10836081 |
| 8529 | I | 51 | 55752 | 117286 |
| 8530 | I | 21 | 1314 | 1429 |
| 8531 | I | 45 | 79479 | 163363 |
| 8532 | I | 49 | 64829 | 39743 |
| 8533 | I | 15 | 45342 | 32933 |
| 8534 | I | 65 | 10768 | 7531 |
| 8535 | I | 66 | 76817 | 94361 |
| 8536 | I | 44 | 636814 | 593175 |
| 8538 | I | 16 | 5871 | 2234 |
| 8539 | I | 48 | 19176 | 13796 |
| 8571 | I | 34 | 434518 | 438536 |
| 8572 | I | 41 | 36011 | 19661 |
| 8573 | I | 7 | 1388 | 1738 |
| 8574 | I | 40 | 1237407 | 794197 |
| 8575 | I | 63 | 2920 | 2840 |
| 8576 | I | 52 | 212953 | 179674 |
| 8577 | I | 22 | 134807 | 191998 |
| 8578 | I | 37 | 59247 | 70484 |
| 8579 | I | 43 | 80660 | 78237 |
| 8531 | B | 72 | 25059 | 10827 |
| 8533 | B | 17 | 658 | 567 |
| 8535 | B | 30 | 1237 | 439 |
| 8539 | B | 36 | 1534 | 974 |
| 8577 | B | 2 | 7615 | 6124 |
| 8578 | B | 46 | 358 | 412 |

Table C.73 Raw Data Returned By Laboratory 31 – screening of unirradiated samples

| SP number | Status | Code number | Lab 31 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 253 | 294 |
| 8513 | U | 31 | 1266 | 1206 |
| 8514 | U | 56 | 1010 | 952 |
| 8515 | U | 42 | 393 | 306 |
| 8516 | U | 25 | 3224 | 3491 |
| 8517 | U | 58 | 571 | 18 |
| 8518 | U | 54 | 322 | 103 |
| 8519 | U | 39 | 227 | 244 |
| 8520 | U | 68 | 375 | 295 |
| 8521 | U | 3 | 444 | 194 |
| 8523 | U | 12 | 381 | 224 |
| 8525 | U | 1 | 611 | 484 |
| 8526 | U | 59 | 423 | 365 |
| 8528 | U | 10 | 295 | 268 |
| 8529 | U | 11 | 308 | 364 |
| 8530 | U | 23 | 289 | 282 |
| 8531 | U | 29 | 2471 | 1063 |
| 8532 | U | 27 | 629 | 378 |
| 8533 | U | 14 | 431 | 248 |
| 8534 | U | 4 | 290 | 329 |
| 8535 | U | 53 | 597 | 666 |
| 8536 | U | 19 | 813 | 883 |
| 8538 | U | 61 | 344 | 353 |
| 8539 | U | 20 | 437 | 485 |
| 8571 | U | 9 | 401 | 218 |
| 8572 | U | 6 | 433 | 389 |
| 8573 | U | 60 | 360 | 294 |
| 8574 | U | 70 | 664 | 791 |
| 8575 | U | 71 | 342 | 23 |
| 8576 | U | 13 | 400 | 416 |
| 8577 | U | 55 | 2178 | 1743 |
| 8578 | U | 28 | 3221 | 3491 |
| 8579 | U | 57 | 883 | 633 |

Non-negative results shaded

Table C.74 Raw Data Returned By Laboratory 31 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 31 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 15957 | 13246 |
| 8513 | I | 26 | 12626482 | 9240995 |
| 8514 | I | 64 | 7886958 | 9224106 |
| 8515 | I | 35 | 4209 | 3653 |
| 8516 | I | 18 | 447490 | 293894 |
| 8517 | I | 5 | 268133 | 200267 |
| 8518 | I | 33 | 169472 | 139027 |
| 8519 | I | 24 | 28446 | 24628 |
| 8520 | I | 69 | 11852 | 13271 |
| 8521 | I | 38 | 14797 | 12644 |
| 8523 | I | 32 | 165124 | 130861 |
| 8525 | I | 8 | 11605 | 9176 |
| 8526 | I | 50 | 9641 | 7367 |
| 8528 | I | 62 | 8802033 | 13377612 |
| 8529 | I | 51 | 34990 | 28603 |
| 8530 | I | 21 | 1582 | 1351 |
| 8531 | I | 45 | 150004 | 126904 |
| 8532 | I | 49 | 39100 | 31920 |
| 8533 | I | 15 | 32517 | 24680 |
| 8534 | I | 65 | 7268 | 5456 |
| 8535 | I | 66 | 56525 | 46074 |
| 8536 | I | 44 | 575704 | 453872 |
| 8538 | I | 16 | 2364 | 2284 |
| 8539 | I | 48 | 12330 | 9698 |
| 8571 | I | 34 | 293134 | 229461 |
| 8572 | I | 41 | 18608 | 14508 |
| 8573 | I | 7 | 1438 | 1215 |
| 8574 | I | 40 | 630333 | 493351 |
| 8575 | I | 63 | 2167 | 1848 |
| 8576 | I | 52 | 121036 | 94883 |
| 8577 | I | 22 | 108293 | 84556 |
| 8578 | I | 37 | 46223 | 36797 |
| 8579 | I | 43 | 74121 | 56997 |
| 8531 | B | 72 | 9120 | 7457 |
| 8533 | B | 17 | 514 | 540 |
| 8535 | B | 30 | 1887 | 1154 |
| 8539 | B | 36 | 863 | 864 |
| 8577 | B | 2 | 3611 | 2978 |
| 8578 | B | 46 | 462 | 468 |

Table C.75 Raw Data Returned By Laboratory 33 – screening of unirradiated samples

| SP number | Status | Code number | Lab 33 | |
|-----------|--------|-------------|-----------|-----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 445 | 326 |
| 8513 | U | 31 | 670 | 810 |
| 8514 | U | 56 | 1392 | 1116 |
| 8515 | U | 42 | 320 | 308 |
| 8516 | U | 25 | 806 | 946 |
| 8517 | U | 58 | 126 | *16777011 |
| 8518 | U | 54 | 374 | 221 |
| 8519 | U | 39 | 433 | 442 |
| 8520 | U | 68 | 426 | 370 |
| 8521 | U | 3 | 323 | 347 |
| 8523 | U | 12 | 60 | 477 |
| 8525 | U | 1 | 453 | 643 |
| 8526 | U | 59 | *16777016 | 357 |
| 8528 | U | 10 | 517 | 272 |
| 8529 | U | 11 | 414 | 472 |
| 8530 | U | 23 | 366 | 361 |
| 8531 | U | 29 | 1837 | 1063 |
| 8532 | U | 27 | 3816 | 1542 |
| 8533 | U | 14 | 263 | 461 |
| 8534 | U | 4 | 852 | 357 |
| 8535 | U | 53 | 578 | 482 |
| 8536 | U | 19 | 1161 | 2280 |
| 8538 | U | 61 | 309 | 264 |
| 8539 | U | 20 | 348 | 236 |
| 8571 | U | 9 | 493 | 484 |
| 8572 | U | 6 | 467 | 235 |
| 8573 | U | 60 | 442 | 204 |
| 8574 | U | 70 | 894 | 991 |
| 8575 | U | 71 | 360 | 333 |
| 8576 | U | 13 | 279 | 354 |
| 8577 | U | 55 | 1899 | 1638 |
| 8578 | U | 28 | 1167 | 425 |
| 8579 | U | 57 | 429 | 289 |

Non-negative results shaded

* zero-crossing error – result negative

Table C.76 Raw Data Returned By Laboratory 33 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 33 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 20224 | 20671 |
| 8513 | I | 26 | 16168898 | 13067819 |
| 8514 | I | 64 | 13841231 | 2816721 |
| 8515 | I | 35 | 13234 | 17734 |
| 8516 | I | 18 | 773410 | 588560 |
| 8517 | I | 5 | 314380 | 324030 |
| 8518 | I | 33 | 214261 | 207000 |
| 8519 | I | 24 | 26398 | 49275 |
| 8520 | I | 69 | 14332 | 8056 |
| 8521 | I | 38 | 22734 | 41586 |
| 8523 | I | 32 | 223017 | 246611 |
| 8525 | I | 8 | 21154 | 73380 |
| 8526 | I | 50 | 39809 | 28018 |
| 8528 | I | 62 | 14824633 | 10787122 |
| 8529 | I | 51 | 69829 | 59587 |
| 8530 | I | 21 | 1226 | 1161 |
| 8531 | I | 45 | 171898 | 182307 |
| 8532 | I | 49 | 63790 | 78516 |
| 8533 | I | 15 | 47548 | 39384 |
| 8534 | I | 65 | 7887 | 6848 |
| 8535 | I | 66 | 89169 | 90837 |
| 8536 | I | 44 | 563570 | 552029 |
| 8538 | I | 16 | 3100 | 4552 |
| 8539 | I | 48 | 27824 | 14840 |
| 8571 | I | 34 | 500627 | 478735 |
| 8572 | I | 41 | 20060 | 31755 |
| 8573 | I | 7 | 2183 | 2897 |
| 8574 | I | 40 | 1093100 | 866587 |
| 8575 | I | 63 | 1976 | 939 |
| 8576 | I | 52 | 182756 | 163922 |
| 8577 | I | 22 | 189073 | 133944 |
| 8578 | I | 37 | 60051 | 69323 |
| 8579 | I | 43 | 137354 | 163131 |
| 8531 | B | 72 | 9081 | 19519 |
| 8533 | B | 17 | 761 | 685 |
| 8535 | B | 30 | 1667 | 671 |
| 8539 | B | 36 | 1295 | 1235 |
| 8577 | B | 2 | 8639 | 8566 |
| 8578 | B | 46 | 373 | 418 |

Table C.77 Raw Data Returned By Laboratory 34 – screening of unirradiated samples

| SP number | Status | Code number | Lab 34 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 640 | 512 |
| 8513 | U | 31 | 507 | 418 |
| 8514 | U | 56 | 3591 | 904 |
| 8515 | U | 42 | 246 | 240 |
| 8516 | U | 25 | 818 | 903 |
| 8517 | U | 58 | 320 | 458 |
| 8518 | U | 54 | 372 | 405 |
| 8519 | U | 39 | 293 | 324 |
| 8520 | U | 68 | 356 | 484 |
| 8521 | U | 3 | 341 | 314 |
| 8523 | U | 12 | 194 | 299 |
| 8525 | U | 1 | 426 | 523 |
| 8526 | U | 59 | 283 | 282 |
| 8528 | U | 10 | 532 | 427 |
| 8529 | U | 11 | 597 | 447 |
| 8530 | U | 23 | 263 | 250 |
| 8531 | U | 29 | 432 | 450 |
| 8532 | U | 27 | 447 | 326 |
| 8533 | U | 14 | 415 | 321 |
| 8534 | U | 4 | 292 | 263 |
| 8535 | U | 53 | 701 | 773 |
| 8536 | U | 19 | 313 | 414 |
| 8538 | U | 61 | 359 | 353 |
| 8539 | U | 20 | 464 | 340 |
| 8571 | U | 9 | 497 | 412 |
| 8572 | U | 6 | 276 | 417 |
| 8573 | U | 60 | 284 | 276 |
| 8574 | U | 70 | 3049 | 1325 |
| 8575 | U | 71 | 436 | 380 |
| 8576 | U | 13 | 447 | 559 |
| 8577 | U | 55 | 4319 | 2193 |
| 8578 | U | 28 | 402 | 609 |
| 8579 | U | 57 | 435 | 397 |

Non-negative results shaded

Table C.78 Raw Data Returned By Laboratory 34 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 34 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 34678 | 34693 |
| 8513 | I | 26 | 7803080 | 6066954 |
| 8514 | I | 64 | 16028790 | 15137986 |
| 8515 | I | 35 | 6883 | 6891 |
| 8516 | I | 18 | 342196 | 199866 |
| 8517 | I | 5 | 175984 | 307263 |
| 8518 | I | 33 | 18202 | 20360 |
| 8519 | I | 24 | 24951 | 29127 |
| 8520 | I | 69 | 14472 | 13318 |
| 8521 | I | 38 | 13481 | 14413 |
| 8523 | I | 32 | 27410 | 250732 |
| 8525 | I | 8 | 27005 | 22388 |
| 8526 | I | 50 | 24649 | 23545 |
| 8528 | I | 62 | 7464843 | 15034157 |
| 8529 | I | 51 | 92559 | 92600 |
| 8530 | I | 21 | 474 | 663 |
| 8531 | I | 45 | 185828 | 48241 |
| 8532 | I | 49 | 74165 | 93443 |
| 8533 | I | 15 | 27504 | 13633 |
| 8534 | I | 65 | 2923 | 4238 |
| 8535 | I | 66 | 56864 | 59995 |
| 8536 | I | 44 | 610369 | 123612 |
| 8538 | I | 16 | 2049 | 1548 |
| 8539 | I | 48 | 9850 | 17132 |
| 8571 | I | 34 | 418119 | 356221 |
| 8572 | I | 41 | 9849 | 12448 |
| 8573 | I | 7 | 1950 | 2878 |
| 8574 | I | 40 | 405238 | 359324 |
| 8575 | I | 63 | 2426 | 3748 |
| 8576 | I | 52 | 338152 | 411458 |
| 8577 | I | 22 | 91503 | 103111 |
| 8578 | I | 37 | 54204 | 41290 |
| 8579 | I | 43 | 41669 | 67018 |
| 8531 | B | 72 | 12259 | 7529 |
| 8533 | B | 17 | 443 | 505 |
| 8535 | B | 30 | 424 | 628 |
| 8539 | B | 36 | 929 | 825 |
| 8577 | B | 2 | 9964 | 9275 |
| 8578 | B | 46 | 544 | 365 |

Table C.79 Raw Data Returned By Laboratory 35 – screening of unirradiated samples

| SP number | Status | Code number | Lab 35 | |
|-----------|--------|-------------|--------|-------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 374 | 441 |
| 8513 | U | 31 | 490 | 467 |
| 8514 | U | 56 | 1419 | 1290 |
| 8515 | U | 42 | 690 | 299 |
| 8516 | U | 25 | 1394 | 4138 |
| 8517 | U | 58 | 677 | 463 |
| 8518 | U | 54 | 394 | 399 |
| 8519 | U | 39 | 454 | 307 |
| 8520 | U | 68 | 525 | 473 |
| 8521 | U | 3 | 477 | 369 |
| 8523 | U | 12 | 439 | 409 |
| 8525 | U | 1 | 808 | 1939 |
| 8526 | U | 59 | 277 | 477 |
| 8528 | U | 10 | 342 | 381 |
| 8529 | U | 11 | 496 | 503 |
| 8530 | U | 23 | 248 | 330 |
| 8531 | U | 29 | 1786 | 1396 |
| 8532 | U | 27 | 338 | 410 |
| 8533 | U | 14 | 383 | 370 |
| 8534 | U | 4 | 405 | 366 |
| 8535 | U | 53 | 773 | 1674 |
| 8536 | U | 19 | 1991 | 3310 |
| 8538 | U | 61 | 347 | 577 |
| 8539 | U | 20 | 647 | 569 |
| 8571 | U | 9 | 1580 | 699 |
| 8572 | U | 6 | 535 | 392 |
| 8573 | U | 60 | 369 | 323 |
| 8574 | U | 70 | 2444 | 1547 |
| 8575 | U | 71 | 248 | 297 |
| 8576 | U | 13 | 3391 | 722 |
| 8577 | U | 55 | 3117 | 3553 |
| 8578 | U | 28 | 566 | 550 |
| 8579 | U | 57 | 484 | 545 |

Non-negative results shaded

Table C.80 Raw Data Returned By Laboratory 35 – screening of irradiated samples and blends

| SP number | Status | Code number | Lab 35 | |
|-----------|--------|-------------|---------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 55652 | 38435 |
| 8513 | I | 26 | 6687907 | 13520225 |
| 8514 | I | 64 | 4855752 | 4691697 |
| 8515 | I | 35 | 40321 | 9952 |
| 8516 | I | 18 | 360637 | 328750 |
| 8517 | I | 5 | 457542 | 346754 |
| 8518 | I | 33 | 328567 | 373917 |
| 8519 | I | 24 | 142006 | 36820 |
| 8520 | I | 69 | 17373 | 17919 |
| 8521 | I | 38 | 32483 | 27417 |
| 8523 | I | 32 | 434999 | 285849 |
| 8525 | I | 8 | 34262 | 48209 |
| 8526 | I | 50 | 21189 | 14474 |
| 8528 | I | 62 | 2481252 | 4706400 |
| 8529 | I | 51 | 160651 | 92441 |
| 8530 | I | 21 | 1517 | 1461 |
| 8531 | I | 45 | 331744 | 268128 |
| 8532 | I | 49 | 105551 | 105819 |
| 8533 | I | 15 | 76534 | 65869 |
| 8534 | I | 65 | 14622 | 9549 |
| 8535 | I | 66 | 158956 | 166736 |
| 8536 | I | 44 | 1018542 | 1082550 |
| 8538 | I | 16 | 2347 | 3681 |
| 8539 | I | 48 | 31027 | 33223 |
| 8571 | I | 34 | 770119 | 762155 |
| 8572 | I | 41 | 31204 | 40733 |
| 8573 | I | 7 | 2361 | 3713 |
| 8574 | I | 40 | 1289735 | 1293027 |
| 8575 | I | 63 | 5885 | 7550 |
| 8576 | I | 52 | 301427 | 238782 |
| 8577 | I | 22 | 220781 | 155387 |
| 8578 | I | 37 | 90358 | 77964 |
| 8579 | I | 43 | 169695 | 225004 |
| 8531 | B | 72 | 68513 | 25644 |
| 8533 | B | 17 | 983 | 874 |
| 8535 | B | 30 | 1054 | 1031 |
| 8539 | B | 36 | 7384 | 3637 |
| 8577 | B | 2 | 11890 | 6550 |
| 8578 | B | 46 | 527 | 503 |

Table C.81 Raw Data Returned By Laboratory 35 – calibrated PSL for unirradiated samples

| SP number | Status | Code number | Lab 35 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | U | 67 | 28255 | 34630 |
| 8513 | U | 31 | 11869210 | 9413035 |
| 8514 | U | 56 | 9453898 | 10651708 |
| 8515 | U | 42 | 7588 | 10981 |
| 8516 | U | 25 | 227744 | 140528 |
| 8517 | U | 58 | 410011 | 305164 |
| 8518 | U | 54 | 219644 | 221910 |
| 8519 | U | 39 | 24673 | 19165 |
| 8520 | U | 68 | 15135 | 7578 |
| 8521 | U | 3 | 20469 | 28468 |
| 8523 | U | 12 | 203771 | 321662 |
| 8525 | U | 1 | 34885 | 31967 |
| 8526 | U | 59 | 14855 | 13056 |
| 8528 | U | 10 | 13145418 | 13028800 |
| 8529 | U | 11 | 67459 | 43159 |
| 8530 | U | 23 | 1416 | 1616 |
| 8531 | U | 29 | 161358 | 166620 |
| 8532 | U | 27 | 58543 | 62625 |
| 8533 | U | 14 | 41800 | 39007 |
| 8534 | U | 4 | 9052 | 6360 |
| 8535 | U | 53 | 102496 | 90155 |
| 8536 | U | 19 | 781618 | 741704 |
| 8538 | U | 61 | 3039 | 4926 |
| 8539 | U | 20 | 22836 | 13477 |
| 8571 | U | 9 | 140340 | 102332 |
| 8572 | U | 6 | 14369 | 18451 |
| 8573 | U | 60 | 1608 | 2603 |
| 8574 | U | 70 | 1012745 | 1226859 |
| 8575 | U | 71 | 5303 | 4687 |
| 8576 | U | 13 | 181198 | 161513 |
| 8577 | U | 55 | 579590 | 673549 |
| 8578 | U | 28 | 55067 | 98299 |
| 8579 | U | 57 | 82637 | 130128 |

Table C.82 Raw Data Returned By Laboratory 35 – calibrated PSL for irradiated samples and blends

| SP number | Status | Code number | Lab 35 | |
|-----------|--------|-------------|----------|----------|
| | | | Pot 1 | Pot 2 |
| 8512 | I | 47 | 48576 | 72308 |
| 8513 | I | 26 | 965942 | 11089817 |
| 8514 | I | 64 | 10908597 | 156770 |
| 8515 | I | 35 | 13135 | 39916 |
| 8516 | I | 18 | 905861 | 469724 |
| 8517 | I | 5 | 632458 | 577179 |
| 8518 | I | 33 | 332197 | 373193 |
| 8519 | I | 24 | 49778 | 145808 |
| 8520 | I | 69 | 15806 | 20547 |
| 8521 | I | 38 | 36864 | 39136 |
| 8523 | I | 32 | 283921 | 337326 |
| 8525 | I | 8 | 48120 | 26110 |
| 8526 | I | 50 | 34434 | 31289 |
| 8528 | I | 62 | 8600844 | 14840439 |
| 8529 | I | 51 | 146832 | 161489 |
| 8530 | I | 21 | 2364 | 2624 |
| 8531 | I | 45 | 629832 | 353883 |
| 8532 | I | 49 | 129985 | 109207 |
| 8533 | I | 15 | 80269 | 71538 |
| 8534 | I | 65 | 17763 | 12250 |
| 8535 | I | 66 | 143069 | 133099 |
| 8536 | I | 44 | 1446309 | 1233774 |
| 8538 | I | 16 | 5752 | 6274 |
| 8539 | I | 48 | 37396 | 102911 |
| 8571 | I | 34 | 1147262 | 1028591 |
| 8572 | I | 41 | 49358 | 54279 |
| 8573 | I | 7 | 2557 | 3368 |
| 8574 | I | 40 | 1771630 | 1572129 |
| 8575 | I | 63 | 7013 | 7631 |
| 8576 | I | 52 | 288663 | 299463 |
| 8577 | I | 22 | 201707 | 206819 |
| 8578 | I | 37 | 129926 | 94605 |
| 8579 | I | 43 | 1898375 | 295421 |
| 8531 | B | 72 | 218681 | 179140 |
| 8533 | B | 17 | 38061 | |
| 8535 | B | 30 | 82009 | 71075 |
| 8539 | B | 36 | 27824 | 16070 |
| 8577 | B | 2 | 633541 | 509575 |
| 8578 | B | 46 | 50535 | 63693 |

| | | | | | | | | | |
|------------------------|--------------|--------------|--------|---------------------------|--------------|--------------|--------|---------------------------|--|
| Lab No. | 1 | | | | | | | | |
| TL reader type | TLD 3500 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 0,37 nC | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 235 | | | | | | | | |
| Chip 2 | 233 | | | | | | | | |
| Mean | 234 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 6.49 | 5810 | 0.001 | | 16.38 | 7830 | 0.002 | | untreated |
| 2 | 13230 | 10750 | 1.23 | | 6231 | 4208 | 1.48 | | irradiated |
| 3 | 73.16 | 673 | 0.109 | 2 peaks by channel 60+100 | 160.3 | 1881 | 0.0085 | 2 peaks by channel 60+100 | mixtures containing irradiated materials |
| 4 | 27.51 | 26.7 | 1.03 | | 1.46 | 10.51 | 0.14 | | irradiated |
| 5 | 0.22 | 15.91 | 0.014 | | 0.46 | 18.46 | 0.025 | | untreated |
| 6 | 0.18 | 46.99 | 0.004 | | 0.41 | 32.91 | 0.012 | | untreated |
| 7 | 2.836 | 299.9 | 0.01 | | 1.344 | 2435 | 0.0006 | | untreated |
| 8 | 1.12 | 4600 | 0.0002 | | 1.226 | 2528 | 0.0005 | | untreated |
| 9 | 2322 | 1917 | 1.21 | | 916.5 | 1029 | 0.89 | | irradiated |
| 10 | 0.2 | 52.05 | 0.004 | | 0.16 | 53.39 | 0.003 | | untreated |
| 11 | 164.4 | 83.49 | 1.97 | | 43.76 | 55.02 | 0.8 | | irradiated |
| 12 | 2.44 | 131.7 | 0.019 | 2 peaks by channel 59+100 | 2.4 | 191.3 | 0.013 | 2 peaks by channel 61+100 | mixtures containing irradiated materials |
| 13 | 11560 | 6990 | 1.65 | | 9760 | 12860 | 0.76 | | irradiated |
| 14 | 130.5 | 3583 | 0.036 | peak-max. by channel 75 | 27.17 | 1660 | 0.016 | peak-max.by channel 80 | mixtures containing irradiated materials |
| 15 | 31.89 | 6571 | 0.005 | | 10.1 | 2653 | 0.004 | | untreated |
| 16 | 2.41 | 2073 | 0.001 | | 10.19 | 3061 | 0.003 | | untreated |
| 17 | 4.496 | 3742 | 0.001 | | 3.389 | 2185 | 0.002 | | untreated |
| 18 | 4695 | 3045 | 1.54 | | 2542 | 2301 | 1.1 | | irradiated |

Table C.83 TL data for TL laboratory 1

| | | | | | | | | | |
|------------------------|------------------|--------------|-------|--|--------------|--------------|-------|--|--|
| Lab No. | 2.00 | | | | | | | | |
| TL reader type | Harshaw TLD 4000 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 7.50 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 273.50 | | | | | | | | |
| Chip 2 | 277.00 | | | | | | | | |
| Mean | 275.25 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 28.03 | 8648.23 | 0.00 | maximum 369°C | 2.52 | 1915.77 | 0.00 | maximum 413°C | untreated |
| 2 | 8671.23 | 3920.89 | 2.21 | maximum at 245°C with shoulder in front of the peak; lower, broad maximum 382°C | 7933.50 | 3743.26 | 2.12 | maximum at 245°C with shoulder in front of the peak; lower, broad maximum 360°C | irradiated |
| 3 | 296.52 | 1494.26 | 0.20 | maximum at 249°C with shoulder in front of the peak; lower, broad maximum 403°C | 122.78 | 1528.71 | 0.08 | maximum at 242°C with shoulder in front of the peak; lower, broad maximum 374°C | irradiated, only a part of a mixture, because |
| 4 | | | | | | | | | no possible |
| 5 | | | | | | | | | no possible |
| 6 | | | | | | | | | no possible |
| 7 | 0.84 | 408.11 | 0.00 | maximum at 327°C | 0.35 | 660.51 | 0.00 | no maximum | untreated |
| 8 | 0.54 | 433.42 | 0.00 | broad maximum at 400°C | 0.55 | 821.97 | 0.00 | broad plateau 313°C - 400°C | untreated |
| 9 | 553.50 | 409.04 | 1.35 | round peak maximum at 229°C | 848.94 | 791.58 | 1.07 | round peak maximum at 220°C | irradiated |
| 10 | 0.30 | 151.13 | 0.00 | maximum at 400°C, like blank | 0.31 | 135.02 | 0.00 | maximum at 400°C, like blank | untreated |
| 11 | 983.36 | 708.65 | 1.39 | maximum at 221°C, lower peak 393°C | 20.79 | 11.93 | 1.74 | maximum at 241°C, lower peak 391°C | irradiated |
| 12 | 0.36 | 13.25 | 0.03 | maximum at 400°C, like blank not enough minerals | 5.56 | 27.33 | 0.20 | sharp maximum at 244°C, lower peak 391°C | irradiated, only a part of a mixture, because |
| 13 | 1316.69 | 1220.52 | 1.08 | maximum at 213°C | 748.55 | 566.32 | 1.32 | maximum at 248°C with fronting | irradiated |
| 14 | 3.49 | 261.08 | 0.01 | plateau from 295-390°C | 4.58 | 640.37 | 0.01 | broad maximum at 361°C | untreated |
| 15 | 4.69 | 1204.83 | 0.00 | plateau from 295-390°C | 7.91 | 954.11 | 0.01 | broad peak, maximum at 296°C | untreated |
| 16 | 17.41 | 1651.68 | 0.01 | lower maximum 209°C, broad peak maximum 391°C | 3.64 | 1078.42 | 0.00 | broad peak, maximum at 362°C | mixture irradiated material |
| 17 | 9.34 | 3143.09 | 0.00 | plateau from 329-383°C | 4.78 | 2849.53 | 0.00 | plateau from 317-370°C | untreated |
| 18 | 2285.56 | 1531.60 | 1.49 | maximum 235°C, fronting | 3150.55 | 2151.98 | 1.46 | maximum 248°C, fronting | irradiated |

Table C.84 TL data for TL laboratory 2

| | | | | | | | | | |
|------------------------|--------------|--------------|-------|---------------------------|--------------|--------------|-------|---------------------------|--|
| Lab No. | 3 | | | | | | | | |
| TL reader type | TLD 3500 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 3,12nC | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 241 | | | | | | | | |
| Chip 2 | 246 | | | | | | | | |
| | | | | | | | | | |
| Mean | 243.5 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 4.375 | 897 | 0.00 | no peak | 6.159 | 47.93 | 0.13 | no peak | untreated |
| 2 | 207.5 | 17.08 | 12.15 | breadth peak in ROI1 | 5429 | 920 | 5.90 | peak in ROI1 | irradiated |
| 3 | 384.3 | 71.59 | 5.37 | peak in ROI1 | 87.35 | 164.5 | 0.53 | small peak | irradiated |
| 4 | 1.448 | 0.7514 | 1.93 | little peak in ROI1 | 1.137 | 7.495 | 0.15 | little peak in ROI1 | (irradiated) |
| 5 | 0.444 | 6.541 | 0.07 | no peak | 0.3103 | 5.749 | 0.05 | no peak | (untreated) |
| 6 | 0.5627 | 3.813 | 0.15 | no peak | 0.2581 | 3.911 | 0.07 | no peak | (untreated) |
| 7 | 4.413 | 72.58 | 0.06 | no peak | 1.143 | 113.2 | 0.01 | no peak | untreated |
| 8 | 3.542 | 160.8 | 0.02 | little breadth peak | 11.43 | 302.5 | 0.04 | little breadth peak | untreated |
| 9 | 220.1 | 108.5 | 2.03 | peak in ROI1 | 553.9 | 179.3 | 3.09 | peak in ROI1 | irradiated |
| 10 | 0.2809 | 0.9046 | 0.31 | no peak | 0.5128 | 14 | 0.04 | no peak | (untreated) |
| 11 | 1.437 | 0.8806 | 1.63 | little breadth peak | 1.223 | 0.5657 | 2.16 | little breadth peak | (irradiated) |
| 12 | 1.039 | 1.049 | 0.99 | little peak in ROI1 | 4.038 | 15.67 | 0.26 | little peak in ROI1 | (mixtures containing irradiated materials) |
| 13 | 662.5 | 286.1 | 2.32 | peak in ROI1 | 1471 | 527 | 2.79 | peak in ROI1 | irradiated |
| 14 | 53.65 | 1291 | 0.04 | breadth peak higher temp. | 10.07 | 416.3 | 0.02 | no peak | untreated |
| 15 | 36.9 | 1921 | 0.02 | breadth peak higher temp. | 16.19 | 544.5 | 0.03 | breadth peak higher temp. | untreated |
| 16 | 6.962 | 159.8 | 0.04 | little peak in ROI1 | 0.5932 | 54.05 | 0.01 | no peak | mixtures containing irradiated materials |
| 17 | 0.4123 | 53.25 | 0.01 | no peak | 3.489 | 815.1 | 0.00 | no peak | untreated |
| 18 | 1105 | 243.3 | 4.54 | peak in ROI1 | 28.47 | 5.773 | 4.93 | peak in ROI1 | irradiated |

Table C.85 TL data for TL laboratory 3

| Lab No. | 4 | | | | | | | | |
|------------------------|-------------------------|--------------|--------|------------------|--------------|--------------|------------|------------------|------------|
| TL reader type | Harshaw TLD-Reader 3500 | | | | | | | | |
| Units for intensity | nano Coulomb (nC) | | | | | | | | |
| MDL | 20 nC | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 245 | | | | | | | | |
| Chip 2 | 245 | | | | | | | | |
| Mean | 245 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 48 | 4856 | 0.0099 | no peak | | | | | untreated |
| 2 | 5253 | 7200 | 0.7296 | gaussian, normal | | | | | irradiated |
| 3 | 645 | 3523 | 0.1831 | gaussian, normal | | | | | irradiated |
| 4 | 98 | 85 | 1.1529 | gaussian, normal | | | | | irradiated |
| 5 | 0 | 48 | 0.0000 | no peak | | | | | untreated |
| 6 | 5 | 10 | 0.5000 | no peak | | | | | untreated |
| 7 | 51 | 2687 | 0.0190 | gaussian, normal | | | | | irradiated |
| 8 | 0 | 1877 | 0.0000 | no peak | 0 | 2447 | 0 | no peak | untreated |
| 9 | 2282 | 1401 | 1.6288 | gaussian, normal | 1798 | 1968 | 0.91361789 | gaussian, normal | irradiated |
| 10 | 169 | 119 | 1.4202 | gaussian, normal | | | | | irradiated |
| 11 | 72 | 36 | 2.0000 | gaussian, normal | | | | | irradiated |
| 12 | 9 | 2 | 4.5000 | very small peak | | | | | ???? |
| 13 | 3078 | 1711 | 1.7989 | gaussian, normal | | | | | irradiated |
| 14 | 74 | 1489 | 0.0497 | no peak | | | | | untreated |
| 15 | 27 | 487 | 0.0554 | no gaussian | | | | | irradiated |
| 16 | 119 | 824 | 0.1444 | no peak | | | | | untreated |
| 17 | 28 | 715 | 0.0392 | gaussian, normal | | | | | irradiated |
| 18 | 8 | 1374 | 0.0058 | gaussian, normal | | | | | irradiated |

Table C.86 TL data for TL laboratory 4

| | | | | | | | | | |
|---------------------|--------------|--------------|--------------------------|----------------------|--------------|------------------------|-------|----------------------|------------|
| Lab No. | 5 | | | | | | | | |
| TL reader type | Harshaw 3500 | | MDL | 0.0238 | | LiF data | | | |
| Units for intensity | nC | | | 0.0153 | | Peak V temperature /°C | | | |
| | | | | 0.0285 | | Chip 1 | 91 | | |
| | | | | 0.0311 | | Chip 2 | 90 | | |
| | | | | 0.0321 | | Mean | 90.5 | | |
| | | | Average x | 0.0262 | | | | | |
| | | | Standard deviation .xs n | 0.0069 | | | | | |
| | | | MDL (x + 3s) | 0.0468 | | | | | |
| | | | 10*MDL | 0.468 | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 6.80 | 2256 | 0.00 | untreated | 3.79 | 1108 | 0.00 | untreated | negativ |
| 2 | 16710 | 4755 | 3.5 | irradiated | 5320 | 1489 | 3.6 | irradiated | positiv |
| 3 | 348 | 1693 | 0.21 | irradiated (mixture) | 352 | 924 | 0.38 | irradiated (mixture) | positiv |
| 4 | 150 | 58.7 | 2.6 | irradiated | | | | | positiv |
| 5 | 5.25 | 28.7 | 0.18 | irradiated (mixture) | | | | | positiv |
| 6 | 0.05 | 0.29 | 0.17 | | | | | | <10 MDL |
| 7 | 6.61 | 1956 | 0.00 | untreated | 6.4 | 2503 | 0.00 | untreated | negativ |
| 8 | 3.77 | 4111 | 0.00 | untreated | | | | | negativ |
| 9 | 11260 | 5510 | 2.0 | irradiated | 10740 | 4387 | 2.4 | irradiated | positiv |
| 10 | 4.72 | 48.6 | 0.1 | irradiated (mixture) | | | | | positiv |
| 11 | 338 | 122 | 2.8 | irradiated | | | | | positiv |
| 12 | 19.1 | 10.8 | 1.8 | irradiated | | | | | positiv |
| 13 | 14930 | 4003 | 3.7 | irradiated | 14690 | 2673 | 5.5 | irradiated | positiv |
| 14 | 310 | 3819 | 0.1 | irradiated (mixture) | 106 | 1247 | 0.1 | irradiated (mixture) | positiv |
| 15 | 50.8 | 5839 | 0.01 | untreated | 62.5 | 3169 | 0.02 | untreated | negativ |
| 16 | 12.4 | 3201 | 0.00 | untreated | | | | | negativ |
| 17 | 14.3 | 2592 | 0.01 | untreated | | | | | negativ |
| 18 | 8475 | 2618 | 3.2 | irradiated | 7572 | failed | | irradiated | positiv |

Table C.87 TL data for TL laboratory 5

| | | | | | | | | | |
|------------------------|---------------------|---------------------|--------------|----------------------|---------------------|---------------------|--------------|----------------------|-------------------|
| Lab No. | 6 | | | | | | | | |
| TL reader type | HARSHAW TLD 3500 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 0,098 nC | | | | | | | | |
| Heating rate | 6 oC/s | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 235.6 | | | | | | | | |
| Chip 2 | 233.8 | | | | | | | | |
| Mean | 234.7 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 3.98 | 1685 | 0.002 | untreated | 5.5 | 2836 | 0.002 | untreated | untreated |
| 2 | 5584 | 3001 | 1.861 | irradiated | 3389 | 1392 | 2.435 | irradiated | irradiated |
| 3 | 749 | 1683 | 0.445 | irradiated | 118 | 936 | 0.126 | irradiated | mixture |
| 4 | | | | | | | | | |
| 5 | | | | | | | | | |
| 6 | | | | | | | | | |
| 7 | 0.81 | 1278 | 0.001 | untreated | | | | | |
| 8 | 1.01 | 1532 | 0.001 | untreated | 1.16 | 1830 | 0.001 | untreated | untreated |
| 9 | 1756 | 1304 | 1.347 | irradiated | | | | | irradiated |
| 10 | 0.15 | 45 | 0.003 | untreated | | | | | untreated |
| 11 | 387 | 382 | 1.013 | irradiated | | | | | irradiated |
| 12 | 8.31 | 194 | 0.043 | irradiated | | | | | mixture |
| 13 | 3698 | 2857 | 1.294 | irradiated | 3363 | 2476 | 1.358 | irradiated | irradiated |
| 14 | 33.56 | 2096 | 0.016 | untreated | 3.28 | 581 | 0.006 | untreated | untreated |
| 15 | 7.66 | 1890 | 0.004 | untreated | 18.72 | 2194 | 0.009 | untreated | untreated |
| 16 | 2.21 | 741 | 0.003 | untreated | 1.28 | 736 | 0.002 | untreated | untreated |
| 17 | 0.92 | 802 | 0.001 | untreated | 1.4 | 597 | 0.002 | untreated | untreated |
| 18 | 1064 | 718 | 1.482 | irradiated | 1840 | 1025 | 1.795 | irradiated | irradiated |

Table C.88 TL data for TL laboratory 6

| | | | | | | | | | |
|------------------------|---------------|--------------|-------|--|--|--------------|--------------|-------|--|
| Lab No. | 7 | | | | | | | | |
| TL reader type | Harshaw M3500 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 0 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 260 | | | | | | | | |
| Chip 2 | 250 | | | | | | | | |
| Mean | 255 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | | Aliquot B | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape |
| 1 | 3 | 42 | 0.06 | no distinct peak in relevant temperature range | | 5 | 101 | 0.05 | no distinct peak in relevant temperature range |
| 2 | 518 | 24 | 21.75 | distinct peak at 225°C | | 817 | 86 | 9.47 | distinct peak at 225°C |
| 3 | 38 | 8 | 4.64 | distinct peak at 215°C | | 113 | 13 | 8.61 | distinct peak at 215°C |
| 4 | 48 | 34 | 1.42 | distinct peak at 200°C | | 10 | 8 | 1.15 | distinct peak at 200°C |
| 5 | 0 | 3 | 0.05 | No signal in the relevant temperature range | | 0 | 2 | 0.07 | No signal in the relevant temperature range |
| 6 | 0 | 7 | 0.03 | No signal in the relevant temperature range | | 0 | 1 | 0.11 | No signal in the relevant temperature range |
| 7 | 5 | 42 | 0.13 | flat signal with peak at 200°C | | 2 | 58 | 0.03 | flat signal with peak at 200°C |
| 8 | 1 | 42 | 0.03 | no distinct peak in relevant temperature range | | 2 | 74 | 0.02 | no distinct peak in relevant temperature range |
| 9 | 200 | 56 | 3.58 | distinct peak at 220°C | | 590 | 251 | 2.35 | distinct peak at 220°C |
| 10 | 0 | 7 | 0.03 | No signal in the relevant temperature range | | 0 | 30 | 0.01 | No signal in the relevant temperature range |
| 11 | 18 | 12 | 1.50 | small peak at 200°C | | 5 | 4 | 1.41 | small peak at 200°C |
| 12 | 12 | 22 | 0.54 | small peak at 210°C | | 0 | 2 | 0.12 | small peak at 210°C |
| 13 | 8455 | 1594 | 5.30 | very large peak at 200°C | | 7801 | 2622 | 2.98 | very large peak at 200°C |
| 14 | 96 | 5270 | 0.02 | signal starts at 140°C and rises from there on up to 400°C. No distinct peak in relevant temperature range | | 99 | 2680 | 0.04 | signal starts at 140°C and rises from there on up to 400°C. No distinct peak in relevant temperature range |
| 15 | 2 | 320 | 0.01 | signal starts in the relevant temperature range, but the small peak is at 300°C and therefore not within the relevant temperature range. | | 2 | 332 | 0.00 | signal starts in the relevant temperature range, but the small peak is at 300°C and therefore not within the relevant temperature range. |
| 16 | 35 | 1691 | 0.02 | signal starts at 90°C and rises from there on up to 400°C. No distinct peak in relevant temperature range | | 7 | 1755 | 0.00 | signal starts at 90°C and rises from there on up to 400°C. No distinct peak in relevant temperature range |
| 17 | 4 | 785 | 0.00 | Signal starts before the relevant temperature range and rises through it up to 400°C without a distinct peak | | 9 | 1845 | 0.00 | Signal starts before the relevant temperature range and rises through it up to 400°C without a distinct peak |
| 18 | 13700 | 1572 | 8.72 | Very high signal with distinct peak at 210°C | | 17550 | 3203 | 5.48 | Very high signal with distinct peak at 210°C |

Table C.89 TL data for TL laboratory 7

| | | | | | | | | | |
|----------------------------------|------------------|--------------|--------|--|--------------|--------------|--------|--|--|
| Lab No. | 8 | | | | | | | | |
| TL reader type | Risø TL/OSL-DA15 | | | | | | | | |
| Units for intensity | counts | | | | | | | | |
| MDL (full Temperaturintervall I) | 502 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 254 | | | | | | | | |
| Chip 2 | 254 | | | | | | | | |
| Mean | 254 | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 118 135 | 28 316 457 | 0.004 | no Max. in temperatur interval I | 357 977 | 36 675 782 | 0.01 | no Max. in temperatur interval I | not irradiated |
| 2 | 48 669 589 | 33 114 212 | 1.47 | wide glow-curve, Max. in temperatur interval I | 23 724 372 | 14 433 661 | 1.64 | wide glow-curve, Max. in temperatur interval I | irradiated |
| 3 | 7 843 999 | 30 551 003 | 0.26 | wide glow-curve, Max. in temperatur interval I | 41 753 048 | 63 846 699 | 0.65 | wide glow-curve, Max. in temperatur interval I | irradiated (probably mixture of irradiated |
| 4 | 181 043 | 166 066 | 1.09 | Max. in temperatur interval I | 3 080 | 4 087* | | small Max. in temperatur interval I | irradiated (conclusion from mea |
| 5 | 327 | 663 191 | 0.0005 | no Max. in temperatur interval I | 366 | 1 321 061 | 0.0003 | no Max. in temperatur interval I | not irradiated |
| 6 | 348 | 7 404 | 0.05 | no Max. in temperatur interval I | 353 | 843* | | | probably (conclusion from measurement |
| 7 | 10 819 | 21 683 979 | 0.0005 | no Max. in temperatur interval I | 3 093 | 9 332 116 | 0.0003 | no Max. in temperatur interval I | not irradiated |
| 8 | 11 483 | 19 475 343 | 0.0006 | no Max. in temperatur interval I | 18 343 | 35 595 880 | 0.0005 | no Max. in temperatur interval I | not irradiated |
| 9 | 41 580 251 | 37 680 886 | 1.1 | Max. in temperatur interval I | 25 775 572 | 22 005 892 | 1.17 | Max. in temperatur interval I | irradiated |
| 10 | 704 | 4 883 669 | 0.0001 | no Max. in temperatur interval I | 2 067 | 4 049 836 | 0.0005 | no Max. in temperatur interval I | not irradiated |
| 11 | 1 580 084 | 1 139 931 | 1.39 | Max. in temperatur interval I | 2 891 780 | 2 324 858 | 1.24 | Max. in temperatur interval I | irradiated |
| 12 | 980 358 | 11 733 319 | 0.08 | Max. in temperatur interval I | 656 250 | 7 717 984 | 0.09 | Max. in temperatur interval I | irradiated (probably mixture of irradiated |
| 13 | 29 353 880 | 26 282 946 | 1.12 | Max. in temperatur interval I | 18 959 289 | 18 616 101 | 1.02 | Max. in temperatur interval I | irradiated |
| 14 | 771 205 | 15 365 051 | 0.05 | Max. in temperatur interval I | 77 171 | 19 865 724 | 0.004 | no Max. in temperatur interval I | irradiated, conclusion (probably mixture of a small am |
| 15 | 116 508 | 12 889 750 | 0.009 | no Max. in temperatur interval I | 341 299 | 49 398 361 | 0.007 | no Max. in temperatur interval I | not irradiated |
| 16 | 22 018 | 10 331 106 | 0.002 | no Max. in temperatur interval I | 13 278 | 6 083 037 | 0.002 | no Max. in temperatur interval I | not irradiated |
| 17 | 9 715 | 2 617 519 | 0.004 | no Max. in temperatur interval I | 23 389 | 7 876 940 | 0.003 | no Max. in temperatur interval I | not irradiated |
| 18 | 26 625 387 | 19 113 548 | 1.39 | Max. in temperatur interval I | 15 056 059 | 10 622 814 | 1.42 | Max. in temperatur interval I | irradiated |

Table C.90 TL data for TL laboratory 8

| | | | | | | | | |
|---------------------|------------------------|--------------|---------|--------------------------------|--------------|--------------|---------|--------------------------------|
| Lab No. | 10 | | | | | | | |
| TL reader type | Harshaw 3500 | | | | | | | |
| Units for intensity | nC | | | | | | | |
| MDL | 0.54 | | | | | | | |
| | | | | | | | | |
| LiF data | Peak V temperature /°C | | | | | | | |
| Chip 1 | 248.00 | | | | | | | |
| Chip 2 | 244.00 | | | | | | | |
| | | | | | | | | |
| Mean | 246.00 | | | | | | | |
| | | | | | | | | |
| Sample data | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape |
| 1 | 3.41 | 1483.00 | 0.00230 | weak peak at about 300°C | 4.32 | 2703.00 | 0.00160 | weak peak at about 300°C |
| 2 | 13030.00 | 1681.00 | 7.75134 | peak at about 212°C | 7341.00 | 2743.00 | 2.67627 | peak at about 229°C |
| 3 | 773.20 | 3868.00 | 0.19990 | peak at about 220°C | 797.60 | 4162.00 | 0.19164 | peak at about 220°C |
| 4 | 2.56 | 13.17 | 0.19438 | peak at about 200°C | 1.10 | 14.28 | 0.07675 | peak at about 200°C |
| 5 | 0.20 | 7.01 | 0.02853 | no peaks are visible | 0.21 | 8.37 | 0.02557 | no peaks are visible |
| 6 | 0.26 | 62.71 | 0.00415 | no peaks are visible | 0.26 | 64.56 | 0.00397 | no peaks are visible |
| 7 | 31.79 | 15848.00 | 0.00201 | peaks at about 200 and 300°C | 26.12 | 12478.00 | 0.00209 | peaks at 200 and 300°C |
| 8 | 11.26 | 10234.00 | 0.00110 | no peaks are visible | 4.50 | 6656.00 | 0.00068 | no peaks are visible |
| 9 | 6443.00 | 6441.00 | 1.00031 | peak at about 200°C | 1265.00 | 1707.00 | 0.74107 | peak at about 200°C |
| 10 | 0.28 | 15.53 | 0.01777 | no peaks are visible | 0.25 | 165.50 | 0.00153 | no peaks are visible |
| 11 | 55.05 | 58.51 | 0.94086 | peak at about 200°C | 0.76 | 10.77 | 0.07057 | peak at about 200°C |
| 12 | 0.43 | 10.02 | 0.04281 | very weak peak at about 200 °C | 0.41 | 7.11 | 0.05753 | very weak peak at about 200 °C |
| 13 | 34700.00 | 23273.00 | 1.49100 | peak at about 200°C | 35300.00 | 24330.00 | 1.45088 | peak at about 200°C |
| 14 | 240.20 | 20469.00 | 0.01173 | weak peak at about 200°C | 243.30 | 16973.00 | 0.01433 | weak peak at about 200°C |
| 15 | 50.58 | 13900.00 | 0.00364 | weak peak at about 300°C | 27.70 | 8778.00 | 0.00316 | weak peak at 300°C |
| 16 | 9.93 | 5302.00 | 0.00187 | peaks at about 200 and 300°C | 6.68 | 3287.00 | 0.00203 | peaks at 200 and 300°C |
| 17 | 4.14 | 10880.00 | 0.00038 | weak peak at about 300 °C | 10.74 | 8670.00 | 0.00124 | no peaks are visible |
| 18 | 10397.00 | 10781.00 | 0.96438 | peak at about 225 °C | 13840.00 | 13907.00 | 0.99518 | peak at about 225 °C |

Table C.91 TL data for TL laboratory 10

| | | | | | | | | | |
|------------------------|----------------|--------------|-------|------------------|--------------|--------------|-------|------------------|--------------|
| Lab No. | TL 11 | | | | | | | | |
| TL reader type | TL-DA-15; RISO | | | | | | | | |
| Units for intensity | | | | | | | | | |
| MDL | 358.49 | | | | | | | | |
| LiF data | | | | | | | | | |
| | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 255 | | | | | | | | |
| Chip 2 | 259 | | | | | | | | |
| Mean | 257 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 65845 | 14792736 | 0.004 | max. 340 | 28386 | 14514844 | 0.002 | max. 340 | unirradiated |
| 2 | 49826572 | 27014391 | 1.844 | max. 200 and 340 | 50477721 | 23558922 | 2.143 | max. 200 and 340 | mixture |
| 3 | 4207445 | 23531350 | 0.179 | max. 200 and 340 | 2801563 | 9169738 | 0.306 | max. 200 and 340 | mixture |
| 4 | 1023013 | 567270 | 1.803 | max. 200 | 803912 | 493360 | 1.629 | max. 200 | irradiated |
| 5 | -13 | 28382 | 0.000 | bckgrd. | 241 | 125206 | 0.002 | bckgrd. | unirradiated |
| 6 | 128 | 56752 | 0.002 | bckgrd. | 59 | 146180 | 0.000 | bckgrd. | unirradiated |
| 7 | 1143 | 2156883 | 0.001 | max. 340 | 8907 | 6385980 | 0.001 | max. 230 and 340 | mixture |
| 8 | 50735 | 19832392 | 0.003 | max. 300 | 17927 | 14962105 | 0.001 | max. 300 | unirradiated |
| 9 | 22523123 | 10090741 | 2.232 | max. 200 | 47781830 | 25135827 | 1.901 | max. 200 | irradiated |
| 10 | 386 | 1307815 | 0.000 | bckgrd. | 349 | 204658 | 0.002 | bckgrd. | unirradiated |
| 11 | 254682 | 98053 | 2.597 | max. 200 | 895046 | 451928 | 1.981 | max. 200 | irradiated |
| 12 | 37323 | 114953 | 0.325 | max. 200 | 1002 | 102507 | 0.010 | bckgrd. | irradiated |
| 13 | 129255206 | 81273060 | 1.590 | max. 200 | 117042721 | 84043022 | 1.393 | max. 200 | irradiated |
| 14 | 1339775 | 47970540 | 0.028 | max. 200 and 300 | 957453 | 48487761 | 0.020 | max. 300 | mixture |
| 15 | 505550 | 53504181 | 0.009 | max. 270 - 340 | 364121 | 24918540 | 0.015 | max. 270 - 340 | unirradiated |
| 16 | 473832 | 5382455 | 0.088 | max. 200 and 300 | 13077 | 3923832 | 0.003 | max. 300 | mixture |
| 17 | 31727 | 7395225 | 0.004 | max. 340 | 44960 | 9492378 | 0.005 | max. 300-340 | unirradiated |
| 18 | 25774335 | 14059233 | 1.833 | max. 200 | 18346005 | 11413995 | 1.607 | max. 200 and 340 | irradiated |

Table C.92 TL data for TL laboratory 11

| | | | | | | | | |
|----------------|----------------|------------------------|---------------------|--|--------------|--------------|-------|--|
| Lab No. | 13 | | | | | | | |
| TL reader type | Risoe TL-DA-15 | | Units for intensity | counts | MDL | 2492 | | |
| LiF data | | Peak V temperature /°C | | | | | | |
| Chip 1 | | 263 | | | | | | |
| Chip 2 | xxx | | | | | | | |
| | Mean | 263 | | | | | | |
| | Aliquot A | | | | Aliquot B | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape |
| | 186-245°C | | | | | | | |
| 1 | 75016 | 26675090 | 0.003 | Max at 350°C | 119775 | 2.5E+07 | 0.005 | Max at 350°C |
| 2 | 41434815 | 16419165 | 2.524 | Max at 215°C; shoulder at 320°C | 46272445 | 1.9E+07 | 2.405 | Max at 330°C |
| 3 | 5032240 | 21032474 | 0.239 | Max at 200°C; shoulder at 330°C | 8817061 | 2.7E+07 | 0.321 | Max at 205°C; shoulder at 330°C |
| 4 | 2950165 | 1363530 | 2.164 | Max at 210°C | 2982827 | 1513771 | 1.97 | Max at 200°C |
| 5 | 45822 | 1293515 | 0.035 | first peak at 220°C;second peak at 350°C very small signals | 17697 | 1193987 | 0.015 | rise from 160°C; shoulder at 200°C;Max at 360°C |
| 6 | 56249 | 2118013 | 0.027 | Max at 190°C new rise from 280°C to 360°C; very small signals | 5086 | 1148071 | 0.004 | Max at 360°C |
| 7 | 47884 | 15195815 | 0.003 | rise from 200°CMax at 320°C;very small signal | 59028 | 1.8E+07 | 0.003 | rise from 160°CMax at 340°C;very small signal |
| 8 | 35879 | 21265834 | 0.002 | rise from 200°C,Max between 330°C and 340°C;very small signal | 39561 | 1.4E+07 | 0.003 | rise from 200°C,Max between 280°C and 360°C;very small signal |
| 9 | 8548668 | 5831549 | 1.466 | Max at 190°C | 8590994 | 5744781 | 1.495 | Max at 200°C |
| 10 | 1996 | 853513 | 0.002 | very small signal no maximum | 3015 | 839956 | 0.004 | Max at 340°C;small signal |
| 11 | 1310899 | 1127275 | 1.163 | Max at 195°C | 5750470 | 2849275 | 2.018 | Max at 200°C |
| 12 | 438296 | 1099691 | 0.399 | Max at 190°C | 243564 | 2681214 | 0.091 | Max at 180°C |
| 13 | 60718671 | 39137580 | 1.551 | Max at 195°C, shoulder between 230°C and 340°C | 85040694 | 5.6E+07 | 1.506 | Max at 200°C, shoulder at 240°C |
| 14 | 683713 | 30674122 | 0.022 | first rise to 190°C, shoulder between 240°C and 280°C, rise from 320°C | 541249 | 2.8E+07 | 0.019 | first rise to 190°C-200°C,shoulder between 250°C and 280°C,rise from 320°C |
| 15 | 205570 | 26733480 | 0.008 | Max at 280°C | 368335 | 5.3E+07 | 0.007 | Max between 270°C and 280°C |
| 16 | 135128 | 22797343 | 0.006 | Max at 320°C | 352618 | 2.9E+07 | 0.012 | first peak at 195°C;new rise with second max at 320°C |
| 17 | 70853 | 18044642 | 0.004 | Max at 320°C | 99460 | 2.6E+07 | 0.004 | Max at 320°C |
| 18 | 43462457 | 26101173 | 1.665 | Max at 200°C,shoulder at 320°C | 49746066 | 3.1E+07 | 1.596 | Max at 200°C, shoulder at 320°C |

Table C.93 TL data for TL laboratory 13

| | | | | | | | | | |
|------------------------|--------------|--------------|-------------------|-----------------|--------------|--------------|-------------------|-----------------|----------------|
| Lab No. | 14 | | | | | | | | |
| TL reader type | Harshaw 3500 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 0.354 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 220 | | | | | | | | |
| Chip 2 | 220 | | | | | | | | |
| Mean | 220 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 2.62 | 1993 | 0.0013146 | absent | 4.705 | 844.9 | 0.00556871 | dominant | mixture |
| 2 | 4272 | 3085 | 1.38476499 | dominant | 3708 | 3247 | 1.14197721 | dominant | positive |
| 3 | 283 | 2115 | 0.13380615 | dominant | 172.1 | 1029 | 0.16724976 | dominant | positive |
| 4 | 8.56 | 9.527 | 0.898499 | dominant | 3.01 | 3.847 | 0.78242787 | dominant | positive |
| 5 | 0.14 | 6.195 | 0.02259887 | absent | 0.361 | 20.22 | 0.01785361 | dominant | mixture |
| 6 | 0.229 | 13.11 | 0.01746758 | minor | 12.46 | 33.2 | 0.3753012 | dominant | mixture |
| 7 | 0.902 | 1010 | 0.00089307 | dominant | 0.713 | 967.3 | 0.0007371 | absent | mixture |
| 8 | 0.692 | 1605 | 0.00043115 | absent | 0.614 | 2934 | 0.00020927 | absent | negative |
| 9 | 4572 | 4827 | 0.94717216 | dominant | 3898 | 4184 | 0.93164436 | dominant | positive |
| 10 | 0.281 | 24.55 | 0.01144603 | absent | 0.218 | 40.92 | 0.00532747 | absent | negative |
| 11 | 33.79 | 24.93 | 1.35539511 | dominant | 51.76 | 100.5 | 0.51502488 | dominant | positive |
| 12 | 13.06 | 67.1 | 0.19463487 | dominant | 7.819 | 17.97 | 0.43511408 | dominant | positive |
| 13 | 26830 | 18960 | 1.41508439 | dominant | 25550 | 20410 | 1.25183733 | dominant | positive |
| 14 | 221.8 | 18250 | 0.01215342 | dominant | 287.7 | 15930 | 0.01806026 | dominant | mixture |
| 15 | 4.582 | 5256 | 0.00087177 | absent | 6.123 | 12160 | 0.00050354 | absent | negative |
| 16 | 2.53 | 2453 | 0.00103139 | absent | 4.292 | 1189 | 0.00360976 | absent | negative |
| 17 | 1.616 | 1931 | 0.00083687 | absent | 0.847 | 1005 | 0.00084279 | absent | negative |
| 18 | 3781 | 3196 | 1.1830413 | dominant | 404.8 | 347.4 | 1.1652274 | dominant | positive |

Table C.94 TL data for TL laboratory 14

| | | | | | | | | | |
|------------------------|--------------|--------------|---------|--------------------|--------------|--------------|---------|--------------------|------------|
| Lab No. | 15 | | | | | | | | |
| TL reader type | TL-DA-10 | | | | | | | | |
| Units for intensity | counts | | | | | | | | |
| MDL | 636 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 298 | | | | | | | | |
| Chip 2 | 288 | | | | | | | | |
| Mean | 293 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 1710 | 5526361 | 3.1E-04 | max. at 356 deg. C | 922 | 3058422 | 3.0E-04 | max. at 484 deg. C | untreated |
| 2 | 369843751 | 20192801 | 1.8E+01 | max. at 224 deg. C | 5587022 | 3115646 | 1.8E+00 | max. at 202 deg. C | irradiated |
| 3 | 6906284 | 13843432 | 5.0E-01 | max. at 202 deg. C | 109639 | 5901970 | 1.9E-02 | max. at 362 deg. C | mixture |
| 4 | 284351 | 167551 | 1.7E+00 | max. at 166 deg. C | - | - | - | - | - |
| 5 | 22951 | 8134351 | 2.8E-03 | max. at 434 deg. C | - | - | - | - | - |
| 6 | 10851 | 1077751 | 1.0E-02 | max. at 484 deg. C | - | - | - | - | - |
| 7 | 2255 | 14684444 | 1.5E-04 | max. at 492 deg. C | 1086 | 11815093 | 9.2E-05 | max. at 372 deg. C | untreated |
| 8 | 735 | 7183306 | 1.0E-04 | max. at 392 deg. C | 955 | 5240851 | 1.8E-04 | max. at 402 deg. C | untreated |
| 9 | 4105132 | 6968942 | 5.9E-01 | max. at 208 deg. C | 13300929 | 20234384 | 6.6E-01 | max. at 212 deg. C | irradiated |
| 10 | 10351 | 639151 | 1.6E-02 | max. at 498 deg. C | - | - | - | - | - |
| 11 | 338651 | 339551 | 1.0E+00 | max. at 198 deg. C | - | - | - | - | - |
| 12 | 41651 | 522151 | 8.0E-02 | max. at 458 deg. C | - | - | - | - | - |
| 13 | 7455246 | 8837855 | 8.4E-01 | max. at 214 deg. C | 10957351 | 15935211 | 6.9E-01 | max. at 200 deg. C | irradiated |
| 14 | 770276 | 20983301 | 3.7E-02 | max. at 188 deg. C | 58375 | 19059194 | 3.1E-03 | max. at 290 deg. C | mixture |
| 15 | 136351 | 110821451 | 1.2E-03 | max. at 298 deg. C | 44091 | 30605051 | 1.4E-03 | max. at 306 deg. C | untreated |
| 16 | 9512 | 8152399 | 1.2E-03 | max. at 330 deg. C | 2283 | 9540472 | 2.4E-04 | max. at 368 deg. C | untreated |
| 17 | 6627 | 10318861 | 6.4E-04 | max. at 360 deg. C | 4347 | 13892094 | 3.1E-04 | max. at 360 deg. C | untreated |
| 18 | 12641168 | 13396818 | 9.4E-01 | max. at 212 deg. C | 7273061 | 10125609 | 7.2E-01 | max. at 208 deg. C | irradiated |

Table C.95 TL data for TL laboratory 15

| | | | | | | | | | |
|------------------------|-----------------|--------------|-------|-------------------------|--------------|--------------|-------|-----------------------------|---------------------|
| Lab No. | No18 | | | | | | | | |
| TL reader type | Harshow TL 3500 | | | | | | | | |
| Units for intensity | n C | | | | | | | | |
| MDL | 2.56 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 229 | | | | | | | | |
| Chip 2 | 235 | | | | | | | | |
| Mean | 232 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 | 183.86 | 5814.26 | 0.032 | Peak at 329°C | 232.40 | 5644.12 | 0.041 | Peak at 335°C | Un-irradiated |
| 2 | 21978.15 | 10832.01 | 2.029 | Peak at 198°C | 16345.85 | 8614.27 | 1.898 | Peak at 198°C | Irradiated |
| 3 | 4362.90 | 16132.80 | 0.270 | Peak at 186°C | 3524.93 | 16228.32 | 0.217 | Peak at 187°C | Irradiated(mixture) |
| 4 | 2555.66 | 1293.03 | 1.976 | Peak at 179°C | 2095.25 | 1010.35 | 2.074 | Peak at 183°C | Irradiated |
| 5 | 8.07 | 622.83 | 0.013 | Peak at 186°C | 12.40 | 736.26 | 0.017 | Peak at 189°C | Irradiated(mixture) |
| 6 | 7.17 | 930.38 | 0.008 | no peak was observed | 8.19 | 1003.16 | 0.008 | no peak was observed | Un-irradiated |
| 7 | 106.75 | 13944.44 | 0.008 | small shoulder at 189°C | 187.41 | 25724.32 | 0.007 | small shoulder at 183°C | Irradiated(mixture) |
| 8 | 88.17 | 14580.86 | 0.006 | Peak at 322°C | 100.71 | 18419.95 | 0.005 | Peak at 317°C | Un-irradiated |
| 9 | 23297.25 | 17080.90 | 1.364 | Peak at 183°C | 24074.51 | 17478.48 | 1.377 | Peak at 183°C | Irradiated |
| 10 | 13.56 | 1075.97 | 0.013 | no peak was observed | 16.38 | 1946.71 | 0.008 | no peak was observed | Un-irradiated |
| 11 | 4857.31 | 1936.39 | 2.508 | Peak at 181°C | 5729.54 | 2400.65 | 2.387 | Peak at 183°C | Irradiated |
| 12 | 172.90 | 675.73 | 0.256 | Peak at 188°C | 62.10 | 626.72 | 0.099 | Peak at 183°C | Irradiated(mixture) |
| 13 | 6691.95 | 2362.35 | 2.833 | Peak at 185°C | 4191.67 | 1746.66 | 2.400 | Peak at 186°C | Irradiated |
| 14 | 432.21 | 10390.84 | 0.042 | Peak at 190°C | 645.84 | 14423.25 | 0.045 | Peak around 200-250C(broad) | Irradiated(mixture) |
| 15 | 330.76 | 14384.48 | 0.023 | Peak at 266°C | 175.75 | 8417.15 | 0.021 | Peak at 263°C | Un-irradiated |
| 16 | 206.18 | 9187.53 | 0.022 | small shoulder at 189°C | 100.41 | 5519.04 | 0.018 | small shoulder at 175°C | Irradiated(mixture) |
| 17 | 92.65 | 5713.64 | 0.016 | Peak at 348°C | 137.20 | 9145.05 | 0.015 | Peak at 348°C | Un-irradiated |
| 18 | 21255.62 | 12226.83 | 1.738 | Peak at 181°C | 17872.33 | 10408.60 | 1.717 | Peak at 183°C | Irradiated |

Table C.96 TL data for TL laboratory 18

| | | | | | | | | | |
|------------------------|------------------|--------------|--------|---------------|--------------|--------------|-------|---------------|------------|
| Lab No. | 24 | | | | | | | | |
| TL reader type | Harshaw TLD 3500 | | | | | | | | |
| Units for intensity | nC | | | | | | | | |
| MDL | 1.9 | | | | | | | | |
| LiF data | | | | | | | | | |
| Peak V temperature /°C | | | | | | | | | |
| Chip 1 | 228 | | | | | | | | |
| Chip 2 | 228 | | | | | | | | |
| Mean | 228 | | | | | | | | |
| Sample data | | | | | | | | | |
| | Aliquot A | | | | Aliquot B | | | | |
| Sample | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | G1 Intensity | G2 Intensity | G1/G2 | G1 peak shape | Evaluation |
| 1 (29) | 7.497 | 1145 | 0.007 | Pmax 202°C | 10.44 | 3946 | 0.003 | no Peak | negative |
| 2(45) | 10340 | 6752 | 1.531 | Pmax 211°C | 5836 | 4793 | 1.218 | Pmax 228°C | positive |
| 3(72) | 619 | 3536 | 0.175 | Pmax 225°C | 3536 | 4339 | 0.815 | Pmax 222°C | positive |
| 4(15) | 56.43 | 6787 | 0.008 | Pmax 215°C | 67.87 | 109.5 | 0.620 | Pmax 222°C | positive |
| 5(17) | 4.191 | 118.1 | 0.035 | Pmax 225°C | 3.142 | 141.9 | 0.022 | Pmax | negative |
| 6(14) | 0.968 | 283.8 | 0.003 | no Peak | 1.06 | 117.3 | 0.009 | no Peak | negative |
| 7(30) | 22.61 | 3122 | 0.007 | no Peak | 32.27 | 7340 | 0.004 | | negative |
| 8(53) | 7.733 | 1660 | 0.005 | Pmax 208°C | 4.698 | 1593 | 0.003 | no Peak | negative |
| 9(66) | 12880 | 7857 | 1.639 | Pmax 205°C | 1334 | 5060 | 0.264 | Pmax 284°C | positive |
| 10(20) | 5.56 | 245.8 | 0.023 | Pmax 234°C | 4.06 | 511.8 | 0.008 | Pmax 185°C | negative |
| 11(48) | 428.2 | 335.9 | 1.275 | Pmax 205°C | 472.2 | 389.1 | 1.214 | Pmax 208°C | positive |
| 12(36) | 6.956 | 140.6 | 0.049 | | 25.55 | 584.9 | 0.044 | | positive |
| 13(22) | 16220 | 996.4 | 16.279 | Pmax 208°C | 2173 | 1494 | 1.454 | Pmax 208°C | positive |
| 14(2) | 123.6 | 3399 | 0.036 | Pmax 222°C | 298.6 | 6733 | 0.044 | Pmax 222°C | negative |
| 15(55) | 62.41 | 6309 | 0.010 | no Peak | 36 | 7573 | 0.005 | no Peak | negative |
| 16(46) | 51.09 | 5981 | 0.009 | no Peak | 44.75 | 4962 | 0.009 | no Peak | negative |
| 17(28) | 38.22 | 6894 | 0.006 | no Peak | 35.44 | 2878 | 0.012 | no Peak | negative |
| 18(37) | 3286 | 2747 | 1.196 | Pmax 228°C | 6710 | 5400 | 1.243 | Pmax 218°C | positive |

Table C.97 TL data for TL laboratory 24

APPENDIX D

SAMPLE ALLOCATION

| Pot order | Lab no |
|-----------|--------|
| 1 | 27 |
| 2 | 9 |
| 3 | 7 |
| 4 | 21 |
| 5 | 19 |
| 6 | 20 |
| 7 | 14 |
| 8 | 2 |
| 9 | 15 |
| 10 | 3 |
| 11 | 11 |
| 12 | 5 |
| 13 | 35 |
| 14 | 10 |
| 15 | 8 |
| 16 | 4 |
| 17 | 28 |
| 18 | 34 |
| 19 | 22 |
| 20 | 6 |
| 21 | 16 |
| 22 | 25 |
| 23 | 24 |
| 24 | 26 |
| 25 | 17 |
| 26 | 23 |
| 27 | 29 |
| 28 | 12 |
| 29 | 33 |
| 30 | 32 |
| 31 | 1 |
| 32 | 30 |
| 33 | 18 |
| 34 | 31 |
| 35 | 13 |

Table D1. Allocation of pots to laboratories for round 2 – 35 sets of materials were prepared although only 30 laboratories returned data

| SP Number | Serial Number | Sample no |
|-----------|---------------|-----------|
| 8512 u | 1 | 67 |
| 8512 i | 2 | 47 |
| 8513 u | 3 | 31 |
| 8513 i | 4 | 26 |
| 8514 u | 5 | 56 |
| 8514 i | 6 | 64 |
| 8515 u | 7 | 42 |
| 8515 i | 8 | 35 |
| 8516 u | 9 | 25 |
| 8516 i | 10 | 18 |
| 8517 u | 11 | 58 |
| 8517 i | 12 | 5 |
| 8518 u | 13 | 54 |
| 8518 i | 14 | 33 |
| 8519 u | 15 | 39 |
| 8519 i | 16 | 24 |
| 8520 u | 17 | 68 |
| 8520 i | 18 | 69 |
| 8521 u | 19 | 3 |
| 8521 i | 20 | 38 |
| 8523 u | 21 | 12 |
| 8523 i | 22 | 32 |
| 8525 u | 23 | 1 |
| 8525 i | 24 | 8 |
| 8526 u | 25 | 59 |
| 8526 i | 26 | 50 |
| 8528 u | 27 | 10 |
| 8528 i | 28 | 62 |
| 8529 u | 29 | 11 |
| 8529 i | 30 | 51 |
| 8530 u | 31 | 23 |
| 8530 i | 32 | 21 |
| 8531 u | 33 | 29 |
| 8531 i | 34 | 45 |
| 8531 b | 35 | 72 |
| 8532 u | 36 | 27 |
| 8532 i | 37 | 49 |
| 8533 u | 38 | 14 |
| 8533 i | 39 | 15 |
| 8533 b | 40 | 17 |
| 8534 u | 41 | 4 |
| 8534 i | 42 | 65 |
| 8535 u | 43 | 53 |
| 8535 i | 44 | 66 |
| 8535 b | 45 | 30 |
| 8536 u | 46 | 19 |
| 8536 i | 47 | 44 |
| 8538 u | 48 | 61 |
| 8538 i | 49 | 16 |
| 8539 u | 50 | 20 |
| 8539 i | 51 | 48 |

| | | |
|--------|----|----|
| 8539 b | 52 | 36 |
| 8571 u | 53 | 9 |
| 8571 i | 54 | 34 |
| 8572 u | 55 | 6 |
| 8572 i | 56 | 41 |
| 8573 u | 57 | 60 |
| 8573 i | 58 | 7 |
| 8574 u | 59 | 70 |
| 8574 i | 60 | 40 |
| 8575 u | 61 | 71 |
| 8575 i | 62 | 63 |
| 8576 u | 63 | 13 |
| 8576 i | 64 | 52 |
| 8577 u | 65 | 55 |
| 8577 i | 66 | 22 |
| 8577 b | 67 | 2 |
| 8578 u | 68 | 28 |
| 8578 i | 69 | 37 |
| 8578 b | 70 | 46 |
| 8579 u | 71 | 57 |
| 8579 i | 72 | 43 |

Table D2. Correspondence between SP number and sample number as presented to participants

APPENDIX E

MIXING THE BLENDS

Mixing the Green Tea

On 25 April 2006 499.5g of unirradiated green tea powder was taken from pots 181-200 and supplementary pots 95-97 and part of 98. 0.5014g was taken evenly from irradiated pots 101-110. The material was placed in a 3l plastic mixing bowl and mixed for 60s with Braun Quattro Multimix cakebeaters at the lowest speed to minimise dust generation. After 60s, approximately 3g of material was placed in each of 10 petri dishes, taking the material from different parts of the bowl. The aliquots were screened following the normal PSL procedure as run MIXGT. After measurement, mean SD and CV were tabulated and the cycle repeated for a total of 6 cycles, tabulated below.

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|-----|------|-----|-----|-----|-----|
| Aliquot | | | | | | |
| 1 | 487 | 395 | 386 | 607 | 541 | 555 |
| 2 | 480 | 506 | 380 | 508 | 761 | 389 |
| 3 | 391 | 529 | 677 | 427 | 391 | 549 |
| 4 | 350 | 515 | 376 | 586 | 605 | 553 |
| 5 | 502 | 608 | 466 | 578 | 549 | 473 |
| 6 | 634 | 520 | 406 | 322 | 495 | 886 |
| 7 | 580 | 509 | 478 | 338 | 623 | 335 |
| 8 | 509 | 1473 | 483 | 605 | 662 | 638 |
| 9 | 375 | 485 | 454 | 441 | 577 | 555 |
| 10 | 492 | 448 | 490 | 506 | 708 | 505 |
| Mean | 480 | 599 | 460 | 492 | 591 | 544 |
| SD | 89 | 312 | 89 | 107 | 107 | 149 |
| CV(%) | 19 | 52 | 19 | 22 | 18 | 27 |

Table E.1 PSL measurements during mixing for SP8578 Green Tea

It can be seen that the signals recorded are smaller than predicted, but that the criterion for CV <50% has been met.

On 28 April 2006 the material was potted for the participants, using 8g per pot.

Mixing the Chilli Powder

On 26 April 2006 the same process was repeated for the chilli powder using 499.5g unirradiated and 0.5329g irradiated. Pots 181-195 were used for the unirradiated material and 101-110 for the spike.

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|------|------|------|------|------|------|
| Aliquot | | | | | | |
| 1 | 587 | 839 | 810 | 725 | 908 | 1345 |
| 2 | 831 | 690 | 595 | 1135 | 1401 | 585 |
| 3 | 1165 | 1036 | 747 | 664 | 932 | 1480 |
| 4 | 810 | 710 | 2563 | 779 | 530 | 772 |
| 5 | 1039 | 1093 | 764 | 666 | 957 | 552 |
| 6 | 615 | 837 | 521 | 1496 | 852 | 797 |

| | | | | | | |
|-------|------|------|------|------|------|------|
| 7 | 735 | 818 | 1464 | 1069 | 1249 | 581 |
| 8 | 4209 | 1121 | 1850 | 1251 | 810 | 1390 |
| 9 | 1189 | 496 | 671 | 599 | 2245 | 2310 |
| 10 | 1296 | 843 | 758 | 595 | 897 | 695 |
| Mean | 1248 | 848 | 1074 | 898 | 1078 | 1051 |
| SD | 1014 | 184 | 636 | 300 | 449 | 541 |
| CV(%) | 81 | 22 | 59 | 33 | 42 | 51 |

Table E.2 PSL measurements during mixing for SP8535 Chilli Powder

Mixing of Ginger

On 26 April 2006 the same process was repeated for the ginger powder using 385g unirradiated and 3.9003g irradiated. Pots 181-195 and 93-100 were used for the unirradiated material and 101-110 for the spike. It was realised that this is the incorrect concentration for this product. Therefore on 27 April 2006 a further 38.9g of irradiated ginger were added, making a total of 42.8g of spike in 427.8g of blended material, i.e. 10%.

Results of 1% mixing were as follows:

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|------|------|------|------|-------|------|
| Aliquot | | | | | | |
| 1 | 4467 | 3011 | 4422 | 1980 | 2872 | 3217 |
| 2 | 1838 | 2299 | 1237 | 1959 | 1093 | 1459 |
| 3 | 5257 | 1506 | 1342 | 1441 | 1601 | 1179 |
| 4 | 1077 | 987 | 1309 | 1760 | 894 | 1606 |
| 5 | 1577 | 1704 | 3162 | 1417 | 1632 | 1090 |
| 6 | 6277 | 6662 | 4130 | 4664 | 2237 | 742 |
| 7 | 2007 | 1550 | 7523 | 1887 | 4414 | 2558 |
| 8 | 1644 | 1625 | 1796 | 1607 | 2761 | 1311 |
| 9 | 1507 | 4566 | 3929 | 2811 | 2785 | 4570 |
| 10 | 2533 | 1246 | 2020 | 3755 | 17750 | 1405 |
| Mean | 2818 | 2516 | 3087 | 2328 | 3804 | 1914 |
| SD | 1826 | 1796 | 1992 | 1087 | 5008 | 1187 |
| CV(%) | 65 | 71 | 65 | 47 | 132 | 62 |

Table E.3 PSL measurements during mixing for SP8531 Ginger 1%

Results of 10% mixing were as follows (from cycle 5 onwards a spoon was used to mix the contents of the petri dishes back into the mix):

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|--------|-------|-------|-------|-------|--------|
| Aliquot | | | | | | |
| 1 | 132317 | 18631 | 29131 | 25373 | 20460 | 29978 |
| 2 | 46022 | 23256 | 32118 | 36264 | 34870 | 28036 |
| 3 | 29889 | 40510 | 30521 | 24243 | 27006 | 27769 |
| 4 | 21053 | 35030 | 49456 | 15367 | 55175 | 20575 |
| 5 | 24754 | 26736 | 20096 | 28848 | 20363 | 23216 |
| 6 | 49485 | 29598 | 41336 | 16371 | 30578 | 45810 |
| 7 | 18082 | 21613 | 36341 | 41747 | 35186 | 111396 |
| 8 | 33840 | 23609 | 17719 | 39367 | 37860 | 24200 |

| | | | | | | |
|-------|-------|-------|-------|-------|-------|-------|
| 9 | 23358 | 15398 | 45279 | 91130 | 34371 | 24200 |
| 10 | 17763 | 29798 | 19610 | 25919 | 63566 | 31277 |
| Mean | 39656 | 26418 | 32161 | 34463 | 35944 | 36646 |
| SD | 34369 | 7582 | 11037 | 21795 | 13872 | 27175 |
| CV(%) | 87 | 29 | 34 | 63 | 39 | 74 |

Table E.4 PSL measurements during mixing for SP8531 Ginger 10%

After a further 60s mixing, approximately 7g of blended ginger (calculated from the total required for participants and for homogeneity testing) was dispensed into pots.

Mixing the Siberian Ginseng

On 27 April 2006 the above process was used to mix and check a 1% blend of Siberian Ginseng. Pots 181-200 plus 93-95 were used to provide 495.1g of unirradiated material which was spiked with 5.0734g of irradiated material from pots 101-110. During the measurement cycle the decision was taken to perform 6 cycles of mixing and screening regardless of the CV obtained for statistical comparability.

PSL results were as follows:

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|-------|-------|-------|-------|-------|-------|
| Aliquot | | | | | | |
| 1 | 13987 | 11029 | 12201 | 10701 | 9176 | 23521 |
| 2 | 8723 | 9404 | 10947 | 8897 | 8115 | 17382 |
| 3 | 9143 | 10991 | 17336 | 15057 | 19073 | 13751 |
| 4 | 8166 | 10827 | 11663 | 11359 | 10624 | 32622 |
| 5 | 16160 | 11439 | 12326 | 9979 | 8527 | 12156 |
| 6 | 10562 | 13188 | 7280 | 13235 | 17283 | 10343 |
| 7 | 8793 | 19298 | 39537 | 12460 | 10901 | 8482 |
| 8 | 10170 | 8515 | 11398 | 9932 | 15523 | 14210 |
| 9 | 15319 | 14326 | 16974 | 7769 | 8834 | 9218 |
| 10 | 10312 | 12158 | 14070 | 12908 | 14864 | 23494 |
| Mean | 11134 | 12118 | 15373 | 11230 | 12292 | 16518 |
| SD | 2923 | 3028 | 8982 | 2210 | 4027 | 7787 |
| CV(%) | 26 | 25 | 58 | 20 | 33 | 47 |

Table E.5 PSL measurements during mixing for SP8577 Ginseng

On 28 April 2006 the blended material was dispensed into 50 pots for TL and PSL, with 8g per pot. Allocation of pots for PSL followed the same sequence as the ginger.

Mixing the rosemary

The same process was followed on 28 April with the rosemary, to make a 1% blend. 409g of unirradiated material was taken from pots 93-100 and 181-200, and mixed with 4.1791g spike from pots 101-110 to give 413g blend. PSL results are tabulated below:

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|------|------|------|-------|------|------|
| Aliquot | | | | | | |
| 1 | 540 | 2299 | 1800 | 66252 | 846 | 657 |
| 2 | 1516 | 2533 | 1256 | 403 | 857 | 529 |
| 3 | 1072 | 1604 | 1047 | 817 | 882 | 1517 |
| 4 | 892 | 1901 | 1191 | 983 | 967 | 1329 |
| 5 | 1499 | 769 | 766 | 2751 | 1319 | 1968 |
| 6 | 557 | 636 | 340 | 500 | 785 | 1375 |
| 7 | 887 | 964 | 1523 | 825 | 1447 | 1188 |
| 8 | 461 | 885 | 1251 | 1139 | 1889 | 1135 |
| 9 | 1214 | 581 | 792 | 2198 | 922 | 844 |
| 10 | 1979 | 637 | 710 | 762 | 959 | 396 |
| Mean | 1062 | 1281 | 1068 | 7663 | 1087 | 1094 |
| SD | 495 | 740 | 429 | 20600 | 354 | 488 |
| CV(%) | 47 | 58 | 40 | 269 | 33 | 45 |

Table E.6 PSL measurements during mixing for SP8533 Rosemary

6g of blended material was put into each of 50 tubs as before.

Mixing the tarragon

The same process was followed on 28 April 2006 with tarragon, using 324g from pots 93-100 and 181-200 and 36.2g of spike from pots 101-110. to give a blended weight of 360.2g.

PSL results are tabulated below:

| Cycle | 1 | 2 | 3 | 4 | 5 | 6 |
|---------|------|------|------|------|-------|------|
| Aliquot | | | | | | |
| 1 | 1585 | 6275 | 3232 | 2294 | 1565 | 5185 |
| 2 | 725 | 1616 | 2258 | 1958 | 943 | 2697 |
| 3 | 6505 | 4374 | 8500 | 2596 | 2538 | 1943 |
| 4 | 2111 | 2410 | 1644 | 4576 | 1636 | 2966 |
| 5 | 2643 | 6820 | 3393 | 5427 | 1703 | 4320 |
| 6 | 6120 | 3403 | 2008 | 2202 | 2077 | 1495 |
| 7 | 1311 | 6345 | 5141 | 1434 | 10127 | 1852 |
| 8 | 1926 | 3583 | 1911 | 3411 | 1786 | 2164 |
| 9 | 4635 | 3002 | 6089 | 2059 | 5335 | 1345 |
| 10 | 1459 | 4553 | 2364 | 2268 | 2329 | 1444 |
| Mean | 2902 | 4238 | 3654 | 2823 | 3004 | 2541 |
| SD | 2084 | 1771 | 2243 | 1269 | 2770 | 1294 |
| CV(%) | 72 | 42 | 61 | 45 | 92 | 51 |

Table E.7 PSL measurements during mixing for SP8539 Tarragon

On 2 May 2006 the chilli powder and tarragon were potted (8g and 6g per pot respectively).